- (GB)** Instructions for use/Technical description
AdTec mini – instrument system for single use
- (USA)** AdTec mini – instrument system for single use
- (D)** Gebrauchsanweisung/Technische Beschreibung
AdTec mini – Instrumentensystem zum Einmalgebrauch
- (F)** Mode d'emploi/Description technique
AdTec mini – système d'instruments à usage unique
- (E)** Instrucciones de manejo/Descripción técnica
AdTec mini instrumental de un solo uso
- (I)** Istruzioni per l'uso/Descrizione tecnica
AdTec mini – sistema strumenti monouso
- (P)** Instruções de utilização/Descrição técnica
AdTec mini – sistema de instrumental para utilização única
- (NL)** Gebruiksaanwijzing/Technische beschrijving
AdTec mini – instrumentsysteem voor eenmalig gebruik
- (S)** Bruksanvisning/Teknisk beskrivning
AdTec mini – instrumentsystem för engångsbruk
- (RUS)** Инструкция по применению/Техническое описание
АдТек мини инструменты для одноразового использования
- (CZ)** Návod k použití/Technický popis
Systém nástrojů k jednorázovému použití AdTec mini
- (PL)** Instrukcja użytkowania/Opis techniczny
AdTec mini – system instrumentów do jednorazowego użytku
- (SK)** Návod na používanie/Technický opis
AdTec mini – systém nástrojov na jedno použitie
- (TR)** Kullanım Kılavuzu/Teknik açıklama
AdTec mini – tek seferlik kullanım için ekipman sistemi


Fig. 1


Fig. 2

Fig. 3

B | BRAUN
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 014296 04/15 V6

CE 0123 – DIR 93/42/EEC

Technical alterations reserved

Aesculap®

AdTec mini – instrument system for single use

Legend

- 1 Actuating lever
- 2 Switch
- 3 Handle part

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture

Intended use

All professional disciplines that use endoscopy: Cutting, dissection, mobilization and coagulation of tissue.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product.

The product is gamma-sterilized and supplied in sterile packaging.

The product must not be reused.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Do not use products from open or damaged sterile packaging.
- Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Do not use the product after its use-by date.

Risk of injury from ignition or explosion of flammable gases! Sparks may occur when using the HF device as directed.

- Observe the safety guidelines in the instructions for use of the HF device.

Thermal injuries to patients/users due to insufficient insulation of leads in active accessories!

- Adjust the HF device to an appropriate setting to ensure that the peak output voltage does match or not exceed the accessory voltage rating specified for the product.

- Adjust the HF power output to the intended surgical intervention. Take into account clinical experience or reference values.
- Select the lowest possible HF power output.
- Keep the product's contact surfaces clean during surgery. Remove encrusted tissue residues or body fluids with a moistened swab.

The product is furnished with the following connection on the plug side: 4 mm pin.

Refer to our brochures to find a compatible cable.

The accessory voltage rating of the product is 1 000 Vp.

The accessory voltage rating must exceed or match the peak output voltage with which the product is operated in combination with a suitable HF device at an appropriate operating mode/setting (see IEC/DIN EN 60601-2-2).

To avoid HF burns:

- Push the HF cable plug onto the HF pin to the positive stop.
- Check for secure seating of the plug before activating the HF current.
- Always keep the working end of the product in the user's field of vision whenever the HF power is activated.
- Prior to activating the HF device, check that the working end of the product is not touching any electrically conductive accessories.
- Prior to each use, visually inspect the product for: damage or surface changes to the insulation.
- Never place the product on or next to the patient.
- When using accessories for endoscopy or laparoscopy, deactivate the automatic switch-on mode of the HF device.
- Follow the instructions for use of the HF device.

Safe operation

Handle part with lock:

- To release the ratchet mechanism, press actuating lever 1 and handle part 3 towards each other, see Fig. 1.
- To activate the ratchet mechanism, release actuating lever 1.

Switch 2 can be used for permanently deactivating the ratchet mechanism.

- To deactivate the ratchet mechanism on the handle part 3: Position switch 2, see Fig. 2.
- To activate the ratchet mechanism on the handle part 3: Position switch 2, see Fig. 3.

Storage

- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Technical Service

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042

Aesculap Repair Hotline

Phone: +1 (800) 214-3392

Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

For accessories and spare parts, see brochure C304, C280, C285 and C951 or visit <http://www.endoscopy-catalog.com>

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

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