Aesculap[®] Excia[®]

10 Years Excellence in Arthroplasty

Aesculap Orthopaedics

Excia®

Excellence in Arthroplasty

The Excia Hip System is based on many years of experience with straight stem implants in France, where the first Excia was implanted in April 2000.

Today, the Excia system is used successfully throughout Europe, the US and Japan.

The Excia system concept of using one instrument set for implantation with or without bone cement is time tested and has been subject to ongoing enhancements.

Excia is one of the most important products in the Aesculap hip endoprosthesis range.

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Excia 10 Years

4th Generation Hip Replacement

Successful hip replacement surgery has benefited from various developments over the past decades.

At the time Excia was developed, less invasive and minimally invasive surgical techniques, such as straight stem systems for universal use and computer-aided applications and navigation, were coming to the fore.

These 4th generation stems were accompanied by significant improvements in articulation materials and further developments in acetabular systems for use with larger prosthesis heads. From 2010 to 2020, we can expect to see sustained clinical results from enhanced systems, an optimized endoprosthesis hip replacement process and an individualized, patient-focused approach to hip replacement.

Excia has made a major contribution to these developments over the course of its 10 years.

2nd

1980 ...

- Porous coatings
- Proximal stem fixation
- Reduced thigh pain

Hip stems
cementless

3rd

1990 ...

- Optimized cementing technique
- Use of distal centralizers
- CoCr cemented stems

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4th

2000 ...

- Minimally invasive surgery
- Universal hip systems
- CAS & Navigation

Hip systems
and surgery

2010 ...

- Long term results
- Patient matched solutions
- Health economics impact

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Excia[®]

Design Concept

The Excia straight stem concept is based on the characteristic design features of the proximal and distal stem components.

This allows it to meet the Excia system requirements for a cementless or cemented implant in a unique way.

The heart of the Excia system is its modern straight stem instrument set featuring high precision rasps and instruments for all minimally invasive procedures.

For cementless Excia implants, proximal pressfit is facilitated by a sub-proximal distal guide function of the tapered stem component. The proximal-lateral wing rasp increases primary rotational stability.

The cemented Excia stem is smaller than the stem rasp.

A special feature is the Excia 8/10 conical trunnion, which can achieve a significantly higher range of motion, even with small head diameters, collarless XL head lengths and bipolar hip replacements.

Proximal

Small Taper
High Range of Motion


Trochanter Wing
Rotational Stability

Plasmapore μ -CaP
Secondary Stability

Bilateral Flanges
Primary Stability

Distal

High Precision Rasps
Pressfit and Formfit

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Excia®

Hip Technology

The modern Excia implant system embodies the use of new technologies in hip endoprosthetics. As early as the year 2000, the Excia system supported the use of less invasive implantation techniques and computer-assisted hip navigation.

The Excia stem rasp, the modular wing rasp, curved handle and MIOS implantation instruments have been continually enhanced for use in less invasive procedures.

Based on the kinematic OrthoPilot acetabular navigation, the Excia hip stem could be used with a combined acetabular and stem implantation, which made possible the first image-free total hip navigation procedure.

The most recent version of OrthoPilot hip navigation has also been documented by Excia clinical trial findings on the use of ultrasound technology for hip navigation.

The cementless Excia straight stem and Plasmacup press-fit system were among the first hip endoprosthetic implantation components to be offered with materials such as combined metal and dicalcium phosphate surface coatings.

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
Excia®


Hip System 2010

The Excia Hip System has benefited from the successes of and further developments in hip endoprosthesis, and today features all the necessary components for a universal straight stem system for primary joint replacement and fracture repair.

The Excia stem is available in all sizes in a standard and lateralized version, in 135° and 128°, and as either cementless or cemented.

The Excia implant product line has been expanded over the course of clinical usage and also through implant components with a 12/14 neck taper to accommodate standardization within the Aesculap hip endoprosthesis offerings for revision and short stem endoprosthesis.



The system is completed by head components for the 8/10 or 12/14 conical trunnion, 22.2 to 36 mm, or bipolar, metal or ceramic as well as cementless or cemented Aesculap acetabular components.

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Excia®

Experiences

Excia hip stems have been successfully implanted in Europe and the US for the past ten years. The X-ray study over time shows the first Excia replacement from the year 2000. The results of follow-up investigations are documented in ongoing studies of both cementless and cemented Excia stems. Approximately 65% of all Excia hip endoprosthesis stems are cementless. At present, 20% use the lateralized stem design. Excia has also been implanted in Japan since the beginning of 2010 where it is frequently used in the treatment of dysplastically malformed and smaller femoral morphologies. The Excia stem design has thus proven itself as a universal implant system over its first ten years, and with its enhancements is ready for the future.

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- Tissue sparing implantation
- Line-to-line stability
- Proximal bone remodelling
- No thigh pain

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Plasmapore µ-CaP, Surface Coating.

Aesculap AG, Scientific Information 2001; 051002

Bobyn JD, Pilliar RM, Cameron HU, Weatherly GC

The optimum pore size for the fixation of porous surfaced metal implants by ingrowth of bone.

Clin Orthop Relat Res 1980; Jul-Aug, (150): p. 263-270

Proximal stem fixation

Aldinger PR, Jung AW, Breusch SJ, Ewerbeck V, Parsch D

Survival of the cementless Spotorno stem in the second decade.

Clin Orthop Relat Res 2009; Sep, 467(9): p. 2297-2304

McLaughlin JR, Lee KR

Total hip arthroplasty with an uncemented tapered femoral component.

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Prospective long-term follow-up of the cementless bicontact hip stem with plasmapore coating.

Z Orthop Unfall 2007; Sep-Oct, 145 Suppl 1: p. S3-8

Reduced thigh pain

Suckel A, Geiger F, Kinzl L, Wulker N, Garbrecht M

Long-term results for the uncemented Zweymüller/Alloclassic hip endoprosthesis.

A 15-year minimum follow-up of 320 hip operations.

J Arthroplasty 2009; Sep, 24(6): p. 846-853

Hallan G, Lie SA, Furnes O, Engesaeter LB, Vollset SE, Havelin LI

Medium- and long-term performance of 11,516 uncemented primary femoral stems from the Norwegian arthroplasty register.

J Bone Joint Surg Br 2007; Dec, 89(12): p. 1574-1580

Optimized cementing technique

El Masri F, Kerboul L, Kerboul M, Courpied JP, Hamadouche M

Is the so-called 'French paradox' a reality?: long-term survival and migration of the Charnley-Kerboul stem cemented line-to-line.

J Bone Joint Surg Br 2010; Mar, 92(3): p. 342-348

Sierra RJ, Timperley JA, Gie GA

Contemporary cementing technique and mortality during and after Exeter total hip arthroplasty.

J Arthroplasty 2009; Apr, 24(3): p. 325-332

Rice J, Prenderville T, Murray P, McCormack B, Quinlan W

Femoral cementing techniques in total hip replacement.

Int Orthop 1998; 22(5): p. 308-311

Use of distal centralizers

Breusch SJ, Lukoschek M, Kreutzer J, Brocai D, Gruen TA

Dependency of cement mantle thickness on femoral stem design and centralizer. J Arthroplasty 2001; Aug, 16(5): p. 648-657

Goldberg BA, al-Habbal G, Noble PC, Paravic M, Liebs TR, Tullos HS

Proximal and distal femoral centralizers in modern cemented hip arthroplasty.

Clin Orthop Relat Res 1998; Apr, (349): p. 163-173

Hanson PB, Walker RH

Total hip arthroplasty cemented femoral component distal stem centralizer.

Effect on stem centralization and cement mantle.

J Arthroplasty 1995; Oct, 10(5): p. 683-688

CoCr cemented stems

Kovac S, Trebse R, Milosev I, Pavlovic V, Pisot V
Long-term survival of a cemented titanium-aluminium-vanadium alloy straight-stem femoral component.
J Bone Joint Surg Br 2006; Dec, 88(12): p. 1567-1573

Jergesen HE, Karlen JW
Clinical outcome in total hip arthroplasty using a cemented titanium femoral prosthesis.
J Arthroplasty 2002; Aug, 17(5): p. 592-599

Tompkins GS, Lachiewicz PF, DeMasi R
A prospective study of a titanium femoral component for cemented total hip arthroplasty.
J Arthroplasty 1994; Dec, 9(6): p. 623-630

Minimally invasive surgery

Oinuma K, Eingartner C, Saito Y, Shiratsuchi H
Total hip arthroplasty by a minimally invasive, direct anterior approach.
Oper Orthop Traumatol 2007; Aug, 19(3): p. 310-326

Judet H
Five years of Experience in hip navigation using a mini-invasive anterior approach.
Orthopedics 2007; Vol. 30, No.10: p. S141-143

Ogonda L, Wilson R, Archbold P, Lawlor M, Humphreys P, O'Brien S
minimal-incision technique in total hip arthroplasty does not improve early postoperative outcomes. A prospective, randomized, controlled trial.
J Bone Joint Surg Am 2005; Apr, 87(4): p. 701-710

Kennon RE, Keggi JM, Wetmore RS, Zatorski LE, Huo MH, Keggi KJ
Total hip arthroplasty through a minimally invasive anterior surgical approach.
J Bone Joint Surg Am 2003; 85-A Suppl 4: p. 39-48

Universal hip systems

D'Angelo F, Murena L, Vulcano E, Zatti G, Cherubino P
Seven to twelve year results with versys et cementless stem. A retrospective study of 225 cases. Hip Int 2010; Jan-Mar, 20(1): p. 81-86

Mainard D
Results of a bone preserving cementless straight hip stem A prospective study with a minimum follow-up of five years.
EFORT Vienna 2009; E-Poster P434

Ajmal M, Ranawat AS, Ranawat CS
A new cemented femoral stem: a prospective study of the Stryker accolade C with 2- to 5-year follow-up.
J Arthroplasty 2008; Jan, 23(1): p. 118-122

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Mainard D
Navigated and Nonnavigated Total Hip Arthroplasty: Results of Two Consecutive Series Using a Cementless Straight Hip Stem.
Orthopedics 2008; Vol. 31, No. 10: 22-26

Hasart O, Perka C, Tohtz S
Comparison between pointer-based and ultrasound-based navigation technique in THA using a minimally invasive approach.
Orthopedics 2008; Oct, 31 (10 Suppl): p. 31-34

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B. Braun Austria GmbH | Aesculap Division | Otto Braun-Straße 3-5 | 2344 Maria Enzersdorf | Österreich
Tel. +43 2236 4 65 41-0 | Fax +43 2236 4 65 41-177 | www.bbraun.at

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B. Braun Medical AG | Seesatz | 6204 Sempach | Schweiz
Tel. +41 58 258 50 00 | Fax +41 58 258 60 00 | www.bbraun.ch

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Deutschland
Tel. (0 74 61) 95-0 | Fax (0 74 61) 95-26 00 | www.aesculap.de

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Die Hauptproduktmarke Aesculap und die Produktmarke Excia sind eingetragene Marken der Aesculap AG.