



SPINE SURGERY

AESCULAP® TSPACE® TITANIUM

TRANSFORAMINAL LUMBAR INTERBODY FUSION SYSTEM SURGICAL MANUAL

AESCULAP® THORACOLUMBAR SPINE

PROTECTING AND PRESERVING SPINAL STABILITY

Modern lifestyle has resulted in increasing physical inactivity among people all over the world. Of the many medical problems associated with this, spinal disorders are among the most critical. This is even more significant as the spinal column is one of the most important structures in the human body.

It supports and stabilizes the upper body and is the center of our musculoskeletal system, which gives the body movement. Our work in the field of spine surgery is dedicated to protecting the spinal column and preserving its stability. We support spine surgeons with durable, reliable products and partner services for reliable procedures and good clinical outcomes (1-7).

Our philosophy of sharing expertise with healthcare professionals and patients allows us to develop innovative implant and instrument systems that help to preserve stability and stabilize the cervical and thoracolumbar spine.

RELIABLE PARTNER IN SPINE SURGERY

Discover our comprehensive product portfolio by clicking on this area.



www.bbraun.com/spine-surgery

AESCULAP® TSPACE® TITANIUM CONTENT

A B



A	GENE	RAL INFORMATION		
В	SURGICAL MANUAL			
	B.1.	POSITIONING OF THE PATIENT AND INCISION MARKING		
	B.2.	FASCIAL INCISION		
	В.З.	EXPOSURE AND BLUNT DISSECTION OF THE PARASPINAL MUSCLES		
	B.4.	INTRODUCTION OF THE SPINE CLASSICS RETRACTOR SYSTEM		
	B.5.	INSERTION OF PEDICLE SCREW		
	B.6.	REMOVAL OF FACET JOINT		
	B.7.	OPENING OF THE DISC AND REMOVAL OF DISC MATERIAL		
	B.8.	RESTORATION OF DISC HEIGHT		
	B.9.	CLEANING OF THE INTERVERTEBRAL SPACE		
	B.10.	DETERMINATION OF IMPLANT SIZE USING TRIAL IMPLANTS		
	B.11.	IMPLANT INSERTION		
	B.12.	FINAL IMPLANT POSITIONING		
	B.13.	APPLICATION OF ROD AND SET SCREW		
	B.14.	PEDICLE SCREW POSITIONING ON THE CONTRA-LATERAL SIDE		
С	IMPL	ANT & INSTRUMENT OVERVIEW		

A | GENERAL INFORMATION

PHILOSOPHY

TSPACE[®] TITANIUM implants are used for stabilization of the lumbar and thoracic spine through transforaminal approach, monosegmental and multisegmental.

Always use TSPACE[®] TITANIUM implants in conjunction with an internal fixator.

It is indicated for the treatment of thoracolumbar instabilities.

TSPACE® TITANIUM IS DESIGNED TO DELIVER

- PRIMARY STABILITY (11).
- RESTORATION OF THE NATURAL DISC HEIGHT AND LORDOSIS.



A | GENERAL INFORMATION

IMPLANT MATERIAL



The heart of this implant is a solid titanium alloy core (Ti6Al4V acc. to ISO 5832-3). The core is mantled with the proven PLASMAPORE[®] coating to increase the contact area between implant and endplate. PLASMAPORE[®] is a pure titanium coating (Ti/ISO 5832-2) which allows ingrowth of bone due to its balanced relationship between pore depth, porosity and roughness (2, 4, 8, 10).

Using a special manufacturing procedure, the implant surface is sprayed with pure titanium powder. Molten titanium particles settle on the core of the implant where they cool rapidly, building a firm form-lock between coating and core (Fig. 1).

In this way, each layer of the coating is built up and a favorable surface for bone ingrowth is created (2, 4, 8, 10) (Fig. 1).

AIM OF THE PLASMAPORE® COATING

I Primary Stability

The increased surface roughness of the PLASMAPORE[®] coating contributes to the primary stability of the motion segment (2).

I Secondary Stability

Bone growth into the coating is enabled due to the supportive features of PLASMAPORE[®], which leads to bone fusion between the adjacent vertebras with the implant (2, 4, 8, 10).

The coating concept, which has been proven as a result of many years of use in the field of hip prosthetics, is widely used in spinal surgery (2, 9, 12, 13).

IMPLANT FEATURES













PLASMAPORE® COATING

- High primary stability due to roughened surface which increases migration resistance and mechanical strength (2).
- High secondary stability due to bone ongrowth and ingrowth into the PLASMAPORE[®] structure (2, 4, 8, 10).

IMPLANT DESIGN

- Anatomical shape for a precise implant fit.
- Lordotic implant design.
- Serrated profile for a high primary stability (11).
- Load sharing on a large contact area.

IMPLANT VARIETY

Adequate range of sizes to enable the choice of implant size to fit the patient.

INSTRUMENT DESIGN

- Specifically designed and clearly arranged instruments.
- Straight and curved interbody inserter available.
- Inserters attach to implant via screw thread connection.

B | SURGICAL MANUAL



B.1. POSITIONING OF THE PATIENT AND INCISION MARKING

- A minimally invasive approach requires the patient to be placed on a radiolucent table which allows for AP views of the various anatomic structures (Fig. 1).
- The appropriate position of the longitudinal incision (4-5 cm in length) is determined by using a C-arm. The intended skin incision is marked paraspinally on the right and respectively on the left side (Fig. 2).



B.2. FASCIAL INCISION

A slightly curved fascial incision 1.5 cm from the mid-line is performed. This allows a firm hold of the speculum and counter retractor, facilitating the exposure of the individual segment (Fig. 3/4).



B.3. EXPOSURE AND BLUNT DISSECTION OF THE PARASPINAL MUSCLES

After splitting of the thoracolumbar fascia a blunt dissection of the paraspinal muscles is performed with the fingertip. In accordance with the palpatory finding, a correction of the skin incision is still possible, as the muscle retractor should be introduced as vertically as possible and in the direction of the interlaminar space. The length of the retractor is selected by using the index finger (Fig. 5/6).



B.4. INTRODUCTION OF THE SPINE CLASSICS RETRACTOR SYSTEM

- Spine Classics retractor system: see brochure 011402
- I The muscle retractor is introduced with closed blades and with the handle in the longitudinal direction. It is then turned 90° and expanded (Fig. 7/8).

B | SURGICAL MANUAL





Using the standard technique, the pedicle screws are inserted, e. g. AESCULAP[®] Ennovate[®] (Fig. 1).

For more information please visit www.bbraun.com/ennovate

INFORMATION

Thoracolumbar pedicle screw system shown in images is the AESCULAP $^{\circ}$ S $^{4^{\circ}}$ Spinal System.



B.6. REMOVAL OF FACET JOINT

A complete unilateral facetectomy should be considered on the side targeted for the implant insertion, for which an osteotome may be used. The inferior articular process of the facet joint is resected first, then the subjacent superior articular process is resected (Fig 2/3).



B.7. OPENING OF THE DISC AND REMOVAL OF DISC MATERIAL

- To open the disc a small window is cut into the annulus (Fig. 4).
- Rongeurs are used to remove the opened annulus (Fig. 5).
- Posterior osteophytes are removed by using Kerrisons.



B.8. RESTORATION OF DISC HEIGHT

- The desired restoration of the natural disc height can be set using the distractors. They are available in heights from 7-17 mm in 2 mm increments and can be attached to the t-handle in the set.
- I The distractor must be inserted horizontally and then rotated. Rotating clockwise the distractors are blunt and allow for a blunt height restoration maneuver. Rotating counterclockwise the distractors have a specially designed sharp rim to allow for removal of disc material (Fig. 6).

B | SURGICAL MANUAL





- I The disc space is cleared using rongeurs, bone curettes and box curettes (Fig. 1).
- I The bone rasps are used for endplate preparation.
- Alternatively, the box curettes can be used.

INFORMATION

Make certain that the endplates of the neighboring vertebral bodies are not weakened, in order to minimize the risk of migration.

Make certain that the implant bed is properly prepared to avoid damage to the implant when it is driven in.

Use the nerve root retractors to protect the dura during insertion.

B.10. DETERMINATION OF IMPLANT SIZE USING TRIAL IMPLANTS

- The TSPACE[®] TITANIUM trial implants are available in heights from 7 17 mm in 2 mm increments.
- Using the slap hammer the desired trial implant is inserted (Fig. 2).

INFORMATION

The trials are essential to ensure the correct implant size to be used.





B.11. IMPLANT INSERTION

- I It is recommended to place bone material harvested from the facet joint around the TSPACE[®] TITANIUM implant.
- The insertion guide supports the implantation of the TSPACE[®] TITANIUM implant (Fig. 3).
- The TSPACE[®] TITANIUM implant is inserted gradually into the disc space using the straight or the curved implant inserter (Fig. 4).
- The TSPACE[®] TITANIUM inserters attach to implant via a screw thread connection.



B | SURGICAL MANUAL



B.12. FINAL IMPLANT POSITIONING

- Carefully pull the inserter out of the implant. Avoid tilting of the instrument.
- Using the straight or angled impactor the implant is rotated 90° to achieve the final positioning (Fig. 1).
- **I** X-ray control verify the implant positioning.
- I It is recommended to put bone material harvested from the facet joint around the TSPACE[®] TITANIUM implant (Fig. 2).







- Final assembly of the thoracolumbar pedicle screw system (Fig. 3).
- Compression is applied to the pedicle screws to support the contact area between the TSPACE[®] TITANIUM implant and the endplates.
- I Final tightening of the pedicle screws and closure.



B.14. PEDICLE SCREW POSITIONING ON THE CONTRA-LATERAL SIDE

I Thoracolumbar pedicle screw system is applied on the contralateral side (Fig. 4).

To view our full portfolio please visit www.bbraun.com

INFORMATION

Thoracolumbar pedicle screw system shown in images is the AESCULAP* $\mathsf{S}^{\scriptscriptstyle 4^{\scriptscriptstyle 0}}$ Spinal System.

C | IMPLANT & INSTRUMENT OVERVIEW



Article No.	Size (Length x Width x Height)	Angle
SJ327T	26 x 12 x 7 mm	5°
SJ328T	26 x 12 x 8 mm	5°
SJ329T	26 x 12 x 9 mm	5°
SJ330T	26 x 12 x 10 mm	5°
SJ331T	26 x 12 x 11 mm	5°
SJ332T	26 x 12 x 12 mm	5°
SJ333T	26 x 12 x 13 mm	5°
SJ335T	26 x 12 x 15 mm	5°
SJ337T	26 x 12 x 17 mm	5°
SJ349T	26 x 12 x 9 mm	8°
SJ351T	26 x 12 x 11 mm	8°
SJ353T	26 x 12 x 13 mm	8°
SJ355T	26 x 12 x 15 mm	8°
SJ357T	26 x 12 x 17 mm	8°
SJ367T	30 x 12 x 7 mm	5°
SJ368T	30 x 12 x 8 mm	5°
SJ369T	30 x 12 x 9 mm	5°
SJ370T	30 x 12 x 10 mm	5°
SJ371T	30 x 12 x 11 mm	5°
SJ372T	30 x 12 x 12 mm	5°
SJ373T	30 x 12 x 13 mm	5°
SJ375T	30 x 12 x 15 mm	5°
SJ377T	30 x 12 x 17 mm	5°
SJ389T	30 x 12 x 9 mm	8°
SJ391T	30 x 12 x 11 mm	8°
SJ393T	30 x 12 x 13 mm	8°
SJ395T	30 x 12 x 15 mm	8°
SJ397T	30 x 12 x 17 mm	8°

FJ630 – TSPACE[®] TITANIUM INSTRUMENTATION

FJ633R TSPACE® TITANIUM – Preparation Trav

INSTRUMENTS	Article No.	Description	Quantity
	FJ679R	Left angled bone curette, 45°	1*
	FJ680R	Right angled bone curette, 45°	1*
	FJ698R	Left angled bone curette, 20°	1
	FJ699R	Right angled bone curette, 20°	1
	FJ681R	Curette, straight	1
	FJ682R	Left angled curette, 45°	1*
	FJ683R	Right angled curette, 45°	1*
	FJ702R	Left angled curette, 20°	1
	FJ703R	Right angled curette, 20°	1
	FJ658R	Straight osteotome, 8 mm	1
	FJ685R	Left angled bone rasp, 45°	1*
	FJ686R	Right angled bone rasp, 45°	1*
	FJ704R	Left angled bone rasp, 20°	1
	FJ705R	Right angled bone rasp, 20°	1
	FJ633R	Tray for preparation instruments	1
	JH217R	Wide Perf. basket lid	1

C | IMPLANT & INSTRUMENT OVERVIEW

FJ630 – TSPACE[®] TITANIUM INSTRUMENTATION

INSTRUMENTS	Article No.	Description	Quantity
-	FJ646R	T-Handle for distractors	1
	FJ647R	Distractor, 7 mm	1
	FJ649R	Distractor, 9 mm	1
5 minator F (40)	FJ651R	Distractor, 11 mm	1
	FJ653R	Distractor, 13 mm	1
	FJ655R	Distractor, 15 mm	1*
	FJ657R	Distractor, 17 mm	1*
	FJ666R	Slap hammer for trials	1
	FJ667R	TSPACE [®] trial, 7 mm	1*
	FJ669R	TSPACE [®] trial, 9 mm	1*
2 BADALAP ANTE AND Day	FJ671R	TSPACE [®] trial, 11 mm	1*
	FJ673R	TSPACE [®] trial, 13 mm	1*
	FJ675R	TSPACE [®] trial, 15 mm	1*
	FJ677R	TSPACE [®] trial, 17 mm	1*
	FJ619R	TSPACE [®] curved trial, 7 mm	1
	FJ621R	TSPACE [®] curved trial, 9 mm	1
	FJ623R	TSPACE [®] curved trial, 11 mm	1
	FJ625R	TSPACE [®] curved trial, 13 mm	1
	FJ627R	TSPACE [®] curved trial, 15 mm	1*
	FJ629R	TSPACE [®] curved trial, 17 mm	1*

FJ635R TSPACE® TITANIUM – Implantation Trav

INSTRUMENTS	Article No.	Description	Quantity
	FJ051R	Retractor S	1
	FJ052R	Retractor M	1
1994 - 2	FJ053R	Retractor L	1
	FJ054R	Retractor XL	1
	FJ700R	TSPACE® TITANIUM straight insertion instrument	1
	FJ701R	TSPACE [®] TITANIUM curved insertion instrument	1
	FJ661R	TSPACE [®] insertion guide	1*
	FJ662R	TSPACE [®] impactor, straight	1
	FJ663R	TSPACE [®] impactor, angled	1
	FJ635R	Tray for implantation instr.	1
	JH217R	Wide Perf. basket lid	1

PLEASE NOTE

Either the straight or the curved trials and inserter fit into the tray.

AESCULAP® TSPACE® TITANIUM REFERENCES

- (1) Tropiano P, Bronsard JJ, Louis C, Tallet JM, Sauget Y. Three Column Stabilisation through a Posterior Approach with a Titanium PLASMAPORE[®] Intervertebral Block (PROSPACE[®]) Radiological and Clinical Study after 4 Years. Rivista di Neuroradiologia. 1999;12(1 Suppl):89–94.
- (2) Kroppenstedt S, Gulde M, Schönmayr R. Radiological comparison of instrumented posterior lumbar interbody fusion with one or two closed-box PLASMAPORE[®] coated titanium cages: follow-up study over more than seven years. SPINE. 2008;33(19):2083-8.
- (3) Kreinest M, Schmahl D, Grützner PA, Matschke S. Radiological Results and Clinical Patient Outcome After Implantation of a Hydraulic Expandable Vertebral Body Replacement following Traumatic Vertebral Fractures in the Thoracic and Lumbar Spine. SPINE. 2017;42(8):E482-E489.
- (4) Takeuchi M, Yasuda M, Niwa A, Wakao N, Nakura T, Osuka K et al. PLASMAPORE[®]-coated titanium cervical cages induce more rapid and complete bone fusion after anterior cervical discectomy and fusion as compared to noncoated titanium cages. World Neurosurg. 2014;82(3-4):519-22.
- (5) Vanek P, Bradac O, Konopkova R, Lacy P de, Lacman J, Benes V. Treatment of thoracolumbar trauma by short-segment percutaneous transpedicular screw instrumentation: prospective comparative study with a minimum 2-year follow-up. J Neurosurg Spine. 2014;20(2):150-6.
- (6) Beisse R. Endoscopic surgery on the thoracolumbar junction of the spine. Eur Spine J. 2006;15(6):687-704.
- (7) Finger T, Bayerl S, Onken J, Czabanka M, Woitzik J, Vajkoczy P. Sacropelvic fixation versus fusion to the sacrum for spondylodesis in multilevel degenerative spine disease. Eur Spine J. 2014;23(5):1013-20.
- (8) Fink U. PLASMAPORE[®]: A plasma-sprayed microporous titanium coating to improve the long-term stability 1996; Acutalités en. Biomatériaux. (III:97-104).

- (9) Swamy G, Pace A, Quah C, Howard P. The Bicontact cementless primary total hip arthroplasty: long-term results. Int Orthop. 2012;36(5):915-20.
- (10) Kashii M, Kitaguchi K, Makino T, Kaito T. Comparison in the same intervertebral space between titanium-coated and uncoated PEEK cages in lumbar interbody fusion surgery. J Orthop Sci. 2019.
- (11) Mechanical testing of the AESCULAP® TSPACE® Titanium implant evaluating shear resistance. Tuttlingen, 2009. The AESCULAP® TSPACE® Titanium implant was tested in a shear force test following internal standards to evaluate the risk of implant dislocation related to shear forces occurring in situ. For this test, a synthetic bone-equivalent polyurethane test block was used. A quasi-static shear force was applied to the implant. The implant was axially loaded during the whole test with a constant force to simulate a clinically relevant loading situation. Based on the maximum shear load obtained in the test, the risk of implant dislocation due to shear forces at a physiological load situation is unlikely to occur in the case of the AESCULAP® TSPACE® Titanium implant.
- (12) Arregui R, Aso J, Martinez-Quinones JV, Consolini F, Lamban N, Dominguez M. Cespace. Cervical interbody fusion system. Preliminary retrospective study in 104 cases (120 implants). Neurocirugia (Astur) 2011; 22(6):542-53; discussion 553.
- (13) Cabraja M, Oezdemir S, Koeppen D, Kroppenstedt S. Anterior cervical discectomy and fusion: Comparison of titanium and polyetheretherketone cages. BMC Musculoskelet Disord 2012; 13(1):1-9.

AESCULAP[®] – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "AESCULAP" and the product trademarks "Ennovate", "PLASMAPORE", "S⁴" and "TSPACE" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.