en **Instructions for use/Technical description**

USA Caiman 5 Articulating, Sealing & Cutting Instruments

Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
Aesculap®
Caiman 5 Articulating, Sealing & Cutting Instruments

Legend

- Instrument**
- 1 RF activation button
 - 2 Rotation wheel
 - 3 Jaw articulation
 - 4 Moveable jaw part with marking
 - 5 Articulation activation lever
 - 6 Blade activation lever
 - 7 Jaw part activation lever
 - 8 RF cable with connector
- Generator**
- 9 SEE LEKTRAFUSE GENERATOR IFU
 - 10 Foot switch connection socket
 - 11 Regrasp signal lamp
 - 12 Display
 - 13 Error signal lamp
 - 14 Instrument connection socket
 - 15 "Power on" indicator light
 - 16 Foot control switch

Symbols on product and packages

	Sterilization with irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product
	Follow the instructions for use
	Date of manufacture

Scope

Aesculap Caiman 5 Articulating, Sealing & Cutting Instruments are sterile, disposable instruments. Caiman 5 instruments with the LEKTRAFUSE GENERATOR form an RF system for applying RF energy. Caiman 5 instruments can only be used in conjunction with the LEKTRAFUSE GENERATOR.

► For item-specific instructions for use and information regarding material compatibility, see Aesculap Extranet at www.extranet.bbraun.com

Intended use

Indications

Caiman Seal & Cut is a bipolar RF sealing system, which consists of the LEKTRAFUSE RF Generator and Caiman instruments. This system can be used for grasping, preparation, sealing and cutting of tissue during open and minimally invasive surgical procedures.

Caiman Seal & Cut can be used on vessels and vessel bundles with diameters up to and including 7 mm as well as soft tissue in general surgery and also surgical specialties such as gynecology, urology and bariatric, colorectal and thoracic surgery.

Absolute contraindications

- Caiman Seal & Cut is not suitable for grasping, preparation, sealing and cutting of hard tissue such as bone or cartilage.
- Caiman Seal & Cut is not approved for tubal sterilization, coagulation of the fallopian tubes, etc.
- Caiman Seal & Cut is not approved for application on parenchyma tissue of organs.

Relative contraindications

- Particular caution must be exercised in case of pathological changes of vessels such as atherosclerosis.

Since the safe and effective use of bipolar RF sealing systems strongly depends on factors such as tissue type, pathology and surgical procedure, the above information can only be used as a general guidance. The clinically successful application is dependent upon the knowledge and experience of the surgeon, who is responsible for deciding which structures can be reasonably treated and whether the patient's condition, e.g. the coagulation status, allows a successful treatment while taking into account the safety and warning instructions given in these instructions for use.

Safe handling and preparation

DANGER

Risk to the patient and/or user when using the instruments with RF sources other than the LEKTRAFUSE GENERATOR!
Non-observance can result in injury or death and restriction in the function of the instruments.

► Use the instruments only with the LEKTRAFUSE GENERATOR.

DANGER

Risk of death or injury to the patient through non-observance of the following directions!

► Use the instruments with caution on patients with cardiac pacemakers or defibrillators.

WARNING

Risk of injury to patient and/or user by using damaged or defective instrument!

► Before the instrument is used for the first time, lock the instrument jaws and then activate the blade trigger to verify that the instrument is functional.

► Prior to each use, inspect the instrument for loose, bent, broken, cracked or fractured components.

► Do not use the instrument if it is damaged or defective.

► If any unusual instrument performance or damage is observed, discontinue using the instrument immediately. Replace the instrument with a new one.

WARNING

Risk of injury due to evaporation of tissue fluid!
During the sealing process, tissue fluid can be converted into steam by the RF energy. This can result in collateral damage to the directly neighboring tissue.

► Absorb any fluids before activating the instrument.

WARNING

Risk of injury or damage to the instrument caused by incorrect operation!

► Do not activate the RF current before the instrument has been completely closed and locked.

► Do not open the instruments during the RF process as this will automatically interrupt the sealing process.

► Do not advance the blade during the RF process. Insufficient sealing or damage to the instrument may occur.

WARNING

Risk of injury from ignition or explosion of flammable gases such as oxygen and nitrous oxide or solvents such as alcohol and ether!
Sparks may occur when using the LEKTRAFUSE GENERATOR as directed.

► Do not use flammable gases for anesthesia or insufflation and wait until solvents volatilized.

WARNING

Risk of insufficient sealing due to contaminated electrode surfaces!

► Keep the instrument's contact surfaces clean during surgery, see chapter Intraoperative cleaning.

WARNING

Risk of injury to the patient and/or damage to the instrument when using incompatible or damaged trocars!

► When using the instrument through a trocar:

- Check the compatibility of the trocar and the instrument.
- Check the trocar for sharp edges and damage.

► Do not use the instrument with metal trocars.

► Do not use the instrument with damaged trocars.

WARNING

Risk of injury to patient and/or damage to the instrument when removing the instrument through a trocar!

► Before inserting or removing the instrument, place the articulation mechanism in the neutral position and lock the instrument jaws.

► Insert the instrument into the trocar and remove it only when locked, and carefully.

- The Caiman system may only be used by surgeons who have been trained in the usage of the products and who have the required knowledge and experience for the specific surgical procedures.
 - To prevent injury or burn to patients and medical personnel, fundamental understanding of the principles and techniques of RF surgery is required.
- Read, follow, and keep the instructions for use.
- Only use product after consulting the instructions for use of the LEKTRAFUSE GENERATOR.
- Use the product only in accordance with its intended use, see Intended use.
- The product is gamma-sterilized and supplied in sterile packaging.
- The product must not be reused.
- Do not use products from open or damaged sterile packaging.
- Do not use the product after its use-by date.
- The accessory voltage rating of the product is 240 Vp.
- To avoid RF burns:
- Always keep the working end of the product in the user's field of vision whenever the RF power is activated.
- Prior to activating the LEKTRAFUSE GENERATOR, check that the working end of the product is not touching any electrically conductive accessories.
- Never place the product on or next to the patient.

Safe operation

WARNING

- Risk of injury when using the product beyond the field of view!**
- Apply and activate the product only under visual control.
 - Before activating product, ensure that no electrically conductive objects are in contact.
 - Do not place the product on or next to the patient.

WARNING

- Risk of damage to the tissue or the instrument due to incorrect use of the instrument with the jaws closed and locked!**
- Avoid the following actions when the jaws are locked:
 - Rotation of the shaft
 - Articulation of the jaw parts
 - Ensure that the jaws are open before articulating or rotating the instrument.

CAUTION

- Risk of insufficient sealing or transection of the tissue!**
- Do not bundle or fold the tissue which is to be sealed and cut, or clamp more than one layer between the jaw parts.
 - Grip the tissue that is to be sealed centrally between the electrodes.
 - Do not start the process if the instrument jaws are in conductive liquids (e.g. blood or saline solution).
 - Do not start the process if there are conductive objects (e.g. clamps, vessel clamps, clips, etc.) between the instrument jaws.

The software of the LEKTRAFUSE GENERATOR recognizes the instruments and applies the corresponding device settings. The duration of the sealing process is controlled by the generator software.

The sealing process is started by tapping on an RF button once. The sealing process is interrupted by again tapping an RF activation button.

Note

In case of unsuccessful or incomplete sealing, as well as insufficient functionality of the instruments, a visual and acoustic alarm will be generated by the generator.

Note

In the event of visibly insufficient sealing of the tissue, the sealing quality can be improved by again activating the RF process.

Tissue manipulation, sealing of vessels and tissue bundles and tissue separation

- Remove the instrument from the sterile packaging.
- Insert the plug **8** of the instrument into the connecting socket **14** of the LEKTRAFUSE GENERATOR.
- In order to bring the jaw parts into the desired position, turn the rotation wheel **2**. The indication feature of the rotation wheel **2** indicates the central position.
- Bring the tip of the instrument into the desired position for surgery by rotating the shaft and applying an articulation movement to the articulation activation lever **5**.
- Grip the tissue that is to be sealed and/or cut between the jaw parts.
- Lock the moveable jaw part **4** using the activation lever **7**.
- In order to seal the clamped tissue, actuate the RF activation button **1** on the instrument or foot switch **16**.

Note

For additional options see the instructions for use of the LEKTRAFUSE GENERATOR.

- In order to cut clamped and sealed tissue, use the blade activation lever **6**.

WARNING

- Risk of injury to patient and/or damage to the instrument due to clamping of metallic objects and/or strong, rigid or very thick tissue!**
- Do not grasp, seal or cut metallic objects such as clips or clamps.
 - Do not grasp, seal or cut strong, rigid or very thick tissue such as bone or cartilage.

- If an error report appears:
 - Do not cut clamped tissue.
 - Unlock and open the instrument jaws.
 - When sealing quality is sufficient, the sealed tissue can be cut using surgical scissors or the instrument.
 - When sealing quality is insufficient, the sealing process must be restarted. If the additional RF cycle leads to a normal seal, the tissue can be cut, and the instrument can be opened and removed.
- If repeated error messages occur, discontinue using the instrument immediately. Replace the instrument with a new one.
- Do not leave the jaw actuation lever in the locked position when the instrument is not in use.

Intraoperative cleaning

WARNING

- Risk of injury to patient and/or user caused by incorrect cleaning!**
- Do not activate the RF current when cleaning the instrument.

CAUTION

- Risk of damage to and/or failure of the instrument caused by incorrect cleaning! This can result in consistent REGRASP INDICATOR – SHORT error messages (SEE LEKTRAFUSE GENERATOR IFU).**
- When cleaning the instrument:
 - Do not immerse the instrument jaws in liquid.
 - Do not use abrasives.
 - Do not apply excessive force on the jaws.

- Clean the instrument only as necessary:
 - If the blade becomes difficult to activate.
 - If the instrument sticks to the tissue.
 - If electrode surfaces become encrusted with tissue residues or body fluids.

When cleaning is required:

- Lock the instrument jaws and then cycle the blade trigger five to six times. This action may remove residual tissue debris from the proximal part of the jaws.
- If further cleaning is necessary, clean jaws with a moistened swab. Remove excess liquid from swab.

Note

If available, sterile (regular or deionized) water is preferred rather than saline or chemical cleaning solutions.

After cleaning the instrument:

- Inspect both inner and outer features of the jaws for bent or broken parts and for other damage.
- Do not use the instrument if it is damaged or defective.
- If repeated error messages occur, discontinue using the instrument immediately. Replace the instrument with a new one.

Single-use products

WARNING

- Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!**
- Do not reprocess the product!

Storage

- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Troubleshooting

Errors, faults and warnings are displayed on the LEKTRAFUSE GENERATOR as follows:

- An individual sound or a clear sound sequence and interruption of the sound sequence of RF operation
- The yellow regrasp error light or the red error light on the front side of the LEKTRAFUSE GENERATOR lights up
- Notification in the display on the front side of the LEKTRAFUSE GENERATOR
- For details on errors, faults, warnings and remedying such, see the instructions for use of the LEKTRAFUSE GENERATOR.

Technical Service

WARNING

- Risk of injury and/or malfunction!**
- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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