Celsite® PICC–Cel

- en – Peripherally Inserted Central Catheters (PICC).
- fr – Cathéters centraux insérés par voie périphérique (CCVP).
- de – Peripher einführbarer Katheter (PICC).
- es – Catéteres centrales de inserción periférica (PICC).
- it – Cateteri centrali inseriti perifericamente (PICC).
- sv – Perifert inlagda centrala katetrar (PICC).
- pt – Catéteres Venosos Centrais Inseridos Perifericamente (PICC).
- zh – 经外周中心静脉导管 (PICC)。
- ru – Периферически вводимый центральный катетер (ПВЦК).
- nl – Perifeer ingebrachte centraal veneuze katheters (PICC).
- no – Perifert innlagte sentrale katetre (PICC).
- da – Perifert anlagte centralvenekatetre (PICC).
- fi – Perifeerisesti asetettavat keskussaskimokatetrit (PICC).
- pl – Cewniki centralne z dostępu obwodowego (PICC).
- el – Περιφερικά εισαγόμενοι κεντρικοί φλεβικοί καθετήρες (PICC).
- tr – Periferal yoldan yerleştirilen Santral Kateterler (PICC).
- cs – Periferně zaváděné centrální katétry (PICC).
- th – อุปกรณ์สอดพวงกับคอมพิวเตอร์สำหรับใส่เข้าในสายสวน (PICC)
- ko – 말초중심정맥 카테터 (PICC).
- vi – Ống thông Tĩnh mạch Trung tâm (PICC).
- ro – Catetere centrale cu inserție periferică (PICC).
- bg – Централни катетри за периферно включване (PICC).
- et – Perifeerselt sisestatud tsentraalveeni kateetrid (PICC).
- hu – Perifériásan bevezetett centrális katéterek (PICC).
- lt – Periferiškai įvedami centrinės venos kateteriai.
- sk – Periférne vložené centrálné katétre (PICC).

B | BRAUN

ENGLISH	P 1
FRANÇAIS	P 3
DEUTSCH	P 5
ESPAÑOL	P 7
ITALIANO	P 9
SVENSKA	P 11
PORTUGUÊS	P 13
中文	P 15
РУССКИЙ	P 17
NEDERLANDS	P 19
NORSK	P 21
DANSK	P 23
SUOMI	P 25
POLSKI	P 27
ΕΛΛΗΝΙΚΑ	P 29
TÜRKÇE	P 31
CESKY JAZYK	P 33
ภาษาไทย	P 35
한국어	P 37
TIẾNG VIỆT	P 39
ROMÂNĂ	P 41
БЪЛГАРСКИ	P 43
EESTI KEEL	P 45
MAGYAR	P 47
LIETUVIŲ KALBA	P 49
SLOVENSKÝ	P 51

en - INSTRUCTIONS FOR USE
fr - MODE D'EMPLOI
de - BEDIENUNGSANLEITUNG
es - INSTRUCCIONES DE USO
it - ISTRUZIONI PER L'USO
sv - BRUKSANVISNING
pt - INSTRUÇÕES DE USO
zh - 使用说明
ru - ИНСТРУКЦИЯ ПО ИСПОЛЬЗОВАНИЮ
nl - GEBRUIKSIINSTRUCTIES
no - BRUKSANVISNING
da - BRUGERVEJLEDNING
fi - KÄYTTÖOHJEET
pl - INSTRUKCJA UŻYCIA
el - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
tr - KULLANIM KILAVUZU
cs - NÁVOD K POUŽITÍ
th - สำหรับใช้ในการเอ็กซ์เรย์
ko - 사용설명서
vi - HƯỚNG DẪN SỬ DỤNG
ro - INSTRUȚIUNI DE UTILIZARE
bg - УКАЗАНИЯ ЗА УПОТРЕБА
et - KASUTUSJUHEND
hu - HASZNÁLATI UTASÍTÁS
lt - NAUDOJIMO INSTRUKCIJOS
sk - NÁVOD NA POUŽITIE

en

Complementary Information Celsite® PICC-Cel IFU

II. INDICATIONS

Celsite® PICC-Cel is indicated for intravenous therapies which includes drugs such as chemotherapy, antibiotics, antivirals and contrast media. The Celsite® PICC-Cel may be used for up to 3 months.

V. WARNING

- Patient movement can cause catheter tip movement.
- The catheter must be secured in place to minimise the risk of catheter breakage and subsequent embolisation.
- To prevent shearing of the guide-wire never withdraw the J guide-wire through the Seldinger needle. Remove the guide-wire and the dilator together. Do not remove the guide-wire through the dilator as this may result in the guide-wire unravelling.

Central Venous Pressure Monitoring:

Central venous pressure (CVP) monitoring should always be used in conjunction with other patient assessment metrics when evaluating cardiac function.

Prior to conducting central venous pressure monitoring:

- ensure proper positioning of the catheter tip
- flush catheter vigorously with sodium chloride (NaCl 0.9%)
- ensure the pressure transducer is at the level of the right atrium

It is recommended that a continuous infusion of sodium chloride (NaCl 0.9%) (3 mL/hr) is maintained through the catheter while measuring CVP to improve the accuracy of CVP results.

Use your institutions protocols for central venous pressure monitoring procedures.

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1 – DESCRIPTION OF THE DEVICE

Celsite® PICC-Cel are peripherally inserted central venous catheters (PICC) made of polyurethane (PUR) with a reverse taper design. The catheters are radio-opaque and MR Safe, if used with the Ultrasite needleless valve the PICC-Cel should be considered MR Conditional.

MR Information

General Information

According to IEC Standard 60601-2-33(2008), the scanner must be operated in Normal Operating Mode (defined as the mode of operation of the MR system in which none of the outputs have a value that causes physiological stress to the patient):

- The whole body averaged specific absorption rate (SAR) must be ≤ 2.0 W/kg
- The head SAR must be < 3.2 W/kg

Device Information

MR Conditional

The ULTRASITE valve was determined to be **MR-Conditional** according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08, Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

Non-clinical testing demonstrated that the ULTRASITE valve is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla and 1.5-Tesla
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less
- Maximum whole body averaged specific absorption rate (SAR) of 2.9W/kg for 15 minutes of scanning

MRI-Related Heating

In non-clinical testing, the ULTRASITE valve produced the maximum temperature rise during MRI performed for 15-min (i.e., per pulse sequence) in 3-Tesla (Excite HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems, as:

MRI Condition	MR System Reported, Maximum Whole Body Averaged SAR (W/kg)	Calorimetry Value (W/kg)	Highest Temperature Change	Time for MRI (per pulse sequence)
3-T / 128-MHz	2.9	2.7	1.6 °C	15 min

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the ULTRASITE valve. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

II – INDICATIONS

PICC lines may be used when a central line is needed for repeated short to long-term intra-venous therapy, blood sampling or transfusion, or nutrition.

They may also be used for Contrast Enhanced Computed Tomography (CECT) using high pressure injection (See "IX Recommendation for High Flow rate / High Pressure Injection) and for Central Venous Pressure (CVP) monitoring.

III – CONTRAINDICATIONS

- Known or suspected infection of the puncture site, bacteraemia or septicæmia.
- Known or suspected allergy to any of the materials contained in the device.
- If the patient's anatomy does not allow the insertion or navigation of the catheter to the chosen site or if the patient has had previous radiotherapy to the insertion site.
- Previous venous thrombosis

IV – POTENTIAL COMPLICATIONS

• Air embolism • Cardiac arrhythmia, tamponade, • Catheter rupture, migration or fragmentation • Catheter occlusion • Embolism	• Endocarditis • Extravasation • Fibrin sheath formation • Haematoma/bleeding • Local or general infection/sepsis • Nerve injury	• Risks associated with local or general anaesthesia or with iodine based contrast media • Thrombo-embolism / Thrombophlebitis • Venous thrombosis
--	---	--

Before implantation of the PICC-Cel these, and other, well documented potential complications should be considered.

V – WARNINGS

- During storage avoid freezing or excessively high temperatures, storage at room temperature is recommended.
- Celsite® PICC-Cel are for single use only, do not re-sterilise any of the components, and destroy after use.
- The device and its accessories are not re-usable nor designed to be re-used. Any re-use would definitely compromise the performance and safety of the device.
- Sterilised using Ethylene Oxide. Sterile in unopened, undamaged individual packages.
- Check the expiry date before use and never implant a device past the expiry date.
- Do not open the tray until ready for use and ensure that all handling, maintenance and use and removal is under strict aseptic conditions.
- Always verify the integrity of the device before placement.
- This product should only be implanted, accessed or removed by personnel experienced in the technical and clinical aspects of PICC lines.
- Never use incompatible drugs together or in sequence: rinse the catheter of all residues with 10 mL of sodium chloride (NaCl 0.9%) between each medication and at the end of treatment.
- Avoid repeated and prolonged exposure to alcohol, antiseptics containing alcohol, acetone and ointments containing polyethylene glycol.

VI – SPECIAL ATTENTION

i The Patient

- Ensure that use of the PICC line is indicated taking into consideration the patient's anatomy, needs and proposed use of the PICC line.
- Check the condition of the skin over the puncture/insertion site.
- During insertion place the patient slightly inclined head down.

ii Catheter and Catheter Position

- The catheter tip should always be placed in the Superior Vena Cava at the entrance to, and not inside, the right atrium.
- Ensure that there is no kinking of the catheter along the vein.
- A radiograph should be taken to confirm catheter tip position. Tip position should be monitored routinely according to institutional policy.
- Always verify that the catheter is functional by aspirating 2 mL of blood into a syringe and injecting 5 mL of sodium chloride (NaCl 0.9%) before attempting to start an infusion.
- Stop treatment immediately if any pain or swelling is noted or if blood return is absent.

iii Occlusion

- Always ensure that the system is patent before using.
- If resistance is encountered to aspiration or injection, the lumen of the catheter may be partially or completely occluded. Attempt to inject a few mL of sodium chloride (NaCl 0.9%) into the catheter. If resistance to injection continues, or if swelling occurs along the catheter, device malfunction should be suspected.
- In case of obstruction, never try to clear the blockage using a small syringe or fluid under high pressure, this carries the risk of catheter fracture and migration.

iv Infection

- In case of infection, if appropriate, treat with antibiotic drugs according to local protocols. If this fails, or is not appropriate, the catheter should be removed.

VII – IMPLANTATION TECHNIQUES

General Considerations

- It is recommended that sodium chloride (NaCl 0.9%), with or without heparin according to local protocols, be used to flush the catheter and stylet during implantation.
- When trimming the catheter do not cut the stylet
- Never use force to remove the stylet. Resistance can damage the catheter. Take care not to damage, tear or perforate the catheter when using the stylet.
- During implantation ensure that the catheter is not damaged by ungarded forceps, suture needles or other sharp instruments.
- Do not clamp the lumen of the catheter. Use only the in-line clamps provided.
- It is recommended that the catheter be positioned above the antecubital fossa to reduce the risk of phlebitis and catheter fracture.

Per-cutaneous Insertion Techniques

- Position the arm at 90° to the patient's body; place a tourniquet above the intended insertion site to dilate the veins.
- Select the vein (recommended veins: basilic, cephalic, medial or lateral brachial). The use of ultrasound guidance for vein location is recommended to avoid accidental arterial or nerve puncture. Flush the catheter and stylet with either normal sodium chloride (NaCl 0.9%) or heparinised saline.
- Insert the Seldinger needle, attached to a syringe, into the chosen vein; verify the position by observing blood reflux.
- Insert the guide-wire through the needle 15-20 cms into the vein. Release the tourniquet.
- Fluoroscopic control is recommended when advancing the guide wire past the axilla.
- Leaving the guide-wire in place, withdraw the needle.
- Thread the assembled dilator and peelable sheath over the guide-wire using a twisting motion to pass through the skin planes.
- Determine correct catheter length by advancing the guide-wire, using fluoroscopy control, to the desired catheter tip location (at the entrance to the right atrium).
- Clamp the guide wire at the hub of the dilator with forceps to indicate the desired catheter length.
- Remove the catheter-wire from the vein. To minimise blood loss and reduce the risk of air embolism while the catheter is being trimmed, the patient should perform the Valsalva manoeuvre or a syringe or end cap may be placed on the peelable sheath.
- Using the clamped guide-wire, measure the catheter, withdraw the stylet and trim the catheter to the desired length. Alternatively, the catheter length may be determined using anatomical landmarks and the tape measure provided in the accessories.
- Re-insert the stylet to the distal end of the catheter and insert the catheter and stylet into the peelable sheath.
- Slowly advance the catheter to the desired position. Check catheter tip position using fluoroscopy.
- Remove the stylet from the catheter, the catheter can be stabilised by maintaining pressure on the catheter distal to the peelable sheath.
- Remove the peelable sheath ensuring that the catheter is not displaced, the '0' mark on the catheter should be at the insertion site.
- Aspirate to ensure free blood reflux and flush with sodium chloride (NaCl 0.9%) or heparinised saline.
- Secure the catheter with the catheter securement device and cover with a transparent dressing. A tape may also be placed around the catheter hub for additional security.
- The final catheter tip position should be checked radiographically.

VIII – MAINTENANCE OF THE PICC LINES

i Flushing and patency

- Rinsing of the catheter is essential.
- Rinse the catheter with 10 mL of sodium chloride (NaCl 0.9%) before treatment, between different drugs and after each use. If heparinised saline is used the system should be rinsed first with sodium chloride (NaCl 0.9%) alone, as some drugs may react with heparin and result in blockage of the catheter due to the formation of precipitates.
- Follow local protocols for flushing frequency and heparin concentration, if used.
- Always rinse the catheter between administration of different drugs. Special care should be taken with drugs which have a higher risk of precipitation, with anti-coagulation agents, after blood sampling or after transfusion to reduce the risk of catheter occlusion.
- Failure to maintain the system may result in occlusion of the catheter.

ii Dressing

- Check under the dressing for accumulation of blood, fluid or moisture in the first 24 hours.
- An occlusive dressing should cover the insertion site at all times.
- Check the dressing regularly for cleanliness, change the dressing according to local protocols or if the dressing becomes soiled, wet or non-occlusive.
- Check catheter position during the dressing changes by measuring the external length of the catheter.
- Catheter tip position, placement and patency should be checked regularly.

iii Catheter removal only

- Remove the dressing
- Hold the catheter near the insertion site and pull slowly. Do not use excessive force. Do not pull against resistance.
- If resistance is felt, apply a warm compress, wait 20-30 minutes and attempt removal again.
- Verify that the length of catheter removed corresponds to the length inserted (noted in patient's file).

IX – RECOMMENDATIONS FOR HIGH FLOW RATE/HIGH PRESSURE INJECTION

Warnings:

- Always verify that the catheter is functional by aspirating 2 mL of blood into a syringe and vigorously injecting a 5 mL of sodium chloride (NaCl 0.9%) into the catheter before attempting to start an infusion.
- Do not exceed the recommended pressure (300 psi - 20.7 bars) and flow rate (5 mL/sec) as system failure may occur.
- Contrast media should be warmed to 37°C (98.6°F) according to drug manufacturers recommendations. Failure to follow this recommendation will result in up to 50% lower flow rates and/or catheter or injection system failure.
- Depending on the technical characteristics of the injector system, the target flow rate might not be attained.
- Flush the catheter with 10 mL of sodium chloride (NaCl 0.9%) before, and after, using the catheter for CECT, followed by usual rinsing procedures.

Ref	Size/N°lumens	Length	Material	Max pressure psi/bar	Max flow rate mL/sec	Priming volume mL
04439001	4F* SL	51 cm	PUR	300	5	0.6
04439002	5F* SL	61 cm	PUR	300	5	0.7
04439003	5F* DL	56 cm	PUR	300	5	0.5
04439004	6F* DL	61 cm	PUR	300	5	0.6

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