- (GB)** Instructions for use
(USA) CO₂ cartridge adapter/CO₂ cartridge
(D) Gebrauchsanweisung
 CO₂-Kartuschenadapter/CO₂-Kartusche
(F) Mode d'emploi
 Adaptateur de cartouche de CO₂/cartouche de CO₂
(E) Instrucciones de manejo
 Adaptador para cartuchos de CO₂/Cartuchos de CO₂
(I) Istruzioni per l'uso
 Adattatore per cartucce di CO₂/cartuccia di CO₂
(P) Instruções de utilização
 Adaptador de cartuchos CO₂/Cartuchos CO₂
(NL) Gebruiksaanwijzing
 CO₂-patroonadapter/CO₂-patroon
(S) Bruksanvisning
 CO₂-patronadapter/CO₂-patron
(RUS) Инструкция по применению
 Рукоятка для картриджа CO₂/картридж CO₂
(CZ) Návod k použití
 Adaptér kartuše CO₂/kartuše CO₂
(PL) Instrukcja użytkowania
 Adapter naboju CO₂/nabój CO₂
(SK) Návod na používanie
 Kartušový adaptér CO₂/Kartuša CO₂
(TR) Kullanım Kılavuzu
 CO₂ kartuş adaptörü/CO₂ kartuşu

B | BRAUN
 SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
 Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 013148 04/13 V6 Änd.-Nr. 47095

CE 0123 – DIR 93/42/EEC

Technical alterations reserved


Fig. 1


Fig. 2


Fig. 3


Fig. 4

Aesculap®

CO₂ cartridge adapter/CO₂ cartridge

Legend

- 1 CO₂ cartridge adapter
- 2 Cartridge sleeve
- 3 Connector piece
- 4 CO₂ cartridge (unsterile)
- 5 Plastic seal
- 6 Sliding sleeve
- 7 Pneumatic holding arm

Symbols on product and packages

	Not for reuse in intended applications as defined by the manufacturer
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The CO₂ cartridge adapter is used for the provision of compressed gas for operating the UNITRAC holding system in unsterile areas.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

CAUTION

Loss of sterility of close-by instruments caused by use of the CO₂ cartridge in the sterile area or by unsterile exhaust air!

- Use the CO₂ cartridge and the CO₂ cartridge adapter only in unsterile areas of the O.R.

DANGER

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the CO₂ cartridge.

Note

Optionally, for joint storage and processing, the CO₂ cartridge adapter can be sterilized. The CO₂ cartridge adapter is intended for use in unsterile areas only.

Note

Use a new CO₂ cartridge for each operation.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product mechanically after removing its transport packaging and prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- Do not use the product after its use-by date.

Assembling

CAUTION

Impact damage to the plastic seal!

- Remove the protective foam cover from the CO₂ cartridge only immediately before use.

- Remove the protective foam cover from the CO₂ cartridge 4.
- Insert CO₂ cartridge 4 into cartridge sleeve 2, see Fig. 2.
- Align connector piece 3 on the top end of cartridge sleeve 2 and screw it down to the positive stop, see Fig. 3. The CO₂ cartridge 4 is pierced.

Note

The CO₂ cartridge sleeve can be difficult to release if tightened too firmly.

Installing the CO₂ cartridge adapter on the pneumatic holding arm

- Pull back sliding sleeve 6 and couple CO₂ cartridge adapter 1 at pneumatic holding arm 7, see Fig. 4.
- Release sliding sleeve 6.
- Sliding sleeve 6 engages audibly, with a click.

Note

Heavy or sluggish movement of the pneumatic holding arm indicates that the CO₂ cartridge is spent.

Changing the CO₂ cartridge

- Dismount the pressure unit and remove the CO₂ cartridge 4, see Disassembling.
- Insert new CO₂ cartridge 4 and couple CO₂ cartridge adapter 1, see Assembling.

Disassembling

WARNING

Risk of injury caused by very cold CO₂ cartridge, due to expansion of residual gas!

- Touch the CO₂ cartridge only at the plastic seal.

CAUTION

Risk of alarm caused by residual gas hissing out of the partially emptied CO₂ cartridge at the moment of removal!

- Inform O.R. personnel in good time of the impending CO₂ cartridge change.

- Pull back sliding sleeve 6 and decouple CO₂ cartridge adapter 1 from pneumatic holding arm 7.
- Unscrew CO₂ cartridge sleeve 2 from connector piece 3. As soon as residual gas escapes audibly, pause unscrewing until all residual gas has escaped.
- Remove CO₂ cartridge 4 from cartridge sleeve 2.

Validated reprocessing procedure

General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Make certain that the CO₂ cartridge 4 is removed from the CO₂ cartridge adapter 1, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

DANGER

Risk to patients!

- The product must only be cleaned mechanically!

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for PEEK plastics,
 - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable temperature of 55 °C.

CAUTION

Contamination or blocking of the CO₂ cartridge adapter with cleaning and disinfecting agents!

- Do not immerse/clean the CO₂ cartridge adapter in an ultrasonic cleaning bath and do not apply immersion disinfection.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Mechanical alkaline cleaning and thermal disinfection	Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).	Chapter Mechanical cleaning/disinfecting and sub-chapter: Chapter Mechanical alkaline cleaning and thermal disinfecting

Mechanical cleaning/disinfecting

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	Concentrate, alkaline: – pH = 13 – <5 % anionic surfactant 0.5 % working solution – pH = 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disin- fecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage to sealing components caused by oil spray on the CO₂ cartridge adapter!

Do not lubricate the CO₂ cartridge adapter.

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Ensure that the CO₂ cartridge has been removed from the CO₂ cartridge adapter, see Disassembling.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

Packaging

- Place the product in its holder or on a suitable tray.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note
The product may only be sterilized when dismantled.

CAUTION

Loss of sterility/damage to the CO₂ cartridge adapter, close-by instruments and the sterilizer caused by sterilization of the CO₂ cartridge adapter with CO₂ cartridge!

Make certain that no CO₂ cartridge is ever put into the sterilizer.

- Make certain that all external and internal surfaces of the product will be exposed to the sterilizing agent.
- Validated sterilization process

–

Disassemble the product

–

Steam sterilization through fractionated vacuum process

–

Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665

–

Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{–6}, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

Technical Service

WARNING

Risk of injury and/or malfunction!

Do not modify the product.

For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Art. no.	Designation
RT040R	UNITRAC holding arm
RT044SU	UNITRAC CO ₂ cartridge (unsterile)
RT043201	Cartridge sleeve for CO ₂ cartridge adapter

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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