
## CRANIOFIX2

|            |                                                                           |
|------------|---------------------------------------------------------------------------|
| <b>GB</b>  | <b>Instructions for use</b>                                               |
| <b>USA</b> | <b>CRANIOFIX2</b> titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm               |
| <b>D</b>   | <b>Gebruiksaanweisung</b>                                                 |
| <b>F</b>   | <b>CRANIOFIX2</b> -Titanklammern Ø 11 mm, Ø 16 mm, Ø 20 mm                |
| <b>E</b>   | <b>Mode d'emploi</b>                                                      |
| <b>I</b>   | <b>Pinces en titane</b> <b>CRANIOFIX2</b> Ø 11 mm, Ø 16 mm, Ø 20 mm       |
| <b>P</b>   | <b>Instrucciones de manejo</b>                                            |
| <b>NL</b>  | <b>Grapas de titânio</b> <b>CRANIOFIX2</b> Ø 11 mm, Ø 16 mm, Ø 20 mm      |
| <b>DK</b>  | <b>Istruzioni per l'uso</b>                                               |
| <b>S</b>   | <b>Clamp in titanio</b> <b>CRANIOFIX2</b> Ø 11 mm, Ø 16 mm, Ø 20 mm       |
| <b>GR</b>  | <b>Instruções de utilização</b>                                           |
| <b>RUS</b> | <b>Grampos de titânio</b> <b>CRANIOFIX2</b> Ø 11 mm, Ø 16 mm, Ø 20 mm     |
| <b>CZ</b>  | <b>Gebruiksaanwijzing</b>                                                 |
| <b>PL</b>  | <b>CRANIOFIX2</b> -Titaniumklemmen Ø 11 mm, Ø 16 mm, Ø 20 mm              |
| <b>SK</b>  | <b>Brugsanvisning</b>                                                     |
| <b>HU</b>  | <b>CRANIOFIX2</b> -Titanklemmer Ø 11 mm, Ø 16 mm, Ø 20 mm                 |
| <b>TR</b>  | <b>Bruksanvisning</b>                                                     |
| <b>CN</b>  | <b>CRANIOFIX2</b> Titanklämmor Ø 11 mm, Ø 16 mm, Ø 20 mm                  |
| <b>KR</b>  | <b>Οδηγίες χρήσης</b>                                                     |
|            | <b>Αγκιτρώων τιτανίου</b> <b>CRANIOFIX2</b> Ø 11 mm, Ø 16 mm, Ø 20 mm     |
|            | <b>Инструкция по применению</b>                                           |
|            | <b>CRANIOFIX2</b> -титановые клипсы Ø 11 mm, Ø 16 mm, Ø 20 mm             |
|            | <b>Návod k použití</b>                                                    |
|            | <b>CRANIOFIX2</b> -titánové svorky Ø 11 mm, Ø 16 mm, Ø 20 mm              |
|            | <b>Instrukcja użytkowania</b>                                             |
|            | <b>Zaciski tytanowe</b> <b>CRANIOFIX2</b> Ø 11 mm, Ø 16 mm, Ø 20 mm       |
|            | <b>Návod na používanie</b>                                                |
|            | <b>CRANIOFIX2</b> titánové svorky Ø 11 mm, Ø 16 mm, Ø 20 mm               |
|            | <b>Használati utasítás</b>                                                |
|            | <b>CRANIOFIX2</b> titán kapcsok Ø 11 mm, Ø 16 mm, Ø 20 mm                 |
|            | <b>Kullanım Kılavuzu</b>                                                  |
|            | <b>CRANIOFIX2</b> titanyum fiksasyon clamp'leri Ø 11 mm, Ø 16 mm, Ø 20 mm |
|            | <b>使用说明书</b>                                                              |
|            | <b>CRANIOFIX2</b> 钛夹 Ø 11 mm, Ø 16 mm, Ø 20 mm                            |
|            | <b>사용방법및사용시주의사항</b>                                                       |
|            | <b>두개성형판 고정제 (CRANIOFIX, 형명별 개별기재 )</b>                                   |

# B | BRAUN

## SHARING EXPERTISE

Aesculap AG

Am Aesculap-Platz  
78532 Tuttlingen  
Germany  
Phone +49 7461 95-0  
Fax +49 7461 95-2600  
www.aesculap.de

TA-Nr. 011200 06/11 Änd.-Nr. 43522

CE marking according to directive 93/42/EEC  
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG  
Marquage CE conforme à la directive 93/42/CEE  
Identificación CE en conformidad con la directriz 93/42/CEE  
Marchio CE conforme alla direttiva 93/42/CEE  
Símbolo CE, em conformidade com a Directiva 93/42/CEE  
CE-certificering conform richtlijn 93/42/EEG  
CE-markering iht. retningslinie 93/42/EEC  
CE-märkning i enlighet med direktiv 93/42/EEG  
Ενδείξη CE σύμφωνα με την Οδηγία 93/42/EEC  
Маркировка CE согласно директиве 93/42/ЕЭС  
CE-označení podle směrnice 93/42/EHS  
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE  
Označenie CE podľa smernice 93/42/EHS  
CE-jelölés a 93/42/EGK irányelv szerint  
93/42/EEG direktívi uyarınca CE işareti  
根据 93/42/EEG 指令的 CE 标志

Technical alterations reserved  
Technische Änderungen vorbehalten  
Sous réserve de modifications techniques  
Sujeto a modificaciones técnicas  
Con riserva di modifichie tecniche  
Salvo alterações técnicas  
Technische wijzigingen voorbehouden  
Retten til tekniske ændringer forbeholdes  
Med reservation för eventuella tekniska ändringar  
Επιφύλασσόμεθα για τεχνικές αλλαγές  
Мы оставляем за собой право вносить технические изменения  
Technické zmeny vyhradeny  
Zmiany techniczne zastrzeżone  
Technické zmeny vyhradené  
A muszaki változtatások joga fenntartva  
Teknik değışiklikler yapma hakkı saklıdır  
保留技术变动的权力

## CRANIOFIX2 titanium clamps

Ø 11 mm, Ø 16 mm, Ø 20 mm

### Figure captions

- Fig. 1**
- 1 **CRANIOFIX2** titanium clamp FF490T/FF490T-UNI, Ø 11 mm
  - CRANIOFIX2** titanium clamp FF491T/FF491T-UNI, Ø 16 mm
  - CRANIOFIX2** titanium clamp FF492T/FF492T-UNI, Ø 20 mm
- Fig. 2**
- 1 **CRANIOFIX2** titanium clamp FF491T/FF491T-UNI, Ø 16 mm/
  - CRANIOFIX2** titanium clamp FF492T/FF492T-UNI, Ø 20 mm
  - CRANIOFIX2** titanium clamp FF490T/FF490T-UNI, Ø 11 mm
- Fig. 3**
- 1 Correct position of **CRANIOFIX2** titanium clamps FF491T/FF491T-UNI/FF492T/FF492T-UNI: Trepanation hole Ø 12/15 mm
  - 2 Incorrect position of **CRANIOFIX2** titanium clamps FF491T/FF491T-UNI/FF492T/FF492T-UNI: Craniotomy gap
- Fig. 4**
- 1 **CRANIOFIX2** titanium clamp
  - 2 New craniotomy line
  - 3 Old craniotomy line

### Symbols on product

|  |                                                                                                  |
|--|--------------------------------------------------------------------------------------------------|
|  | Sterilization using irradiation                                                                  |
|  | This product is not for reuse in terms of its proper applications as defined by the manufacturer |
|  | Use by                                                                                           |
|  | Attention, see instructions for use                                                              |

### Product description

The Aesculap **CRANIOFIX2** titanium clamp system comprises sterile **CRANIOFIX2** titanium clamps for single use and diverse reusable neurosurgical instruments for applying titanium clamps.

### Intended use

The Aesculap **CRANIOFIX2** titanium clamp system is used for fixing craniotomized bone flaps and fractures at the brain skull (neurocranium). The **CRANIOFIX2** titanium clamp FF490T/FF490T-UNI, Ø 11 mm, is applied in the craniotomy gap (see **Fig. 2**) or in the fracture gap. The **CRANIOFIX2** titanium clamp FF491T/FF491T-UNI, Ø 16 mm, is applied in a trepanation hole of Ø 12 mm (see **Fig. 3**) or for bridging the defect in case of a fracture. The **CRANIOFIX2** titanium clamp FF492T/FF492T-UNI, Ø 20 mm, is applied in a trepanation hole of Ø 15 mm (see **Fig. 3**) or for bridging the defect in case of a fracture.

### Materials

The **CRANIOFIX2** titanium clamps are made of titanium alloy (Ti6Al4V) acc. to ISO 5832-3. This titanium alloy is a biocompatible material. The **CRANIOFIX2** titanium clamps are non-ferromagnetic and MR-compatible up to 3.0 Tesla. Consequently, the magnetic fields applied during NMR tomography do not present any risk. Stronger magnetic fields or an enlargement of the imaging field can lead to a dramatically increased incidence of positioning errors and artifacts. The size of the artifacts can vary considerably, depending on the NMR-pulse frequency, and compromise the significance of the NMR image if the area concerned is in the immediate vicinity (within a distance of some millimeters) of the clamp position.

### Contraindications

- The following list shows some, but not all conceivable contraindications.
- Patients with systemic illnesses and metabolic disorders, oversensitivity to metals or allergies against the implant materials, inflammations in the region of the implant site, and patients that show insufficient compliance and co-operation
  - Bone conditions that rule out the application of **CRANIOFIX2** titanium clamps
  - Artificial cranial bone flaps
  - Bone tumors in the area on which the implant is supported
  - Degenerative bone diseases
  - Missing dura mater
  - Application in the facial skull (viscerocranium) and in the orbita or skull-base regions

### Notes

- Improper fixation or combination with other fixation methods (sutures, mini plates or clamp systems from other manufacturers) can lead to step formation, position shifts of the skull calotte and loosening or breakage of implant components. For maximum stability of a craniotomized bone flap, we recommend using three **CRANIOFIX2** titanium clamps in a triangular arrangement. For larger bone flaps or multiple fragments of a skull fracture, additional **CRANIOFIX2** titanium clamps may be required to ensure adequate stability.

- During application of **CRANIOFIX2** titanium clamps, friction between the upper titanium disc and the pin indentation can lead to spark formation. Appropriate safety measures must be taken in order to prevent this. Because of the fire hazard, surgical interventions of this kind must not be carried out in rooms where flammable anesthetics or other flammable gases, liquids, objects or oxidizing agents are stored.
- After the operation, exceptionally high mechanical loads on the fixated bone flap must be avoided.

### Safety measures

The operating surgeon should be thoroughly familiar with the surgical techniques required for using this product. The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. The operating surgeon is also responsible for complications resulting from misdiagnosis and flawed operation techniques.

### Possible side effects

Infections, damage to nervous tissue, hematomas and wound dehiscence are among the general perioperative risks. The patient should be made aware of these risks as well as other risks of neurosurgery, general surgery and the administration of anesthetics.

### Delivery format

Aesculap **CRANIOFIX2** titanium clamps are delivered in sterile condition and are designed for single use only. The **CRANIOFIX2** titanium clamps must not be resterilized.

**Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!**

- Do not reprocess the product.

### Provision and handling of the implants

#### CAUTION

Federal law restricts this device to sale by or on order of a physician

**Damage to/destruction of the titanium clamps caused by using a non-compatible applier!**

- Use the **CRANIOFIX2** applier (FF494R) only for **CRANIOFIX2**-titanium clamps (FF490T/FF490T-UNI–FF492T/FF492T-UNI).
- Use the **CRANIOFIX** applier (FF102R, FF106R and FF107R) only for the **CRANIOFIX** titanium clamps (FF099T–FF101T).

Prior to the surgical intervention, the operating surgeon must prepare an operation plan in order to ensure that a sufficient number of sterile **CRANIOFIX2** titanium clamps are available, aseptic operation conditions are in place, the **CRANIOFIX2** applier required and other **CRANIOFIX** instruments are complete and ready for operation, and the operating surgeon and the O.R. team are familiar with the instruments and operation technique. The craniotomy incisions must be performed in an osteoplastic fashion.

Correct handling of the **CRANIOFIX2** titanium clamps before and during the operation is essential for successful bone flap fixation.

The **CRANIOFIX2** titanium clamps are individually wrapped in protective packaging and gamma-ray sterilized.

- The **CRANIOFIX2** titanium clamps have to be kept in their original packaging until immediately before use. Prior to use, check the expiry date and make sure that the sterile packaging is intact.
- Make certain that the surfaces of the **CRANIOFIX2** titanium clamps are not bent or damaged.
- The **CRANIOFIX2** titanium clamps have been designed as a disposable product and must not be reused.
- Damaged or surgically excised **CRANIOFIX2** titanium clamps must not be used again.
- To avoid spark formation, immerse the **CRANIOFIX2** titanium clamps in a sterile liquid prior to application.

### Application of CRANIOFIX2 titanium clamps

#### Caution

The Aesculap **CRANIOFIX2** titanium clamp of Ø 16 mm may only be applied in a trepanation hole of Ø 12 mm (see **Fig. 3**) or for bridging the defect in case of a fracture.

The Aesculap **CRANIOFIX2** titanium clamp of Ø 20 mm may only be applied in a trepanation hole of Ø 15 mm (see **Fig. 3**) or for bridging the defect in case of a fracture.

The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. For maximum stability of one bone flap approx. 20 cm<sup>2</sup> in size, we recommend applying three **CRANIOFIX2** titanium clamps in the craniotomy gap or the trepanation hole, respectively, at equal distances to each other. Figure 2 shows a typical application of a craniotomized bone flap.

The following is a brief description of the surgical procedure:

- 1 Position the three **CRANIOFIX2** titanium clamps in equal distances to each other at the bone flap.
- 2 The lower discs are pushed in between the dura and the calvarium.
- 3 Replace the bone flap in its original position.
- 4 Apply the holding forceps on the first **CRANIOFIX2** titanium clamp. To ensure a safe fixation and avoid damage to the pin indentation, the pin of the **CRANIOFIX2** titanium clamp must be held in the recess intended for this purpose.
- 5 Position the **CRANIOFIX2** applier over the pin of the **CRANIOFIX2** titanium clamps.
- 6 Compress the handle of the **CRANIOFIX2** applier so that the upper disk approaches the lower disc along the pin indentation. The spring reeds engage in the pin indentation of the **CRANIOFIX2** titanium clamp, with an audible "click".
- 7 Do not remove the holding forceps until immediately before closing the **CRANIOFIX2** titanium clamp.
- 8 The **CRANIOFIX2** applier is equipped with an automatic strain relief mechanism, which ensures that the **CRANIOFIX2** applier is disengaged as soon as the maximum allowable force is applied. Once the **CRANIOFIX2** applier has disengaged, the handle should be fully compressed several times until there is no more audible clicking of the spring reeds as they engage in the indentation of the pin of the **CRANIOFIX2** titanium clamp.

The top disk can loosen after application if the pin is not cut off properly!

- After tightening the top disk, cut the excess pin off at a right angle to the pin's axis using the **CRANIOFIX** cutting forceps.
  - Be sure that the indentation above the top disk remains completely intact.
- 9 Use the cutting forceps to remove the remaining, outward projecting section of the pin. When doing so, the protruding pin of the **CRANIOFIX2** titanium clamp must be held securely to prevent the pin from snapping off into the operation field in an uncontrolled manner.
- 10 The same procedure is carried out with the other **CRANIOFIX2** titanium clamps.

The surgical procedure for applying **CRANIOFIX2** titanium clamps is described in the **CRANIOFIX2** brochure and in the literature list ([www.aesculap.de](http://www.aesculap.de)). Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact us directly:

#### Aesculap AG

Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 7461 95-0  
Fax: +49 7461 95-2600  
Website: [www.aesculap.de](http://www.aesculap.de)

Or in the US:

Aesculap Inc.  
Attn: Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood, MO 63042  
Aesculap Repair Hotline  
Phone: +1 800 214-3392  
Fax: +1 314 895-4420

### Patient ID

The Patient ID serves to keep a written record of important information on the implanted **CRANIOFIX2** titanium clamps and the surgical intervention that was carried out:

- Patient data
- Hospital data
- Date of provision
- Physician in charge and hospital
- Art. no. and lot no. of the **CRANIOFIX2** titanium clamps implanted

To facilitate the postoperative radiological examination, each patient should be issued with a Patient ID, which can be ordered separately from Aesculap. Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact Aesculap directly.

#### Removing the CRANIOFIX2 titanium clamps in a revision operation

In the rare case that the operation has to be repeated, the operating surgeon can leave the **CRANIOFIX2** titanium clamps in place on the bone flap. Aesculap recommends applying a second incision around the original clamp position. If this is not feasible (e.g. for anatomic reasons), Aesculap recommends a second incision within the original clamp site. Under no circumstances must the incision cut through a **CRANIOFIX2** titanium clamp. The operating surgeon is responsible for any decision about how to approach the revision intervention (see **Fig. 4**).

### Accessories

| Designation                                 | Art. no.   |
|---------------------------------------------|------------|
| <b>CRANIOFIX2</b> titanium clamps Ø 11 mm   | FF490T     |
| <b>CRANIOFIX2</b> titanium clamps Ø 11 mm   | FF490T-UNI |
| <b>CRANIOFIX2</b> titanium clamps Ø 16 mm   | FF491T     |
| <b>CRANIOFIX2</b> titanium clamps Ø 16 mm   | FF491T-UNI |
| <b>CRANIOFIX2</b> titanium clamps Ø 20 mm   | FF492T     |
| <b>CRANIOFIX2</b> titanium clamps Ø 20 mm   | FF492T-UNI |
| <b>CRANIOFIX2</b> applier, non-dismountable | FF494R*    |
| <b>CRANIOFIX</b> holding forceps            | FF105R*    |
| <b>CRANIOFIX</b> cutting forceps            | FF103R*    |
| <b>CRANIOFIX</b> extractor forceps          | FF104R*    |
| <b>CRANIOFIX</b> storage tray               | FF094P     |

\* Concerning the handling, functioning and reprocessing of this product, the instructions for use TA011161 must be observed.

### Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034  
USA

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