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
Aesculap®

ELAN 4 shortening sleeve (1-ring or 2-ring) GB861P to GB863P and lateral burr guard (1-ring) GB866R

Legend

- 1 GB861P ELAN 4 shortening sleeve (1-ring)
- 2 GB862P ELAN 4 MIS sleeve (1-ring)
- 3 GB863P ELAN 4 shortening ring (2-ring)
- 4 GB866R ELAN 4 lateral burr guard (1-ring)
- 5 GB727R ELAN 4 ECCOS holder

Symbols on product and packages

Caution
Observe important safety information such as warnings and precautions in the instructions for use.

Contents

1.	Applicable to	2
2.	Intended use	2
2.1	Task/function within system	2
2.2	Application Environment	2
2.3	Indications	3
2.4	Contraindications	3
3.	Safe handling	3
3.1	Non-sterile ELAN 4 shortening sleeves, MIS sleeves and lateral burr guard	4
4.	Product description	4
4.1	Package includes	4
4.2	Components required for operation	4
4.3	Operating principle	4
5.	Preparation	4
6.	Working with the applied part	4
6.1	System set-up	4
6.2	Function check	5
6.3	Operation	5
7.	Validated reprocessing procedure	5
7.1	General safety notes	5
7.2	General information	6
7.3	Preparations at the place of use	6
7.4	Preparation before cleaning	6
7.5	Cleaning/disinfection	6
7.6	Manual cleaning with immersion disinfection – ELAN 4 shortening sleeves and lateral burr guard	8
7.7	Mechanical cleaning/disinfection with manual pre-cleaning – ELAN 4 shortening sleeves and lateral burr guard	8
7.8	Inspection, maintenance and checks	9

7.9	Packaging	9
7.10	Steam sterilization	9
7.11	Sterilization for the US market	10
7.12	Storage	10
8.	Troubleshooting	10
9.	Technical Service	11
10.	Accessories/Spare parts	11
11.	Technical Data	11
12.	Disposal	11
13.	Distributor in the US/Contact in Canada for product information and complaints	11

1. Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

2. Intended use

2.1 Task/function within system

The ELAN 4 standard shortening sleeve for handpieces (1-ring), the ELAN 4 MIS standard sleeve for handpieces (1-ring) and the ELAN 4 HD shortening sleeve for handpieces (2-ring) are accessories of the ELAN 4 motor system. The ELAN 4 standard shortening sleeve for handpieces (1-ring), the ELAN 4 MIS standard sleeve for handpieces (1-ring) and the ELAN 4 HD shortening sleeve for handpieces (2-ring) are only for use in combination with the ELAN 4 air/electro handpieces standard (1-ring) or HD (2-ring).

The ELAN 4 lateral burr guard 1-ring is an accessory of the ELAN 4 motor system.

The ELAN 4 lateral burr guard 1-ring can only be used in combination with the ELAN 4 air/electro handpieces (1-ring). When using the ELAN 4 lateral burr guard 1-ring in combination with the ELAN 4 air/electro handpieces (1-ring), only 1-ring tools with a burr head diameter of max. 4.5 mm must be used.

Note

Use accessories only in combination with the specified handpieces!

2.2 Application Environment

Application in the sterile area

2.3 Indications

Modes of application	Separating, removing and shaping hard tissue, cartilage and similar, as well as for drilling holes in bone and bone substitute materials.
Surgical discipline/areas of application	Neurosurgery, ENT surgery, OMF surgery, orthopedics and trauma surgery

Note

The type and area of application depend on the tools selected.

2.4 Contraindications

The ELAN 4 electro motor system is not authorized for use in the central nervous system or central circulatory system.

Note

The safe and effective use of electrically driven applied parts depends heavily on influences that only the user can control him or herself. Therefore the specifications provided represent framework conditions only.

Note

The clinical success of the use of the ELAN 4 electro motor system is dependent on the knowledge and experience of the surgeon. The surgeon must decide which structures it is sensible to treat and take into account the safety and warning information contained in these instructions for use.

3. Safe handling

CAUTION

Federal law restricts this device to sale by or on order of a physician!

WARNING

Risk of injury and material damage due to inappropriate use of the product!

- Use the product only in accordance with the intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

This product is an accessory of the following products:

- ELAN 4 electro handpieces Standard (1-ring) GA861 – GA864
- ELAN 4 air handpieces Standard (1-ring) GA761 – GA764
- ELAN 4 electro handpieces HD (2-ring) GA865 – GA869
- ELAN 4 air handpieces HD (2-ring) GA765 – GA769
- Follow the instructions for use for ELAN 4 handpieces air/electro.
- Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Prior to use, check that the product is in good working order.
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.
- Ensure that the electrical installations at the premises where the device is to be used meet all requirements acc. to IEC/DIN EN.
- Do not use the control unit and applied part in potentially explosive areas.
- Carry out sterile processing in any reusable applied parts, motor cables and tools prior to use.
- When handling, observe instructions for use that are relevant to Aesculap holder systems TA009721, see Aesculap Extranet at <https://extranet.bbraun.com>

Aesculap®

ELAN 4 shortening sleeve (1-ring or 2-ring) GB861P to GB863P and lateral burr guard (1-ring) GB866R

3.1 Non-sterile ELAN 4 shortening sleeves, MIS sleeves and lateral burr guard

Note

Sterilize the products before each use!

The products are packaged unsterile and must be reprocessed before each use. Steam sterilization is used, see Steam sterilization.

The products may be reused.

- Only have the products operated and used by suitably qualified, trained and experienced personnel.
- Read, follow, and keep the instructions for use.
- Use the products only in accordance with their intended use, see Intended use.
- Prior to each use, visually inspect the products for: loose, bent, broken, cracked, worn, or fractured components.
- Never use damaged or faulty products. Immediately set any damaged products aside.

4. Product description

4.1 Package includes

Art. no.	Designation
GB861P	ELAN 4 shortening sleeve (1-ring) - or -
GB862P	ELAN 4 MIS sleeve (1-ring) - or -
GB863P	ELAN 4 shortening sleeve (2-ring) - or -
GB866R	ELAN 4 lateral burr guard (1-ring) - or -
TA015081	Instructions for use for ELAN 4 shortening sleeves (1-ring or 2-ring) GB861P to GB863P and lateral burr guard (1-ring) GB866R

4.2 Components required for operation

- Operation-ready ELAN 4 electro control unit GA800, see TA014401
- ELAN 4 electro handpieces Standard (1-ring) or HD (2-ring)
 - or -
- ELAN 4 air handpieces Standard (1-ring) or HD (2-ring)
- ELAN 4 tool (1 ring or 2 ring)

4.3 Operating principle

The ELAN 4 shortening sleeves and MIS sleeves are accessories for the ELAN 4 handpieces. They extend the holding area in the direction of the tool and enable a further holding position.

The ELAN 4 lateral burr guard is an accessory for the ELAN 4 handpieces, for the lateral covering of the rotating tool head. This prevents unwanted contact between tissue and the tool head.

5. Preparation

Non-compliance with these rules will result in complete exclusion of liability on the part of Aesculap.

- Prior to use, inspect the product and its accessories for any visible damage.
- Use the products and their accessories only if they are in perfect technical condition.

6. Working with the applied part

WARNING

Damage to the product if dropped!

- Use the products only if they are in perfect technical condition, see Function check.

6.1 System set-up

Note

Control elements on the system components of the ELAN 4 air/electro motor system are labeled with a gold marking.

Connecting the accessories

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

All configurations must comply with basic standard IEC/DIN EN 60601-1. Any individual connecting devices with one another is responsible for such configuration and must ensure compliance with basic standard IEC/DIN EN 60601-1 or applicable national standards.

- Follow the instructions for use of individual accessories.
- Please contact your B. Braun/Aesculap Partner or Aesculap Technical Service, address see Technical Service with any inquiries in this respect.

ELAN 4 Assembling the MIS sleeve and shortening sleeves

Note

First connect the holding sleeve and then the tool to the handpiece.

- ▶ Slide the holding sleeve onto the shaft from the front as far as the holding area of the handpiece. Turn slightly until it locks in place. Then the holding sleeve automatically locks with the holding area of the handpiece.
- ▶ Pull gently on the holding sleeve to check that it is securely coupled.
- ▶ To disassemble, remove the holding sleeve from the handpiece.

Assembling the ELAN 4 lateral burr guard

Note

When using the ELAN 4 lateral burr guard 1-ring in combination with the ELAN 4 air/electro handpieces (1-ring), only 1-ring tools with a maximum burr head diameter of $\varnothing$ 4.5 mm must be used.

- ▶ Select the position of the burr guard 4 on the handpiece shaft.
- ▶ Lock the burr guard on the handpiece shaft.
- ▶ If necessary, disassemble the burr guard, select a new position and re-assemble. The burr guard should be disassembled without a tool in the shaft.
- ▶ Pull on the lateral burr guard to check that it is securely coupled.


Fig. 1 Assembling the lateral burr guard.


Fig. 2 Disassembling the lateral burr guard.

6.2 Function check

- ▶ Prior to each use, check that all products to be used are in good working order.
- ▶ Check the secure connection of all products to be used.
- ▶ Briefly run the applied parts at maximum speed.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.

6.3 Operation

WARNING

Risk of injury and/or malfunction!

- ▶ Always carry out a function test prior to each use of the product.

WARNING

Risk of injury due to application of the product outside the field of view!

- ▶ Only use the product under vision.

Note

For further information, see instructions for use of ELAN 4 electro control unit GA800 (TA014401).

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

The recommended chemical was used for validation.

Aesculap®

ELAN 4 shortening sleeve (1-ring or 2-ring) GB861P to GB863P and lateral burr guard (1-ring) GB866R

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion:

- Do not exceed a period of 6 hours between use and reprocessing.
- Do not use fixating pre-cleaning temperatures >45 °C.
- Do not use fixating disinfecting agents (active ingredient: aldehyde, alcohol).

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result damage to metallic products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Material damage, e.g. corrosion, cracks, fracturing, premature aging, or swelling
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion.
- ▶ For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

7.3 Preparations at the place of use

- ▶ Separate the products from each other immediately after use.
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

7.4 Preparation before cleaning

ELAN 4 shortening sleeves and lateral burr guard

- ▶ Prior to first mechanical cleaning/disinfection process: Mount the bottom part of the ECCOS holder GB727R in a suitable tray.
- ▶ Place shortening sleeves with the larger open side facing upward in the tool holder of the ECCOS holder GB727R 5, see foldout page.
- ▶ Place the burr guard in a suitable tray.

7.5 Cleaning/disinfection

Product-specific safety guidelines on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning and disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents which
 - are approved for plastic material and high-grade steel,
 - do not attack softeners (e.g., in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C during chemical cleaning and/or disinfection.
- ▶ Do not exceed the maximum temperature of 90 °C during thermal disinfection with DI water.
- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

Note

The indicated drying time is a guide time only. It must be checked taking into account the specific conditions (e.g. load) and if applicable adjusted.

Validated cleaning and disinfection procedure

Product	Validated procedure	Reference
Shortening sleeves and lateral burr guard	Manual cleaning with immersion disinfection	see Manual cleaning with immersion disinfection – ELAN 4 shortening sleeves and lateral burr guard
	Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection	see Mechanical cleaning/disinfection with manual pre-cleaning – ELAN 4 shortening sleeves and lateral burr guard
Shortening sleeves and lateral burr guard	Steam sterilization	see Steam sterilization

Aesculap®

ELAN 4 shortening sleeve (1-ring or 2-ring) GB861P to GB863P and lateral burr guard (1-ring) GB866R

7.6 Manual cleaning with immersion disinfection – ELAN 4 shortening sleeves and lateral burr guard

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Pre-cleaning	RT (cold)	-	-	D-W	-
II	Cleaning	RT (cold)	>5	1	D-W	pH-neutral, pH ~ 9*
III	Intermediate rinse	RT (cold)	-	-	D-W	-
IV	Disinfection	RT (cold)	>15	2	D-W	Concentrate-free, aldehyde-free, phenol-free, QUAT-free, pH ~ 9**
V	Final rinse	RT (cold)	-	-	D-W	-
VI	Drying	RT	-	-	-	-

D-W: Drinking water
 RT: Room temperature
 *Recommended: BBraun Helizyme
 **Recommended: BBraun Stabimed

Phase I

Note

Only clean the outside of the shortening sleeves with a suitable cleaning brush.

Do not clean the inside of the shortening sleeves with a brush.

- Clean the ELAN 4 shortening sleeves and ELAN 4 lateral burr guard under running water until all visible residues have been removed from the surfaces.
- If necessary, use a suitable cleaning brush for outer surfaces.

Phase II

- Fully immerse the applied part in enzymatic cleaning solution for at least 5 minutes.

Phase III

- Rinse the entire product (all accessible surfaces) under running water.

Phase IV

- Prior to manual disinfection, allow rinsing water to drip off the product.
- Fully immerse the applied part in enzymatic cleaning solution for at least 15 minutes.

Phase V

- Rinse the entire product (all accessible surfaces) under running water.

Phase VI

- Dry the product in the drying phase with suitable equipment (e.g. cloth,

7.7 Mechanical cleaning/disinfection with manual pre-cleaning – ELAN 4 shortening sleeves and lateral burr guard

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Pre-cleaning	RT (cold)	-	-	D-W	-

D-W: Drinking water
 RT: Room temperature

Note

Only clean the outside of the shortening sleeves with a suitable cleaning brush.

Do not clean the inside of the shortening sleeves with a brush.

- ▶ Clean the ELAN 4 shortening sleeves and ELAN 4 lateral burr guard under running water until all visible residues have been removed from the surfaces.
- ▶ If necessary, use a suitable cleaning brush for outer surfaces.

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH ~ 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Place ELAN 4 shortening sleeves with the larger open side facing upward on the ECCOS holder GB727R 5, see fold-out page.
- ▶ Connect the ECCOS holder GB727R 5 to the rinsing connection of the rinsing cart.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

7.8 Inspection, maintenance and checks

- ▶ Allow the ELAN 4 shortening sleeves and ELAN 4 lateral burr guard to cool to room temperature.
- ▶ Check the ELAN 4 shortening sleeves and ELAN 4 lateral burr protection after every cleaning and disinfection process to ensure they are: clean, functional, and undamaged.
- ▶ Set aside the product if it is damaged.

7.9 Packaging

- ▶ Always follow the instructions for use of the respective packaging and storage devices (e.g. instructions for use TA009721 for the Aesculap ECCOS holding system).
- ▶ Place ELAN 4 shortening sleeves with the larger open side facing upward on the ECCOS holder GB727R 5, see figure on front page.

7.10 Steam sterilization

- ▶ Make certain that all external and internal surfaces of the product will be exposed to the sterilizing agent.
- ▶ Validated sterilization process:
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer: ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Aesculap®

ELAN 4 shortening sleeve (1-ring or 2-ring) GB861P to GB863P and lateral burr guard (1-ring) GB866R

7.11 Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

7.12 Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

8. Troubleshooting

Malfunction	Detection	Cause	Remedy
Shortening sleeve cannot be slid onto the shaft or handpiece from the front.	Shortening sleeve cannot be assembled with manual force.	Shortening sleeve defective	Use a new shortening sleeve.
Shortening sleeve cannot be securely coupled with the shaft or handpiece.	No clicking in place in the end position	Shortening sleeve defective	Use a new shortening sleeve.
Shortening sleeve cannot be securely coupled with the shaft or handpiece.	Shortening sleeve cannot be gently slid up and down the shaft.	Shortening sleeve defective	Use a new shortening sleeve.
Lateral burr guard cannot be attached to the shaft of the handpiece.	Lateral burr guard cannot be assembled with manual force.	Lateral burr guard defective/bent	Use a new lateral burr guard.
Lateral burr guard cannot be locked on the shaft.	Lateral burr guard cannot be gently slid up and down the shaft.	Lateral burr guard defective/bent	Use a new lateral burr guard.

9. Technical Service

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- ▶ **Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.**
- ▶ **Do not modify the product.**

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen /Germany

Phone: +49 7461 95-1601

Fax: +49 7461 14-939

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042

Aesculap Repair Hotline

Phone: +1 800 214-3392

Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

10. Accessories/Spare parts

Art. no.	Designation
TA015080	Instructions for use for ELAN 4 shortening sleeves (1-ring or 2-ring) GB861P to GB863P and ELAN 4 lateral burr guard GB866R (Extranet version)
TA015081	Instructions for use for ELAN 4 shortening sleeves (1-ring or 2-ring) GB861P to GB863P and ELAN 4 lateral burr guard GB866R (Leaflet)

11. Technical Data

Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
GB861P	ELAN 4 shortening ring (1-ring)	I
GB862P	ELAN 4 MIS sleeve (1-ring)	I
GB863P	ELAN 4 shortening ring (2-ring)	I
GB866R	ELAN 4 lateral burr guard (1-ring)	Ila

12. Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

13. Distributor in the US/Contact in Canada for product information and complaints

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