




Aesculap Endoscopic Technology

- (GB)** Instructions for use/Technical description
- (USA)** Full HD CMOS camera system PV470
- (D)** Gebrauchsanweisung/Technische Beschreibung
Full HD CMOS Kamerasystem PV470
- (F)** Mode d'emploi/Description technique
Système de caméra Full HD CMOS PV470
- (E)** Instrucciones de manejo/Descripción técnica
Sistema de cámara CMOS Full HD PV470
- (I)** Istruzioni per l'uso/Descrizione tecnica
Sistema videocamera Full HD CMOS PV470
- (P)** Instruções de utilização/Descrição técnica
Sistema de câmara Full HD CMOS PV470
- (NL)** Gebruiksaanwijzing/Technische beschrijving
Full HD CMOS-camerasysteem PV470


- (S)** Bruksanvisning/Teknisk beskrivning
Full HD CMOS kamerasystem PV470
- (RUS)** Инструкция по применению/Техническое описание
Система камеры Full HD с КМОП-матрицей PV470
- (CZ)** Návod k použití/Technický popis
Systém kamery CMOS full HD PV470
- (PL)** Instrukcja użytkowania/Opis techniczny
System kamery Full HD CMOS PV470
- (SK)** Návod na použitie/Technický opis
Full HD CMOS Kamerový systém PV470
- (TR)** Kullanım Kılavuzu/Teknik açıklama
Full HD CMOS kamera sistemi PV470


1.

2.

3.

Aesculap®

Full HD CMOS camera system PV470

Legend

- 1 "Power off" switch
- 2 "Power on" switch
- 3 LED display on/off
- 4 White balance button
- 5 Camera socket
- 6 Mains power
- 7 Fuse holder
- 8 DVI outputs
- 9 HD-SDI outputs
- 10 Equipotentialization connector
- 11 Communication interface (master)
- 12 S-Video output
- 13 Remote connection
- 14 USB connector for service purposes
- 15 Interface for MAQUET TEGRIS system (RS232 socket)
- 16 Safety lock
- 17 Endocoupler latch
- 18 Focusing ring
- 19 Zoom ring
- 20 Button A (up)
- 21 Button B (right)
- 22 Button D (left)
- 23 Button C (down)
- 24 Cap
- 25 Connector plug
- 26 Camera cable
- 27 Pendulum-locking lever

Symbols on product and packages

	Caution, general warning symbol Caution, see documentation supplied with the product
	Meaning dependent on the standard that applies in each case: 2nd edition IEC 60601-1: Please take note of the accompanying document! 3rd edition IEC 60601-1: Caution
	Follow the instructions for use
	Equipotentialization connection
	Unlock
	Lock
	Hazardous electrical voltage
	Fuse
	Alternating current
	Power ON switch
	Power OFF switch
	Applied part type CF
	Marking of electric and electronic devices according to directive 2002/96/EC (WEEE). see Disposal
	Date of manufacture

Contents

1.	Applicable to	3
2.	Safe handling	3
3.	Product description	4
3.1	Scope of supply	4
3.2	Components required for operation	4
3.3	Intended use	4
3.4	Operating principle	4
4.	Preparation and setup	6
4.1	Preparation	6
5.	Working with the camera control unit Full HD PV470 and the camera head PV471/PV472/PV473	7
5.1	System set-up	7
5.2	Function checks	8
5.3	Safe operation	8
6.	Validated reprocessing procedure	14
6.1	General safety instructions	14
6.2	Single-use products	14
6.3	General information	14
6.4	Preparations at the place of use	15
6.5	Cleaning/disinfection	15
6.6	Wipe disinfection for electrical devices without sterilization ..	16
6.7	Manual cleaning/disinfecting	16
6.8	Sterilization	17
6.9	Inspection, maintenance and checks	17
6.10	Packaging	17
6.11	Storage	17
7.	Maintenance	18
8.	Troubleshooting list	18
8.1	Fuse replacement	19
9.	Technical Service	20
10.	Accessories/Spare parts	20
11.	Technical data	20
11.1	Ambient conditions	21
12.	Disposal	21
13.	Distributor in the US/Contact in Canada for product information and complaints	21

1. Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

2. Safe handling

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

DANGER

Risk of death by electric shock!

- Do not open the product.
- Connect the product only to a grounded power supply.

WARNING

Risk to the patient if the product is not used gently/a damaged product is used!

- Handle the product with appropriate care.
- Do not use the product if it has been subjected to strong mechanical stress or if it has been dropped, and send the product to the manufacturer for inspection.

CAUTION

Malfunction in connection with magnetic resonance!

- Do not use the product in a magnetic resonance environment.

- Manually clean the new product after removing its transport packaging and prior to first use, and disinfect if necessary.
- Prior to use, check that the product is in good working order.
- Observe "Notes on Electromagnetic Compatibility (EMC)", see TA022130.
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Only allow products and accessories to be operated and used by individuals who have been trained in the use of endoscopic procedures.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.
- Ensure that the electrical facilities in the room in which the device is being used conform to IEC norms.
- Set up the devices in such a way that the ventilation slots at the back of the device are clear.
- Connect the devices to a suitable power source.
- Unplug the device by pulling on the plug, and never on the power cord.

Aesculap®

Full HD CMOS camera system PV470

- In order to avoid damage to the electrical components and inoperability:
 - Do not bend the product.
 - Avoid bending, torsional, tensile or compression loading.
- Operate the product with original manufacturer's parts only: System cable, video cable and power cord.
- The Full HD CMOS camera system is intended for daily to weekly use.

Note

The connection of electrical devices to a multiple socket creates a system and can lead to a reduced level of safety.

Note

The camera case is only transport packaging.

3. Product description

3.1 Scope of supply

Camera control unit

Art. no.	Designation
PV470	Full HD CMOS camera control unit
–	DVI-cable
TA013748	Instructions for use
TA022130	Notes on electromagnetic compatibility (EMC)

Camera head

Art. no.	Designation
PV471	Pendulum camera head
–or–	
PV472	Zoom camera head
–or–	
PV473	Fixed focal length camera head

3.2 Components required for operation

In addition to the camera system, the following components are required for an endoscopic intervention:

- Power cord
- Endoscope
- Light source
- Light cable
- Monitor

3.3 Intended use

The Full HD CMOS camera system PV470/PV471/PV472/PV473 is used in medical endoscopy, together with an endoscope, a light source and a monitor.

It serves to visualize the endoscopic image for display on a color monitor.

Indication

The product is indicated for supporting diagnosis and treatment in endoscopic imaging.

Contraindication

Due to the universal suitability of the product for endoscopy, no specific contraindications for its use can be given. The use of the product is contraindicated if the endoscopic procedure are contraindicated or if its use could pose a risk to the patient.

3.4 Operating principle

The Endoscope visualizes the operation site and the endoscopic lens focuses this image onto the CMOS sensor, where it is converted into electrical signals that can be further amplified. The image brightness is automatically adjusted to the light intensity (autosshutter).

Through camera cable **26**, the image signals are transferred to camera control unit, where they are further amplified for display on the monitor.

The video signal processed in this way can then be picked up at the video output sockets for connecting monitors, video printers or any other digital archiving systems.

For areas that are particularly difficult to illuminate, the camera system offers a dynamic electronic brightness enhancement function (GAIN), which automatically adapts to different light conditions, and a choice of mode settings for different situations during application.

The video camera system is also equipped with a manual white balance system for optimal color rendition.

Connections on the front panel of the camera control unit

Control elements / indicators	Description
"Power off" switch 1 "Power on" switch 2	To switch the device on and off. LED display 3 above the "Power on" switch 2 lights up when the device is switched on.
Pushbutton white balance 4	For precise white-balancing, ensuring realistic rendition of all other colors
Camera socket 5	For connecting the camera plug at the end of the camera cable

Connections on the rear panel of the camera control unit

Connector	Description
Interface for MAQUET TEGRIS system 15	Interface for combination with the MAQUET TEGRIS system
USB connector for service purposes 14	Is only used by the manufacturer for service purposes
Remote connection 13	Connection for remote operation of peripheral devices (e.g. video printers) via cables with a 3.5 mm mono-plug
S-Video output 12	S-Video signal output
Communication interface 11	Connector to connect further communication interface devices
Equipotentialization connector 10	Connector for equipotentialization cables acc. to DIN 42801
HD-SDI outputs 9	HD SDI video signal output (2 connectors)
DVI outputs 8	DVI video signal output (2 connectors)
Fuse holder 7	For standard fuse 1.6 A (T)
Mains connector 6	Connection for power cord IEC 320, see accessories / spare parts

Control element on the camera head with endocoupler

Connector	Description
Safety lock 16	Prevents inadvertent opening of endocoupler
Locking mechanism 17	Spring-loaded locking for the endoscope ocular
Focusing ring 18	Turnable focusing ring
Zoom ring 19	Turnable zoom ring (only PV472)
Button A (up) 20	Function key
Button B (right) 21	Function key
Button D (left) 22	Function key
Button C (down) 23	Function key
Cap 24	Cap to prevent moisture from entering the plug
Connector plug 25	Connector plug to the camera control unit
Camera cable 26	Connection cable between the camera head and camera control unit
Endocoupler lens	Coupler for endoscope ocular

4. Preparation and setup

Non-compliance with the following instructions will preclude all responsibility and liability in this respect on the part of Aesculap.

DANGER

Danger of explosion from a device that is improperly set up!

- ▶ Make certain that the power plug is connected to the power supply outside of any areas where there is a risk of explosion.
- ▶ Have ready replacement products/camera system.

- ▶ When setting up and operating the product, adhere to
 - national regulations for installation and operation,
 - national regulations on fire and explosion protection.

Note

For the safety of patients and users it is essential that the mains power cord and, especially, the protective earth connection are intact. In many cases defective or missing protective earth connections are not registered immediately.

Note

The potential equalization plugs of the device must be connected to the potential equalization band (see IEC/DIN 60601-1-1 or corresponding national standards).

Note

The connection of electrical devices to a multiple socket creates a system and can lead to a reduced level of safety.

Note

If the equipment cart is set up outside of the sterile area, measure the cable length to give enough room for movement so that the sterility is not compromised.

- ▶ Connect the device via the potential equalization terminal at the rear panel of the device to the potential equalization system of the room used for medical purposes.
- ▶ Ensure that all devices operated in the vicinity meet their relevant EMC requirements.

Note

The potential equalization lead can be ordered from the manufacturer as art. no. GK535 (4 m length) or TA008205 (0.8 m length).

4.1 Preparation

Non-compliance with the following instructions will preclude all responsibility and liability in this respect on the part of Aesculap.

First use

WARNING

Risk of injury and/or product malfunction due to incorrect operation of the electromedical system!

- ▶ Do not use the device in areas where there is a risk of explosion or where there is a high level of oxygen.
- ▶ Adhere to the instructions for use of any medical device.
- ▶ Do not use the device in a magnetic resonance environment.

There are three different connection variants, depending on system requirements.

- ▶ If a parallel video image from the documentation device is not required, always connect the components as shown under variant 1, see fold-out page.
- ▶ If a parallel video image from the documentation device is required on a monitor, connect the components as shown under variant 2 or 3, see fold-out page.
- ▶ Always connect the video output (Video out) of one device to the video input (Video in) of the next device.

1. variant: Connecting the camera, monitor and one documentation device to the camera control unit (no parallel video images)

- ▶ Select signal or type of video cable according to the common signal standard of the camera and monitor, always choosing the highest-quality signal option (video signals in descending order of quality: DVI, HD-SDI, S-Video).
- ▶ If necessary, connect remote cable.

2. variant: Connecting the camera and monitor with one documentation device daisy-chained between them (for parallel video images)

- ▶ If video images from the documentation device are intended to be overlaid on the monitor display, connect the device (e.g. digital documentation system) between the camera and the monitor.
- ▶ Select signal or type of video cable according to the common signal standard of the camera and monitor, always choosing the highest-quality signal option (video signals in descending order of quality: DVI, HD-SDI, S-Video).
- ▶ If necessary, connect remote cable.

Note

The optional direct video connection between the camera and the monitor is recommended as an additional signal line if the video signal can be looped through the documentation device only with the device switched on (see instructions for use provided by the manufacturer of the respective device). This ensures that, in case of a documentation device failure, the monitor can be immediately switched to the direct video signal.

3. variant: Connecting the camera and monitor with several documentation devices daisy-chained between them (for parallel video images on the monitor)

- ▶ If video images from the documentation devices are intended to be overlaid on the monitor display, connect Device 1 and Device 2 (e.g. digital documentation system) between the camera and Monitor 1.
- ▶ Select signal or type of video cable according to the common signal standard of the camera and monitor, always choosing the highest-quality signal option (video signals in descending order of quality: DVI, HD-SDI, S-Video).
- ▶ Connect Monitor 2 to Monitor 1 if the video images from the documentation devices are to be displayed on both monitors.
- ▶ If necessary, connect remote cable.

Note

The optional direct video connection between the camera and the monitor is recommended as an additional signal line if the video signal can be looped through device 1 and device 2 only with the devices switched on (see instructions for use of the devices). This ensures that, in case of a failure of Device 1 or Device 2, the monitor can be immediately switched to the direct video signal.

The alternative direct video connection between Monitor 2 and the camera is recommended if Monitor 2 is only to display the live image.

Stacking of units

- ▶ Make certain the system is set up on a sufficiently stable support (e.g. a table, ceiling support, equipment cart, etc.).
- ▶ Do not exceed a stack height of 300 mm.

Moving stacks of units

- ▶ Always lift stacks from underneath the bottommost unit.

Presetting the correct voltage

The mains voltage must correspond to the voltage indicated on the type plate at the back of the unit.

5. Working with the camera control unit Full HD PV470 and the camera head PV471/PV472/PV473

5.1 System set-up

Connecting the accessories

DANGER

Risk of injury due to unapproved configuration using additional components!

- ▶ For all applied components, ensure that their classification matches that of the application component (Type CF) of the respective device.

DANGER

Risk of injury when using devices or components that are not of type CF!

- ▶ Do not touch the device and patient at the same time.
- ▶ Do not place the patient in direct contact with the device.

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

Also note that any equipment connected at the interfaces must demonstrably meet the respective IEC standards (e.g. IEC 60950 for data processing equipment, IEC/DIN EN 60601-1 for electromedical devices).

All configurations must comply with basic standard IEC/DIN EN 60601-1. Any individual connecting devices with one another is responsible for such configuration and must ensure compliance with basic standard IEC/DIN EN 60601-1 or applicable national standards.

- ▶ Please address your B. Braun/Aesculap partner or Aesculap Technical Service with any inquiries in this respect; for a contact address, see Technical Service.

Connect video cable

Connect video cables according to diagrams, see First use.

Communication interface connection (optional accessory)

With the communication interface connection, the LED light source OP940 can be switched on and to standby via the FULL HD CMOS camera system – see connection diagrams 1, 2 or 3. In order to use this remote control of the light source, the communication interface cable OP942 must be connected to the two devices.

- ▶ Connect the communication interface cable OP942 between the LED light source OP940 and the camera control unit 11.

Aesculap®

Full HD CMOS camera system PV470

Use the MAQUET TEGRIS interface

The MAQUET TEGRIS interface **15** is used to connect the device to the MAQUET TEGRIS system for central device control.

When the device is connected, all the main functions and settings of the camera system:

- can be controlled via the TEGRIS control and
- they are displayed on the TEGRIS screen.

When the MAQUET TEGRIS system is connected and active, all the functions and settings of the camera system can be set and displayed via the device and via the MAQUET TEGRIS system.

If both are being operated at the same time, entries on the device take precedence. Doctors and the nursing or service personnel must be thoroughly familiar with how to set and operate the MAQUET TEGRIS system, see TEGRIS documentation for information on settings, use and troubleshooting.

Connecting the power supply

DANGER

Risk of death by electric shock!

- ▶ **Connect the product only to a grounded power supply.**

- ▶ Insert the power cord on the rear panel of the camera control unit into the mains connector **6**.
- ▶ Plug in the mains plug at a socket of the building mains.
- ▶ Plug in the equipotentialization cable into the equipotentialization connector **10** of the camera control unit.
- ▶ Connect the plug of the equipotentialization cable to the building mains equipotentialization connector.
- ▶ Make certain that the rear of the device is accessible at all times.

5.2 Function checks

Camera control unit

- ▶ Check camera control unit for external damage and signs of impact or other mechanical impacts.

Power cord

- ▶ Check that the power cord is not broken, kinked or twisted.

Video cable

- ▶ Check that the video cable is not broken, kinked or twisted.

Camera head with endoscopic lens and camera cable

Note

The endoscopic lens is firmly attached to the camera head and therefore cannot be removed!

WARNING

Risk of infection due to soiled or contaminated products!

- ▶ **Disinfect camera head and camera cable before each use, or use with sterile cover.**

- ▶ Check that the front glass of endocoupler on the camera head is clean. If necessary, clean the glass with isopropyl alcohol (70 %).
- ▶ Make certain that adjusting ring for the focusing ring **18** turns freely without sticking.
- ▶ Check that camera cable **26** is not broken, kinked or twisted. Do not use the camera head in the event of damage.
- ▶ Check that pins in connection plug **25** are not bent and that there is no moisture in the plug.
- ▶ Ensure that the housing of the camera head is not damaged or broken. Do not use the camera head in the event of damage.
- ▶ Aim camera head to an object in the room, focus the image and check for good image quality.

Sterile cover

Note

Aesculap recommends using disposable sterile covers, (e.g. art. no. JG904). Use of the sterile covers is described in the relevant instructions for use.

- ▶ Check that the packaging of the sterile covers is unopened, undamaged and dry before using the sterile covers in the operating room.

5.3 Safe operation

WARNING

Risk of injury and/or malfunction!

- ▶ **Always carry out a function check prior to using the product.**

WARNING

Risk of injury when using the product beyond the field of view!

- ▶ **Apply the product only under visual control.**

Note

The camera control unit requires approx. 30 s to boot up.

Connecting the camera cable

- ▶ Withdraw cap **24** on the connection plug **25**.
- ▶ Connect the connection plug **25** to the camera socket **5** on the front panel of the camera control unit. When doing this take care that the red marking on the plug is aligned with the red marking on camera socket **5**.

Switching on the system devices

- ▶ Switch on the monitor and set it to the required channel (e.g. DVI, HD-SDI).
- ▶ Press the "Power on" switch **2** on the camera control unit. LED display **3** lights up.

Connecting the endoscope

PV471/PV472/PV473:

- ▶ Press both endocoupler latches **17** at the same time in order to spread apart the holding clamps.
- ▶ Insert the endoscope and release the endocoupler latches **17**.
- ▶ In order to avoid inadvertent opening of the endocoupler: Rotate safety lock **16** to the locked position (**LOCK**).

Connecting the light cable

- ▶ Connect the light cable to the endoscope.
- ▶ Switch on the light source at a low light intensity setting.

Adjusting the image focus and enlargement

- ▶ Adjust the intensity of the light source until sufficient illumination is achieved.
- ▶ In order to sharpen the image, turn the focusing ring **18** on the camera head.
- ▶ In order to set the desired image size, turn the black zoom ring **19** on the camera head (only PV472).

Pendulum function (only PV471)

This camera head is self-aligning. The image always remains upright as the camera head remains aligned perpendicularly to the endoscope axis due to the force of gravity.

- ▶ In order to inhibit rotatability around the endoscope axis: Actuate pendulum-locking lever **27**.

Performing white-balancing

Note

The system will request a white balance after a Factory Reset.

- ▶ Point endoscope at a white object (e.g. a sterile swab).
- ▶ Press the white balance button **4** on the camera control unit or the white balance button on the camera head (depending on the button programming), and hold the endoscope pointing at the white object until the camera has completed white balancing.

The image will be displayed in pure white. This procedure ensures realistic rendition of all colors.

The following messages will be displayed on the monitor during white-adjustment:

During white-balancing	After white-balancing	White balance error
White balance in progress	White balance OK	White balance failed

- ▶ If white-adjustment fails, restart white-adjustment, see Performing white balance.

Buttons on the camera head

Buttons **20–23** on the camera head can each be programmed with up to two functions. The first function is triggered by pressing the button briefly. The second function is triggered by pressing the button and holding (for at least 1 second). The buttons are freely allocable apart from the "Menu" function, see Assigning the buttons on the camera head. The menu can always be called up by pressing and holding button **20**.

Default setting

Pressing the button	Button A 20	Button B 21	Button C 23	Button D 22
Short	free	free	free	Photo release
Long	Menu (fixed)	White balance	LED light source on/standby	Video start/stop

Aesculap®

Full HD CMOS camera system PV470

Calling up the menu

The camera has a menu that appears as text on the screen ("on-screen display" or OSD). The menu appears on the screen by pressing and holding (>1 s) button A **20** (up).

Navigation in the menu

By actuating buttons A **20** (up) and C **23** (down) it is possible to navigate within the menu. As soon as a menu point is marked, this can be activated by actuating button B **21** (right) and then modified by pressing button B **21** (right) or button D **22** (left).

Overview of the main menu

Menu	
Profile	Preinstalled camera settings, optimized to the use. The choices are: Laparoscopy Arthroscopy Hysteroscopy Neuro 1 Neuro 2 Fiberscope Users 1, 2, 3 (programmable by the user)
Brightness	Setting brightness manually
Saturation	Setting color saturation manually
Digital zoom	Setting digital zoom up to factor 2.5
Sharpness	Setting edge enhancement manually
Contrast	Setting contrast manually
System setup	Opens a submenu with other settings
Save and exit	Save setting and exit menu
Exit	Changes will not be saved, exit menu

Overview of the "profile" menu

Profile	
Laparoscopy	Camera settings optimized for this use
Arthroscopy	Camera settings optimized for this use
Hysteroscopy	Camera settings optimized for this use
Neuro 1	Camera settings optimized for this use
Neuro 2	Camera settings optimized for this use
Fiberscope	Camera settings optimized for this use
User 1	Camera settings saved by the user; disk space 1
User 2	Camera settings saved by the user; disk space 2
User 3	Camera settings saved by the user; disk space 3
User settings	Opens submenu. The "User 1/2/3" settings can be defined there.
Return	Returns the user to the main menu

Note

Once a profile has been chosen, the menu will close automatically and the user will be back in the operation mode (live image).

Overview of the "system setup" menu

System setup	
Language	Setting the menu language
Button layout	Opens a submenu in which the camera head buttons can be assigned
Image format	Sets the image refresh rate Selectable options: 50 Hz/60 Hz
Light source	Switch the light source on or to standby (only in conjunction with Aesculap LED light source OP940) Selectable options: On/standby
System information	Opens a submenu with other settings
Return to default settings	Restore all settings to the original values
Return	Takes the user back to the previous menu

Overview of the "user settings" menu

User settings	
Ref. profile	Here, a profile can be selected as a basis for the other settings. The choices are: Laparoscopy Arthroscopy Hysteroscopy Neuro 1 Neuro 2 Fiberscope
Brightness	Setting brightness manually
Saturation	Setting color saturation manually
Digital zoom	Setting digital zoom manually
Sharpness	Setting edge enhancement manually
Contrast	Setting contrast manually
Save as	Here, disk space can be established for the settings. The choices are: User 1 User 2 User 3
Return and save	Saves the settings and takes the user back to the previous menu
Return	Takes the user back to the previous menu without saving the changes

Aesculap®

Full HD CMOS camera system PV470

Overview of the "system information" menu

System information	
Software status	
Camera head	Camera head software status
Control unit	Control unit software status
Serial number	
Camera head	Camera head serial number
Control unit	Control unit serial number
Light source	Light source serial number
Operating hours	
Camera head	Camera head operating hours
Control unit	Control unit operating hours
Return	Takes the user back to the previous menu

Overview of the "button layout" menu

Note

"Top long" is set as "Menu" and cannot be changed. All other buttons can be assigned freely with different functions.

Button layout (default settings)	
Top short	-
Up long	Menu
Bottom short	-
Down long	LED on/standby
Left short	Photo
Left long	Video
Right short	-
Right long	White balance
Return and save	Saves the settings and takes the user back to the previous menu
Return	Takes the user back to the previous menu without saving the changes

Apply user-defined camera settings and save

Here, up to three different camera settings can be set and saved according to the user's preference. These settings can then be called up again later by selecting "User 1", "User 2" or "User 3" in the main menu.

- ▶ Start menu (button A 20 (up), press and hold).
- ▶ Select "User settings".
- ▶ Press button B 21 (right) briefly. A submenu appears.
- ▶ Under "Ref. profile", select the desired entry, e.g. "Laparoscopy".
- ▶ Now apply the desired image settings.
- ▶ Under "Save as", select "User 1", "User 2", or "User 3".
- ▶ Select "Return and save" or "Return" (changes are not saved).
- ▶ These user-defined settings become active immediately if the profile "User 1/2/3" is selected as a profile in the main menu.

Assigning the buttons on the camera head

By pressing and holding the button, button A 20 (up) is assigned to the "Menu" function and cannot be changed. All other buttons can be assigned freely with different functions. In order to allocate functions to the buttons, proceed as follows:

- ▶ Access the start menu (button A 20 (top), press and hold).
- ▶ Select "System Settings".
- ▶ Press button B 21 (right) briefly. A submenu appears.
- ▶ Select "Button layout".
- ▶ Press button B 21 (right) briefly.
- ▶ Select the desired button.
- ▶ A function can be selected by pressing briefly on button B 21 (right) or button D 22 (left).
- ▶ Define all buttons as required with the same procedure.
- ▶ If the button is not intended to have any function, select "empty" as the function.
- ▶ Finally, select "Return and save" or just "Return" (changes are not saved).

Note

For safety reasons, not all the camera functions are available for key assignment.

The followings table shows the preset key assignment and which combinations are possible:

Function/push-button	Up short	Top long	Down short	Bottom long	Left short	Left long	Right short	Right long
Menu		■ (default)						
White balance				■		■		■ (default)
Light source on/off				■ (default)		■		■
Photo	■		■	■	■ (default)	■	■	■
Video	■		■	■	■	■ (default)	■	■
Brightness +	■		■	■	■	■	■	■
Brightness -	■		■	■	■	■	■	■
Saturation +	■		■	■	■	■	■	■
Saturation -	■		■	■	■	■	■	■
Digital zoom +	■		■	■	■	■	■	■
Digital zoom -	■		■	■	■	■	■	■
Sharpness +	■		■	■	■	■	■	■
Sharpness -	■		■	■	■	■	■	■
Contrast +	■		■	■	■	■	■	■
Contrast -	■		■	■	■	■	■	■
Empty	■ (default and fixed)		■ (default)	■	■	■	■ (default)	■

Return to factory settings

In order to return to factory settings:

- ▶ Start menu (button A 20 (up), press and hold).
- ▶ Mark "System setting".
- ▶ Actuate button B 21 (right). A submenu appears.
- ▶ Select "Return to default settings".

End application

- ▶ Press the "Power off" switch 1 on the camera control unit.
- ▶ Disconnect light guide cable.
- ▶ Disconnect endoscope.
- ▶ To completely disconnect the device from the power supply, remove the plug from the mains connector 6 on the rear of the device.

6. Validated reprocessing procedure

6.1 General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob-disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

6.2 Single-use products

WARNING

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

► **Do not reprocess the product!**

Designation	Art. no.
Disposable sterile cover	JG904

6.3 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

6.4 Preparations at the place of use

- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

CAUTION

Damage to the product due to ultrasound or inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Only damp-wipe the camera control unit (disinfection by wiping).
- ▶ Do not immerse the camera control unit in any fluid.
- ▶ Do not clean the camera control unit with endo-coupler in an ultrasound-bath.
- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for plastic material and high-grade steel,
 - do not attack softeners (e.g. in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 80 °C.

6.5 Cleaning/disinfection

Product-specific safety instructions for the reprocessing procedure

DANGER

Risk of electric shock and fire hazard!

- ▶ Unplug the device before cleaning.
- ▶ Do not use flammable or explosive cleaning or disinfecting solutions.
- ▶ Ensure that no fluids will penetrate the product.

DANGER

Damage to, or destruction of the product caused by mechanical cleaning/disinfection!

- ▶ Only clean and disinfect the product manually.

Validated cleaning and disinfection procedure

Validated procedure	Special features	Reference
Wipe disinfection for electrical devices without sterilization <ul style="list-style-type: none"> ■ Camera head ■ Camera control unit 	<ul style="list-style-type: none"> ■ Keep working ends open for cleaning. ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. 	Chapter Wipe disinfection for electrical devices without sterilization
Manual cleaning with immersion disinfection <ul style="list-style-type: none"> ■ Camera head 	<ul style="list-style-type: none"> ■ Firmly slide open cap 24 on the connection plug 25. ■ Keep working ends open for cleaning. ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or compressed air for medical purposes 	Chapter Manual cleaning/disinfecting and sub-section: <ul style="list-style-type: none"> ■ Chapter Manual cleaning with immersion disinfection

Aesculap®

Full HD CMOS camera system PV470

6.6 Wipe disinfection for electrical devices without sterilization

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Wipe disinfection	RT	≥1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol

RT: Room temperature

Phase I

- ▶ Remove any visible residues with a disposable disinfectant wipe.
- ▶ Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- ▶ Observe the specified application time (1 min minimum).

6.7 Manual cleaning/disinfecting

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	34–45/ 95–113	5	0,8	D–W	Enzymatic Cleaner Cidezyme/Enzol
II	Intermediate rinse	RT (cold)	3x1	-	D–W	-
III	Disinfection	20–25/ 68–77	12	0,55	D–W	Cidex OPA (process solution 0.55 % ortho-phthalaldehyde)
IV	Final rinse	RT (cold)	3x2	-	FD–W sterile	-
V	Drying	RT	-	-	-	-

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

Phase I

- ▶ Fully immerse the product for at least 5 min in the cleaning solution. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product immersed in the solution, using a soft cloth or, if necessary, a soft cleaning brush, until all visible residues have been removed from the surfaces.
- ▶ If applicable brush through all surfaces that are not accessible to visual inspection for at least 1 minute.
- ▶ During cleaning, mobilize non-rigid components, such as set screws, joints, slides, etc., 3 times in each direction to the positive stop.
- ▶ After cleaning, use a suitable disposable syringe to thoroughly rinse these parts of the product that are difficult to access at least 3 times.

Phase II

- ▶ Rinse/flush the product (all accessible surfaces) thoroughly for at least 1 minute under running water.
- ▶ During rinsing, mobilize non-rigid components, such as set screws, joints, slides, etc., 3 times in each direction to the positive stop.
- ▶ Perform these steps at least 3 times.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product for at least 12 minutes in the disinfectant solution.
- ▶ During rinsing, mobilize non-rigid components, such as set screws, joints, slides, etc., 3 times in each direction to the positive stop.
- ▶ Rinse lumens at least 3 times at the beginning of the exposure time, using a disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Following disinfection, rinse/flush the product completely (all accessible surfaces) for at least 2 minutes.
- ▶ During rinsing, mobilize non-rigid components, such as set screws, joints, slides, etc., 3 times in each direction to the positive stop.
- ▶ Rinse lumens at least 3 times at the beginning of the exposure time, using a disposable syringe.
- ▶ Perform these steps at least 3 times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product with suitable tools.

6.8 Sterilization

CAUTION

CAUTION

Damage to the camera control unit by sterilization!

- ▶ Do not sterilize the camera control unit under any circumstances.
- ▶ Clean the camera control unit by surface disinfection only.

Damage to the camera head and the adapter caused by inappropriate sterilization processes!

- ▶ Never sterilize the camera head and the adapter with steam or in an autoclave.
- ▶ Do not carry out chemical sterilization.

Working with disposable sterile covers

- ▶ Use suitable Aesculap disposable sterile covers. Follow the instructions for use of the sterile covers.

Approval of compatibility of materials exists for:

ETO-sterilization

6.9 Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Set aside the product if it is damaged.

6.10 Packaging

- ▶ Always follow the instructions for use of the respective packaging and storage devices (e.g. instructions for use TA009721 for the Aesculap Eccos® storage system).
- ▶ Insert the product in its proper position in the Eccos® holder, or put it on a tray in such a way that the product is protected against damage. Ensure that any cutting edges present are protected.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

6.11 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- ▶ Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Aesculap®

Full HD CMOS camera system PV470

7. Maintenance

To ensure reliable operation, Aesculap recommends maintaining the product at least once a year. For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

Note

If a component in the visualization system is replaced (e.g. service or upgrade), the system must be started up again.

Note

If this is not possible for urgent reasons, prepare product as far as is possible and label it correspondingly.

Note

The specialized repair company may refuse to repair soiled or contaminated products for safety reasons. The manufacturer reserves the right to return contaminated products to the sender.

WARNING

Risk of infection due to soiled or contaminated products!

- ▶ Thoroughly clean, disinfect and sterilize the product before sending. Otherwise prepare the product in as far as is possible, and label it correspondingly.
- ▶ Select suitable and safe packaging (ideally the original packaging). Package the product in such a way that the packaging does not become contaminated.

8. Troubleshooting list

Malfunction	Cause	Remedy
No image on the monitor	Device not powered	Verify that all units are connected to mains power and switched on
	Fuse is broken	Replace the fuse in the camera
	Video cable defective	Replace video cable
	Video cable not connected	Connect video cable
Image too dark	Camera cable not connected or defective	Ensure that the camera head is connected to the camera control unit and the camera plug is dry
	Light guide cable defective	Connect a new light guide cable. Have repaired by manufacturer
	Light source dimmed down	Turn up the light source intensity
	Camera brightness regulated down	Regulate camera brightness up
	Endoscope optical system defective	If the image is too dark even without the camera: Use a different endoscope and have the defective endoscope repaired by the manufacturer
Image fuzzy or poorly defined	Object out of focus	Turn focusing ring 18 on the endoscopic lens and adjust the focus until the image is sharp
	Debris on the endoscope tip	Clean tip with a sterile cloth and isopropyl alcohol (70 %)
	Moisture on the camera cable plug	Dry plug with gauze before plugging it in at the camera socket

Malfunction	Cause	Remedy
Image fogged	-	Check if the front and rear windows of the endoscope and the front window of the endoscopic lens are clean; if necessary, clean with isopropyl alcohol (70 %)
Noise or excessive grain in the picture	-	Increase intensity of the light source
Image too bright or glary	-	Reduce the intensity of the light source
Color bar diagram on the monitor	Camera cable not connected to the camera control unit	Connect camera cable at the camera control unit
Unrealistic color representation on the monitor	White balance failed	Start white balance-adjustment
	Monitor color setting misadjusted	Reset to factory setting
	Video cable defective	Replace video cable
The LED light source is no longer controlled via the camera head	Problems with the communication interface between the light source and camera control unit	Connect the communication interface cable correctly Remove the communication interface cable and operate the light source manually
Camera system is not switching the light source on/off properly when the communication interface cable is connected	Problems with the communication interface between the light source and camera control unit	Remove the communication interface cable and manually operate the light source
MAQUET TEGRIS system is not operating the camera system or light source connected to the camera system via the communication interface properly.	Problems with the communication interface between the MAQUET TEGRIS system and camera control unit	Check the cable connections Operate the camera system and, if applicable, the light source manually, and if necessary remove the TEGRIS interface cable and communication interface cable.
"System error" appears on the monitor OSD	Internal product error	Restart camera
"System error" appears on the monitor OSD even after restart	Internal product error	Send the camera to the Service department, see Technical Service

8.1 Fuse replacement

Risk of death by electric shock!

► **Unplug the device before changing the fuses!**

- Use a small screwdriver to release the clip on the fuse holder.
- Pull out the fuse holder.
- Replace both fuses.
- Reinsert the fuse holder in such a way that it audibly snaps into place.

Note

If the fuses burn out frequently, the device is faulty and should be repaired, see Technical Service.

9. Technical Service

WARNING

Risk of injury and/or malfunction!

► **Do not modify the product.**

► For service and repairs, please contact your national B. Braun/Aesculap agency.

CAUTION

Malfunction of the product due to damage in transport!

► **Select suitable and safe packaging (ideally the original packaging).**

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 14-939

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042 USA

Aesculap Repair Hotline

Phone: +1 (800) 214 -3392

Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

10. Accessories/Spare parts

Art. no.	Designation
PV470	Camera control unit
PV471	Camera head with 90° endoscopic lens and camera cable
PV472	Camera head with zoom endoscopic lens and camera cable
PV473	Camera head with fixed focal length Endoscopic lens and camera cable
JG904	Disposable sterile covers, pack of 25
A085442	Fuse T1.6AH 250 V
TA013748	Instructions for use
TA022130	Notes on electromagnetic compatibility (EMC)
PV436	DVI-cable
OP942	Communication interface cable

11. Technical data

Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
PV470	Camera control unit	I
PV471	Camera head with fixed focal length, 90° endoscopic lens	I
PV472	Camera head with zoom endoscopic lens	I
PV473	Camera head with fixed focal length endoscopic lens	I







Camera control unit PV470

Mains supply	100–240 V~, 0.63–0.35 A, 50/60 Hz
Power consumption	max. 75 VA
Protection class (acc. to IEC/DIN EN 60601-1)	I
Device fuse	T1,6AH 250 V
Weight	5.2 kg
Dimensions (W x H x D)	305 mm x 75 mm x 305 mm
TV system	PAL (50 Hz), NTSC (60 Hz) for standard definition
Shutter range (automatic)	1/50 – 1/10,000 seconds
Video outputs	DVI 1,080p (x2), HD-SDI (x2), S video
Interfaces	Communication interface, remote output, RS232
EMC	IEC/DIN EN 60601-1-2 CISPR11 Class A
Conforming to standard	IEC/DIN EN 60601-1

Camera head PV471/PV472/PV473

Image recording	CMOS sensor 1/3"
Resolution	1,920 x 1,080 pixels
Focal length of endoscopic lens	PV471: 14.25 mm PV472: 14.25 mm–28 mm PV473: 14.25 mm
Camera head dimensions incl. endoscopic lens: (W x H x D)	PV471: approx. 52 mm x 63 mm x 130 mm PV472: approx. 52 mm x 52 mm x 130 mm PV473: approx. 52 mm x 52 mm x 110 mm
Camera head weight incl. endoscopic lens:	PV471: 610 g without cable: approx. 200 g PV472: 730 g without cable: approx. 335 g PV473: 670 g without cable: approx. 300 g
Camera cable length	4 m
EMC	IEC/DIN EN 60601-1-2
Conforming to standard	IEC/DIN EN 60601-1
IP classification	IPX7
Applied part in accordance with IEC/DIN EN 60601-1	Type CF


11.1 Ambient conditions

	Operation	Storage and transport
Temperature	 10 °C – 40 °C	 -10 °C – 50 °C
Relative humidity	 30 % – 90 %	 10 % – 90 %
Atmospheric pressure	 770 hPa – 1 060 hPa	 500 hPa – 1 060 hPa

12. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

Adhere to national regulations when disposing of or recycling the product, its components and its packaging! The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.) Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

13. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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