GAV 2.0[®] GAV 2.0[®] LP

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Ⓔ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
Abb. 10: GAV 2.0 LP (U-Form)

TABLE OF CONTENTS

INDICATION	12
TECHNICAL DESCRIPTION	12
FUNCTION OF THE VALVE	12
SELECTION OF THE APPROPRIATE PRESSURE LEVEL	13
PRESSURE LEVELS RECOGNITION IN X-RAY IMAGES	13
POSSIBLE SHUNT COMPONENTS	14
TUBE SYSTEMS	14
IMPLANTATION	14
VALVE TEST	15
PRESSURE-FLOW CHARACTERISTICS	16
PRECAUTIONS AND CONTRAINDICATIONS	17
FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	17
ADVERSE REACTIONS AND INTERACTIONS	17
STERILISATION	17
REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)	17
MEDICAL PRODUCTS CONSULTANTS	18
VARIANTS	18

- 1 Micro spiral spring
- 2 Tantalum ball
- 3 Sapphire ball
- 4 Coding


Fig. 1: GAV 2.0 in cross section

INDICATION

The GAV 2.0 is used for cerebrospinal fluid (CSF) drainage in the treatment of hydrocephalus

TECHNICAL DESCRIPTION

GAV 2.0 is a valve made from titanium. It consists of a ball-cone unit and a gravitational unit. Thus a physiological intraventricular pressure (IVP) can be reached in any body position.

In the proximal part of the GAV 2.0, micro spiral spring (1) enables the opening pressure of the ball-cone unit. The gravitational unit in the distal part of the valve consists of a tantalum ball (2) which defines the opening pressure of this unit, and a sapphire ball (3) which enables its precise closure. Coding (4) permits the identification of pressure levels in an X-ray.

FUNCTION OF THE VALVE

The GAV 2.0 is a position-dependent valve.

Horizontal position

In the horizontal position, the gravitational unit is always open and does not present any resistance. For that reason, the GAV 2.0 opening pressure in this position is only defined by the ball-cone unit. Fig. 2a shows the ball-cone unit in closed status. If the patient's intraventricular pressure (IVP) exceeds the opening pressure of the micro spiral spring, the ball from the cone, thus opening a gap for drainage (Fig. 2b).


Fig. 2: GAV 2.0 in horizontal position
a) closed b) open

Vertical position

When the patient moves into an upright position, the gravitational unit is activated and the opening pressure of the GAV 2.0 is significantly increased (Fig. 3a). In addition to the opening pressure of the ball-cone unit, the weight of the tantalum ball (opening pressure of the gravitational unit) has to be exceeded. Drainage is only possible once the sum of intraventricular pressure (IVP) and hydrostatic pressure is greater than the opening pressure of both units (Fig. 3b).


Fig. 3: GAV 2.0 in vertical position
a) closed b) open

SELECTION OF THE APPROPRIATE PRESSURE LEVEL

Horizontal position

Depending on the clinical picture and age of the patient, the opening pressure for this position can be selected between pressure levels 5 and 10 cmH₂O.

Vertical position

Patient size, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure for this position (see pressure level recommendations at <https://www.miethke.com/en/products/downloads/>).

PRESSURE LEVELS RECOGNITION IN X-RAY IMAGES

The selected pressure levels can be detected radiographically after implantation.

Valve opening pressure horizontal vertical		Coding
5 cmH ₂ O	20 cmH ₂ O	
5 cmH ₂ O	25 cmH ₂ O	
5 cmH ₂ O	30 cmH ₂ O	
5 cmH ₂ O	35 cmH ₂ O	
10 cmH ₂ O	25 cmH ₂ O	
10 cmH ₂ O	30 cmH ₂ O	

Fig. 4: GAV 2.0 pressure levels encodings on X-ray

POSSIBLE SHUNT COMPONENTS

The GAV 2.0 can be ordered as a shunt system in a range of configurations. The configurations can be combined with the accessories presented in brief below. In each case, versions for paediatric hydrocephalus and for normal pressure hydrocephalus (NPH) in adults are available.

Reservoirs

The use of a reservoir in combination with shunt systems provides options for the withdrawal of cerebrospinal fluid, administration of drugs and pressure control.

Thanks to the one-way flux (flow) system of the *SPRUNG RESERVOIR* and the *CONTROL RESERVOIR*, cerebrospinal fluid can be pumped towards the valve, thus making it possible to check the distal part of the drainage system as well as (proximal) ventricular catheter. During the pump action, access to the ventricular catheter is closed. The use of reservoirs does not increase the opening pressure of the shunt system. A puncture should be performed as perpendicular as possible to the reservoir surface with a maximum cannula diameter of 0.9 mm. 30 punctures are possible without any restrictions.

Warning notice: Frequent pumping can result in excessive drainage and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.

Burrhole deflector

Because of the tight fit on the ventricular catheter, the burrhole deflector makes it possible to choose the length of catheter penetrating into the skull prior to implantation. The ventricular catheter is deflected at a right angle in the burrhole (see chapter "Implantation").

TUBE SYSTEMS

The GAV 2.0 can be ordered as an individual valve unit or as a shunt system with integrated catheters (interior diameter 1.2 mm, exterior diameter 2.5 mm). The supplied catheters do not fundamentally change the pressure-flow characteristics. If catheters by other manufacturers are used, a tight fit must be ensured. In any

case, catheters have to be carefully fixed with a ligature to the valve's titanium connectors.

IMPLANTATION

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The required skin incision should be made in form of a lobule pedicled towards the draining catheter. If a burrhole deflector is used, the skin incision should not be located right above the reservoir. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burrhole.

The GAV 2.0 is available in a range of different configurations: If a burrhole reservoir is used, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the reservoir is connected, whereby the connection is secured with a ligature. A shunt system with prechamber comes with a burrhole deflector. The deflector is used for adjusting the length of catheter to be implanted and for its positioning inside the ventricle. The ventricular catheter is deflected, connected to the prechamber, and the prechamber is put into place. The position of the ventricular catheter should be inspected after the procedure by imaging (such as CT or MRI).

Positioning the valve

The operating principle of the GAV 2.0 is posture-dependent. For that reason, care must be taken to implant the valve parallel to the body axis. For VP drainage, a suitable position for implantation is behind the ear. After skin incision and tunnelling under the skin, the catheter is pushed forward from the burrhole to the intended valve implantation site, shortened where necessary and secured with a ligature. For LP drainage, the valve is placed in a subcutaneous skin pocket in the abdominal or back region. The valve should not be located directly under the skin incision. The valve is marked with an arrow indicating the distal direction of flow.

Warning notice: The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. For example, it can be applied para-umbilically in a horizontal direction or transrectally at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. The recommendation is to pull the peritoneal catheter using a subcutaneous tunneling tool from the valve to the intended position, if necessary with the aid of an auxiliary incision. The peritoneal catheter which is usually securely attached to the valve has an open distal end and no wall slits. Following the exposure of the peritoneum or with the aid of a trocar, the peritoneal catheter (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

Reimplantation

Products that have previously implanted must not subsequently be reimplanted into the same or another patient.

VALVE TEST

Preoperative valve test

The most careful way of filling the valve is by aspiration through a sterile single-use syringe attached to the distal end of the catheter. The distal end of the valve is connected and immersed in a sterile physiological saline solution. The valve is patent if saline solution can be extracted (Fig. 5).

Warning notice: Contamination in the solution used for testing can impair the product's performance.


Fig. 5: Patency test

Warning notice: Pressurisation by the single-use syringe should be avoided both at the proximal and the distal end (fig. 6).


Fig. 6: Avoidance of pressurisation

Postoperative valve test

The GAV 2.0 has been constructed as a reliably functioning unit without pump or test function. The valve test can be performed by flushing, pressure measurement or pumping.

PRESSURE-FLOW CHARACTERISTICS

The following diagrams show the pressure-flow characteristics for the available pressure settings of the GAV 2.0 (Fig. 7). The total opening pressure refers to a reference flow of 5 ml/h. For flow rates of 20 ml/h, stated pressures are approx. 1-2 cmH₂O higher.


Fig. 7: Pressure-flow characteristics for the available pressure settings of the GAV 2.0

PRECAUTIONS AND CONTRAINDICATIONS

Patients must be carefully monitored after implantation. Reddening of skin or tightness in the area of the drained tissue may be indications of infections at the shunt system. Symptoms such as headache, dizziness, confusion or vomiting often occur in conjunction with shunt dysfunction. These symptoms and a leakage within the shunt system require the immediate replacement of the affected shunt component or the entire shunt system.

The implantation of medical devices is contraindicated if the patient has an infection or suspected infection (e.g. meningitis, ventriculitis, peritonitis, bacteraemia, septicaemia) in the region affected by the implantation.

FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. However, no guarantee can be given that these medical devices may not require replacement for medical or technical reasons. These medical devices are able to resist positive and negative pressures up to 200 cmH₂O during and after implantation. These medical devices have to be stored in a clean and dry environment at all times.

Nuclear magnetic resonance examinations up to a field strength of 3 Tesla or computed tomography examinations can be performed without risk or impairment to the valve function. The valve is MR Conditional. Supplied catheters are MR Safe. Reservoirs, deflectors and connectors are MR Conditional.

ADVERSE REACTIONS AND INTERACTIONS

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases noise development. Violent shocks for the outside (accident, fall) may put the integrity of the shunt system at risk.

The GAV 2.0 must not be used in conjunction with hydrostatic valves as this may result in increased ventricular pressure outside of the physiological range. In case of doubt, please contact the medical device consultants at Christoph Miethke GmbH & Co. KG.

STERILISATION

The products are sterilised with steam under strictly controlled conditions. The double wrapping in sterile bags ensures sterility for a five-year period. The expiry date is printed on the wrapping of each individual product. If the packaging is damaged, the product must not be used in any circumstances. No guarantee can be given for the functional safety and reliability of resterilised products.

REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the implanted valve should therefore be recorded in the patient's medical records and patient data card to ensure complete traceability.

Translations of these instructions for use into additional languages can be found on our website (<https://www.miethke.com/en/products/downloads/>).

MEDICAL PRODUCTS CONSULTANTS

In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical products consultants as contacts for all product-related questions:

Dipl.-Ing. Christoph Miethke
 Dipl.-Ing. Roland Schulz
 Michaela Funk-Neubarth
 Dipl.-Ing. Thoralf Knitter
 Dr. Andreas Bunge
 Jan Mgel

Contact details can be found on the reverse of these instructions for use.

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany