# Lyoplast®

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# Lyoplast®

## Description

Lyoplast® is an implant of pure collagen that has been obtained from bovine pericardium. The special preparation process ensures that Lyoplast® is purified from noncollagenous components, for example, enzymes, lipids and noncollagenous proteins.

The gentle freeze-drying guarantees that Lyoplast® retains its loose fibrous structure, which provides ideal healing conditions after implantation. After implantation Lyoplast® is slowly degraded enzymatically and is successively replaced by endogenous connective tissue.

Lyoplast® does not possess any independent pharmacological effect. The collagen, which forms the implant, is degraded to amino acids in the body. Lyoplast® is sterilized with ethylene oxide.

## Indications

Replacement and extension of connective tissue structure in neurosurgery:

- for covering cerebral and cerebellar dura defects
  - after extirpation of tumours
  - for sealing cerebrospinal fluid fistulae
- for cerebral decompression surgery when there is elevated intracranial pressure
- for covering spinal dura defects
  - after removal of spinal tumours
  - after spinal traumata
- for spinal decompression surgery

## Mode of action

After implantation the Lyoplast® implant is colonized by endogenous connective tissue cells, i. e. revitalization takes place. The loosely fibrous structure, the low material thickness and the good connection of Lyoplast® with the host tissue are all of importance in this respect.

The colonization of the implant by connective tissue cells begins after only a few days. The complete revitalization takes place within a period of 1 to 3 months.

## Contra-indications

Lyoplast® should not be applied

- in infected regions
- to replace connective tissue structures that are subject to mechanical stress
- to replace parts of the arterial system or the cardiac wall
- in case of observed hypertensive reactions against proteins of bovine origin

## Warning notes

None

## Precautionary measures in application

See under mode of application

## Interactions with other agents

None known

## Mode of application

If not otherwise prescribed the suitable Lyoplast® size is chosen to fit the application area and is cut to fit the size of the defect.

Before the start of the operation it is advisable to place the Lyoplast® in sterile physiological saline solution or another isotonic solution in order to obtain better suppleness and flexibility. The implant should be cut to fit the defect as closely as possible in order to obtain bedding, which is as stress-free as possible. Lyoplast® can be sutured and additionally fastened with fibrin adhesive.

Before implantation of the Lyoplast® the person using it must be familiar with the surgical technique required, the specific applications and the properties of Lyoplast® in vivo.

## Side effects

- Tissue adhesion to the bio-implant can occur occasionally.
- It is not possible to exclude immune reactions with certainty.

## Notes and details of stability

Lyoplast® should not be used after the expiry date given.

Lyoplast® should be used immediately once the package is opened.

Lyoplast® should be stored at 25 + 5 °C.

## Other information

Lyoplast® may only be used if the packaging is undamaged. Opened packs of Lyoplast® and pieces of implant that are no longer required, must not be used later. Resterilization must not be carried out, because the structure of the Lyoplast® implant and, hence, its behaviour in vivo can be detrimentally affected.

Lyoplast® may not be resterilized!

Safety with regard to the transmission of zoo-anthropozoonoses In view of the fact that bovine material from New Zealand is regarded as safe by the European authorities with respect to BSE (bovine spongiform encephalopathy), the raw material is imported from there. Furthermore, Lyoplast® is subjected to treatment with NaOH during processing, in order to further reduce any theoretical risk, by means of this recognized decontamination method.

## Product range

Packs of 1 piece

| Size      | REF      |
|-----------|----------|
| 6 x 14 cm | 106 6021 |
| 4 x 10 cm | 106 6048 |
| 6 x 8 cm  | 106 6242 |
| 8 x 9 cm  | 106 6030 |
| 5 x 6 cm  | 106 6050 |

Packs of 2 pieces

| Size       | REF      |
|------------|----------|
| 4 x 5 cm   | 106 6064 |
| 2 x 10 cm  | 106 6080 |
| 1,5 x 3 cm | 106 6102 |

## Symbols used on labelling

Do not reuse!

Use until Year + Month

Temperature limitation

Sterile unless package is opened or damaged.  
Method of Sterilization - Ethylene oxide

0123 CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.

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# Lyoplant®

## Product description

Lyoplant® is a lyophilized dura substitute of pure collagen derived from bovine pericardium. The product is provided sterile, for single use only.

## Composition

Lyoplant® consists of specially purified collagen. It is EO sterilized.

## Intended use

Lyoplant® is intended as dura mater substitute in neurological procedures for soft tissue reconstruction of damaged, impaired, or missing tissue, such as for:

- Closing cerebral and cerebellar dural defects, such as after excision of tumors, operations on the posterior cranial fossa, after the removal of basal meningiomas and to seal CSF fistulas such as after frontobasal fractures;  
and
- Closure of spinal dural defects after removal of spinal tumors and after spinal trauma.

## Contraindications

Lyoplant® should not be:

- implanted in infected areas, open cranial trauma or open spina bifida.
- used as a replacement for mechanically stressed connective tissue structures, or as substitute for parts in the arterial system or the cardiac wall.

## How supplied

The product is supplied sterile, non-pyrogenic and packaged as single or double pieces. The products are available in various dimensions ranging from 1.5 x 3 cm to 6 x 14 cm pieces. Specific sizes are noted by catalog number below.

## Precautions

Before using Lyoplant®, the surgeon should be familiar with the surgical technique, special applications and the properties of Lyoplant® in vivo.

## Warnings

Lyoplant® should only be used if the package is intact. Visually inspect the packaging; do not use product if the package is found opened, punctured, torn or tampered with as sterility may be compromised.

Do not resterilize, as the sterilization process may alter the structure of Lyoplant® and thus its properties in vivo.

Do not use left-over implant pieces.

Do not use product beyond the expiration date.

U.S.A. Federal Law restricts this device to sale by or on the order of a physician.

Applied in case of observed hypertensive reactions against proteins of bovine origin.

## Possible adverse reactions

Hypersensitivity or an immune response to Lyoplant® may occur. The physician should closely monitor the patient for such a reaction and treat accordingly.

Other possible adverse reactions may include, but are not limited to infection, postoperative adhesions, leakage of cerebrospinal fluid and hematoma.

## Directions for use

For improved flexibility and easier shaping, it is recommended that the graft be adequately rehydrated in sterile physiological saline or another isotonic solution prior to implantation. To minimize tension after implantation, the correct size of Lyoplant® should be chosen and cut as closely as possible to fit the defect. Lyoplant® can be cut before or after rehydration. Proper sizing of Lyoplant® is important for optimal results.

It is recommended that Lyoplant® be sutured in place using non-absorbable (polyester, polypropylene) suture material to improve suture security. The use of atraumatic round-bodied needles will allow suturing without significant damage to the implant.

## Storage

Store the Lyoplant® box in a clean, dry and protected area with protection from extreme changes in temperature and humidity.

## Product range

Packs of 1 piece

| Dimensions | REF      |
|------------|----------|
| 6 x 14 cm  | 106 6021 |
| 4 x 10 cm  | 106 6048 |
| 6 x 8 cm   | 106 6242 |
| 8 x 9 cm   | 106 6030 |
| 5 x 6 cm   | 106 6050 |

Packs of 2 pieces

| Dimensions | REF      |
|------------|----------|
| 4 x 5 cm   | 106 6064 |
| 2 x 10 cm  | 106 6080 |
| 1,5 x 3 cm | 106 6102 |

## Distributor in the US/Contact in Canada for product information and complaints

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CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.

Batch Number

See Instructions for Use

Cat. No.

**DIM**

Size

**B | BRAUN**  
SHARING EXPERTISE

**Aesculap AG**  
Am Aesculap-Platz  
78532 Tuttlingen  
Phone +49 7461 95-0  
Fax +49 7461 95-2600  
[www.aesculap.de](http://www.aesculap.de)

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Aesculap Inc.