Lyoplast® Onlay

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Lyoplant® Onlay

Product Description

Lyoplant® Onlay is a lyophilized dura substitute of pure collagen derived from bovine pericardium and bovine split hide. The product is provided sterile, for single use only.

Federal law restricts this device to sale by or on the order of a physician.

Composition

Lyoplant® Onlay consists of specially purified collagen. It is EtO sterilized.

Indications for Use

Lyoplant® Onlay is indicated as a dura substitute for the repair of the dura mater.

Contraindications

- Do not use in infected areas, open cranial trauma or open spina bifida.
- Do not use as a replacement for mechanically stressed connective tissue structures.
- Do not use as a substitute for parts in the arterial system or the cardiac wall.
- Do not use if there are known allergies to proteins of bovine origin.

Precautions

Before using Lyoplant® Onlay, the surgeon should be familiar with the surgical technique, special applications and the properties of Lyoplant® Onlay in vivo.

Warnings

- Lyoplant® Onlay should only be used if the package is intact. Visually inspect the packaging; do not use product if the package is found opened, punctured, torn or tampered with as sterility may be compromised.
- Do not re-sterilize, as the sterilization process may alter the structure of Lyoplant® Onlay and thus its properties in vivo.
- Do not use left-over implant pieces.
- Do not use product beyond expiration date.

How Supplied

The product is supplied sterile and packaged as a single piece. The product is available in various dimensions ranging from 2.5 x 2.5 cm to 10.0 x 12.5 cm pieces. Specific sizes are listed by catalog number below:

Catalog Number	Dimensions
1067010	2.5 x 2.5 cm
1067020	5.0 x 5.0 cm
1067030	2.5 x 7.5 cm
1067040	7.5 x 7.5 cm
1067050	10.0 x 12.5 cm

Possible Adverse Reactions

Hypersensitivity or an immune response to Lyoplant® Onlay may occur. The physician should closely monitor the patient for such a reaction and treat accordingly. Other possible adverse reactions may include, but are not limited to infection, postoperative adhesions, leakage of cerebrospinal fluid and hematoma.

Recommendations for Use

To minimize tension after implantation, the correct size of Lyoplant® Onlay should be chosen and cut as closely as possible to fit the defect. Proper sizing of Lyoplant® Onlay is important for optimal results.

For improved flexibility and easier shaping, it is recommended that the graft be adequately rehydrated before use. Soak Lyoplant® Onlay in sterile 0.9% saline for a minimum of 3 minutes. Hydration is indicated by a color change from white (dry) to grey (wet).

Ensure that the fleece-like, porous surface is facing the dura (stamped "DURA SIDE" on graft). Identify the direction in which each surface should be facing before hydration.

Lyoplant® Onlay can be used as an onlay or sutured into place. When suturing, non-absorbable (polyester, polypropylene) suture material is recommended to improve suture security. The use of atraumatic round-bodied needles will allow suturing without significant damage to the implant.

Storage

Store unopened product in a clean, dry and protected area between 20–30°C (68–86°F).

Lyoplast® Onlay

Description

Lyoplast® Onlay is an implant made from pure collagen obtained from bovine pericardium and bovine split-skin. In a special processing procedure, non-collagen components such as enzymes, fats and non-collagen proteins are removed from the Lyoplast® Onlay.

The gentle freeze-drying process ensures that the loose fiber architecture of Lyoplast® Onlay is preserved, which provides optimum integration and healing conditions after implantation.

After implantation, Lyoplast® Onlay is slowly broken down by enzymes and replaced by the body's own connective tissue. Lyoplast® Onlay has no pharmacological action of its own. The collagen implant is broken down by the body into amino acids.

Composition

Lyoplast® Onlay is a biphasic three-dimensional collagen type I/III product. Lyoplast® Onlay is composed of 12 ± 4 mg/cm² pericardium and 10 ± 2 mg/cm² sponge components. Lyoplast® Onlay is not chemically cross-linked.

Indication

Replacement and extension of connective tissue structures in neurosurgery:

- to close holes in the cerebral and cerebellar dura
 - after extirpation of tumors
 - to close CSF fistulas
- for procedures to relieve increased intracranial pressure
- to close holes in the spinal dura
 - after removal of spinal tumors
 - after spinal traumas
- for procedures to relieve spinal pressure

Mode of application

The appropriate size of the Lyoplast® Onlay is chosen according to the area of application and cut according to the size of the defect, leaving an excess edge on the implant of approximately one centimeter (1 cm). Prior to implantation, Lyoplast® Onlay is placed in sterile saline solution or in another isotonic solution to ensure better suppleness and flexibility. The implant should be fixed without being pulled tautly. Lyoplast® Onlay consists of two different sides. Ensure that the fleece-like, porous side (labeled „DURA SIDE“) is facing the dura. Identify the direction in which each side should be facing before hydration. Lyoplast® Onlay can be used as an Onlay graft. Ensure that the implant is lying flat against the defect edges and is not taut. In all other cases it is recommended to fixate the implant as well.

Lyoplast® Onlay can simply be applied, or, if required and if considered necessary by the user, it can also be sutured in place and sealed with fibrin glue.

Prior to implanting Lyoplast® Onlay, the user must be acquainted with the surgical technique, special applications and the characteristics of Lyoplast® Onlay in vivo.

Mode of operation

After implantation, the Lyoplast® Onlay implant integrates with the body's own connective tissue cells, i.e. revitalization takes place. The fiber architecture of the compact side of the implant, the narrow thickness of the material, and a good connection with the recipient tissue are particularly important for this. The fleece-like quality of the porous side allows the Lyoplast® Onlay to adhere to the edges of the defect.

The implant starts to integrate with the connective tissue cells after just a few days. Full revitalization takes place within 1 to 3 months.

Contraindications

Lyoplast® Onlay should not be used

- in infected regions
- as a replacement for connective tissue structures that undergo mechanical stress (tendons/ligaments)
- as a replacement for parts of the arterial system or walls of the heart
- in case of known hypersensitivity to proteins of bovine origin.

Warnings

WARNING: Risk of changing the material properties of Lyoplast® Onlay if resterilized!

Lyoplast® Onlay must not be resterilized as this can adversely alter the structure of the implant and therefore its behavior in vivo.

Lyoplast® Onlay may only be used if the packaging is undamaged.

Do not reuse already opened packages of Lyoplast® Onlay or implant pieces that were not required.

WARNING: Do not reuse the product!

Reuse of the implant can cause infections in the patient or user and functional loss of the implant. Risk of injury, illness or death through contamination and/or limited functionality of the product.

Safety and precautions

Safety against the transmission of zoonoses.

Due to the fact that bovine material from New Zealand and Australia is considered BSE (bovine spongiform encephalopathy) -safe by the European authorities, the raw material is imported from these countries. Lyoplast® Onlay is also treated with sodium hydroxide solution during its further processing, and this recognized depletion process additionally reduces the theoretical risk of transmission.

Side effects

- Occasionally, tissue adhesions can occur.
- Immune reactions cannot be ruled out with absolute certainty.
- The occurrence of CSF leaks cannot be fully ruled out in some individual cases.

Sterilization

Lyoplast® Onlay is EO-sterilized.

Other instructions

Do not use Lyoplast® Onlay after the expiry date indicated.

Use Lyoplast® Onlay immediately after opening the packaging.

Store Lyoplast® Onlay at 25 ± 5 °C.

Delivery forms and packaging sizes

Pack with 1 unit

DIM	REF
2,5 × 2,5 cm	1067010
5,0 × 5,0 cm	1067020
2,5 × 7,5 cm	1067030
7,5 × 7,5 cm	1067040
10,0 × 12,5 cm	1067050

Explanation of the symbols on the packaging

	Follow the instructions for use!
	Do not reuse!
	Date of manufacture
	Use until year, month, day
	Maximum storage temperature
	Sterile, as long as the packaging is undamaged and unopened. Sterilization method: ethylene oxide.
	CE mark and identification number of the notified body. The product complies with the essential requirements of Council Directive 93/42/EEC concerning medical devices
	Batch designation
	Art. no.
DIM	Size

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Lyoplant® Onlay

Description

Lyoplant® Onlay is an implant made from pure collagen obtained from bovine pericardium and bovine split-skin. In a special processing procedure, non-collagen components such as enzymes, fats and non-collagen proteins are removed from the Lyoplant® Onlay.

The gentle freeze-drying process ensures that the loose fiber architecture of Lyoplant® Onlay is preserved, which provides optimum integration and healing conditions after implantation.

After implantation, Lyoplant® Onlay is slowly broken down by enzymes and replaced by the body's own connective tissue. Lyoplant® Onlay has no pharmacological action of its own. The collagen implant is broken down by the body into amino acids.

Composition

Lyoplant® Onlay is a biphasic three-dimensional collagen type I/III product. Lyoplant® Onlay is composed of 12 ± 4 mg/cm² pericardium and 10 ± 2 mg/cm² sponge components. Lyoplant® Onlay is not chemically cross-linked.

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Mode of application

The appropriate size of the Lyoplant® Onlay is chosen according to the area of application and cut according to the size of the defect, leaving an excess edge on the implant of approximately one centimeter (1 cm). Prior to implantation, Lyoplant® Onlay is placed in sterile saline solution or in another isotonic solution to ensure better suppleness and flexibility. The implant should be fixed without being pulled tautly. Lyoplant® Onlay consists of two different sides. Ensure that the fleece-like, porous side (labeled "DURA SIDE") is facing the dura. Identify the direction in which each side should be facing before hydration. Lyoplant® Onlay can be used as an Onlay graft. Ensure that the implant is lying flat against the defect edges and is not taut. In all other cases it is recommended to fixate the implant as well.

Lyoplant® Onlay can simply be applied, or, if required and if considered necessary by the user, it can also be sutured in place and sealed with fibrin glue.

Prior to implanting Lyoplant® Onlay, the user must be acquainted with the surgical technique, special applications and the characteristics of Lyoplant® Onlay in vivo.

Mode of operation

After implantation, the Lyoplant® Onlay implant integrates with the body's own connective tissue cells, i.e. revitalization takes place. The fiber architecture of the compact side of the implant, the narrow thickness of the material, and a good connection with the recipient tissue are particularly important for this. The fleece-like quality of the porous side allows the Lyoplant® Onlay to adhere to the edges of the defect.

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Warnings

WARNING: Risk of changing the material properties of Lyoplant® Onlay if resterilized!

Lyoplant® Onlay must not be resterilized as this can adversely alter the structure of the implant and therefore its behavior in vivo.

Lyoplant® Onlay may only be used if the packaging is undamaged.

Do not reuse already opened packages of Lyoplant® Onlay or implant pieces that were not required.

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Due to the fact that bovine material from New Zealand and Australia is considered BSE (bovine spongiform encephalopathy) -safe by the European authorities, the raw material is imported from these countries. Lyoplant® Onlay is also treated with sodium hydroxide solution during its further processing, and this recognized depletion process additionally reduces the theoretical risk of transmission.

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Lyoplant® Onlay is EO-sterilized.

Other instructions

Do not use Lyoplant® Onlay after the expiry date indicated.

Use Lyoplant® Onlay immediately after opening the packaging.

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Use by year, month, day

Temperature limitation

Sterile, as long as the packaging is undamaged and unopened.

Sterilization method: ethylene oxide.

CE mark and identification number of the notified body. The product complies with the essential requirements of Council Directive 93/42/EEC concerning medical devices

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Aesculap AG, Am Aesculap-Platz, 78532 Tuttlingen,
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany