

Aesculap Spine


- (GB)** Instructions for use/Technical description
MACS II
- (USA)** Gebrauchsanweisung/Technische Beschreibung
MACS II
- (D)** Mode d'emploi/Description technique
MACS II
- (F)** Instrucciones de manejo/Descripción técnica
MACS II
- (E)** Istruzioni per l'uso/Descrizione tecnica
MACS II
- (I)** Instruções de utilização/Descrição técnica
MACS II
- (P)** Gebruiksaanwijzing/Technische beschrijving
MACS II

- (S)** Bruksanvisning/Teknisk beskrivning
MACS II
- (RUS)** Инструкция по применению/Техническое описание
MACS II
- (CZ)** Návod k použití/Technický popis
MACS II
- (PL)** Instrukcja użytkowania/Opis techniczny
MACS II
- (SK)** Návod na použitie/Technický opis
MACS II
- (TR)** Kullanım Kılavuzu/Teknik açıklama
MACS II


A


C


D



F


G

H

Aesculap®

MACS II

Legend

A Mounting instruments

- 1 Torque wrench for mounting the retaining bolt
- 2 Mounting tool for SX804T

B K-wire instruments

- 3 K-wire protection
- 4 K-wire set, sterile
 - a K-wire (x2)
 - b Removal instrument (x1)
- 5 K-wire insertion instrument
 - a Aiming tube
 - b Handle rod
 - c Impactor
 - d Counter holder

C Insertion instruments for polyaxial screw

- 6 Cannulated screwdriver
- 7 Guide tube
- 8 Counter torque (cleaning rod, counter torque with-out cleaning rod)
- 9 Ball-end screwdriver

D Insertion instruments for stabilization plate

- 7 Guide tube
- 8 Counter torque (cleaning rod, counter torque with-out cleaning rod)
- 10 Holding forceps for rods/plates
- 11 Socket wrench for nut
- 12 Torque wrench for nut
- 13 Screwdriver for retaining bolt (optional)

E Insertion instruments for SSE connecting rods

- 7 Guide tube
- 8 Counter torque (cleaning rod, counter torque with-out cleaning rod)
- 10 Holding forceps for rods/plates
- 11 Socket wrench for nut
- 12 Torque wrench for nut
- 13 Screwdriver for retaining bolt (optional)
- 14 SSE rod bending pliers

F Insertion instruments for stabilization screw

- 6 Cannulated screwdriver
- 15 Screw revision instrument
- 16 Lumbar centering tool
- 17 Aiming sleeve for lumbar stabilization screw
- 18 Aiming sleeve for thoracic stabilization screw
- 19 Centering tool for purple stabilization screw
- 20 Screwdriver for purple stabilization screw, see TA011177

G Insertion instruments for stabilization screw

- 7 Guide tube
- 8 Counter torque (cleaning rod, counter torque with-out cleaning rod)
- 13 Screwdriver for retaining bolt (optional)
- 21 Torque wrench for retaining bolt

H Revision instruments

- 6 Cannulated screwdriver
- 9 Ball-end screwdriver
- 10 Holding forceps for rods/plates
- 11 Socket wrench for nut
- 13 Screwdriver for retaining bolt (optional)
- 15 Screw revision instrument
- 20 Screwdriver for purple stabilization screw, see TA011177
- 22 T-handle with 1/4" coupling

I Instruments for intraoperative disassembly/ mounting of the pre-assembled retaining bolt

- 1 Torque wrench for mounting the retaining bolt
- 7 Guide tube
- 8 Counter torque (cleaning rod, counter torque without cleaning rod)
- 13 Screwdriver for retaining bolt (optional)
- 22 T-handle with 1/4" coupling

J Implant components

- 23 SSE connecting rod, 5.2 mm
- 24 Stabilization plate
- 25 Clamping element T
- 26 Clamping element TL
- 27 Pre-assembled retaining bolt
 - a Centralizer
 - b Retaining bolt
- 28 Locknut
- 29 Stabilization screw, 4.5 mm, purple
- 30 Stabilization screw, 6.5 mm, yellow
- 31 Polyaxial screw, 7.0 mm, blue
- 32 Polyaxial screw XL, 10 mm, blue

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbBraun.com>

Intended use

The MACS II system is used for ventral monosegmental and multisegmental stabilization of the lumbar and thoracic vertebrae.

Aesculap®

MACS II

System components and available sizes

Item no.	Art. no.	Designation
1	FW419R	Torque wrench for mounting the retaining bolt
2	FW429	Mounting tool for SX804T
3	FW394R	K-wire protection
4	FW406S	K-wire set, sterile (K-wire x2, removal instrument x1)
5	FW408R	K-wire insertion instrument (aiming tube, handle rod, impactor, counter holder)
6	FW398R	Cannulated screwdriver
7	FW397R	Guide tube
8	FW399R	Counter torque (cleaning rod, counter torque without cleaning rod)
9	FW396R	Ball-end screwdriver
10	FW326R	Holding forceps for rods/plates
11	FW417R	Socket wrench for nut
12	FW416R	Torque wrench for nut (casing, handle, wrench)
13	FW411R	Screwdriver for retaining bolt
14	FW013R	SSE rod bending pliers
15	FW336R	Screw revision instrument
16	FW407R	Lumbar centering tool
17	FW409R	Aiming sleeve for lumbar stabilization screw

Item no.	Art. no.	Designation
18	FW413R	Aiming sleeve for thoracic stabilization screw
19	FW405R	Centering tool for purple stabilization screw
20	FW342R	Screwdriver for purple stabilization screw (see TA011177 instructions for use)
21	FW412R	Torque wrench for retaining bolt
22	FW395R	T-handle, 1/4" coupling
23	SX110T	100 mm SSE connecting rod, 5.2 mm
	SX112T	120 mm SSE connecting rod, 5.2 mm
	SX115T	150 mm SSE connecting rod, 5.2 mm
	SX118T	180 mm SSE connecting rod, 5.2 mm
	SX120T	200 mm SSE connecting rod, 5.2 mm
24	SX810T	Stabilization plate, 40 mm
	SX811T	Stabilization plate, 45 mm
	SX812T	Stabilization plate, 50 mm
	SX813T	Stabilization plate, 55 mm
	SX814T	Stabilization plate, 60 mm
	SX834T	Stabilization plate, 65 mm
	SX835T	Stabilization plate, 70 mm
	SX839T	Stabilization plate, 75 mm
	SX816T	Stabilization plate, 80 mm
	SX840T	Stabilization plate, 85 mm
25	SX817T	Stabilization plate, 90 mm
	SX841T	Stabilization plate, 95 mm
26	SX818T	Stabilization plate, 100 mm
	SX800T	Clamping element TL

Item no.	Art. no.	Designation
27	SX804T	Pre-assembled retaining bolt (retaining bolt, centralizer)
28	SX802T	Locknut
29	SX826T	Stabilization screw, 4.5 mm x 15 mm, purple
	SX827T	Stabilization screw, 4.5 mm x 20 mm, purple
	SX828T	Stabilization screw, 4.5 mm x 25 mm, purple
	SX829T	Stabilization screw, 4.5 mm x 30 mm, purple
30	SX782T	Stabilization screw, 6.5 mm x 25 mm, yellow
	SX783T	Stabilization screw, 6.5 mm x 30 mm, yellow
	SX784T	Stabilization screw, 6.5 mm x 35 mm, yellow
	SX785T	Stabilization screw, 6.5 mm x 40 mm, yellow
	SX786T	Stabilization screw, 6.5 mm x 45 mm, yellow
	SX787T	Stabilization screw, 6.5 mm x 50 mm, yellow

Item no.	Art. no.	Designation
31	SX790T	Polyaxial screw, 7.0 mm x 15 mm, blue
	SX791T	Polyaxial screw, 7.0 mm x 20 mm, blue
	SX792T	Polyaxial screw, 7.0 mm x 25 mm, blue
	SX793T	Polyaxial screw, 7.0 mm x 30 mm, blue
	SX794T	Polyaxial screw, 7.0 mm x 35 mm, blue
	SX795T	Polyaxial screw, 7.0 mm x 40 mm, blue
	SX796T	Polyaxial screw, 7.0 mm x 45 mm, blue
	SX797T	Polyaxial screw, 7.0 mm x 50 mm, blue
	SX798T	Polyaxial screw, 7.0 mm x 55 mm, blue
32	SX821T	Polyaxial screw XL, 10 mm x 30 mm, blue
	SX823T	Polyaxial screw XL, 10 mm x 40 mm, blue

Aesculap®

MACS II

System configurations

System incompatibility may result in increased operating time or abortion of the procedure!

- ▶ Follow only these instructions for use and the TA013373 instructions for use when using the pre-assembled retaining bolt SX804T.
- ▶ Do not use MACS II in conjunction with the HMA monoaxial or polyaxial screws of MACS TL®.

- Polyaxial Twin Screw system TL **26, 30, 31** with stabilization plate **24** or two SSE 5.2 mm connecting rods **23**
- Polyaxial Twin Screw system T **25, 29, 31** with stabilization plate **24** or two SSE 5.2 mm connecting rods **23**
- Polyaxial Twin Screw system XL **30, 32** with stabilization plate **24** or two SSE 5.2 mm connecting rods **23**

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

Risk of injury caused by incorrect operation of the product!

- ▶ Attend appropriate product training before using the product.
- ▶ For information about product training, please contact your national B. Braun/Aesculap agency.

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.
- ▶ Do not use products that are in open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.

Pre-assembled retaining bolt and sterile K-wire set

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- ▶ Do not reprocess the product.

The product is gamma-sterilized and supplied in sterile packaging.

The product must not be reused.

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Do not use the product after its use-by date.

Safe operation

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

Risk of injury when using the product beyond the field of view!

- ▶ Apply the product only under visual control.

Note

For a detailed description of the operating procedure, see *MACS II operating technique*, brochure number 042702.

Mounting the clamping element and 7.0 mm blue polyaxial screw or 10 mm blue polyaxial screw XL with pre-assembled retaining bolt

- ▶ Select either clamping element TL, T **25**, **26** or 10 mm blue polyaxial screw XL **32**.
- ▶ Select 7.0 mm blue polyaxial screw **31** and stabilization screw **29**, **30** (this step is missed out if the 10 mm blue polyaxial screw XL is used **32**).

Aesculap®

MACS II

WARNING

Risk of injury if incorrect anchoring screws are selected!

- ▶ To ensure safe implantation of the anchoring screws, ensure that the diameter of the vertebral body is carefully measured in the planned screw position and direction.
- ▶ Choose a screw length shorter than the diameter measured.

- ▶ Screw the 7.0 mm blue polyaxial screw **31** into the clamping element **25, 26**.
- ▶ Insert the assembled 7.0 mm blue polyaxial screw **31** and clamping element **25, 26** or XL screw **32** into the mounting tool for SX804T **2**.

WARNING

Risk of injury if centralizer comes loose during the procedure!

- ▶ If the 10 mm blue polyaxial screw XL is used, insert this in the appropriate side of the mounting tool for SX804T (labeled with XL).
- ▶ If the 7.0 mm blue polyaxial screw is used, insert this in the appropriate side of the mounting tool for SX804T (labeled with POLY).

- ▶ Remove the pre-assembled retaining bolt **27** from the packaging using the torque wrench for mounting the retaining bolt **1**.

WARNING

Take care not to drop the pre-assembled retaining bolt!

- ▶ Open the packaging correctly and do not turn it to prevent the retaining bolt from falling out.
- ▶ Apply the torque wrench for mounting the retaining bolt as far as it will go onto the pre-assembled retaining bolt.

- ▶ Screw the pre-assembled retaining bolt **27** onto the clamping element **25, 26** with 1.8 Nm.

WARNING

Risk of injury if pre-assembled retaining bolt comes loose during the procedure!

- ▶ To mount the pre-assembled retaining bolt, use the mounting tool for SX804T and the torque wrench for mounting the retaining bolt and tighten with 1.8 Nm.
- ▶ Prior to applying the torque wrench, make certain the display arrow points to zero.
- ▶ Replace the torque wrench if it is damaged.

- Check that the retaining lug of the centralizer **27a** is correctly positioned in the clamping element **25**, **26**.

Risk of injury if the necessary counter torque is not applied when tightening the locknut or pre-assembled retaining bolt SX804T!

- Make sure that the retaining lugs of the centralizer engage fully in the appropriate grooves in the clamping element.

Applying the K-wires

- Assemble the K-wire insertion instrument **5**.
- Insert the threaded end of the K-wire **4a** into the end of the aiming tube **5a** and screw tight.

Damage to the K-wire insertion instrument if the thread of a K-wire is damaged!

- Do not re-use a sterile K-wire set.

- Position the K-wire **4a** correctly and drive it in under radiographic control.

Risk of injury if K-wires are positioned incorrectly!

- Always use radiographic control when positioning K-wires.
- Use radiographic control to constantly monitor the position of the K-wires throughout the operation.
- Apply K-wires in such a way that neither nerve roots, the spinal column, nor blood vessels or organs are injured in any way!

- Unscrew the impactor **5c** counterclockwise.
- Open the cortical bone.

The K-wire will twist if the aiming tube is tilted!

- Do not tilt the K-wire insertion instrument and K-wire assembly.

Aesculap®

MACS II

- Insert the counter holder **5d** into the K-wire insertion instrument **5**.

WARNING

Risk of injury if incorrect instruments are used!

- Do not use the screw revision instrument to apply counterforce at the K-wires.
- Only use the counter holder to apply counterforce at the K-wires.

- Remove the K-wire insertion instrument **5** from the K-wire **4a**. While doing so, use the counter holder **5d** to apply counterforce to the K-wire **4a**.

WARNING

Risk of injury if K-wire is driven forwards with counter holder!

- Remove the K-wire insertion instrument from the K-wire under radiographic control.

- Attach K-wire protection **3** if desired, using the holding forceps for rods/plates **10** in order to protect the adjacent soft tissue.

WARNING

Risk of injury if K-wire is driven forwards with K-wire protection!

- When attaching the K-wire protection, ensure that the K-wire is not pushed forwards.
- Attach the K-wire protection under radiographic control.

Implanting the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL

- Remove K-wire protection **3** using the holding forceps for rods/plates **10**.

WARNING

Risk of injury caused by the tip of the K-wire!

- Ensure that, when the K-wire protection is being removed with the holding forceps for rods/plates, the K-wire is not withdrawn from the vertebral body.

- Attach the counter torque **8** onto the guide tube **7**.

WARNING

Damage to the guide tube caused by incorrect operation!

- Fit the counter torque onto the guide tube until it clicks into position.
- To detach the counter torque from the guide tube, release the snap connection by pulling the handle of the cleaning rod, and pull the counter torque up and out.
- Fit the counter torque onto the correct side of the guide tube. Look for the gold marking.

- Attach the guide tube **7** onto the centralizer **27a** ex situ. Ensure that this is done in the correct direction.

Risk of the cannulated screwdriver getting jammed in the guide tube caused by incorrect operation!

- **Ensure that the guide tube snaps fully onto the centralizer and a click is heard.**
- Insert the cannulated screwdriver **6** into the guide tube **7**.
- Guide the 7.0 mm blue polyaxial screw **31** over the K-wire **4a** and use the cannulated screwdriver **6** to screw it in, until the first turns of the thread are held in the vertebral body.

Risk of injury if the K-wire is driven forwards with the 7.0 mm blue polyaxial screw, the 10 mm blue polyaxial screw XL, or the cannulated screwdriver!

- Insert the cannulated screwdriver fully into the guide tube when threading and when screwing in the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL.
- Never withdraw the cannulated screwdriver from the guide tube and reinsert it while screwing the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL into the vertebral body.
- Thread the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL over the K-wire in an orthograde direction.
- Screw in the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL a maximum of 2 – 3 turns before withdrawing the K-wire using the removal instrument.

Aesculap®

MACS II

- ▶ Insert the removal instrument **4b** for the K-wire **4a** into the cannulated screwdriver **6** and screw onto K-wire **4a**.
- ▶ Remove the K-wire **4a**.
- ▶ Screw down the 7.0 mm blue polyaxial screw **31** or the 10 mm blue polyaxial screw XL **32**.
The clamping element **25, 26** lies almost on the vertebral body.
- ▶ Ensure that the polyaxiality of the clamping element **25, 26** is not restricted.

WARNING

Risk of injury if posterior wall of the vertebral body is penetrated!

- ▶ Screw in the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL under radiographic control.

If SSE connecting rods or stabilization plates cannot be screwed down:

- ▶ Only screw in the 7.0 mm blue polyaxial screw **31** or the 10 mm blue polyaxial screw XL **32** as far as necessary so that polyaxiality is not restricted.
- ▶ Pull the cannulated screwdriver **6** out of the guide tube **7**.
- ▶ Remove the guide tube **7**, including the counter torque **8**, from the centralizer **27a**.

WARNING

Risk of injury if the 7.0 mm blue polyaxial screw or 10 mm blue polyaxial screw XL is pulled out of the vertebral body through incorrect operation or non-observance of the correct sequence of operating steps!

- ▶ For easy withdrawal of the cannulated screwdriver, apply the insertion instrument orthograde to the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL.
- ▶ Remove the cannulated screwdriver from the guide tube before detaching the guide tube together with counter torque from the centralizer.
- ▶ Apply pulling force over the head part of the counter torque.

Fastening the stabilization plate

- Measure the distance between the clamping elements **25**, **26** using bone graft measuring instrument FG045R.

Note

The correct plate length is derived from this measured value plus 30 mm.

WARNING

Risk of injury to the bordering mobile segment caused by projecting stabilization plate!

- Choose the smallest stabilization plate possible (at 5 mm intervals).

- Thread the stabilization plate **24** onto the centralizer **27a** using the holding forceps for rods/plates **10**. Ensure that the flat, labeled side faces upwards.

WARNING

Risk of injury caused by insufficient system stability!

- Place the stabilization plate onto the clamping element or the 10 mm blue polyaxial screw XL with the flat, labeled side facing upwards.

- Pick up the locknut **28** using the socket wrench for nut **11**.
- Insert the guide tube **7** with the counter torque **8** into the socket wrench for nut **11**.
- Snap the guide tube **7** using the socket wrench for nut **11** onto the centralizer **27a** in situ.

If you cannot snap the guide tube **7** onto the centralizer **27a**:

- To snap the guide tube onto the centralizer: Remove the screwdriver **6**, **13** from the guide tube **7**.
– or –
- Guide the guide tube **7** onto the centralizer **27a** in an orthograde direction.
- Optional: Insert the screwdriver for retaining bolt **13** into the guide tube **7**.
- Screw on the locknut **28** using the T-handle of the socket wrench for nut **11**.

WARNING

Risk of injury if the locknut has insufficient clamping stability!

- Apply locknuts in an orthograde direction.
- Tighten the locknuts smoothly to avoid damaging the external thread of the clamping element or the 10 mm blue polyaxial screw XL.
- Replace the external thread of the clamping element or the 10 mm blue polyaxial screw XL if it is damaged.

Aesculap®

MACS II

If you cannot screw the locknut **28** onto the clamping element **25**, **26** or the 10 mm blue polyaxial screw XL **32**:

- ▶ Apply tensile force via the guide tube **7**. To do this, insert the screwdriver for retaining bolt **13** into the guide tube **7**.
- ▶ Tighten the locknut **28** using the torque wrench for nut **12**. Use the counter torque **8** to apply counterforce.

Risk of injury if the locknut has insufficient clamping stability!

- ▶ Set locknuts in place correctly.
- ▶ Tighten locknuts with the torque wrench for nut.
- ▶ Have the torque wrench for the nut recalibrated by the manufacturer at intervals no longer than every two years.

Risk of injury if no or insufficient counterforce is applied!

- ▶ Use the counter torque provided to apply the counterforce.
- ▶ Ensure that as little force as possible is transferred to the spinal column.

Damage to the implant caused by overtightening of the locknut!

- ▶ Tighten the locknut to the exact torque specified for this purpose.
- ▶ Do not block the release mechanism of the torque wrench for nut.
- ▶ Have the torque wrench for the nut recalibrated by the manufacturer at intervals no longer than every two years.

Injuries to blood vessels and/or soft tissue caused by the torque wrench for the nut slipping off from the socket wrench for the nut!

- ▶ Fit the wrench part of the torque wrench for the nut fully onto the socket wrench for the nut.

- ▶ Optional for the transfer of tensile forces: Pull the screwdriver for retaining bolt **13** out of the guide tube **7**.

- Remove the guide tube **7** with the counter torque **8** from the centralizer **27a**. Use the socket wrench for nut **11** to apply counterforce.

Risk of injury if the 7.0 mm blue polyaxial screw or 10 mm blue polyaxial screw XL is pulled out of the vertebral body through incorrect operation or non-observance of the correct sequence of operating steps!

- For easy detachment from the centralizer, apply the guide tube with counter torque in an orthograde direction to the centralizer.
- Remove the screwdriver for retaining bolt from the guide tube before detaching the guide tube together with counter torque from the centralizer.
- Apply pulling force over the head part of the counter torque.

- Remove the socket wrench for nut **11** from the locknut **28**.

Fastening the SSE 5.2 mm connecting rods

- Measure the distance between the clamping elements **25**, **26** or the 10 mm blue polyaxial screws XL **32** using bone graft measuring instrument FG045R.

Note

The correct rod length is derived from this measured value plus 30 mm.

Risk of injury caused by projecting SSE connecting rods or sharp edges!

- Choose the shortest SSE connecting rod possible.
- If necessary, shorten the SSE connecting rod with a suitable instrument. Also round off sharp edges with a suitable tool if necessary.

- Adjust the SSE connecting rods **23** to the required curvature of the spine using SSE rod bending pliers **14**.

Risk of injury caused by insufficient system stability!

- Do not bend or bend back the SSE connecting rods excessively.
- Use SSE rod bending pliers to adjust the curvature of the SSE connecting rods.
- Do not bend the SSE connecting rods in the region of the support on the clamping elements or the 10 mm blue polyaxial screws XL.

Aesculap®

MACS II

- ▶ Pick up the locknut **28** using the socket wrench for nut **11**.
- ▶ Insert the guide tube **7** with the counter torque **8** into the socket wrench for nut **11**.
- ▶ Snap the guide tube **7** using the socket wrench for nut **11** onto the centralizer **27a** in situ.

If you cannot snap the guide tube **7** onto the centralizer **27a**:

- ▶ To snap the guide tube onto the centralizer: Remove the screwdriver **6**, **13** from the guide tube **7**.
– or –
- ▶ Guide the guide tube **7** onto the centralizer **27a** in an orthograde direction.
- ▶ Fit the first SSE 5.2 mm connecting rod **23** into the posterior hollow of the first clamping element **25**, **26** or the 10 mm blue polyaxial screws XL **32** using the holding forceps for rods/plates **10**.
- ▶ Ensure that the polyaxiality of the clamping element or the 10 mm blue polyaxial screws XL is not restricted.
- ▶ Screw on the locknut **28** using the T-handle of the socket wrench for nut **11** with a maximum of two turns.

The posterior rod is secured.

Risk of injury if the locknut has insufficient clamping stability!

- ▶ Apply locknuts in an orthograde direction.
- ▶ Tighten the locknuts smoothly to avoid damaging the external thread of the clamping element or the 10 mm blue polyaxial screw XL.
- ▶ Replace the external thread of the clamping element or the 10 mm blue polyaxial screw XL if it is damaged.

If you cannot screw the locknut **28** onto the clamping element **25**, **26** or the 10 mm blue polyaxial screw XL **32**:

- ▶ Apply tensile force via the guide tube **7**. To do this, insert the screwdriver for retaining bolt **13** into the guide tube **7**.
- ▶ Fit the first SSE connecting rod **23** into the posterior hollow of the second clamping element **25**, **26** or the 10 mm blue polyaxial screw XL **32** using the holding forceps for rods/plates **10**.
- ▶ Ensure that the polyaxiality of the clamping element is not restricted.
- ▶ Screw on the locknut **28**, using the T-handle of the socket wrench for nut **11**.
The posterior rod is secured.
- ▶ Fit the second SSE connecting rod **23** into the anterior hollow of the first and second clamping element **25**, **26** or the 10 mm blue polyaxial screw XL **32** using the holding forceps for rods/plates **10**.
- ▶ Tighten the locknut **28** using the torque wrench for nut **12**. Use the counter torque **8** to apply counterforce.

WARNING

Risk of injury if the locknut has insufficient clamping stability!

- ▶ Tighten locknuts with the torque wrench for nut.
- ▶ Have the torque wrench for the nut recalibrated by the manufacturer at intervals no longer than every two years.

WARNING

Risk of injury if no or insufficient counterforce is applied!

- ▶ Use the counter torque provided to apply the counterforce.
- ▶ Ensure that as little force as possible is transferred to the spinal column.

WARNING

Damage to the implant caused by overtightening of the locknut!

- ▶ Tighten the locknut to the exact torque specified.
- ▶ Do not block the release mechanism of the torque wrench for nut.
- ▶ Have the torque wrench for the nut recalibrated by the manufacturer at intervals no longer than every two years.

WARNING

Injuries to blood vessels and/or soft tissue caused by the torque wrench for the nut slipping off from the socket wrench for the nut!

- ▶ Fit the wrench part of the torque wrench for the nut fully onto the socket wrench for the nut.

- ▶ Optional for the transfer of tensile forces: Pull the screwdriver for retaining bolt 13 out of the guide tube 7.
- ▶ Remove the guide tube 7 with the counter torque 8 from the centralizer 27a. Use the socket wrench for nut 11 to apply counterforce.

WARNING

Risk of injury if the 7.0 mm blue polyaxial screw or 10 mm blue polyaxial screw XL is pulled out of the vertebral body through incorrect operation or non-observance of the correct sequence of operating steps!

- ▶ For easy detachment from the centralizer, position the guide tube with counter torque in an orthograde direction to the centralizer.
- ▶ Remove the screwdriver for retaining bolt from the guide tube before detaching the guide tube together with counter torque from the centralizer.
- ▶ Apply pulling force over the head part of the counter torque.

- ▶ Remove the socket wrench for nut 11 from the locknut 28.

Aesculap®

MACS II

Final application of the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL

- Screw down the polyaxial screw or XL screw with the ball-end screwdriver.

Risk of pulling the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL out of the vertebral body through incorrect operation or non-observance of the correct sequence of operating steps!

- Tighten the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL a maximum of one turn per side.
- Remove the ball-end screwdriver from the rear section of the 7.0 mm blue polyaxial screw or the 10 mm blue polyaxial screw XL by gently rotating it left and right.

Implantation of the stabilization screw

- Snap on the aiming sleeve for the stabilization screw 17, 18.

Risk of injury due to use of incorrect aiming sleeve, incorrect placement of the aiming sleeve, or failure to use an aiming sleeve for the stabilization screw!

- To ensure the aiming sleeve for the stabilization screw snaps on easily, guide the aiming sleeve onto the centralizer in an ortho-grade direction.
- Clamping element T is for use with the aiming sleeve for thoracic stabilization screw, see handle geometry T.
- Clamping element TL is for use with the aiming sleeve for lumbar stabilization screw, see handle geometry L.
- The 10 mm blue polyaxial screw XL is for use with the aiming sleeve for lumbar stabilization screw, see handle geometry L.
- Fit the aiming sleeve for stabilization screw so that it sits above the hole for the stabilization screw in the clamping element or the 10 mm blue polyaxial screw XL.
- Use the appropriate aiming sleeve for the stabilization screw when inserting the stabilization screw.

- Use centering tool **16, 19** to center it.

To avoid incorrect positioning of the stabilization screw **29, 30** or insufficient centering due to use of the incorrect instrument:

- When using clamping element **T 25**, use the centering tool for purple stabilization screw **19**.

– or –

- When using clamping element **TL 26**, use the lumbar centering tool **16**.
- When using clamping element **TL 26** or the 10 mm blue polyaxial screw **XL 32**, insert the screw revision instrument **15** into the cannulated screwdriver **6**.
- Apply stabilization screw **29, 30** onto the hexagon of the screwdriver **6, 20**.
- When using clamping element **TL 26** or the 10 mm blue polyaxial screw **XL 32**, screw the screw revision instrument **15** into the 6.5 mm yellow stabilization screw **30**.
- When using clamping element **T 25**, turn the knurled nut clockwise to fix the 4.5 mm purple stabilization screw **29** onto the screwdriver for purple stabilization screw **20**, see TA011177 instructions for use.
- Screw the stabilization screw **29, 30** into the vertebral body until it is fully seated in the clamping element **25, 26** or the 10 mm blue polyaxial screws **XL 32**.

WARNING

Risk of injury if stabilization screw comes out!

- Screw the stabilization screw in fully so that the self-locking mechanism engages.

WARNING

Risk of injury if posterior wall of the vertebral body is penetrated!

- Screw in the stabilization screw under radiographic control.
- Use the appropriate aiming sleeve for the stabilization screw when inserting the stabilization screw.

- When using clamping element **TL 26** or the 10 mm blue polyaxial screw **XL 32**, unscrew the screw revision instrument **15** from the 6.5 mm yellow stabilization screw **30**.
- When using clamping element **T 25**, turn the knurled nut counterclockwise to release the 4.5 mm purple stabilization screw **29**, see TA011177 instructions for use.
- When using clamping element **TL 26** or the 10 mm blue polyaxial screw **XL 32**, pull the cannulated screwdriver **6** out of the 6.5 mm yellow stabilization screw **30**.

WARNING

Risk of injury if the 7.0 mm blue polyaxial screw or 10 mm blue polyaxial screw **XL is pulled out of the vertebral body through incorrect operation or non-observance of the correct sequence of operating steps!**

- For easy detachment from the stabilization screw, apply the cannulated screwdriver in an ortho-grade direction to the stabilization screw.
- Before removing the cannulated screwdriver from the stabilization screw, unscrew the screw revision instrument from the stabilization screw.

Aesculap®

MACS II

- Pull the aiming sleeve for the stabilization screw **17, 18** off the centralizer **27a**.

Risk of injury if the aiming sleeve for the stabilization screw becomes jammed in the centralizer as a result of severe angular misalignment between the approach direction and centralizer axis!

- For the controlled removal of a jammed aiming sleeve for a stabilization screw, unscrew it with a 6.5 mm yellow stabilization screw at least 40 mm long.

Risk of injury and/or damage to the implant components through non-observance of the correct sequence of operating steps!

- Only fit the torque wrench for retaining bolt onto the screwdriver for retaining bolt when this is fully inserted in the retaining bolt.
- Do not slide the screwdriver for retaining bolt in the guide tube together with the torque wrench for retaining bolt, and do not let it slide down too quickly.

Locking the polyaxiality

- Snap the guide tube **7** with the counter torque **8** onto the centralizer **27a** in situ.

If you cannot snap the guide tube **7** onto the centralizer **27a**:

- To snap the guide tube onto the centralizer: Remove the screwdriver **6, 13** from the guide tube **7**.
- or -
- Guide the guide tube **7** onto the centralizer **27a** in an orthograde direction.
- Insert the screwdriver for retaining bolt **13** into the guide tube **7**.
- Attach the torque wrench for retaining bolt **21** to the screwdriver for retaining bolt **13**.

- Tighten the retaining bolt **27b**.

The retaining bolt is tightened down when two clicks are heard.

WARNING

Risk of injury if the retaining bolt is insufficiently tightened!

- ▶ Do not remove instruments from the retaining bolt after the first click.
- ▶ Only remove instruments from the retaining bolt after the second click has been heard.
- ▶ Ensure that the centralizer does not slip off the clamping element or the 10 mm blue polyaxial screws XL after the first click.
- ▶ Use only the torque wrench for retaining bolt to tighten the retaining bolt.
- ▶ Have the torque wrench for the retaining bolt recalibrated by the manufacturer at intervals no longer than every two years.
- ▶ Ensure that the screwdriver for retaining bolt is fully engaged in the retaining bolt to prevent it slipping off.
- ▶ Ensure that the polished groove on the screwdriver for retaining bolt is fully covered by the guide tube.

WARNING

Damage to the implant caused by overtightening of the retaining bolt!

- ▶ Tighten the retaining bolt to the torque specified for this purpose.
- ▶ Have the torque wrench for the retaining bolt recalibrated by the manufacturer at intervals no longer than every two years.
- ▶ Use only the torque wrench for retaining bolt to tighten the retaining bolt.

WARNING

Risk of injury if no or insufficient counterforce is applied!

- ▶ To apply a counterforce, always use the counter holder provided, ensuring that as little force as possible is transferred to the spinal column.

WARNING

Risk of injury if the centralizer is left in the body!

- ▶ If the torque wrench for retaining bolt is not available, the retaining bolt can be secured using the T-handle.

Aesculap®

MACS II

- ▶ Remove instruments **7, 8, 13, 21**.
- ▶ Pull the torque wrench for retaining bolt **21** with the screwdriver for retaining bolt **13** out of the guide tube **7**.
- ▶ Remove the centralizer **27a** from the guide tube **7**.

WARNING

Risk of cross-contamination from inadequate reprocessing and non-observance of instructions for use!

- ▶ **Before reprocessing, remove the centralizer from the guide tube.**
- ▶ **To detach the centralizer from the guide tube, remove the screwdriver for retaining bolt from the guide tube.**

Revision

- ▶ Insert the socket wrench for nut **11** into the locknut **28**.
- ▶ Insert the T-handle **22** into the screwdriver for retaining bolt **13**.
- ▶ Insert the screwdriver for retaining bolt **13** with the T-handle **22** into the retaining bolt **27b**.
- ▶ Unscrew the retaining bolt **27b** using the screwdriver for retaining bolt with T-handle **22**.
- ▶ Apply counterforce over the socket wrench for nut **11**.
- ▶ When using clamping element **TL 26** or the 10 mm blue polyaxial screws **XL 32**, insert the cannulated screwdriver **6** into the 6.5 mm yellow stabilization screw **30**.
- ▶ When using clamping element **TL 26** or the 10 mm blue polyaxial screws **XL 32**, insert the screw revision instrument **15** into the cannulated screwdriver **6** and screw it into the 6.5 mm yellow stabilization screw **30**.
- ▶ When using clamping element **T 25**, insert the screwdriver for purple stabilization screw **20** into the 4.5 mm purple stabilization screw **29** and tighten the knurled nut clockwise, see TA011177 instructions for use.
- ▶ Unscrew the stabilization screw **29, 30** while applying tensile force. This will overcome the secure fastening of the screw.
- ▶ Unscrew the 7.0 mm blue polyaxial screw **31** or the 10 mm blue polyaxial screw **XL 32** from the vertebral body using the ball-end screwdriver **9**.
- ▶ Remove the entire implant construction using the holding forceps for rods/plates **10**.

Intraoperative disassembly of pre-assembled retaining bolt

- ▶ Snap the guide tube **7** with the counter torque **8** onto the centralizer **27a** in situ.

If you cannot snap the guide tube **7** onto the centralizer **27a**:

- ▶ To snap the guide tube onto the centralizer: Remove the screwdriver **6, 13** from the guide tube **7**.
– or –
- ▶ Guide the guide tube **7** onto the centralizer **27a** in an orthograde direction.
- ▶ Insert the screwdriver for retaining bolt **13** into the guide tube **7**.
- ▶ Insert the screwdriver for retaining bolt **13** into the retaining bolt **27b**.
- ▶ Insert the T-handle **22** into the screwdriver for retaining bolt **13**.

WARNING

Risk of injury and/or damage to the implant components through non-observance of the correct sequence of operating steps!

- ▶ Only fit the T-handle onto the screwdriver for retaining bolt when this is fully inserted in the retaining bolt.
- ▶ Do not let the screwdriver for retaining bolt slide down too quickly in the guide tube together with the T-handle.

- ▶ Unscrew the pre-assembled retaining bolt **27b** from the clamping element **25, 26** or the 10 mm blue polyaxial screws XL **32** counterclockwise.

Intraoperative mounting of pre-assembled retaining bolt

- ▶ Screw the pre-assembled retaining bolt **27** with the mounted guide tube **7**, counter torque **8**, and screwdriver for retaining bolt **13** onto the clamping element **25, 26** or the 10 mm blue polyaxial screw XL **32** by hand.

WARNING

The pre-assembled retaining bolt may break off if instructions for use are not observed!

- ▶ Never use the T-handle or torque wrench for retaining bolt to screw on the pre-assembled retaining bolt.

WARNING

Risk of injury if the necessary counter torque is not applied when tightening the locknut or retaining bolt!

- ▶ Ensure that the retaining lugs of the centralizer are fully engaged in the appropriate grooves in the clamping element or the 10 mm blue polyaxial screw XL.

WARNING

Risk of injury if the retaining bolt has insufficient clamping stability!

- ▶ Apply the pre-assembled retaining bolt in an ortho-grade direction.
- ▶ Ensure that the pre-assembled retaining bolt is screwed on smoothly, otherwise the external thread of the pre-assembled retaining bolt and/or the internal thread of the clamping element or the 10 mm blue polyaxial screw XL can be damaged.
- ▶ If the internal thread of the clamping element or the 10 mm blue polyaxial screw XL is damaged, replace the part.

Aesculap®

MACS II

- ▶ Insert the torque wrench for mounting the retaining bolt 1 into the screwdriver for retaining bolt 13.
- ▶ Tighten the torque wrench for mounting the retaining bolt 1 with 1.8 Nm.

WARNING

Risk of injury if pre-assembled retaining bolt comes loose during the procedure!

- ▶ To mount the pre-assembled retaining bolt, use the torque wrench for mounting the retaining bolt and tighten with 1.8 Nm.
- ▶ Prior to applying the torque wrench, make certain the display arrow points to zero.
- ▶ Replace the torque wrench if it is damaged.

- ▶ Remove the torque wrench for mounting the retaining bolt 1 from the screwdriver for retaining bolt 13.
- ▶ Remove the screwdriver for retaining bolt 13 from the guide tube 7.
- ▶ Remove the guide tube 7, including the counter torque 8, from the centralizer 27a.

WARNING

Risk of injury if the 7.0 mm blue polyaxial screw or 10 mm blue polyaxial screw XL is pulled out of the vertebral body through incorrect operation or non-observance of the correct sequence of operating steps!

- ▶ For easy removal from the centralizer, ensure that the guide tube with counter torque is positioned ortho-grade to the centralizer.
- ▶ First remove the screwdriver for retaining bolt from the guide tube before detaching the guide tube together with counter torque from the centralizer.
- ▶ Apply pulling force over the head part of the counter torque.

Disassembling

Screwdriver for stabilization screw, purple, FW342R

Note

Assembling and disassembling are described in the TA011177 instructions for use.

K-wire insertion instrument FW408R

Note

To disassemble the product, follow the steps for assembling in reverse order.

Counter torque FW399R

Note

To disassemble the product, follow the steps for assembling in reverse order.

Torque wrench for nut FW416R

Note

To disassemble the product, follow the steps for assembling in reverse order.

Assembling

Screwdriver for stabilization screw, purple, FW342R

Note

Assembling and disassembling are described in the TA011177 instructions for use.

K-wire insertion instrument FW408R

- Screw in the handle rod 5b.
- Insert the impactor 5c into the aiming tube 5a, see Fig. 1.


Fig. 1

Counter torque FW399R

- Screw the cleaning rod 8b into the counter torque 8a, see Fig. 2.


Fig. 2

Torque wrench for nut FW416R

- Tilt wrench 12c into the casing, until it clicks into position.
- Push on the handle 12a in the opposite direction to the arrow and screw on, see Fig. 3.


Fig. 3

Aesculap®

MACS II

Validated reprocessing procedure

General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Single-use products

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

► Do not reprocess the product!

Applies to:

Art. no.	Designation
SX804T	Pre-assembled retaining bolt (retaining bolt, centralizer)
FW406S	K-wire set, sterile (K-wire x2, removal instrument x1)

- Before reprocessing the guide tube **7**, make sure that the centralizer **27a** has been removed from the working end and discarded.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- ▶ Disassemble the product immediately after use, as described in the respective instructions for use.
- ▶ Open up instruments with hinges.

Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.

- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- ▶ Dismantle the product prior to cleaning, see Disassembling.
- ▶ Disassemble the screwdriver for the purple stabilization screw **20** prior to cleaning, see TA011177.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for (e.g. aluminum, plastics, high-grade steel),
 - do not attack softeners (e.g. in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 60 °C.

- ▶ Do not use oxidizing chemicals (e.g. H₂O₂), which could cause bleaching/layer loss of the product.

Aesculap®

MACS II

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FW013R ■ FW326R ■ FW336R ■ FW342R ■ FW394R ■ FW396R-FW399R ■ FW405R ■ FW407R-FW409R ■ FW411R-FW413R ■ FW417R ■ FW419R ■ FW429	■ Cleaning brush: – e.g., TA007747 for FW394R, FW397R, FW409R, FW413R, FW417R, FW429R – e.g., F0764200 for FW394 – e.g., TE654202 for FW408R – e.g., GK469200 for FW398R – e.g., TA011944 for FW399R ■ Disposable syringe, 20 ml ■ Keep working ends open for cleaning. ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ FW395R ■ FW416R	■ Cleaning brush: – e.g., TA007747 for FW416R – e.g., TA006874 for FW395R ■ Disposable syringe, 20 ml ■ Keep working ends open for cleaning. ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with ultrasound and immersion disinfection

Validated procedure	Specific requirements	Reference
Mechanical alkaline cleaning and thermal disinfection ■ FW013R ■ FW396R ■ FW405R ■ FW407R ■ FW411R ■ FW419R	■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfecting and sub-chapter: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW326R ■ FW336R ■ FW342R ■ FW394R-FW395R ■ FW397R-FW399R ■ FW408R-FW409R ■ FW412R-FW413R ■ FW417R ■ FW429	■ Cleaning brush: – e.g., TA007747 for FW397R, FW413R, FW417R, FW419R, FW429 – e.g., F0764200 for FW394R – e.g., TE654202 for FW408R – e.g., GK469200 for FW398R – e.g., TA011944 for FW399R – e.g., TA006874 for FW395R ■ Disposable syringe, 20 ml ■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Aesculap®

MACS II

Validated procedure	Specific requirements	Reference
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection FW416R	- Reinigungsbürste: z. B. TA007747 - Disposable syringe, 20 ml - Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots). - Connect components with lumens and channels directly to the rinsing port of the injector carriage. - Keep working ends open for cleaning. - Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: - Chapter Manual pre-cleaning with ultrasound and brush - Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed

Aesculap®

MACS II

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

Aesculap®

MACS II

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfecting

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Aesculap®

MACS II

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Aesculap®

MACS II

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- **Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).**

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note

The guide tube 7 must be disconnected from the counter torque 8 prior to sterilization.

Note

The K-wire insertion instrument 5, counter torque 8, and torque wrench for nut 12 must only be sterilized after they have been dismantled.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Disassemble the instrument
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Aesculap®

MACS II

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom)			
Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/ 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Maintenance

The torque wrench FW419R is factory-calibrated to 1.8 Nm.

If the display arrow does not point to zero with the torque wrench under no load, the wrench is defective and must not be used under any circumstances.

- Replace the damaged torque wrench.
- Have the FW412R and FW416R torque wrench recalibrated by the manufacturer at intervals no longer than every two years.

Technical Service

Risk of injury and/or malfunction!

► **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems LLC

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042

Aesculap Repair Hotline

Phone: +1 (800) 214-3392

Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Incineration of the product, its components and its packaging produces harmless residues.

- Comply with national regulations for incineration under all circumstances!

To avoid the release of toxic incineration products, the incineration of the product, its components and its packaging must be carried out in an appropriate incinerator (>1 000 °C).

- Comply with national regulations for incineration under all circumstances!

Distributor in the US/Contact in Canada for product information and complaints

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany