AESCULAP®


- D** **Gebrauchsanweisung**
MACST^{TL}-Polyaxialschraube XL
- GB** **Instructions for use**
- USA** **MACST^{TL} polyaxial screw XL**
- F** **Mode d'emploi**
Vis polyaxiale MACST^{TL} XL
- E** **Instrucciones de manejo**
Tornillo poliaxial MACST^{TL} XL
- I** **Istruzioni per l'uso**
Vite poliassiale MACST^{TL} XL
- P** **Instruções de utilização**
Parafuso poliaxial MACST^{TL} XL
- NL** **Gebruiksaanwijzing**
MACST^{TL}-polyaxiale schroef XL

AESCULAP AG & CO. KG
Am Aesculap-Platz
78532 Tuttlingen/Germany

Phone +49 (74 61) 95-0
Fax +49 (74 61) 95-26 00
E-mail information@aesculap.de

B | BRAUN

CE 0123

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– Not applicable for the USA – See version for the USA –

Intended use

The MACST^{TL} polyaxial screw XL, in combination with the other MACST^{TL} implants of the MACST^{TL} system, is used for the ventral mono- and multisegmental stabilization of the lumbar and thoracic spine.

The MACST^{TL} system comprises:

- Anchoring screws for the vertebral bodies
- Tensioning elements
- Fixation nuts and clamping screws
- Connecting plates and, if required, connecting rods of the Aesculap Spine System Evolution, diameter 5.2 mm, as connection elements

Materials

The materials used in the implant are listed on the package.

- ISOTAN^{®F} titanium forged alloy Ti6Al4V according to ISO 5832-3

MACST^{TL}[®] and ISOTAN[®] are registered trademarks of AESCULAP AG & CO. KG, 78532 Tuttlingen/Germany.

System configurations

- Polyaxial twin-screw system with plate or 2 rods (see Fig.)

– Not applicable for the USA – See version for the USA –

Indications

Surgically installed implants serve to support normal healing processes. They are not supposed to replace normal body structures or to support permanently loads occurring in cases where healing did not take place.

Use for the following indications:

- Fractures
- Spinal tumor
- Degenerative disc disorder (spondylolisthesis, spondylolysis, spinal stenosis)
- Post-trauma instability

Contraindications

Do not apply in the presence of:

- Infection
 - systemic
 - in the spine
 - local
- Pregnancy

– Not applicable for the USA – See version for the USA –

- Acute osteopenia

In cases of osteopenia, the operating surgeon has to weigh the risks that may endanger the success of a therapy with the MACST^{TL} system.

If it must be expected that the standard polyaxial screw does not provide sufficient stability, the polyaxial screw XL can be used. If even more anchoring stability is required, the polyaxial screw XL can be fixated with an Aesculap cement suitable for screw augmentation. It is absolutely imperative to observe the safety advice concerning the Aesculap cement application system.

- Medical or surgical conditions that could negatively impact the success of the implantation
- Foreign body sensitivity to the implant materials
- Inadequate patient compliance
- Fever
- Cases not listed under indications

– Not applicable for the USA – See version for the USA –

Side-effects or adverse interactions

- Implant failure resulting from excess load
 - bending
 - loosening
 - breakage
- Inadequate fixation
- Absence of, or delayed, bony union
- Cement leakages (cement emerging from the vertebral body)
- Infection
- Fractured vertebral body or bodies
- Injuries to
 - nerve roots
 - spinal cord
 - vessels
 - organs

– Not applicable for the USA – See version for the USA –

Safety information

This system is not licensed to be screwed onto, or fastened to, the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

- It is the operating surgeon's responsibility to ensure that implant components are applied correctly during operative procedures.
- General risk factors associated with surgical procedures are not described in the present instructions for use.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Under no circumstances may modular implant components from different suppliers be combined.

– Not applicable for the USA – See version for the USA –

- Under no circumstances may damaged components or surgically excised components be used.
- Implants once used must not be reused.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keep the patient well informed.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.
In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures.
- Delayed healing can cause implant components to fracture as a result of metal fatigue.
- The physician in charge decides whether an implant component implanted should be removed, taking into account the potential risks of another operation and the difficulties involved in the removal of an implant.

– Not applicable for the USA – See version for the USA –

Sterility

- The implant components are delivered unsterile.
- The implant components are delivered in individual packaging.
- Store the implant components in their original packaging and only remove them from such packaging immediately before use.
- Pre-clean and steam-sterilize the implants (according to the respective hospital guidelines for the provision of sterile materials) before use.

Sterilization method and parameters

- Use the implant system storage unit for sterilization and sterile preparation.
- Sterilize with steam, taking note of the following:
The sterilization has to be done according to a validated steam sterilization procedure (e.g. in a sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, and validated in conformance with EN 554/ISO 13683). In case of application of the fractionated vacuum procedure the sterilization has to be carried out for a minimum of 5 minutes at 134 °C and at 2 bar pressure.

– Not applicable for the USA – See version for the USA –

Application

The MACST^{TL} system is described in detail in the O.R. manual. Detailed user instructions can be found in the product documentation.

The O.R. manual can be ordered from the manufacturer at any time.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bones
- Location of intraoperative landmarks

Caution:

For a secure implantation of the anchoring screws, careful measurement of the vertebral body's diameter at the planned screw position and direction is crucial. A screw length smaller than the measured diameter has to be chosen.

The following conditions must be fulfilled prior to application:

- All requisite implant components must be ready to hand
- Operating conditions must be highly aseptic
- Implantation instruments including the special Aesculap implant system instruments are complete and in working condition
- The special instrument set for implanting the elements and for the distraction, compression and reposition of the lumbar and thoracic spine is ready to hand

– Not applicable for the USA – See version for the USA –

- The operating surgeon and operating room team must be thoroughly conversant with the operating technique, as well as the range of implants and instruments to be applied; complete information on these subjects must be readily available at the workplace
- Those performing operations must be thoroughly conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant articles by medical specialists from the professional literature
- Information was obtained from the manufacturer in case of an ambiguous preoperative situation and if implants are present in the area to be provided for

The operative procedure has been explained to the patient, and the latter's understanding of the following information has been documented:

- If the bone union is delayed or does not take place, the implants may break or come loose as a result of the considerable strains arising in such cases.
- The life-span of the implant is determined by body weight.
- The implant components must not be subjected to overload due to extreme strains or through work-related or sports activities.
- Revision surgery may become necessary in cases of an implant coming loose, implant fracture or correction loss.
- The patient must undergo regular check-ups of the implant components, performed by a physician.

– Not applicable for the USA – See version for the USA –

The implantation of the MACS^{TL} system has to be carried out in the following steps:

- Select the appropriate MACS^{TL} system variant and anchoring screws according to indication, preoperative planning and bone situation found intraoperatively.
- To determine the optimum implant position of the MACS^{TL} system:
 - Preposition the implant with Kirschner wires and
 - prepare the bone correctly, using cannulated instruments.
- Apply the implantation instruments in the correct way.
- To avoid internal stress on, and weakening of, the implant: Avoid scoring or scratching any implant components.
- Do not bend or warp any metal implants with the exception of the rods of the Aesculap Spine System Evolution.
- Do not bend or re-bend the rods excessively.
- For bending the rods, do not use any other tool than the bending instruments that are part of the MACS^{TL} instrument set.
- Mount the centering sleeve on the clamping element, using the special mounting instrument.
- Insert the implant components, using the special insertion instruments.
- Do not tighten or screw in the polyaxial screw XL completely, so that the clamping element can still be moved.
- Prior to unlocking the insertion instrument, remove the Kirschner wire through the cannulated insertion element, using the appropriate removal instrument.

– Not applicable for the USA – See version for the USA –

- Remove the insertion instrument.
- Ensure that the centering sleeves remain mounted, temporarily, on the clamping elements.
- Distract and insert placeholder, if necessary.

WARNING

Risk of injury caused by screws loosening or tearing out during reposition, distraction or compression, especially in cases of poor bone quality!

- Apply force carefully and well-dosed.
 - Maintain permanent visual control (direct or through imaging techniques).
- Insert stabilizing plate or rods.
- Implant stabilizing plate and fixation nut in such a way that the inscribed side (top side) is visible for the user.

WARNING

Risk of injury caused by nuts with insufficient clamping stability!

- Set nuts in place correctly.
 - Use a torque wrench to tighten the nuts at 15 Nm.
- Remove centering devices.
- Tighten the polyaxial screws XL evenly.
- Insert the stabilizing screw through the appropriate guiding sleeve.

– Not applicable for the USA – See version for the USA –

- If necessary, cement the polyaxial screw XL with an Aesculap cement suitable for screw augmentation.
- Always follow the safety advice for the Aesculap cement application system.

WARNING

Risk of injury from cement leakages!

- Never prepare and position the implants without the aid of imaging processes.
- Make certain that the slits of the polyaxial screw XL are fully seated in the bone.
- Do not penetrate the opposite corticalis, and verify that the opposite corticalis is not damaged (e.g. by a previously installed, longer screw).
- Do not cement the polyaxial screw XL in a fractured or partly resected vertebral body. Carry out preoperative thin-section CT control.
- Apply cement only under high-quality radiographic control. Max. amount of cement per polyaxial screw XL: 2 ml.

– Not applicable for the USA – See version for the USA –

- Insert clamping screw.

WARNING

Risk of injury caused by clamping screws that cannot be inserted or tightened because the cement is already hardening!

- Put in place and tighten the clamping screws immediately after cementing the polyaxial screw XL.

WARNING

Risk of injury caused by insufficiently tightened clamping screws of the polyaxial elements!

- Set clamping screws in place correctly.
- Use a torque wrench to tighten the clamping screws at 10 Nm.

Note

Once the polyaxial screw XL has been cemented, the explantation of the screws can become extremely difficult, especially in cases of osteopenia.

Further information on Aesculap implant systems is available from Aesculap or the Aesculap office responsible.

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Intended use

This anterolateral/anterior system consists of several vertebral screws, locking nuts, spine plates and rods. The points of attachment are screw fixation to the anterolateral vertebral bodies of the lumbar and thoracic spine (T1 to L5). This system is intended to provide stabilization during development of a solid spinal fusion. When used as an anterolateral/anterior spine plate MACS^{TL} Modular Anterior Construct System is indicated for patients with:

- Degenerative disc disease defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies
- Spondylolisthesis
- Spondylolysis
- Fracture
- Spinal stenosis
- Deformities (i.e., scoliosis, kyphosis, lordosis, whether neuromuscular or related to deficient posterior elements)
- Tumors (neoplastic disease)
- Pseudarthrosis
- Revision of previous surgery

Materials

The materials used in the implant are listed on the package.

- ISOTAN®_F titanium forged alloy Ti6Al4V according to ISO 5832-3

MACSTL® and ISOTAN® are registered trademarks of AESCULAP AG & CO. KG, 78532 Tuttlingen/Germany.

System configurations

- Polyaxial twin-screw system with plate or 2 rods (see Fig.)

Indications

Surgically installed implants serve to support normal healing processes. They are not supposed to replace normal body structures or to support permanently loads occurring in cases where healing did not take place.

Use for the following indications:

- Fractures
- Spinal tumor
- Degenerative disc disorder (spondylolisthesis, spondylolysis, spinal stenosis)
- Post-trauma instability

Contraindications

Do not apply in the presence of:

- Infection
 - systemic
 - in the spine
 - local

- Pregnancy

- Acute osteopenia

In cases of osteopenia, the operating surgeon has to weigh the risks that may endanger the success of a therapy with the MACST^{TL} system.

If it must be expected that the standard polyaxial screw does not provide sufficient stability, the polyaxial screw XL can be used.

- Medical or surgical conditions that could negatively impact the success of the implantation
- Foreign body sensitivity to the implant materials
- Inadequate patient compliance
- Fever
- Cases not listed under indications

Side-effects or adverse interactions

- Implant failure resulting from excess load
 - bending
 - loosening
 - breakage
- Inadequate fixation
- Absence of, or delayed, bony union
- Infection
- Fractured vertebral body or bodies
- Injuries to
 - nerve roots
 - spinal cord
 - vessels
 - organs

Safety information

This system is not licensed to be screwed onto, or fastened to, the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

- It is the operating surgeon's responsibility to ensure that implant components are applied correctly during operative procedures.
- General risk factors associated with surgical procedures are not described in the present instructions for use.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Under no circumstances may modular implant components from different suppliers be combined.

- Under no circumstances may damaged components or surgically excised components be used.
- Implants once used must not be reused.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keep the patient well informed.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.

In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures.

- Delayed healing can cause implant components to fracture as a result of metal fatigue.
- The physician in charge decides whether an implant component implanted should be removed, taking into account the potential risks of another operation and the difficulties involved in the removal of an implant.

Sterility

- The implant components are delivered unsterile.
- The implant components are delivered in individual packaging.
- Store the implant components in their original packaging and only remove them from such packaging immediately before use.
- Pre-clean and steam-sterilize the implants (according to the respective hospital guidelines for the provision of sterile materials) before use.

Sterilization method and parameters

- Use the implant system storage unit for sterilization and sterile preparation.
- Sterilize with steam, taking note of the following:
The sterilization has to be done according to a validated steam sterilization procedure (e.g. in a sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, and validated in conformance with EN 554/ISO 13683). In case of application of the fractionated vacuum procedure the sterilization has to be carried out for a minimum of 5 minutes at 134 °C and at 2 bar pressure.

Sterilization for the US market:

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270–275 °F	4 min	4 min

Application

The MACST^{TL} system is described in detail in the O.R. manual. Detailed user instructions can be found in the product documentation.

The O.R. manual can be ordered from the manufacturer at any time.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bones
- Location of intraoperative landmarks

Caution:

For a secure implantation of the anchoring screws, careful measurement of the vertebral body's diameter at the planned screw position and direction is crucial. A screw length smaller than the measured diameter has to be chosen.

The following conditions must be fulfilled prior to application:

- All requisite implant components must be ready to hand
- Operating conditions must be highly aseptic
- Implantation instruments including the special Aesculap implant system instruments are complete and in working condition
- The special instrument set for implanting the elements and for the distraction, compression and reposition of the lumbar and thoracic spine is ready to hand

- The operating surgeon and operating room team must be thoroughly conversant with the operating technique, as well as the range of implants and instruments to be applied; complete information on these subjects must be readily available at the workplace
- Those performing operations must be thoroughly conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant articles by medical specialists from the professional literature
- Information was obtained from the manufacturer in case of an ambiguous preoperative situation and if implants are present in the area to be provided for

The operative procedure has been explained to the patient, and the latter's understanding of the following information has been documented:

- If the bone union is delayed or does not take place, the implants may break or come loose as a result of the considerable strains arising in such cases.
- The life-span of the implant is determined by body weight.
- The implant components must not be subjected to overload due to extreme strains or through work-related or sports activities.
- Revision surgery may become necessary in cases of an implant coming loose, implant fracture or correction loss.
- The patient must undergo regular check-ups of the implant components, performed by a physician.

The implantation of the MACS^{TL} system has to be carried out in the following steps:

- Select the appropriate MACS^{TL} system variant and anchoring screws according to indication, preoperative planning and bone situation found intraoperatively.
- To determine the optimum implant position of the MACS^{TL} system:
 - Preposition the implant with Kirschner wires and
 - prepare the bone correctly, using cannulated instruments.
- Apply the implantation instruments in the correct way.
- To avoid internal stress on, and weakening of, the implant: Avoid scoring or scratching any implant components.
- Do not bend or warp any metal implants with the exception of the rods of the Aesculap Spine System Evolution.
- Do not bend or re-bend the rods excessively.
- For bending the rods, do not use any other tool than the bending instruments that are part of the MACS^{TL} instrument set.
- Mount the centering sleeve on the clamping element, using the special mounting instrument.
- Insert the implant components, using the special insertion instruments.
- Do not tighten or screw in the polyaxial screw XL completely, so that the clamping element can still be moved.
- Prior to unlocking the insertion instrument, remove the Kirschner wire through the cannulated insertion element, using the appropriate removal instrument.

- Remove the insertion instrument.
- Ensure that the centering sleeves remain mounted, temporarily, on the clamping elements.
- Distract and insert placeholder, if necessary.

**WARNING**

Risk of injury caused by screws loosening or tearing out during reposition, distraction or compression, especially in cases of poor bone quality!

- Apply force carefully and well-dosed.
 - Maintain permanent visual control (direct or through imaging techniques).
- Insert stabilizing plate or rods.
- Implant stabilizing plate and fixation nut in such a way that the inscribed side (top side) is visible for the user.

**WARNING**

Risk of injury caused by nuts with insufficient clamping stability!

- Set nuts in place correctly.
 - Use a torque wrench to tighten the nuts at 15 Nm.
- Remove centering devices.
- Tighten the polyaxial screws XL evenly.
- Insert the stabilizing screw through the appropriate guiding sleeve.
- Insert clamping screw.

**WARNING**

Risk of injury caused by insufficiently tightened clamping screws of the polyaxial elements!

- Set clamping screws in place correctly.
- Use a torque wrench to tighten the clamping screws at 10 Nm.

Further information on Aesculap implant systems is available from Aesculap or the Aesculap office responsible.

Distributor in the US/Contact in Canada for product information and complaints

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