GB  
USA

**Instructions for use**  
Micro-instruments

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

Legend

- 1 Locking button
- 2 Handle stems
- 3 Jaw part mounting (on handle)
- 4 Jaw parts with notched rod

Symbols on product and packages

| Symbol                                                                           | Explanation                                                                             |
|----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
|  | Caution: General warning symbol<br>Caution: See documentation supplied with the product |

Intended use

The micro-instruments are intended for use in neurosurgery, and are specially designed for use in tight spaces. They can be used for cutting, grasping, dissecting, etc. of blood vessels and tissues.

Available sizes

The micro-instruments are available as complete instruments or as separate, combinable components. Various jaw parts (e.g. straight scissors, curved scissors or forceps) can be used with the appropriate handles.

*Note*  
The current product range is described in Aesculap brochure C 74511.

Safe handling and preparation

- CAUTION
- Federal law restricts this device to sale by or on order of a physician!
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

- 

Risk of injury and/or malfunction!  
➤ Always carry out a function check before using the product.
- 

Risk of injury when applying the product outside the field of view!  
➤ Apply the product only under visual control.

Disassembling

- Turn the locking button **1** by 180° so that it engages.
- Hold the handle stems **2** with one hand and the jaw parts with notched rod **4** with the other hand.
- Pull to separate the jaw part mounting **3** and the jaw parts with notched rod **4**.
- Turn the locking button **1** back to its original position.

Assembling

- Ensure that the link of the jaw parts with notched rod **4** and the groove of the jaw part mounting **3** are in alignment when the handle stems **2** are closed.
- Push the jaw part mounting **3** and jaw parts with notched rod **4** together until the notched rod engages with an audible click.
- Function check: Press the handle stems **2** together, all the way to the stop.  
Check that the jaw parts operate properly.

Validated processing procedure

*Note*  
Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

*Note*  
Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

*Note*  
Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

*Note*  
Up-to-date information on processing can be found on the Aesculap Extranet at [www.aesculap-extra.net](http://www.aesculap-extra.net)

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Therefore the time interval between application and cleaning should not exceed 6 h, and neither potentially fixating pre-cleaning temperatures >45 °C nor any fixating disinfectants (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfecting and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

• Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process solution can already cause visible surface changes.

➤ Do not use process chemicals that would cause stress cracking or brittleness of plastics.

For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see [www.a-k-i.org](http://www.a-k-i.org)

➤ Use suitable cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open jointed products.
- Rinse surfaces that cannot be visually inspected, e.g. on products with hidden gaps or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.

Preparation prior to cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/Disinfecting

- 

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!  
➤ Use cleaning and disinfecting agents according to the manufacturer's instructions. The cleaning and disinfecting agents must  
– be approved e.g. for aluminum, plastic materials and high-grade steel,  
– not attack softeners (e.g. silicone).  
➤ Observe specifications regarding concentration, temperature and exposure time.  
➤ Do not exceed the maximum allowable cleaning temperature of 55 °C.
- Carry out ultrasound cleaning:  
– as an effective mechanical supplement to manual cleaning/disinfecting.  
– as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/ disinfecting.  
– as an integrated mechanical support measure for mechanical cleaning/disinfecting.  
– for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically, provided they can be securely fixed in machines or storage devices in such a way that they will be thoroughly cleaned.

Manual cleaning/disinfecting

- Keep working ends open during cleaning.
- Clean products with movable links with the links opened or mobilized.
- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfection

| Stage | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                 |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------|
| I     | Cleaning              | RT<br>(cold) | 15         | 2            | D–W              | BBraun Stabimed; aldehyde-free, phenol-free and QAV-free |
| II    | Intermediate<br>rinse | RT<br>(cold) | 1          | –            | D–W              | –                                                        |
| III   | Disinfecting          | RT<br>(cold) | 15         | 2            | D–W              | BBraun Stabimed; aldehyde-free, phenol-free and QAV-free |
| IV    | Final rinse           | RT<br>(cold) | 1          | –            | FD–W             | –                                                        |
| V     | Drying                | RT           | –          | –            | –                | –                                                        |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized)  
RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product under running tap water, using a suitable cleaning brush if necessary, until all visible residues have been removed from the surfaces.
- Brush through surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush (TA no. FM32800).
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, thoroughly (at least five times) rinse through these components with the cleaning/disinfecting solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.

Stage IV

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting with manual pre-cleaning

*Note*  
The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

*Note*  
For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that Ao is >3 000 for the process.

*Note*  
The disinfectant used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with brush

| Stage | Step                     | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                    |
|-------|--------------------------|--------------|------------|--------------|------------------|-------------------------------------------------------------|
| I     | Disinfecting<br>cleaning | RT<br>(cold) | 15         | 2            | D–W              | BBraun Stabimed; aldehyde-free,<br>phenol-free and QAV-free |
| II    | Irrigation               | RT<br>(cold) | 1          | -            | D–W              | -                                                           |

D–W: Drinking water  
RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush (TA no. FM32800).
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfectant without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Keep working ends open during cleaning.
- Ensure that product links and joints are opened when placing the product on the tray.

| Stage | Step                    | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical/Note                                                                     |
|-------|-------------------------|--------------|------------|------------------|-----------------------------------------------------------------------------------|
| I     | Prerinse                | <25/77       | 3          | D–W              | -                                                                                 |
| II    | Cleaning                | 55/131       | 10         | FD–W             | BBRAUN HELIMATIC CLEANER alkaline<br>with tensides; application solution<br>0.5 % |
| III   | Intermediate rinse      | >10/50       | 1          | FD–W             | -                                                                                 |
| IV    | Thermal<br>disinfecting | 90/194       | 5          | FD–W             | -                                                                                 |
| V     | Drying                  | -            | -          | -                | According to disinfectant program                                                 |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized)

Inspection, maintenance and checks

**CAUTION**

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Assemble the separable product, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product into the appropriate storage holders or on an appropriate tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

*Note*  
The product can be sterilized either in disassembled or in assembled condition.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
  - Steam sterilization through fractionated vacuum process
  - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665
  - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
  - Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.
- To achieve a sterility assurance level of 10<sup>-6</sup>, Aesculap recommends the following parameters:

| Aesculap Orga Tray/sterile container (perforated bottom)<br>Minimum cycle parameters* |               |       |                     |
|---------------------------------------------------------------------------------------|---------------|-------|---------------------|
| Sterilization method                                                                  | Temp.         | Time  | Minimum drying time |
| Pre-vacuum                                                                            | 270 °F–275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Technical Service

**WARNING**

**Risk of injury and/or malfunction!**

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 7461 95-1602  
Fax: +49 7461 16-5621  
E-Mail: ats@aesculap.de

Or in the US:  
Aesculap Inc.  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood, MO 63042  
Aesculap Repair Hotline  
Phone: +1 800 214-3392  
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034  
USA

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