GB

Instructions for use/Technical description

MINOP®/PaediScope single-use introducer for neuro-endoscopy

USA

Instructions for use/Technical description

MINOP®/PaediScope single-use introducer for neuro-endoscopy

Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
Aesculap®
MINOP®/PaediScope single-use introducer for neuro-endoscopy

Legend

- 1 Obturator
- 2 Sheath
- 3 Obturator and sheath assembled and interlocked
- 4 Interlock flange
- 5 Release locking knob
- 6 Sheath handles
- 7 Set screws
- 8 Handles

Symbols on product and packages

	Sterilization using ethylene oxide
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

Intended use

The single-use introducer serves the creation and preservation of a temporary access canal to the patient's ventricular system. For the adjustment to different entry depths, the single-use introduce is constructed in such a way that the sheath covering can be cleaved axially and the sheath halves can be hinged away right up to the cranial surface. The single-use introducer consists of two main components: Obturator and sheath, which are interlocked.

Indications

Indications, see Intended use.

Note
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Contraindications

None known.

Risks and side effects

As part of the legal duty to inform, the following typical risks and side effects associated with the use of surgical instruments are referred to. These are predominantly procedure-specific, non product-specific and not limited to unwanted damage to surrounding tissue resulting in e.g. bleeding, infection, material incompatibilities or instrument parts remaining unnoticed in the patient, etc.

Available sizes

Art. no. Single-use introducer	Color of sheath handle	Size (interior part of sheath Ø)	For use in conjunction with (typical combination)
FH603SU	Purple	10 French (3.2 mm)	PaediScope PF010A
FH604SU	Pink	19 French (6.1 mm)	MINOP® 6 mm trocar FF399R

Note
These products are LATEX-FREE.

Safe handling and preparation

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product!

- The product is EO sterilized and sterile packed.
The product must not be reused.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
 - Read, follow, and keep the instructions for use.
 - Use the product only in accordance with its intended use, see Intended use.
 - Do not use products from open or damaged sterile packaging.
 - Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Do not re-sterilize the product.
 - Do not use the product after its use-by date.

Safe operation

- Turn the interlock flange 4 of the obturator, thereby interlocking the obturator into the set screws 7 of the sheath, in order to ensure that the obturator 1 is interlocked into the sheath 2. Compare image 3. The tip of the obturator protrudes from the sheath. At the end of the sheath, there is no visible gap between the obturator and the shaft.
- With the aid of the scale on the side of the obturator, determine the depth to which the introducer should be inserted.
- Insert the single-use introducer to the desired depth and unlock and remove the obturator. When doing this, ensure that the obturator remains sterile.
- Carefully insert the endoscope with pushing forward the sheath. It can be visually confirmed with the endoscope that the sheath is positioned correctly.
- If the sheath is not in the desired position, carefully remove the endoscope and the sheath.
- Insert the obturator into the sheath again and interlock.
- Repeat the ventricular insertion. After removing the obturator, the sheath end should remain in the ventricular system. This reduced the risk of the tip causing hemorrhage when the endoscope is pushed forwards.

The round sheath end may cause hemorrhages!

- Under no circumstances push the sheath forwards into the brain without the obturator being interlocked in the correct position.
- Lower the obturator completely and interlock it in the sheath.

Note
Two sterile individuals are required for fixing the sheath.

- Person 1: Hold the sheath depth steady.
- Person 2: Hold sheath handles 6 securely.
- Separate handles 8. The sheath splits lengthways.
- Split the sheath to just above the cranial bone.
- Fix the split sheath securely to the cranium by sewing or tacking both sheath halves. The access canal to the ventricle is created.
- When introducing objects through the sheath, attention must be paid to ensuring that the sheath does not move. Do not insert any objects that barely fit through the sheath.

Validated reprocessing procedure

Single-use products

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product!

Storage

- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Technical Service

Risk of injury and/or malfunction!

- Do not modify the product.

► For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging! Incineration of the product, its components and its packaging produces harmless residues.
- Comply with national regulations for incineration under all circumstances!

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