

Aesculap Neurosurgery

-
- | | |
|------------|--|
| GB | Instructions for use/Technical description |
| USA | MINOP TREND (TRansnasal ENDoscopy System) |
| D | Gebrauchsanweisung/Technische Beschreibung |
| | MINOP TREND (TRansnasales ENDoskopiesystem) |
| F | Mode d'emploi/Description technique |
| | MINOP TREND (système TRansnasal ENDoscopique) |
| E | Instrucciones de manejo/Descripción técnica |
| | MINOP TREND (sistema TRansnasal de ENDoscopia) |
| I | Istruzioni per l'uso/Descrizione tecnica |
| | MINOP TREND (sistema per ENDoscopia TRansnasale) |

MINOP TREND (TRansnasal ENDoscopy System)

Legend

- 1 MINOP TREND (TRansnasal ENDoscopy System)
- 2 MINOP TREND suction/irrigation sheath, rotatable
- 3 MINOP TREND handle
- 4 Irrigation button
- 5 MINOP TREND endoscope
- 6 Adapter for holding arm
- 7 Cone

Symbols on product

Attention, see instructions for use

Intended use

The MINOP TREND (TRansnasal ENDoscopy System) is used for operating sites that can be reached through a transnasal approach. The MINOP TREND (TRansnasal ENDoscopy System) keeps the field of view of the adapted MINOP TREND endoscope free of blood and secretions and thereby ensures optimal visibility of the operating site in a transnasal approach.

The MINOP TREND (TRansnasal ENDoscopy System) can be used either in free-hand operation or, with adapter for holding arm RT099R, on the appropriate holding arm, combined with an irrigation fluid container or a suction pump and the appropriate instruments for the transnasal approach.

Available sizes

Designation	Art. no.
MINOP TREND handle (with irrigation button)	FH615
MINOP TREND suction/irrigation sheath 0°	FH610R
MINOP TREND suction/irrigation sheath 30°	FH611R
MINOP TREND endoscope 0°, Ø 4 mm	PE487A
MINOP TREND endoscope 30°, Ø 4 mm	PE507A
Adapter for holding arm	RT099R

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to the initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

The MINOP TREND handle can be combined with the following suction/irrigation sheaths:

- MINOP TREND suction/irrigation sheath 0° and MINOP TREND endoscope 0°
- MINOP TREND suction/irrigation sheath 30° and MINOP TREND endoscope 30°

When irrigation button **4** is pressed, the field of view of the adapted MINOP TREND endoscope **5** is rinsed and thus freed from any blood and secretion.

- To rinse MINOP TREND endoscope **5**: Press irrigation button **4** down to the stop. Suction is interrupted while irrigation is active.
- To stop rinsing MINOP TREND endoscope **5**: Release irrigation button **4**. The valve automatically switches to continuous suction and interrupts the irrigation.

Note

Suction is permanently active, except for the short periods when irrigation is activated.

Disassembling

Dismounting the MINOP TREND suction/irrigation sheath from the MINOP TREND handle

- Turn MINOP TREND suction/irrigation sheath **2** with MINOP TREND endoscope **5** inserted while pulling it out rearwards from MINOP TREND handle **3**.

Dismounting the tubes and the irrigation button


Fig. 1

- Remove the tubes from the tube connectors.
- Unscrew irrigation button **4**.

Dismounting adapter for holding arm RT099R

- Undo the cap nut.
- Remove adapter for holding arm RT099R.

Removing the MINOP TREND endoscope from the MINOP TREND suction/irrigation sheath

- Unscrew the light cable from MINOP TREND endoscope **5**.
- Remove the camera from the eyepiece funnel of MINOP TREND endoscope **5**.
- Unlock the clamping lever of MINOP TREND suction/irrigation sheath **2** and remove MINOP TREND endoscope **5**.

MINOP TREND (TRansnasal ENDoscopy System)

Assembling

Mounting the MINOP TREND endoscope in the MINOP TREND suction/irrigation sheath


Fig. 2

CAUTION

Damage to components caused by excessive force!

➤ Handle the components with appropriate care.

- Axially insert MINOP TREND endoscope 5 in MINOP TREND suction/irrigation sheath 2. When doing this make certain that the lock of MINOP TREND suction/irrigation sheath 2 is open and the pin fits into the groove at the bayonet coupling of MINOP TREND endoscope 5, see Fig. 2.
- Turn the clamping lever of MINOP TREND suction/irrigation sheath 2 as far as it will go. MINOP TREND endoscope 5 is now safely locked in position.

Connecting the tubes and the irrigation button


Fig. 3

WARNING

Restricted functionality of the MINOP TREND (TRansnasal ENDoscopy System) if the irrigation button is not screwed in completely!

- Make certain that irrigation button 4 is screwed into MINOP TREND handle 3 down to the stop.

CAUTION

Restricted functionality of the MINOP TREND (TRansnasal ENDoscopy System) due to incorrect connecting direction of the tubes!

- Note the arrow on the MINOP TREND handle.
- Carry out function checks before using the product.
- Verify that irrigation button 4 does not show any damage (valve seals, surfaces, etc.)
- Screw irrigation button 4 into MINOP TREND handle 3 down to the limit stop, see Fig. 3.

- Connect the suction tube to the suction connector, making certain that the arrow points away from MINOP TREND handle 3.
- Connect the irrigation tube to the irrigation connector, making certain that the arrow points towards MINOP TREND handle 3.

Note

To facilitate the connection, a MINOP TREND suction/irrigation tube is available. Its irrigation tube is marked by a continuous, blue line along its entire length. The suction tube is not marked.

Connecting the MINOP TREND suction/irrigation sheath to the MINOP TREND handle


Fig. 4

- Slide MINOP TREND suction/irrigation sheath 2 with MINOP TREND endoscope 5 into the rear of the handle until MINOP TREND suction/irrigation sheath 2 engages, see Fig. 4. MINOP TREND suction/irrigation sheath 2 is rotatable.
- Check the irrigation and suction tubes and the cone seating for any leaks.
- Press irrigation button 4 to check that it functions properly.

Mounting adapter for holding arm RT099R


Fig. 5

Fig. 6

- Undo the cap nut for adapter for holding arm 6.
- Insert adapter for holding arm 6 from either side (left/right) into the borehole until the flange engages in the appropriate groove, see Fig. 5.
- Screw on the cap nut and tighten it until adapter for holding arm 6 is fixated, see Fig. 6.

The MINOP TREND (TRansnasal ENDoscopy System) is now ready for being adapted on the holding arm.

MINOP TREND (TRansnasal ENDoscopy System)

Validated processing procedure

Note

Observe all relevant national regulations and standards with regard to processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, manufacturer's specifications can only serve as guide values for the processing procedures applied by the individual user.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. To avoid this, the time interval between application and processing should not exceed 6 h, and neither fixating pre-cleaning temperatures > 45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading of the laser engraving on stainless steel.

Residues containing chlorine or chlorides, as found e.g. in surgical residues, tinctures, medicines, saline solutions, service water used for cleaning, and in cleaning/disinfecting agents, cause corrosion damage (pitting, stress corrosion) and thereby destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals recommended by the chemicals' manufacturers as effective cleaning/disinfecting agents compatible with the materials used in the respective product. Strictly observe all application instructions, e.g. regarding temperatures, concentrations, exposure times, etc. Failure to do so can result in the following problems:

- optical changes to the material, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH > 8 to cause visible surface changes or
 - Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
 - Do not use process chemicals that cause stress cracking of plastics such as PPSU, or attack and cause brittleness e.g. of silicone softeners.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, Publications, Red Brochure – Proper maintenance of instruments.

- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove visible residues as completely as possible, using a lint-free, single-use cleaning tissue.

Preparation prior to cleaning

- Have the product dry in a disposal container and ready for immediate cleaning and disinfecting within 30 min after use.

Cleaning/Disinfecting

WARNING

Inadequate cleaning result if the irrigation button remains screwed in!

- **Remove the irrigation button prior to cleaning.**

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to manufacturer's instructions. The cleaning and disinfecting agent must
 - be approved for plastics and high-grade steel,
 - not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55 °C.

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting,
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting,
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting,
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically, provided they can be securely fixed in machines or storage devices in such a way that they will be thoroughly cleaned.

MINOP TREND (TRansnasal ENDoscopy System)

Manual cleaning/disinfecting

- Check visible surface for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfecting and cleaning with brush and ultrasound

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting ultrasound cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	1 st Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
IV	2 nd Intermediate rinse	RT	1	-	D-W	-
V	Intermediate drying	-	-	-	-	-
VI	Final rinse	RT (cold)	0.5	-	FD-W	-
VII	Drying	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage III

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush (e.g. PM995200) until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use metal brushes or other abrasives for cleaning the product.

Stage IV

- Rinse the product completely (all accessible surfaces) under running water.

Stage V

- Dry the product with compressed air.

Stage VI

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage VII

- Dry the product completely with compressed air.

MINOP TREND (TRansnasal ENDoscopy System)

Mechanical cleaning/disinfecting

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the injector carriage attachment.

Mechanical alkaline cleaning and thermal disinfecting

MINOP TREND handle

Machine type: Single-chamber washer/disinfector without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	< 25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	- Concentrate, alkaline: pH = 10.9 - < 5 % anionic tensides - 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD-W	- Concentrate, acid: pH = 2.6 Basis: Citric acid - 0.2 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD-W	-
V	Thermal disinfecting	94/201	10	FD-W	-
VI	Drying	90/194	40	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

Mechanical cleaning/disinfecting with manual pre-cleaning

Manual pre-cleaning with brush

MINOP TREND suction/irrigation sheaths

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush (e.g. PM995200) until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives that would damage the surfaces for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

MINOP TREND (TRansnasal ENDoscopy System)

Manual pre-cleaning with ultrasound

Adapter for holding arm/irrigation button

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting ultrasound cleaning	RT (cold)	15	1.5	D-W	Shylke & Meyr Gigasept AF aldehyde-free, phenol-free, QAV- free; pH = 8.5
II	Intermediate rinse	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Adapter for holding arm/irrigation button and MINOP TREND suction/irrigation sheaths

Machine type: Single-chamber washer/disinfector without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the injector carriage attachment.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	< 25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	- Concentrate, alkaline: pH = 10.9 < 5 % anionic tensides - 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD-W	- Concentrate, acid: pH = 2.6 Basis: Citric acid - 0.2 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD-W	-
V	Thermal disinfecting	94/201	10	FD-W	-
VI	Drying	90/194	40	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

MINOP TREND (TRansnasal ENDoscopy System)

Control, care and inspection

- Allow the product to cool down to room temperature.
- Lightly lubricate movable parts (irrigation sheath lock, irrigation button, handle cone) with maintenance oil suitable for the respective sterilization process (e.g. Aesculap STERILIT® Spray JG600 or maintenance oil JG598).
- Assemble the separable product, see Assembling.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Check for compatibility with associated products.
- Set aside the product if it is damaged.

Packaging

- Appropriately protect products with fine working tips.
- Sort the product into its appropriate storage device or put it on a suitable tray.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters

Note

The product can be sterilized either in disassembled or in assembled condition.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer acc. to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST 46-1993, validated acc. to EN 554/ISO 13683

– Sterilization through fractionated vacuum process at 134 °C/holding time 5 min

- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterilcontainer (perforated bottom)			
Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	132 °C (275 °F)	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap Sterilcontainer cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store processed products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

Technical Service

For service, maintenance or repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Designation	Art.-no.
MINOP TREND tray for trocars and optics	FF357R
MINOP TREND suction/irrigation tube	FH605SU
Irrigation button	FH615801
Cleaning brush	PM995200

Disposal

Note

No special measures are required for disposal of the products, although national rules and regulations must be adhered to under all circumstances!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG