- (GB)** Instructions for use/Technical description
- (USA)** MINOP optics and angled optics
- (D)** Gebrauchsanweisung/Technische Beschreibung  
MINOP- und Winkeloptiken
- (F)** Mode d'emploi/Description technique  
Optiques MINOP et optiques coudées
- (E)** Instrucciones de manejo/Descripción técnica  
Ópticas MINOP y ópticas anguladas
- (I)** Istruzioni per l'uso/Descrizione tecnica  
Ottiche MINOP e angolari
- (P)** Instruções de utilização/Descrição técnica  
Endoscópios MINOP e angulares
- (NL)** Gebruiksaanwijzing/Technische beschrijving  
MINOP- en hoekendoscopen
- (S)** Bruksanvisning/Teknisk beskrivning  
MINOP- och vinkelendoskop
- (RUS)** Инструкция по применению/Техническое описание  
MINOP и угловая оптика
- (CZ)** Návod k použití/Technický popis  
MINOP a úhlové optiky
- (PL)** Instrukcja użytkowania/Opis techniczny  
Optyki MINOP i optyki kątowe
- (SK)** Návod na použitie/Technický opis  
Optika MINOP a uhlové optiky
- (TR)** Kullanım Kılavuzu/Teknik açıklama  
MINOP ve açılı optikler


**B | BRAUN**  
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

TA-Nr. 013238 2017-12 V6 Änd.-Nr. 57829

**CE** 0123 – DIR 93/42/EEC

Aesculap®

MINOP optics and angled optics

Legend

- 1 Distal window
- 2 Sheath
- 3 Optical cable connection
- 4 Illumination surface of the optical cable connection
- 5 Eyepiece housing
- 6 Eyepiece window
- 7 Inscription ring
- 8 ACMI adapter (fixed)
- 9 Wolf adapter
- 10 Storz adapter
- 11 Color ring

Symbols on product and packages

|  |                                                                                         |
|--|-----------------------------------------------------------------------------------------|
|  | Caution, general warning symbol<br>Caution, see documentation supplied with the product |
|  | Date of manufacture                                                                     |
|  | Steam sterilization at temperatures up to 134 °C                                        |
|  | Low temperature and plasma sterilization                                                |
|  | High-Definition resolution                                                              |

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

Intended use

MINOP optics and angled optics are used for visualizing the ventricles and other structures in the brain, as well as structures at the spine.

The following MINOP and angled optics are approved for direct contact with the central nervous system:

- PE184A
- PE204A
- PE486A
- PE506A
- PE526A

Color coding

The viewing direction of the endoscope is indicated by color ring 11 at optical cable connection 3.

| Colour | Viewing direction |
|--------|-------------------|
| Green  | 0°                |
| Red    | 30°               |
| Yellow | 70°               |

Indications

Indications, see Intended use.

*Note*  
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Contraindications

None known.

Risks and side effects

As part of the legal duty to inform, the following typical risks and side effects associated with the use of surgical instruments are referred to. These are predominantly procedure-specific, non product-specific and not limited to unwanted damage to surrounding tissue resulting in e.g. bleeding, infection, material incompatibilities or instrument parts remaining unnoticed in the patient, etc.

Available sizes

The rigid, autoclavable neuroendoscopes are available in the following sizes and designs:

- Sheath Ø 2.7 mm and 4.0 mm
- Angled endoscopes

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

This manual contains instructions for the preparation, reprocessing, and disposal of the endoscope. It does not contain information on the actual application of the endoscope.

**Risk of injury from defective endoscopes!**

- Only use an endoscope if it is in perfect condition.

- The products are supplied in an unsterile condition and must be reprocessed and sterilized prior to use.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
  - Read, follow, and keep the instructions for use.
  - Use the product only in accordance with its intended use, see Intended use.
  - When using the endoscope in combination with any electromedical device, make sure that the respective BF conditions are adhered to (insulated, floating patient application piece).
  - Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
  - Store any new or unused products in a dry, clean, and safe place.
  - Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
  - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
  - Replace any damaged components immediately with original spare parts.
  - To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).
  - Make certain a reserve product is available.
  - Do not, under any circumstances, put down the product on the patient or on the surgical drape covering the patient.

The following table shows the approved combinations of endoscopes, sheaths, and working trocars:

| Insert the endoscope | Sheath and working trocar      |
|----------------------|--------------------------------|
| PE184A               | FF397R, FF398R, FF399R, FH601R |
| PE204A               | FF397R, FF398R, FF399R, FH620R |

Safe operation

**Risk of injury and/or malfunction!**

- Always carry out a function check prior to using the product.

Function checks

- Check the optical function of the endoscope. The image has to be clear and distinct.
- Ensure that the distal window 1, the eyepiece window 6 and the illumination surface 4 of the optical cable connection 3 are not cloudy, dirty or scratched.
- Hold optical cable connection 3 of the endoscope against the light and check that the light fibers at the distal end light up evenly.
- Inspect the sheath for dents, bends, and scratches.

Use of the endoscope

**Burns caused by high-frequency current (HF current)!**

- When using an endoscope and HF electrodes at the same time, be cautious to activate the HF current under visual control only.
- Ensure that the active electrode is never in the immediate vicinity of electrically conductive components (e.g. trocar, endoscope).

**Burns suffered by patient or user from exposure to high light intensity!**

- Make certain the distal end of the endoscope or optical cable connection does not touch human tissue or any flammable or heat-sensitive materials while the light source is active.
- Do not put down the endoscope on the patient.
- Do not touch the distal end of the endoscope and the optical cable connection.
- Adjust the light source to the minimum required power for optimal illumination of the endoscopic image.
- Only use light sources of up to 300 W electric power.

**Damage to the endoscope by bending the endoscope sheath!**

- Do not bend the endoscope sheath.
- Use the endoscope with its appropriate sheaths and working trocars only.
- Always hold the endoscope at the eyepiece housing but never at the sheath.

- Use endoscopes only with halogen light sources equipped with a spare lamp or with xenon light sources.
- Observe the respective manufacturer's instructions when using the endoscope in combination with other equipment, instruments, or optical cables.

Disassembling

- Unscrew adapter 9 or 10, respectively, from the endoscope.

Assembling

- Screw adapter 9 or 10, respectively, onto the endoscope.

Validated reprocessing procedure

General safety notes

*Note*  
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

*Note*  
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

*Note*  
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

*Note*  
It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.  
The specified chemistry was used for validation.

*Note*  
If there is no final sterilization, then a virucidal disinfectant must be used.

*Note*  
Only clean medical products can be sterilized safely and effectively. Cleaning is therefore of special importance in the reprocessing procedure.

*Note*  
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at [www.aesculap-extra.net](http://www.aesculap-extra.net).

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.  
Excessive dosages of neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking on stainless steel becoming unreadable either visually or by machine.  
Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.  
Perform additional drying, if necessary.  
Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes in materials, e.g., fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see [www.a-k-i.org](http://www.a-k-i.org) Rubrik Veröffentlichungen Rote Broschüre – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- If applicable, rinse surfaces that are not accessible to visible inspection (preferably with demineralized water), using a disposable syringe, for example.
- Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

**DANGER**

**Risk to patient due to cross contamination!**

- Do not clean contaminated products together with uncontaminated products in a tray.

**CAUTION**

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
  - are approved for rigid endoscopes,
  - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum permitted cleaning temperature of 55 °C.

**CAUTION**

**Damage to the optical system caused by loosening of connections during ultrasound cleaning!**

- Do not clean the endoscope by ultrasonic treatment under any circumstances.

- For prion inactivation, clean/disinfect the endoscope mechanically using a prion-inactivating cleaning agent.
- Do not use oxidizing chemicals (e.g. H<sub>2</sub>O<sub>2</sub>), which could cause bleaching/layer loss of the product.
- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.
- Only use bactericidal, fungicidal, and virucidal disinfectants.
- Preferably use thermal disinfecting processes.
- Following the chemical disinfection, thoroughly rinse the product with clear running water. Always adhere to the manufacturer's instructions.

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

| Phase | step                  | T<br>[°C/°F]    | t<br>[min] | Conc.<br>[%] | Water<br>quality  | Chemicals                                                    |
|-------|-----------------------|-----------------|------------|--------------|-------------------|--------------------------------------------------------------|
| I     | Cleaning              | 35–45<br>95–113 | 5          | 0,8          | D–W               | Cidezyme/Enzol                                               |
| II    | Intermediate<br>rinse | RT (cold)       | 3 x 1      | –            | D–W               | –                                                            |
| III   | Disinfection          | 20–25<br>68–77  | 12         | 0,55         | D–W               | Cidex OPA (process solution,<br>0.55 % ortho-phthalaldehyde) |
| IV    | Final rinse           | RT (cold)       | 3 x 2      | –            | FD–W<br>(sterile) | –                                                            |
| V     | Drying                | RT              | –          | –            | –                 | –                                                            |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination:  
max. 10 microorganisms/ml, low endotoxin: max. 0.25 endotoxin units/ml)  
RT: Room temperature

Phase I

- Completely immerse the instrument in the cleaner/disinfectant solution. Ensure that all accessible surfaces are moistened.
- Clean the product immersed in the solution, using a soft cloth or, if necessary, a suitable cleaning brush, until all visible residues have been removed from the surfaces.
- Use a suitable cleaning brush to clean all surfaces that are not visible, e.g. in products with concealed crevices, lumens, or complex geometries, for at least 1 minute or until no further residues can be removed. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop.
- Thoroughly rinse these components with the cleaning disinfectant solution (at least five times), using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion

Phase II

- Completely rinse/ rinse through the product (all accessible surfaces) 3 times for at least 1 min. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop. Use fresh water for every rinse.
- Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml).
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution. Ensure that all accessible surfaces are moistened.
- During disinfection, mobilize non-rigid components, such as set screws, joints etc., 3 times in each direction to the positive stop.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml). Ensure that all accessible surfaces are moistened.

Phase IV

- Following disinfection, rinse/ rinse through the product completely (all accessible surfaces) 3 times for at least 2 minutes. Use fresh water for every rinse.
- During the final rinse, mobilize non-rigid components, such as set screws, joints etc., 3 times in each direction to the positive stop.
- Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml).
- Drain any remaining water fully.

Phase V

- Dry the product with a soft, lint-free cloth.
- Dry inaccessible areas with medical compressed air (P<sub>max</sub> = 5 bar).

Mechanical cleaning/disinfecting

*Note*  
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).  
*Note*  
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

*Note*  
Only use cleaning agents that are suitable for this product class (rigid endoscope).

- Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots).
- Products with a sheath diameter of ≤4 mm must be cleaned and sterilized in special Aesculap optics trays only.
- Process the product in a cleaning/disinfecting machine. Follow the instructions of the washing machine's manufacturer.
- Avoid abrupt cooling of the product (e.g. in water).

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemicals                                                                   |
|-------|----------------------|--------------|------------|------------------|-----------------------------------------------------------------------------|
| I     | Pre-rinse            | <25/77       | 3          | D–W              | –                                                                           |
| II    | Cleaning             | 55/131       | 10         | FD–W             | Dr. Weigert neodisher® SeptoClean 1 %*<br>working solution                  |
| III   | Intermediate rinse   | >10/50       | 1          | FD–W             | –                                                                           |
| IV    | Thermal disinfection | 90/194       | 5          | FD–W             | –                                                                           |
| V     | Drying               | –            | –          | –                | In accordance with the program for the<br>cleaning and disinfecting machine |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
\*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

**CAUTION**

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) at the marked lubrication points, using maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Inspect surfaces, cavities, lumens and openings for visible debris.
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- For products with locking mechanism (e.g. MINOP): Check the locking mechanism for smooth movement.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Sterilization

*Note*  
Alternating sterilization by different sterilization processes can result in damage to materials and adapters of the product.

- The following processes can be used to sterilize the endoscope:
- Steam sterilization
  - Sterrad® sterilization process: using Sterrad® sterilizers Sterrad® 50, Sterrad® 100S, or Sterrad® 200

**CAUTION**

**Damage to the optical system caused by flash sterilization!**

- Do not flash sterilize the endoscope.
- Do not expose the endoscope to temperatures above 134 °C.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Do not autoclave damaged products.
- Protect the products against mechanical impacts.

Steam sterilization

*Note*  
The product can be steam-sterilized either in disassembled or in assembled condition.

*Note*  
Aesculap guarantees steam sterilization of the product up to 50 autoclave cycles at a holding time of 18 min, or for an unlimited number of cycles at a holding time of 5 min.

- Validated sterilization process
  - Steam sterilization using the fractional vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractional vacuum process at 134 °C, holding time 5 min to 18 min
  - For prion inactivation: Sterilization using fractional vacuum process at 134 °C, holding time 18 minutes
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterrad® sterilization process Sterrad® 50, Sterrad® 100S, Sterrad® 200

*Note*  
The Sterrad® sterilization process can cause cosmetic changes to the product. These changes will, however, not affect its functionality.

*Note*  
The product must be disassembled prior to sterilization through the Sterrad® process.

- Disassemble the product.
- Only sterilize endoscopes carrying the SDS symbol for low-temperature and plasma sterilization on the inscription ring 7 by means of the Sterrad® sterilization process 50, 100S, 200.
- Unscrew dismountable parts of the endoscope, see Disassembling.
- Sterilize using the Sterrad® sterilization process 50, 100S, 200. When doing so, observe the following: Follow the Sterrad® system manufacturer's instructions. Use of a biological indicator is recommended to confirm effective sterilization.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
  - Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.
- To achieve a sterility assurance level of 10<sup>-6</sup>, Aesculap recommends the following parameters:

| Aesculap Orga Tray/Sterile container (perforated bottom)<br>Minimum cycle parameters* |               |       |                     |
|---------------------------------------------------------------------------------------|---------------|-------|---------------------|
| Sterilization method                                                                  | Temp.         | Time  | Minimum drying time |
| Prevacuum                                                                             | 270 °F/275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

**WARNING for the US market**  
If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Maintenance

The product requires no maintenance for the first 50 autoclave cycles at a holding time of 18 min, or for an unlimited number of cycles at a holding time of 5 min.

After that, the image quality of the product can be diminished by normal wear during steam sterilization.

Technical Service

**WARNING**

**Risk of injury and/or malfunction!**

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- Prior to dispatching the product for repairs:
  - Have the product cleaned and disinfected or sterilized, and mark it as "disinfected" or "sterilized", respectively.
  - Pack the endoscope in such a way that it will be protected against transport damage.

**Service addresses**

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 7461 95-1601  
Fax: +49 7461 14-939  
E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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