Aesculap Neurosurgery

(GB)

Instructions for use

MINOP optics and angled optics

(USA)

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(D)

(F)

Mode d'emploi

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
MINOP optics and angled optics

Legend

- 1 Distal window
- 2 Shaft
- 3 Optical cable connection
- 4 Illumination surface of the optical cable connection
- 5 Eyepiece housing
- 6 Eyepiece window
- 7 Inscription ring
- 8 ACMI adapter (fixed)
- 9 Wolf adapter
- 10 Storz adapter
- 11 Color ring

Symbols on product and packaging

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture
	Steam sterilization at up to 134 °C
	Low temperature and plasma sterilization
	High-Definition resolution

Color coding

The viewing direction of the endoscope is indicated by color ring **11** at optical cable connection **3**.

Colour	Viewing direction
Green	0°
Red	30°
Yellow	70°

Intended use

MINOP optics and angled optics are used for visualizing the ventricles and other structures in the brain, as well as structures at the spine.

The following MINOP and angled optics are approved for direct contact with the central nervous system:

- PE184A
- PE204A
- PE486A
- PE506A
- PE526A

Available sizes

The rigid, autoclavable neuroendoscopes are available in the following sizes and designs:

- Sheath diameter 2.7 mm and 4.0 mm
- Angled endoscopes

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

This manual contains instructions for the preparation, reprocessing, and disposal of the endoscope. It does not contain information on the actual application of the endoscope.

WARNING

Risk of injury from defective endoscopes!

- Only use an endoscope if it is in perfect condition.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, and experience.
- Read, follow, and keep the instructions for use.
- Use the product only for its intended purpose, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working tip: Exercise caution when passing products through the working channel (e.g. trocar).
- Make certain a reserve product is available.
- Never, under any circumstances, put down the product on the patient or on the surgical drape covering the patient.

The following table shows the approved combinations of endoscopes, sheaths, and working trocars:

Endoscope	Sheath and working trocar
PE184A	FF397R, FF398R, FF399R, FH601R
PE204A	FF397R, FF398R, FF399R, FH620R

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

Function checks

- Check the optical function of the endoscope. The image has to be clear and distinct.
- Make certain that the window **1** at the distal end, eyepiece window **6** and illumination surface **4** of optical cable connection **3** are not cloudy, dirty, or scratched.
- Hold optical cable connection **3** of the endoscope against the light and check that the light fibers at the distal end light up evenly.
- Inspect the sheath for dents, bends, and scratches.

Use of the endoscope

WARNING

Burns caused by high-frequency current (HF current)!

- When using an endoscope and HF electrodes at the same time, be cautious to activate the HF current under visual control only.
- Ensure that the active electrode is never in the immediate vicinity of electrically conductive components (e.g. trocar, endoscope).

MINOP optics and angled optics

**WARNING**

Burns suffered by patient or user from exposure to high light intensity!

- Make certain the distal end of the endoscope or optical cable connection does not touch human tissue or any flammable or heat-sensitive materials while the light source is active.
- Do not put down the endoscope on the patient.
- Do not touch the distal end of the endoscope and the optical cable connection.
- Adjust the light source to the minimum required power for optimal illumination of the endoscopic image.
- Only use light sources of up to 300 W electric power.

**CAUTION**

Damage to the endoscope by bending the endoscope sheath!

- Do not bend the endoscope sheath.
- Use the endoscope with its appropriate sheaths and working trocars only.
- Always hold the endoscope at the eyepiece housing but never at the sheath.

- Use endoscopes only with halogen light sources equipped with a spare lamp or with xenon light sources.
- Observe the respective manufacturer's instructions when using the endoscope in combination with other equipment, instruments, or optical cables.
- When using the endoscope in combination with any electromedical device, make sure that the respective BF conditions are adhered to (insulated, floating patient application piece).

Disassembly

- Unscrew adapter **9** or **10**, respectively, from the endoscope.

Assembly

- Screw on adapter **9** or **10**, respectively.

Validated reprocessing procedure

Note

National laws, national and international standards and directives, and product-specific hygiene regulations for reprocessing must be observed.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/sterile processing technician is responsible for this.

Note

Up-to-date information on reprocessing can be found on the Aesculap Extranet at www.aesculap-extra.net

Note

Only clean medical products can be sterilized safely and effectively. Cleaning is therefore of special importance in the reprocessing procedure.

Note

Endoscopes are delicate optical devices. Aesculap recommends reprocessing them separately. Aesculap offers special optics trays for the sterilization, storage, and sterile setup of endoscopes.

General information

To prevent increased contamination of loaded instrument trays during use, please ensure that contaminated instruments are collected separately and not returned to the instrument tray.

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion of stainless steel. Consequently, the time interval between application and processing should not exceed 6 hours, pre-cleaning fixating temperatures >45°C should not be used, and no fixating disinfectants (containing aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chloride - e.g. in surgical residues, drugs, saline solutions and water for cleaning, disinfection and sterilization - may cause corrosion damage to stainless steel (pinholing, stress corrosion), thus rendering the products unusable. These must be removed by rinsing thoroughly with demineralized water and then drying.

Only use process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All the chemical manufacturer's application specifications regarding temperature, concentration and contact time should be strictly observed. Failure to do so can result in the following problems:

- Visible changes to materials e.g. fading or changes in the colour of titanium or aluminium. As regards aluminium, visible changes to the surface may already occur at a pH level of >8 in the application/working solution.
- Material damage such as corrosion, cracks, fractures, premature deterioration, or swelling.
- Do not use oxidizing chemicals (e.g. H₂O₂) for mechanical cleaning, as these can cause bleaching or layer loss.

MINOP optics and angled optics

- Do not use process chemicals that cause stress cracks or brittleness in plastics.
- Clean the product immediately after use.

Please see www.a-k-i.org for more detailed information on hygienically safe reprocessing which is protective of materials and retains their value.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in these instructions for use, see Disassembly.
- Rinse non-visible surfaces such as those in instruments with concealed crevices, lumens or complex geometries, preferably with distilled water, using a disposable syringe for instance.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Disassemble the product prior to cleaning, see Disassembly.

Cleaning/disinfection

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions, which
 - are approved for rigid endoscopes,
 - are approved for use e.g. on aluminum, plastic materials, and stainless steel,
 - do not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55 °C.

- Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

CAUTION

Damage to the optical system caused by loosening of connections during ultrasound cleaning!

- Do not clean the endoscope by ultrasonic treatment under any circumstances.

- Only use bactericidal, fungicidal, and virucidal disinfecting agents.
- Preferably use thermal disinfecting processes.
- Following the chemical disinfection, thoroughly rinse the product with clear running water. Always adhere to the manufacturer's instructions.

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Where necessary, repeat the cleaning process.
- To remove residues of cleaning/disinfecting agents, always wipe distal window **1**, illumination window **4** of the light cable connector, and eyepiece glass **6** with a swap moistened with alcohol.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	35-45 95-113	5	0.8	D-W	Cidezyme/Enzol
II	Intermediate rinse	RT (cold)	3 x 1	-	D-W	-
III	Disinfection	20-25 68-77	12	0.55	D-W	Cidex OPA (process solution, 0.55 % ortho-phthalaldehyde)
IV	Final rinse	RT (cold)	3 x 2	-	FD-W (sterile)	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: max. 10 microorganisms/ml, low endotoxin: max. 0.25 endotoxin units/ml)

RT: Room temperature

MINOP optics and angled optics

Phase I

- Completely immerse the instrument in the cleaner/ disinfectant solution. Ensure that all accessible surfaces are moistened.
- Clean the product immersed in the solution, using a soft cloth or, if necessary, a suitable cleaning brush, until all visible residues have been removed from the surfaces.
- Use a suitable cleaning brush to clean all surfaces that are not visible, e.g. in products with concealed crevices, lumens, or complex geometries, for at least 1 minute or until no further residues can be removed. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop.
- Thoroughly rinse these components with the cleaning disinfectant solution (at least five times), using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Phase II

- Completely rinse/ rinse through the product (all accessible surfaces) 3 times for at least 1 min. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop. Use fresh water for every rinse.
- Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml).
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfecting solution. Ensure that all accessible surfaces are moistened.
- During disinfection, mobilize non-rigid components, such as set screws, joints etc., 3 times in each direction to the positive stop.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml). Ensure that all accessible surfaces are moistened.

Phase IV

- Following disinfection, rinse/ rinse through the product completely (all accessible surfaces) 3 times for at least 2 minutes. Use fresh water for every rinse.
- During the final rinse, mobilize non-rigid components, such as set screws, joints etc., 3 times in each direction to the positive stop.
- Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml).
- Drain any remaining water fully.

Phase V

- Dry the product with soft, lint-free cloth.
- Dry inaccessible areas with compressed air ($P_{\max} = 5 \text{ bar}$).

Mechanical cleaning/disinfecting

Note

The disinfectant must be of proven effectiveness (e.g. DGHM or FDA approval or CE mark).

Note

Only use cleaning agents that are suitable for this product class (rigid endoscope).

Note

For thermal disinfection, always use fully desalinated water (demineralized, low microbial contamination: drinking water quality at least) and ensure that Ao is >3 000 for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

- Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots).
- Clean and disinfect endoscopes with a sheath diameter of ≤ 4 mm only in special Aesculap optics trays.
- Process the product in a washer/disinfectant machine. Follow the instructions of the washing machine's manufacturer.
- Avoid abrupt cooling of the product (e.g. in water).

MINOP optics and angled optics

Mechanical alkaline cleaning and thermal disinfection

Machine type: Single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	2	D-W	-
II	Cleaning	55/131	10	FD-W	neodisher MediClean forte 0.5 % (5 ml/L) pH >10
III	Rinse I	>10/50	1	D-W	-
IV	Rinse II	>10/50	1	FD-W	-
V	Thermal disinfection	90/194	5	FD-W	-
VI	Drying	-	-	-	According to disinfectant program

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: max. 10 microorganisms/ml, low endotoxin: max. 0.25 endotoxin units/ml)

Inspection, maintenance, and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Inspect surfaces, cavities, lumens and openings for visible debris.
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- To remove residues of cleaning/disinfecting agents, wipe distal window 1, illumination window 4 of the light cable connector, and eyepiece glass 6 with a swap moistened with alcohol.
- Check that the product functions correctly.
- For products with locking mechanism (e.g. MINOP): Check the locking mechanism for smooth movement.

- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable instruments, see Assembly.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Observe the individual weight limit for each tray.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization

Note

Alternating sterilization by different sterilization processes can result in damage to materials and adapters of the product.

The following processes can be used to sterilize the endoscope:

- Steam sterilization
- Sterrad® sterilization process: in Sterrad® sterilizers Sterrad® 50, Sterrad® 100S, Sterrad® 200

CAUTION

Damage to the optical system caused by flash sterilization!

- Do not flash sterilize the endoscope.
 - Do not expose the endoscope to temperatures above 134 °C.
- Ensure that the sterilant can reach all external and internal surfaces (e.g. by opening valves and stopcocks).

MINOP optics and angled optics

- Do not autoclave damaged products.
- Protect the products against mechanical impacts.

Steam sterilization

Note

The product can be steam-sterilized either in disassembled or in assembled condition.

Note

Only autoclave endoscopes marked with the symbol for steam sterilization on inscription ring 7.

- Validated sterilization process
 - Disassemble the product, if applicable.
 - Steam sterilization using the fractionated vacuum process
 - Steam sterilizer specified in DIN EN 285 and validated in accordance with DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134°C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

Sterrad® sterilization process Sterrad® 50, Sterrad® 100S, Sterrad® 200

Note

The Sterrad® sterilization process can cause cosmetic changes of the product. These changes will, however, not affect its functionality.

Note

The product must be disassembled prior to sterilization through the Sterrad® process.

- Disassemble the instrument.
- Only sterilize endoscopes marked with the SDS symbol for low-temperature and plasma sterilization on inscription ring 7 through the Sterrad® sterilization process 50, 100S or 200.

- Unscrew dismountable parts of the endoscope, see Disassembly.
- Sterilize by the Sterrad® sterilization process 50, 100S or 200, observing the following rules:
Follow the advice of the manufacturers of the Sterrad® system. Use of a biological indicator is recommended to confirm effective sterilization.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F– 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- Prior to dispatching the product for repairs:
 - Have the product cleaned and disinfected or sterilized, and mark it as “disinfected” or “sterilized”, respectively.
 - Pack the endoscope in such a way that it will be protected against transport damage.

Service addresses

AESCLAP FLEXIMED GmbH
Robert-Bosch-Strasse 1
79211 Denzlingen / Germany
Phone: +49 (0) 7666 / 9321-0
Fax: +49 (0) 7666 / 9321-580
E-mail: techserv@aesculap-fleximed.de
www.aesculap-fleximed.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG