
Aesculap Neurosurgery

GB	Instructions for use
USA	MINOP optics and angled optics
D	Gebrauchsanweisung
	MINOP- und Winkeloptiken
F	Mode d'emploi
	Optiques MINOP et optiques coudées
E	Instrucciones de manejo
	Ópticas MINOP y ópticas anguladas
I	Istruzioni per l'uso
	Ottiche MINOP e angolari
P	Instruções de utilização
	Endoscópios MINOP e angulares
NL	Gebruiksaanwijzing
	MINOP- en hoek-endoscopen
DK	Brugsanvisning
	MINOP- og vinkeloptik
N	Bruksanvisning
	MINOP- og vinkeloptikk
S	Bruksanvisning
	MINOP- och vinkelendoskop
FIN	Käyttöohjeet
	MINOP- ja kulmaoptiikat
GR	Οδηγίες χρήσης
	Οπτικά στοιχεία MINOP και γωνίας
CN	使用说明
	刚性耐高压高温内窥镜

MINOP optics and angled optics

Legend

- 1 Distal window
- 2 Sheath
- 3 Optical cable connection
- 4 Illumination surface of the optical cable connection
- 5 Eyepiece housing
- 6 Eyepiece window
- 7 Inscription ring
- 8 ACMI adapter (fixed)
- 9 Wolf adapter
- 10 Storz adapter
- 11 Color ring

Symbols on product and packaging

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture
	Product supplied in unsterile condition
	Steam sterilization at up to 134 °C
	Low temperature and plasma sterilization
	High-Definition resolution

Color coding

The viewing direction of the endoscope is indicated by color ring **11** at optical cable connection **3**.

Color	Viewing direction
Green	0°
Red	30°
Yellow	70°

Intended use

MINOP optics and angled optics are used for visualizing the ventricles and other structures in the brain, as well as structures at the spine.

The following MINOP and angled optics are approved for direct contact with the central nervous system:

- PE184A
- PE204A
- PE486A
- PE506A
- PE526A

Available sizes

The rigid, autoclavable neuroendoscopes are available in the following sizes and designs:

- Sheath diameter 2.7 mm and 4.0 mm
- Angled endoscopes

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

This manual contains instructions for the preparation, reprocessing and disposal of the endoscope. It does not contain information on the actual application of the endoscope.

Risk of injury from defective endoscopes!

- Only use an endoscope if it is in perfect condition.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Exercise caution when passing products through the working channel (e.g. trocar).
- Make certain a reserve product is available.
- Never, under any circumstances, put down the product on the patient or on the surgical drape covering the patient.

The following table shows the approved combinations of endoscopes, sheaths and working trocars:

Endoscope	Sheath and working trocar
PE184A	FF397R, FF398R, FF399R, FH601R
PE204A	FF397R, FF398R, FF399R,

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

Function checks

- Check the optical function of the endoscope. The image has to be clear and distinct.
- Make sure that the window **1** at the distal end, the eyepiece window **6** and the illumination surface **4** of the optical cable connection **3** are not cloudy, dirty or scratched.
- Hold the optical cable connection **3** of the endoscope against the light and check that the light fibers at the distal end light up evenly.
- Inspect the sheath for dents, bends and scratches.

MINOP optics and angled optics

Use of the endoscope

WARNING

Burns caused by high-frequency current (HF current)!

- When using an endoscope and HF electrodes at the same time, be cautious to activate the HF current under visual control only.
- Be sure that the active electrode never is in the immediate vicinity of electrically conductive components (e.g. a trocar, an endoscope).

WARNING

Burns suffered by patient or user from exposure to high light intensity!

- Make certain the distal end of the endoscope or optical cable connection does not touch human tissue or any flammable or heat-sensitive materials while the light source is active.
- Do not put down the endoscope on the patient.
- Do not touch the distal end of the endoscope and the optical cable connection.
- Adjust the light source to the minimum required power for optimal illumination of the endoscopic image.
- Only use light sources of up to 300 W electric power.

CAUTION

Damage to the endoscope by bending the endoscope sheath!

- Do not bend the endoscope sheath.
- Use endoscopes with their appropriate sheaths and working trocars only.
- Always hold the endoscope at the eyepiece housing but never at the sheath.
- Use endoscopes only with halogen light sources equipped with a spare lamp or with xenon light sources.
- Observe the respective manufacturer's instructions when using the endoscope in combination with other equipment, instruments or optical cables.
- When using the endoscope together in combination with any electromedical device, make sure that the respective BF conditions are adhered to (insulated, floating patient application piece).

Disassembling

- Unscrew the adapter 9 or 10, respectively, from the endoscope.

Assembling

- Screw on the adapter 9 or 10 respectively.

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

Note

Only clean medical products can be sterilized safely and effectively. Cleaning is therefore of special importance in the reprocessing procedure.

Rigid autoclavable endoscopes are delivered in an unsterile condition.

- The endoscopes must be cleaned and sterilized before being used.

Endoscopes are delicate optical devices. Therefore, Aesculap recommends processing them separately.

For the sterilization, storage and sterile setup of endoscopes, Aesculap offers special optics trays.

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray. Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Therefore the time interval between application and cleaning should not exceed 6 h, no potentially fixating pre-cleaning temperatures >45 °C be applied, and no fixating disinfectants (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfecting and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

MINOP optics and angled optics

Only use process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process solution can already cause visible surface changes.
 - Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use oxidizing chemicals (e.g. H₂O₂) for mechanical cleaning, as these can cause bleaching or layer loss.
 - Do not use process chemicals that would cause stress cracking or brittleness of plastics.
 - Clean the product immediately after use.
 - Preclean, using non-fixating/NaCl-free method immediately after use.

For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

- Use appropriate cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in these instructions for use, see Disassembling.
- Rinse surfaces that cannot be visually inspected, e.g. on products with hidden crevices or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.

Preparation prior to cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions. The cleaning and disinfecting agents must
 - be approved for rigid endoscopes,
 - not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature, exposure time and use-by date.
- Do not exceed the maximum allowable cleaning temperature of 55 °C.

CAUTION

Damage to the optical system caused by loosening of connections during ultrasound cleaning!

➤ **Do not clean the endoscope with ultrasound.**

- Only use bactericidal, fungicidal and virucidal disinfecting agents.
 - Preferably use thermal disinfecting processes.
 - Following the chemical disinfection, thoroughly rinse the product with clear running water.
- Always adhere to the manufacturer's instructions.

Manual cleaning/disinfecting

- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.
- To remove residues of cleaning/disinfecting agents, always wipe distal window **1**, illumination surface **4** of the light cable connector, and eyepiece window **6** with a swap moistened with alcohol.

MINOP optics and angled optics

Manual cleaning with immersion disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	35–45/ 95–113	5	0.8	D–W	Cidezyme/Enzol
II	Intermediate rinse	RT (cold)	3 x 1	–	D–W	–
III	Disinfecting	20–25/ 68–77	12	0.55	D–W	Cidex OPA (process solution, 0.55 % ortho-phthalaldehyde)
IV	Final rinse	RT (cold)	3 x 2	–	FD–W (sterile)	–
V	Drying	RT	–	–	–	–

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: max. 10 microorganisms/ml, low endotoxin: max. 0.25 endotoxin units/ml)

RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning solution. Make certain that all accessible surfaces are moistened.
- Clean the product immersed in the solution, using a soft cloth or, if necessary, a suitable cleaning brush, until all visible residues have been removed from the surfaces.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 minute or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).

- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Completely rinse/ rinse through the product (all accessible surfaces) 3 times for at least 1 min. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop. Use fresh water for every rinse.
- Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml).
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml). Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop.

Stage IV

- Following disinfection, rinse/ rinse through the product completely (all accessible surfaces) 3 times for at least 2 minutes. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop. Use fresh water for every rinse.
- Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml).
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with soft, lint-free tissue.
- Dry inaccessible areas with compressed air $P_{\max} = 5 \text{ bar}$.

Mechanical cleaning/disinfecting

Note

The disinfectant must be of proven effectiveness (e.g. DGHM or FDA approval or CE mark).

Note

Only use cleaning agents that are suitable for this product class (rigid endoscope).

Note

For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that A_0 is $>3 \text{ 000}$ for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Ensure that the water can drain out through the openings.
- Clean and disinfect endoscopes with a sheath diameter of $\leq 4 \text{ mm}$ only in special Aesculap optics trays.
- Process the product in a washer/disinfectant machine. Follow the instructions of the washing machine's manufacturer.
- Avoid abrupt cooling of the product (e.g. in water).
- After mechanical cleaning/disinfecting, inspect all surfaces, cavities, lumens, and openings for visible debris.
- Clean manually if necessary.

MINOP optics and angled optics

- To remove residues of cleaning/disinfecting agents, always wipe distal window **1**, illumination surface **4** of the light cable connection, and eyepiece window **6** with a swap moistened with alcohol.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfector without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	2	D-W	-
II	Cleaning	55/131	10	D-W	neodischer MediClean forte 0.5 % (5 ml/L) pH >10
III	Rinse I	>10/50	1	D-W	-
IV	Rinse II	>10/50	1	FD-W	-
V	Thermal disinfecting	90/194	5	FD-W	-
VI	Drying	-	-	-	According to disinfector program

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: max. 10 microorganisms/ml, low endotoxin: max. 0.25 endotoxin units/ml)

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e. g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I drip lubricator JG598 or STERILIT® I oil spray JG600).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Check surfaces, cavities, lumens and openings for visible debris.
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- For products with locking mechanism (e.g. MINOP): Check the locking mechanism for smooth movement.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.

- Assemble the separable product, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected. Observe the individual weight limit for each tray.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

Note

The product must be disassembled prior to sterilization through the Sterrad® process.

Note

Alternating sterilization by different sterilization processes can result in damage to product materials and adapters.

The following processes can be used to sterilize the endoscope:

- Steam sterilization
- Sterrad® sterilization process: in Sterrad® sterilizers Sterrad® 50, Sterrad® 100S, Sterrad® 200

MINOP optics and angled optics

CAUTION

Damage to the optical system caused by flash sterilization!

- Do not flash sterilize the endoscope.
- Do not expose the endoscope to temperatures above 134 °C.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Only autoclave endoscopes marked with the symbol for steam sterilization on inscription ring 7.
- Do not autoclave damaged products.
- Protect the product against mechanical impacts.

Steam sterilization

Note

The product can be steam-sterilized either in disassembled or in assembled condition.

- Validated sterilization process
 - Disassemble the product, if applicable.
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/2 bar, holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterrad® sterilization process Sterrad® 50, Sterrad® 100S, Sterrad® 200

Note

The Sterrad® sterilization process can cause cosmetic changes of the product. These changes will, however, not affect its functionality.

- Disassemble the product.
- Only sterilize endoscopes marked with the SDS symbol for low-temperature and plasma sterilization on inscription ring 7 through the Sterrad® sterilization process 50, 100S or 200.
- Unscrew dismountable component parts of the endoscope, see Disassembling.
- Sterilize by the Sterrad® sterilization process 50, 100S or 200, observing the following rules:
Follow the advice of the manufacturers of the Sterrad® system. Use of a biological indicator is recommended to confirm effective sterilization.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization of the product may be accomplished by steam or Sterrad® process.
- Products may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270 °F–275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Technical Service

WARNING

Risk of injury and/or malfunction!

- **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- Prior to dispatching the product for repairs:
 - Have the product cleaned and disinfected or sterilized, and mark it as “disinfected” or “sterilized”, respectively.
 - Pack the endoscope in such a way that it will be protected against transport damage.

Service addresses

AESCULAP FLEXIMED GmbH
Robert-Bosch-Strasse 1
79211 Denzlingen / Germany
Phone: +49 (0) 7666 / 9321-0
Fax: +49 (0) 7666 / 9321-580
E-mail: techserv@aesculap-fleximed.de
www.aesculap-fleximed.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG