- (GB)** Instructions for use/Technical description
- (USA)** Neuropilot IV + EA
- (D)** Gebrauchsanweisung/Technische Beschreibung
Neuropilot IV + EA
- (F)** Mode d'emploi/Description technique
Neuropilot IV + EA
- (E)** Instrucciones de manejo/Descripción técnica
Neuropilot IV + EA
- (I)** Istruzioni per l'uso/Descrizione tecnica
Neuropilot IV + EA
- (P)** Instruções de utilização/Descrição técnica
Neuropilot IV + EA
- (NL)** Gebruiksaanwijzing/Technische beschrijving
Neuropilot IV + EA
- (S)** Bruksanvisning/Teknisk beskrivning
Neuropilot IV + EA
- (RUS)** Инструкция по применению/Техническое описание
Neuropilot IV + EA
- (CZ)** Návod k použití/Technický popis
Neuropilot IV + EA
- (PL)** Instrukcja użytkowania/Opis techniczny
Neuropilot IV + EA
- (SK)** Návod na používanie/Technický opis
Neuropilot IV + EA
- (TR)** Kullanım Kılavuzu/Teknik açıklama
Neuropilot IV + EA

B | BRAUN
SHARING EXPERTISE








Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 010662 03/13 V6 Änd.-Nr. 46931

CE – DIR 93/42/EEC

Technical alterations reserved

Aesculap®
Neuropilot IV + EA

Legend

- 1 Adjustment knob (direction X)
- 2 Adjustment knob (direction Z)
- 3 Adjustment knob (direction Y)
- 4 Reduction sleeve
- 5 Locking nut
- 6 Holding arm connector
- 7 Fixation knob (fixation)
- 8 Holding arm locking/release sleeve

Symbols on product and packages

Caution, general warning symbol

Caution, see documentation supplied with the product

Intended use

The Neuropilot IV + EA is used as a holder in neurosurgery to fix and fine tune the position of endoscopes or trocars. It is designed specially for endoscope-assisted surgery, but can also be used for fully endoscopic intraventricular procedures.

The Neuropilot IV + EA may only be used together with the holding arm (RT040R and FF16B). Adjustment knobs allow the position of the trocar or endoscope to be fine tuned in three different directions in Neuropilot IV + EA.

Safe handling and preparation

- CAUTION**
- Federal law restricts this device to sale by, or on order of a physician!**
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
 - Read, follow, and keep the instructions for use.
 - Use the product only in accordance with its intended use, see Intended use.
 - Clean the new product mechanically after removing its transport packaging and prior to its initial sterilization.
 - Store any new or unused products in a dry, clean, and safe place.
 - Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Replace any damaged components immediately with original spare parts.

Safe operation

- Ensure that the reduction sleeve fits the endoscope or trocar to be inserted.
- Follow the instructions for the endoscope or trocar.
- Rotate the fixation knob 7 in a counterclockwise direction as far as it will go, see Fig. 1.
- Slide the endoscope/trocar into the reduction sleeve 4 in Neuropilot IV + EA from above.
- Insert the endoscope into the reduction sleeve 4 as far as it will go to avoid adverse leverage effects and breakage of the glass lenses in the endoscope.
- Rotate the fixation knob 7 in a clockwise direction until the endoscope/trocar is seated securely.
The endoscope/trocar is now fixed.
- Prior to the operation, rotate the adjustment knobs 1, 2 and 3 so that the trocar spindle of Neuropilot IV + EA is center-positioned in all directions, see Fig. 2.
- During the operation, use the adjustment knobs 1, 2 and 3 to fine position the endoscope/trocar in Neuropilot IV + EA:
 - Use the adjustment knob 1 in direction X.
 - Use the adjustment knob 2 in direction Z.
 - Use the adjustment knob 3 in direction Y.

Disassembling

Removing the Neuropilot IV + EA from the holding arm

- Follow the instructions for use of the holding arm.
- Unscrew locking nut 5, see Fig. 3.
- Turn the sleeve 8 of the holding arm to release the quick-action coupling.
- Slide the sleeve 8 forward and remove the adapter 6 from the quick-action coupling on the support arm.

Removing the reduction sleeve

- Rotate the fixation knob 7 in a counterclockwise direction, see Fig. 4.
- Push the reduction sleeve 4 out of the opening from underneath and pull it out from above.

Assembling

Mounting the Neuropilot IV + EA on the holding arm

- Follow the instructions for use of the holding arm.
- Rotate the adjustment knobs 1, 2 and 3 so that the endoscope or trocar spindle of Neuropilot IV + EA is center-positioned in all directions.
- Unscrew the locking nut 5 in the opposite direction of the arrow.
- Slide the holding arm connector 6 into the quick-action coupling of the holding arm until the sleeve 8 snaps back, see Fig. 5.
- To secure the quick-release coupling against inadvertent opening, turn sleeve 8 of the support arm to the locking position.
- Screw on the locking nut 5 until Neuropilot IV + EA is seated securely.

Inserting the reduction sleeve

WARNING

Accidental pushing of the reduction sleeve into the operating field!

- Always insert the reduction sleeve in a proximal to distal direction.

Validated reprocessing procedure

General safety instructions

- Note*
- Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.
- Note*
- For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.
- Note*
- Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.
- Note*
- Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.
- The recommended chemistry was used for validation.
- Note*
- If there is no final sterilization, then a virucidal disinfectant must be used.
- Note*
- For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com
- The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.

- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.
- Turn the adjustment knobs to their maximum open position.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

DANGER

Risk to patients!

- The product must only be cleaned mechanically!

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for use, for example, on aluminum, plastic materials, and high-grade steel,
 - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable temperature of 60 °C.

- Do not use oxidizing chemicals (e.g. H₂O₂), which could cause bleaching/layer loss of the product.
- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection	Cleaning brush: e.g. TA007747 20 ml disposable syringe Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: Chapter Manual pre-cleaning with ultrasound and brush Chapter Mechanical alkaline cleaning and thermal disinfecting

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	–	D–W	–

D–W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I
► Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
► Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
► If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
► Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
► Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II
► Rinse/flush the product thoroughly (all accessible surfaces) under running water.
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	Concentrate, alkaline: pH = 13 <5 % anionic surfactant 0.5 % working solution pH = 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

Packaging

- Place the product in its holder or on a suitable tray.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note
The product may only be sterilized when dismantled (disassembled reduction sleeve).

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{–6}, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Or in the US:
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Reduction sleeve	Inner diameter	Compatible endoscopes/trocars
RT061R	4.0 mm	Angled neuroendoscopes PE486A, PE506A, PE526A
RT062R	6.2 mm	Ventriculoscope trocars FF370R, FF372R
RT063R	3.2 mm	Minop-trocar FF397R
RT064R	4.6 mm	Minop-trocar FF398R
RT065R	6.0 mm	Minop-trocar FF399R
RT066R	3.0 mm	PaediScope PF010A

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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