Aesculap Spine

GB	Instructions for use PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
USA	Instructions for use PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK <i>Note for U.S. users</i> <i>U.S. cleared device instructions for use can be found at</i> <i>www.aesculapimplantsystems.com "Products & Services".</i> <i>If you do not have internet access please call your local Aesculap</i> <i>representative or customer service at 1-800-258-1946.</i>
D	Gebruiksaanweisung PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
F	Mode d'emploi PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
E	Instrucciones de manejo PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
I	Istruzioni per l'uso PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
P	Instruções de utilização PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
RUS	Инструкция по применению PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
CZ	Návod k použití PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
PL	Instrukcja użytkowania PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
SK	Návod na používanie PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
CN	使用说明书 腰椎后路椎间融合器： A-Space PEEK， CeSpace PEEK， ProSpace PEEK， T-Space PEEK
KR	사용방법및사용시주의사항 추간체고정보형체 CESPAC PEEK CERVICAL IMPLANT, 형명개별기재)

B | BRAUN

SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 012095 08/11 V5 Änd.-Nr. 43592

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Símbolo CE, em conformidade com a Directiva 93/42/CEE
Маркировка CE согласно директиве 93/42/ЕЭС
CE-označení podle směrnice 93/42/EHS
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
Označenie CE podľa smernice 93/42/EHS
根据 93/42/EWG 指令的 CE 标志

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Мы оставляем за собой право вносить технические изменения
Technické zmeny vyhradeny
Zmiany techniczne zastrzeżone
Technické zmeny vyhradené
保留技术变动的权力

PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK

Intended use

PEEK cages are used for the following purposes:

- A-Space PEEK:
 - Stabilization of the lumbar spine through anterior, lateroanterior or lateral approach, mono- and multisegmental
- CeSpace PEEK:
 - Stabilization of the cervical spine C2-T1 through anterior approach mono- and multisegmental
- ProSpace PEEK:
 - Stabilization of the lumbar and thoracic spine through posterior approach, mono- and multisegmental
- T-Space PEEK:
 - Stabilization of the lumbar and thoracic spine through transforaminal approach, mono- and multisegmental

Materials

The materials used in the implant are listed on the packaging. The PEEK cages are made of PEEK-OPTIMA®. They contain marker pins made of tantalum to ensure radiological visibility for inspecting the implant position. PEEK-OPTIMA® is a registered trademark of Invivo Ltd, Lancashire FY5 4QD, UK.

Indications

Surgically installed implants serve to support normal healing processes. They should neither replace normal structures of the body nor permanently bear the loads occurring in the case of incomplete healing.

Use for:

- Degenerative instability
- Spondylolisthesis
- Post-discectomy syndrome
- Post-traumatic instabilities

Contraindications

Do not use in the presence of:

- Acute or chronic infections or severe defects of the osseous structures of the vertebral bodies, which need to be sound for the stable implantation of the PEEK cages
- Bone tumors in the region of implant fixation
- Unwillingness or inability of the patient to follow the instructions for postoperative treatment
- Any medical or surgical condition that could preclude the potential success of the implantation
- Pregnancy
- Osteoporosis or similar bone density loss
- Systemic or metabolic diseases
- Drug abuse or alcoholism
- Generally poor condition of the patient
- Adiposity
- Psychosocial issues; lack of co-operation by the patient
- Any conditions not listed under Indications

Side effects and interactions

The application or improper use of this system entails the following risks:

- Neurological complications caused by overdistraction or trauma of the nerve roots or dura
- Loss of intervertebral disk height due to removal of healthy bone material

Complications that can generally occur in connection with intervertebral surgery:

- Pseudarthrosis
- Incorrect implant position
- Spondylolisthesis
- Loss of fixation; dislocation or migration

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Delayed healing can result in implant rupture due to material fatigue.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.

Damage to implants caused by processing and resterilization!

- Do not process or resterilize the implants.

Application

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information regarding such training, please contact your national B. Braun/Aesculap agency.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating team are aware of information concerning the operating technique and the range of implants and associated instruments; this information is complete and ready to hand.
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific publications by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The intervention has been explained to the patient, whose consent concerning the following information has been documented:

- In the case of delayed or incomplete fusion, the implants can break and loosen due to high loads.
- The life-span of the implant depends on the patient's body weight.
- Corrective surgery may be necessitated by implant loosening, fracture or loss of correction.
- Smokers present an increased risk of bone fusion failure.
- The patient must undergo regular medical follow-up examinations of the implant components.

Implanting the PEEK cages

Increased risk of migration due to over-preparation of the vertebral body endplates!

- When preparing the implant bed, make certain that the base and cover plates of the adjacent vertebral bodies are not weakened.

- Select the appropriate PEEK implant size according to the indication, the preoperative planning and the bone situation found intraoperatively.
- Correctly apply the preparation instruments (rasps, curettes and chisels) for preparing the implant bed, as well as the implantation instrument.
- Apply appropriate care when inserting the implant.

A-Space PEEK

- Check implant height and angle using the trial implants.

Damage to the implant thread!

- Keep to the thread axis when screwing the implant onto the insertion instrument.
- Screw in the implant as far as it will go so that the dihedral of the instrument rests in the groove on the implant.

CeSpace PEEK

- Determine the appropriate implant size, using the trial implants.
- Apply the insertion instrument correctly. Pay attention to any inscriptions/labels.
- Use a cervical plate for additional stabilization.

ProSpace PEEK

Damage to the connector on the implant!

- When inserting the implant into the intervertebral space, avoid canting and take care to maintain an alignment parallel to the cover plate.

- Always implant two implants per layer (PLIF technique).
- Always use ProSpace PEEK in conjunction with an internal fixator.

T-Space PEEK

- Determine the implant height, using the trial impacts.
- Always use T-Space PEEK in conjunction with an internal fixator.

Further information on B. Braun/Aesculap implant systems is always available from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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