PREVISION®

|     |                                                                                                             |
|-----|-------------------------------------------------------------------------------------------------------------|
| GB  | <b>Instructions for use</b><br>PREVISION® revision hip prosthesis                                           |
| D   | <b>Gebrauchsanweisung</b><br>PREVISION®-Revisionshüftprothesenschaft                                        |
| F   | <b>Mode d'emploi</b><br>PREVISION® Prothèse de révision                                                     |
| E   | <b>Instrucciones de manejo</b><br>PREVISION® Vástago de revisión                                            |
| I   | <b>Istruzioni per l'uso</b><br>PREVISION® Stelo da Revisione                                                |
| P   | <b>Instruções de utilização</b><br>PREVISION® Haste de prótese de revisão                                   |
| RUS | <b>Инструкция по применению</b><br>Ножка протеза тазобедренного сустава для ревизионных операций PREVISION® |
| CZ  | <b>Návod k použití</b><br>PREVISION® Revizní dílek endoprotézy kyčelního kloubu                             |
| PL  | <b>Instrukcja użytkowania</b><br>PREVISION®-Rewizyjny trzpień endoprotezy biodrowej                         |
| SK  | <b>Návod na používanie</b><br>PREVISION® drielk revíznej bedrovej protézy                                   |
| KR  | 사용방법및사용시주의사항<br>인공영덩이뼈관절 (Prevision® Modular Revision Hip system, 형명별개별기재 )                                 |

B | BRAUN  
SHARING EXPERTISE

Aesculap AG

Am Aesculap-Platz  
78532 Tuttlingen  
Germany  
Phone +49 7461 95-0  
Fax +49 7461 95-2600  
www.aesculap.de

TA-Nr. 010988 02/11 Änd.-Nr. 40983

CE marking according to directive 93/42/EEC  
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG  
Marquage CE conforme à la directive 93/42/CEE  
Identificación CE en conformidad con la directriz 93/42/CEE  
Marchio CE conforme alla direttiva 93/42/CEE  
Simbolo CE, em conformidade com a Directiva 93/42/CEE  
Маркировка CE согласно директиве 93/42/ЕЭС  
CE-označení podle směrnice 93/42/EHS  
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE  
Označenie CE podľa smernice 93/42/EHS

Technical alterations reserved  
Technische Änderungen vorbehalten  
Sous réserve de modifications techniques  
Sujeto a modificaciones técnicas  
Con riserva di modifiche tecniche  
Salvo alterações técnicas  
Мы оставляем за собой право вносить технические изменения  
Technické zmeny vyhradeny  
Zmiany techniczne zastrzeżone  
Technické zmeny vyhradené

### Intended use

- The implant is used
- as a component of a human hip endoprosthesis: Hip endoprosthesis stem
  - to combine with Aesculap hip endoprosthesis components
  - for implantation without bone cement
  - for revision operations
- Modular implant components:
- Proximal component (PLASMAPORE®  $\mu$ -CaP) with spring nut and locking screw
  - Distal stem component
    - PREVISION® curved
    - or –
    - PREVISION® straight

### Materials

The materials used in the implant are listed on the packaging:

| Materials                                                                                                                   |
|-----------------------------------------------------------------------------------------------------------------------------|
| ISOTAN®; titanium forged alloy Ti6Al4V according to ISO 5832-3                                                              |
| PLASMAPORE® $\mu$ -CaP surface layer made of pure titanium acc. to ISO 5832-2, with an additional calcium phosphate coating |

ISOTAN® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

### Indications

Use in revision surgery to replace hip endoprosthesis shafts:

- Cement-free revision operations
- Periprosthetic fractures
- Subprosthetic fractures
- Treatment of large bone defects on the proximal femur

### Contraindications

The replacement of an endoprosthetic implant is a difficult intervention subject to individual conditions.

To arrive at a decision on whether to provide such replacement, the operating surgeon must assess and, if necessary, weigh the contraindications with regard to the success of the intervention.

Do not use in the presence of:

- Acute or chronic infections near the joint, or systemic infections
- Secondary diseases that could influence joint implant functionality
- Systemic diseases and metabolic disorders
- Acute osteoporosis or osteomalacia
- Severe damage to the bone structures that could prevent the stable implantation, or a sufficient fixation, of the implant components
- Non-regenerative bone conditions with lack of proximal bone support and failure of defect union
- Bone tumors in the region of implant fixation
- Anticipated excessive load on the joint implant
- Anticipated excessive load on the joint implant, especially due to high patient weight or heavy physical strains and activities
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials
- Prosthesis heads with neck length XXL

### Side-effects and adverse interactions

- Changes in position, loosening, wear and tear on, or fracture of implant components
- Secondary operation to remove the distal locking screw (if applicable)
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Venous thrombosis, lung embolism, and cardiac arrest
- Tissue reaction to implant materials
- Injury to blood vessels and nerves
- Hematomas and wound healing disorders
- Periarticular calcification
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise

#### Note

*There are no known side effects or interactions in connection with standard MRI scans up to 3 tesla.*

### Safety notes

- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon is responsible for combining the implant components and implanting them without bone cement.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The instructions for use of the respective Aesculap implant components must be observed.
- Only combine Aesculap modular implant components with each other.

- Do not use damaged or surgically excised components under any circumstances.
- Implants that have been used before must not be reused.
- The operating surgeon is responsible for combining and implanting the implant components, with or without bone cement.
- Delayed or failed healing of the bony implant bed can lead to implant breakage due to excessive strain on the implant material.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Damage to load-transferring bone structures can result in loosening of the components, bone or implant fractures or other severe complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- Damage to the load-bearing structures of the implant may result in loosening of components, their dislocation and migration, and other severe complications.
- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- Avoid damaging the implant, especially in the neck or cone region, by applying instruments (e.g. HF surgical devices) close to the implant.

**CAUTION**

**Risk of implant component breakage due to combination with implant components from other manufacturers!**

- Use Aesculap implant components only.

### Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiry date or whose packaging is damaged.

**WARNING**

**Damage to implants caused by processing and resterilization!**

- Do not process or resterilize the implants.

### Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- Implantation instruments, including the dedicated Aesculap implant system instruments, must be complete and in working condition
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical intervention has been explained to the patient and their agreement concerning the following information has been recorded:

- The revision of a hip endoprosthesis stem is a complex joint-replacement operation; it is generally inferior to the primary joint replacement.
- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
- The prosthetic joint can come loose due to excessive load, wear and tear, or infection.
- The lifespan of the prosthetic joint depends on the body weight and the load put on the joint.
- The prosthetic joint must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- More revision surgery may become necessary if the implant loosens.
- In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
- The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.
- Where distal locking screws have been used, a second operation to remove the metal parts is necessary.
- The patient must be informed of the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The risks of transgressing these rules must be explained to the patient.

The implantation involves the following steps:

- Remove, via a suitable access, the hip endoprosthesis stem to be revised.
- Remove bone cement residues, granulation tissue and, if necessary, avital bone parts.
- A careful operating technique helps preventing further damage to the bone.
- **PREVISION® curved:**  
To prepare the implant bed, use distal intramedullary drills and modular shaping rasps (with proximal and distal components), starting from distal with the marrow cavity drills and proceeding with the modular shaping rasps.

- **PREVISION® straight:**  
To prepare the implant bed, use reamers and proximal shaping rasps (assembled with trial components), starting the preparation with the reamers and proceeding with the shaping rasps.
- Note the depth markings on the instruments when preparing the bone.
- Select implant components (distal and proximal) according to the shaping rasps last inserted in the correct position (PREVISION® curved) or to the proximal shaping rasp with the trial component (PREVISION® straight).
- Implant the distal implant components, taking note of the depth markings (orientate to the intended head center or trochanter major, respectively).
- Carry out distal locking with two locking screws, if necessary.
- Use Aesculap 5-mm locking screws for this step. Observe the instructions for use of the locking screws.
- Perform test reposition with the proximal test implant; check for joint mobility, stability and leg length.
- Remove the test implant; carefully rinse, clean and dry the cone coupling of the distal implant.

**CAUTION**

**Implant breakage or loosening of implant couplings caused by contaminated or damaged coupling surfaces!**

- Check coupling surfaces for cleanness and damage.
- Introduce the proximal implant carefully and in proper axial alignment.

- Insert the modular proximal implant component. Make certain the component is implanted at the correct position (axial and rotational) and implantation depth.
- Insert the spring nut into the proximal implant and screw it on hand-tight, using the socket of the torque wrench.
- Attach the adapter rod of the counterholder to the distal implant and allow it to engage.
- Apply the torque wrench and couple the counterholder to the adapter rod.
- Fixate the spring nut, making sure the counterholder is not twisted in the process.

**CAUTION**

**Insufficient coupling of the implant components due to the spring nut tightened too firmly or not firmly enough!**

- Adhere to the application instructions for the torque wrench.
- Have the torque wrench inspected by the manufacturer every 2 years.

- Tighten the spring nut until the torque wrench engages audibly.
- Remove the torque wrench and the counterholder.
- Screw in the locking screw.

#### Notes on alternative application techniques

For the transfemoral approach, the following points must be observed:

- Proximal bone portions must remain connected to the soft tissue.
- Usually, only distal intramedullary drills (PREVISION® curved) or distal reamers are used for this approach.
- Distal implant components must be locked in position.
- The distal fixation length must exceed 10 cm.

When assembling the implant outside the bone, the following points have to be noted:

- The sizes of the implant components are selected as in the procedure described above.
- The rotational position of the proximal implant is adjusted analog to the in-vivo position of the shaping rasp (PREVISION® curved).
- The procedure for mounting the implant components and using the torque wrench remain unchanged.
- Insert the implant with the insertion tool.

- After implantation, apply the torque wrench again to check that the implant components are properly tightened.
- Carry out another test reduction, checking the joint mobility, the joint stability and the leg length.
- Inspect and document the implant position during the operation (using an image converter).
- Remove the protective cap of the prosthesis cone only immediately before installing the cone.
- Rinse, clean and dry the outside cone of the stem and, if necessary, the inside cone of the prosthesis spheres.
- Install the prosthesis head.
- Check that the cone size of the prosthesis head in use corresponds to the cone plug-in connector of the prosthesis stem (cone size 12/14).
- Couple the prosthesis head and cone at room temperature only.
- If necessary, allow the implant to cool to room temperature.
- In order to forestall abnormal wear and tear on the prosthesis: Remove all open bone and, if necessary, bone cement residues before closing the wound.

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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