- GB

Instructions for use/Technical description
- USA

PREVISION® revision hip endoprosthesis stem
- D

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Intended use

The implant is used

- as a component of a human hip endoprosthesis: Hip endoprosthesis stem
- to combine with Aesculap hip endoprosthesis components
- for implantation without bone cement
- for revision operations

Modular implant components:

- Proximal component (PLASMAPORE®_{μ-CaP} coating) with spring nut and locking screw
- Distal stem component
 - PREVISION® curved
 - PREVISION® straight

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_F titanium forged alloy Ti6Al4V acc. to ISO 5832-3
 - PLASMAPORE®_{μ-CaP} surface coating pure titanium acc. to ISO 5832-2 with additional calcium phosphate coating
- ISOTAN® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for acute disorders of the hip joint that cannot be treated by other therapies:

- Degenerative osteoarthritis
- Rheumatic arthritis
- Fractures of the joint
- Femoral head necrosis

Use in operations to replace hip endoprosthesis stems:

- Cement-free revision operations
- Periprosthetic fractures
- Subprosthetic fractures
- Treatment of large bone defects on the proximal femur

Absolute contraindications

Do not use in the presence of:

- Joint diseases that can be treated by reconstructive surgery (e.g. displacement osteotomy)
- Acute or chronic infections near the joint, or systemic infections
- Secondary diseases that could influence joint implant functionality
- Acute osteoporosis or osteomalacia
- Severe damage to the bone structures that could prevent the stable implantation, or sufficient fixation, of the implant components
- Non-regenerative bone conditions with lack of proximal bone support and failure of defect union
- Bone tumors in the region of implant fixation
- Bone malformations, defective axial positioning, or other bone conditions that rule out implantation of a hip-joint prosthesis
- Foreign body sensitivity to the implant materials
- Prosthesis heads with neck length XXL

Relative contraindications

The following conditions, individual or combined, can compromise the success of the operation:

- Systemic diseases and metabolic disorders
- Expected overload on the joint implant, especially due to high patient weight or major physical strains and activities
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance

Side effects and interactions

- Changes in position, loosening, wear and tear on, or fracture of implant components
- Secondary operation to remove the distal locking screws (if applicable)
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Venous thrombosis, lung embolism, cardiac arrest
- Tissue reaction to implant materials
- Injury to blood vessels and nerves
- Hematomas and wound healing disorders
- Periarticular calcification
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise

Safety notes

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon is responsible for combining the implant components and implanting them without bone cement.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or lack of asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Combine modular implant components only with suitable Aesculap hip endoprostheses.
- Damage to load-bearing bone structures can result in loosening or migration of the components, bone or implant fractures, and other acute complications. To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Delayed or failed healing of the bony implant bed can lead to implant breakage due to excessive strain on the implant material.

WARNING

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

- A precise diagnosis requires X-rays taken in the anterior-posterior and medial-lateral directions.
- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- Avoid damage to the implant, especially in the neck or cone region, caused by instruments (e.g. HF surgical devices) applied close to the implant.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are sterilized by irradiation.
- ▶ Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- ▶ Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- ▶ Do not use implant components that are past their expiration date or whose packaging is damaged.

WARNING

Damage to implants caused by processing and resterilization!

- ▶ Do not reprocess or resterilize the implants.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
 - All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
 - The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
 - The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
 - Information has been obtained from the manufacturer in the event of an ambiguous preoperative situation and if implants are present in the area to be treated
- The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:
- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
 - The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
 - The prosthetic joint can come loose due to excessive load, wear and tear, or infection.
 - The lifespan of the prosthetic joint depends on the patient's body weight and the strain to which the joint is subjected.
 - The prosthetic joint must not be overstrained by extreme loads, hard physical labor or sporting activities.
 - Corrective surgery may become necessary if the implant loosens.
 - In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
 - The revision of a hip endoprosthesis stem is a complex joint replacement operation.
 - Generally, the revision joint replacement will be inferior in its performance to the primary joint replacement.
 - The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.

The implantation and implantation site are prepared in the following way:

- ▶ Remove the hip endoprosthesis stem to be revised via a suitable access point, if necessary.
- ▶ Remove bone cement residues, granulation tissue, and, if necessary, devitalized bone parts.
- ▶ Employ a careful operating technique in order to prevent further damage to the bone.
- ▶ PREVISION® curved: To prepare the implant bed, use distal intramedullary drills and modular form rasps (with proximal and distal components); begin the treatment distally, using the smallest intramedullary drill to start with and getting steadily larger, and proceed with the modular form rasps.
- ▶ PREVISION® straight: To prepare the implant bed, use reamers and proximal form rasps (assembled with trial components); begin the treatment using the smallest reamer, getting steadily larger, and proceed with the form rasps.
- ▶ Note the depth markings on the instruments when preparing the bone.
- ▶ Select implant components (distal and proximal) according to the form rasps (PREVISION® curved) or proximal form rasp and trial component (PREVISION® straight) last inserted in the correct position.
- ▶ Implant the distal implant components, taking note of the depth markings (orientate to the intended head center or trochanter major, respectively).
- ▶ Carry out distal locking of the distal stem component with two locking screws, if necessary.
- ▶ Use Aesculap 5 mm locking screws for this step. Observe the instructions for use of the locking screws.
- ▶ Perform trial reduction with the proximal trial implant; check for joint mobility, stability, and leg length.
- ▶ Remove the trial implant; carefully rinse, clean, and dry the cone coupling of the distal implant.

CAUTION

Implant breakage or loosening of implant couplings caused by contaminated or damaged coupling surfaces!

- ▶ Check coupling surfaces for contamination and damage.
- ▶ Introduce the proximal implant carefully and in proper axial alignment.

- ▶ Insert the modular proximal implant component. Make certain that the component is implanted at the correct position (axial and rotational) and implantation depth.
- ▶ Insert the spring nut into the proximal implant and screw it on hand-tight, using the socket of the torque wrench.
- ▶ Attach the adapter rod of the counterholder to the distal implant and allow it to engage.
- ▶ Apply the torque wrench and couple the counterholder to the adapter rod.
- ▶ Fix the spring nut, making sure that the counterholder is not twisted in the process.

CAUTION

Insufficient coupling of the implant components due to the spring nut being tightened too firmly or not firmly enough!

- ▶ Follow the instructions for use of the torque wrench.
- ▶ Have the torque wrench inspected by the manufacturer every 2 years.

- ▶ Tighten the spring nut until the torque wrench releases audibly.
- ▶ Remove the torque wrench and the counterholder.
- ▶ Screw in the locking screw.

Notes on alternative application techniques

For the transfemoral approach, the following points must be observed:

- Proximal bone portions must remain connected to the soft tissue.
- Usually, only distal intramedullary drills (PREVISION® curved) or distal reamers are used for this approach.
- Distal implant components must be locked in position.
- The distal anchoring length of the implant must be greater than 10 cm.

When assembling the implant outside the bone, the following points must be observed:

- The sizes of the implant components are selected as in the procedure described above.
- The rotational position of the proximal implant is adjusted in the same way as the in-vivo position of the form rasp (PREVISION® curved).
- The procedure for assembling the implant components and using the torque wrench remain unchanged.
- Insert the implant with the insertion tool.
- After implantation, apply the torque wrench again to check that the implant components are properly tightened.
- Carry out another trial reduction, checking the joint mobility, stability, and leg length.
- Inspect and document the implant position during the operation (using an image converter).
- Remove the protective cap from the prosthesis cone only immediately before installing the cone.
- Rinse, clean, and dry the outside cone of the stem and, if necessary, the inside cone of the prosthesis spheres.
- Check that the cone size of the prosthesis head in use corresponds to the cone plug-in connector of the prosthesis stem (cone size 12/14).
- Install the prosthesis head.
- Attach the prosthesis head and prosthesis cone to each other at room temperature only.
- If necessary, allow the implants to cool to room temperature.
- To prevent abnormal wear and tear of the prosthesis: Remove all open bone and, if necessary, bone cement residues before closing the wound.

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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