- (GB) Instructions for use/Technical description
 (USA) Tubular shank instruments
 (D) Gebrauchsanweisung/Technische Beschreibung
 Rohrschaftinstrumente
 (F) Mode d'emploi/Description technique
 Instruments à chemise
 (E) Instrucciones de manejo/Descripción técnica
 Instrumentos de vástago tubular
 (I) Istruzioni per l'uso/Descrizione tecnica
 Strumenti a stelo tubolare
 (P) Instruções de utilização/Descrição técnica
 Instrumentos tubulares
 (NL) Gebruiksaanwijzing/Technische beschrijving
 Buisschaftinstrumenten
 (S) Bruksanvisning/Teknisk beskrivning
 Rörskafinstrument
 (RUS) Инструкция по применению/Техническое описание
 Тубусные инструменты
 (CZ) Návod k použití/Technický popis
 Nástroje s trubkovým dílkem
 (PL) Instrukcja użytkowania/Opis techniczny
 Rozkładane Instrumenty chwytające
 (SK) Návod na používanie/Technický opis
 Nástroje s rúrovou náradou
 (TR) Kullanım Kılavuzu/Teknik açıklama
 Boru shaft ekipmanları

B|BRAUN

SHARING EXPERTISE


Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
 Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 012978 11/12 V6 Änd.-Nr. 45456

Technical alterations reserved

Aesculap®

Tubular shank instruments

Legend

- 1 Rotation wheel
- 2 Handle
- 3 Button
- 4 Movable handle
- 5 Shaft, complete (outer tube, inner tube, working insert)
- 6 Instrument jaw
- 7 Working insert
- 8 Alignment grooves
- 9 Inner tube
- 10 Outer tube
- 11 Spring elements (tactile feedback)

Symbols on product and packages

Caution, general warning symbol
Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

Detachable tubular shaft instruments are used for cutting, dissection and grasping, and removal of biopsies, with different working tips for each intended use. They are introduced into the body by means of access instruments (e.g. trocar).

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).
- If using in combination with a HF instrument within a trocar, insert the tubular shaft instrument first, where possible, to prevent compromising the HF insulation. When removing the instruments, remove the HF instrument first, if possible.

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

WARNING

Risk of injury when using the product beyond the field of view!

- Apply the product only under visual control.

Note

The spring elements give the user some tactile feedback just before the instrument emerges from trocar FF399R into the operating site. The functionality of this tactile feedback mechanism is only ensured when using the product with trocar FF399R.

Disassembling

The tubular shaft instruments are disassembled into the following components:

- Handle 2
- Outer tube 10
- Inner tube 9
- Working insert 7
- Detach handle 2 from working insert 7, see Fig. 1.
 - Fully open moveable handle part 4.
 - Push and hold down button 3.
 - At the same time, pull down moveable handle part 4.
- Detach the shaft 5 from handle 2, see Fig. 2
 - Pull back rotation wheel 1 towards handle 2 to the positive stop.
 - Hold rotation wheel 1 at its positive stop and extract shaft 5.
- To disassemble shaft 5, see Fig. 3:
 - Remove outer tube 10 and inner tube 9 from working insert 7.
 - Remove inner tube 9 from working insert 7, see Fig. 4.

Assembling

- To assemble shaft 5, see Fig. 5:
 - Slide inner tube 9 with alignment grooves 8 in the direction of the working tip 7 until the alignment grooves 8 engage.
 - Slide outer tube 10 over inner tube 9 onto working insert 7 as far as it will go, see Fig. 6.

- Connect the shaft 5 to handle 2:
 - Ensure that outer tube 10 is pushed over inner tube 9 with working insert 7 to the positive stop, see Fig. 7.
 - Slide back and hold rotation wheel 1 at its positive stop, see Fig. 8.
 - Hold the shaft 5 at its closed working tip.
 - Slightly turn shaft 5 and handle 2 during insertion until they connect to each other.
 - Allow the rotation wheel 1 to slide forward, see Fig. 9. Ensure that shaft 5 is securely locked in handle 2 and cannot be removed even by pulling.
- Connect handle 2 entirely to shaft 5:
 - Fully open handle 2.
 - Slide the movable part of the handle 4 upwards whilst simultaneously holding the working insert 7 until the button 3 engages, see Fig. 10.
 - Check that moveable handle part 4 is securely locked.

Validated reprocessing procedure

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open up products with hinges.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which,
 - are approved for the material in question (e.g., aluminum, plastics, high-grade steel).
 - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.

- Do not use oxidizing chemicals (e.g. H₂O₂), which could cause bleaching/layer loss of the product.
- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection	Cleaning brush: e.g., TAO12889 20 ml disposable syringe Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). Connect components with lumens and channels directly to the rinsing port of the injector carriage. Keep working ends open for cleaning. Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: Chapter Manual pre-cleaning with a brush Chapter Mechanical alkaline cleaning and thermal disinfecting

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	–	D–W	–

D–W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.
- Phase I

► Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.

► Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.

► If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

► Mobilize non-rigid components, such as set screws, links, etc. during cleaning.

► Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: pH = 13 <5 % anionic surfactant 0.5 % working solution pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	–
IV	Thermal disinfecting	90/194	5	FD-W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

► Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.

- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note
The product can be sterilized either in disassembled or in assembled condition.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process

– Steam sterilization using fractional vacuum process

– Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665

– Sterilization using fractional vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{–6}, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

► Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Or in the US:
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

For MINOP inVent tubular shaft instruments FH635R, FH636R, FH637R, FH638R, FH639R

Art. no.	Designation
FH633R	MINOP inVent hand grip for tubular shaft instruments
FH635200	MINOP inVent tubular shaft instruments outer tube
FF435R	Jaw part with MINOP micro-scissors Sh/Sh
FF436R	Jaw part with MINOP micro-scissors BI/BI
FF437R	Jaw part with MINOP micro cupped biopsy forceps
FF438R	Jaw part with MINOP micro grasping and preparation forceps
FF439R	Jaw part with MINOP toothed micro forceps

For MINOP tubular shaft instruments FF385R, FF386R, FF387R, FF388R, FF384R

Art. no.	Designation
FF432R	Handle only for instruments FF385R–FF389R
FF433R	Outer tube only for FF385R–FF389R
FF435R	Jaw part with MINOP mirco-scissors Sh/Sh
FF436R	Jaw part with MINOP mirco-scissors BI/BI
FF437R	Jaw part with MINOP micro cupped biopsy forceps
FF438R	Jaw part with MINOP micro grasping and preparation forceps
FF439R	Jaw part with MINOP toothed micro forceps

Disposal

► Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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