- (GB)** Instructions for use/Technical description  
S<sup>4</sup> Element instruments
- (USA)** S<sup>4</sup> Element-Instrumente
- (D)** Gebrauchsanweisung/Technische Beschreibung  
S<sup>4</sup> Element-Instrumente
- (F)** Mode d'emploi/Description technique  
Instruments S<sup>4</sup> Element
- (E)** Instrucciones de manejo/Descripción técnica  
Instrumental S<sup>4</sup> Element
- (I)** Istruzioni per l'uso/Descrizione tecnica  
Strumenti S<sup>4</sup> Element
- (P)** Instruções de utilização/Descrição técnica  
Instrumentos S<sup>4</sup> Element
- (NL)** Gebruiksaanwijzing/Technische beschrijving  
S<sup>4</sup> Element-Instrumenten
- (S)** Bruksanvisning/Teknisk beskrivning  
S<sup>4</sup> Element-instrument
- (RUS)** Инструкция по применению/Техническое описание  
Инструменты S<sup>4</sup> Element
- (CZ)** Návod k použití/Technický popis  
Nástroje S<sup>4</sup> Element
- (PL)** Instrukcja użytkowania/Opis techniczny  
Instrumenty S<sup>4</sup> Element
- (SK)** Návod na používanie/Technický opis  
Nástroje S<sup>4</sup> Element
- (TR)** Kullanım Kılavuzu/Teknik açıklama  
S<sup>4</sup> Element-Araçları

**B | BRAUN**  
SHARING EXPERTISE


Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

TA-Nr. 013272    01/13    V6    Änd.-Nr. 46033

**CE** 0123 – DIR 93/42/EEC

Technical alterations reserved

# Aesculap®

## S<sup>4</sup> Element instruments

### Legend

- 1 S<sup>4</sup> Element monoaxial screwdriver (FW276R)
- 2 S<sup>4</sup> Element polyaxial screwdriver (FW277R)
- 3 S<sup>4</sup> Element counter torque (FW283R)

### Symbols on product and packages

|  |                                                                                         |
|--|-----------------------------------------------------------------------------------------|
|  | Caution, general warning symbol<br>Caution, see documentation supplied with the product |
|  | Date of manufacture                                                                     |

### Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at [www.extranet.bbraun.com](http://www.extranet.bbraun.com)

### Intended use

The S<sup>4</sup> Element instruments are used for the implantation of S<sup>4</sup> Element implants.

### Safe handling and preparation

**CAUTION**  
Federal law restricts this device to sale by, or on order of a physician!

**Risk of injury caused by incorrect operation of the product!**

- Adhere to the instructions for use of any medical device.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).

### Safe operation

**Risk of injury and/or malfunction!**

- Always carry out a function check prior to using the product.

**Damage to the implant through use of the incorrect screwdriver!**

- Use monoaxial screws only with monoaxial screwdriver FW276R.
- Use polyaxial screws only with polyaxial screwdriver FW277R.

- Note**

Do not use screwdriver FW277R or FW276R to apply bending forces or leverage.
- Note**

Do not overtighten the threaded shaft of FW277R or FW276R or tighten with mechanical aids.
- Note**

Do not use screwdriver FW277R/FW276R with torque wrenches or in combination with locking screws.
- Note**

The tips of the screwdrivers FW277R and FW276R are brittle and should not be used as a lever.
- Note**

Instrument parts detached from the instrument body should be removed from the patient because they are not implantable.

### Validated reprocessing procedure

#### General safety instructions

- Note**

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.
- Note**

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.
- Note**

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.
- Note**

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

- Note**

If there is no final sterilization, then a virucidal disinfectant must be used.
- Note**

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at [www.extranet.bbraun.com](http://www.extranet.bbraun.com)

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

#### General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary. Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-l.org](http://www.a-k-l.org), link to Publications, Red Brochure – Proper maintenance of instruments.

#### Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

#### Cleaning/disinfection

##### Product-specific safety notes on the reprocessing procedure

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
  - are approved for plastic materials,
  - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
  - as an effective mechanical supplement to manual cleaning/disinfecting.
  - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
  - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
  - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

#### Validated cleaning and disinfection procedure

| Validated procedure                                                                                                                                                                                                           | Specific requirements                                                                                                                                                                                                                                                                                                                                                                                                  | Reference                                                                                                                                                                                                                                             |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection <ul style="list-style-type: none"><li>■ Counter torque:<ul style="list-style-type: none"><li>– FW283R</li></ul></li></ul>                                                         | <ul style="list-style-type: none"><li>■ Cleaning brush: e.g. TA011944, TE654202 or GK469200</li><li>■ 20 ml disposable syringe</li><li>■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.</li><li>■ Drying phase: Use a lint-free cloth or medical compressed air</li></ul>                                                    | Chapter Manual cleaning/disinfection and sub-chapter: <ul style="list-style-type: none"><li>■ Chapter Manual cleaning with immersion disinfection</li></ul>                                                                                           |
| Manual cleaning with ultrasound and immersion disinfection <ul style="list-style-type: none"><li>■ Screwdriver:<ul style="list-style-type: none"><li>– FW276R</li><li>– FW277R</li></ul></li></ul>                            | <ul style="list-style-type: none"><li>■ Cleaning brush: e.g. TA011944, TE654202 or GK469200</li><li>■ 20 ml disposable syringe</li><li>■ Keep working ends open for cleaning.</li><li>■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.</li><li>■ Drying phase: Use a lint-free cloth or medical compressed air</li></ul>     | Chapter Manual cleaning/disinfection and sub-chapter: <ul style="list-style-type: none"><li>■ Chapter Manual cleaning with ultrasound and immersion disinfection</li></ul>                                                                            |
| Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection <ul style="list-style-type: none"><li>■ Counter torque:<ul style="list-style-type: none"><li>– FW283R</li></ul></li></ul> | <ul style="list-style-type: none"><li>■ Cleaning brush: e.g. TA011944, TE654202 or GK469200</li><li>■ 20 ml disposable syringe</li><li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li><li>■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li><li>■ Place products in the tray with their hinges open.</li></ul> | Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: <ul style="list-style-type: none"><li>■ Chapter Manual pre-cleaning with a brush</li><li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li></ul> |

| Validated procedure                                                                                                 | Specific requirements                                                                                                                                                                                                                                                                                                                                                                                              | Reference                                                                                                                                                                                                                                                         |
|---------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection | <ul style="list-style-type: none"> <li>Cleaning brush: e.g. TA011944, TE654202 or GK469200</li> <li>20 ml disposable syringe</li> <li>Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> <li>Place products in the tray with their hinges open.</li> </ul> | Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: <ul style="list-style-type: none"> <li>Chapter Manual pre-cleaning with ultrasound and brush</li> <li>Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul> |

## Manual cleaning/disinfection

### Manual cleaning with immersion disinfection

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfecting cleaning | RT (cold)    | >15        | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse    | RT (cold)    | 1          | –            | D–W              | –                                                              |
| III   | Disinfection          | RT (cold)    | 15         | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse           | RT (cold)    | 1          | –            | FD–W             | –                                                              |
| V     | Drying                | RT           | –          | –            | –                | –                                                              |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
RT: Room temperature  
\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

#### Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

#### Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

### Manual cleaning with ultrasound and immersion disinfection

| Phase | Step                | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|---------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Ultrasonic cleaning | RT (cold)    | >15        | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse  | RT (cold)    | 1          | –            | D–W              | –                                                              |
| III   | Disinfection        | RT (cold)    | 15         | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse         | RT (cold)    | 1          | –            | FD–W             | –                                                              |
| V     | Drying              | RT           | –          | –            | –                | –                                                              |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
RT: Room temperature  
\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

#### Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

#### Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

## Mechanical cleaning/disinfection with manual pre-cleaning

#### Note

*The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).*

#### Note

*The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.*

### Manual pre-cleaning with a brush

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfectant cleaning | RT (cold)    | >15        | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing               | RT (cold)    | 1          | –            | D–W              | –                                                              |

D–W: Drinking water  
RT: Room temperature  
\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

### Manual pre-cleaning with ultrasound and brush

| Phase | Step                | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|---------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Ultrasonic cleaning | RT (cold)    | >15        | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing             | RT (cold)    | 1          | –            | D–W              | –                                                              |

D–W: Drinking water  
RT: Room temperature  
\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting  
Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                      | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical                                                                                                                                                                                                                                   |
|-------|---------------------------|--------------|------------|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse                  | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                          |
| II    | Cleaning                  | 55/131       | 10         | FD-W             | <div><div></div>Concentrate, alkaline:<ul style="list-style-type: none"><li>pH = 13</li><li>&lt;5 % anionic surfactant</li></ul><div><div></div>0.5 % working solution<ul style="list-style-type: none"><li>pH = 11*</li></ul></div></div> |
| III   | Intermediate<br>rinse     | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                          |
| IV    | Thermal disin-<br>fecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                          |
| V     | Drying                    | -            | -          | -                | According to the program for cleaning<br>and disinfection device                                                                                                                                                                           |

D-W: Drinking water  
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

Packaging

- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10<sup>-6</sup>, Aesculap recommends the following parameters:

| Aesculap Orga Tray/Sterile container (perforated bottom)<br>Minimum cycle parameters* |               |       |                     |
|---------------------------------------------------------------------------------------|---------------|-------|---------------------|
| Sterilization method                                                                  | Temp.         | Time  | Minimum drying time |
| Prevacuum                                                                             | 270 °F/275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

**WARNING for the US market**

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

**Risk of injury and/or malfunction!**

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

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Phone: +49 (7461) 95-1602  
Fax: +49 (7461) 16-5621  
E-Mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems LLC  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042  
Aesculap Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420  
Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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