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
Aesculap®

S⁴® injection cannula for augmentation screws SR146SU, S⁴® percutaneous injection cannula for augmentation screws SR148SU

Legend

- 1 Injection cannulas for open access SR146SU
- 2 Injection cannulas for percutaneous access SR148SU
- 2a Cannula
- 2b Pestle

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

The S⁴ Spinal System implants are used for dorsal monosegmental and multisegmental stabilization of the lumbar and thoracic spine.
The S⁴ Spinal System injection cannula is used for applying bone cement into the S⁴ Spinal System augmentation screw.

Safe handling and preparation

CAUTION
Federal law restricts this device to sale by, or on order of a physician!

Note
Generally, the S⁴ Spinal System–augmentation screw is covered by the instructions for use of the S⁴ Spinal System – Lumbar/Deformity TA011187, as well as the S⁴ Spinal System augmentation screw TA012865!
The present information on the S⁴ Spinal System injection cannula supplement the respective statement in the instructions for use for the S⁴ Spinal System Lumbar/Deformity and augmentation screw.

Note
The S⁴ Spinal System injection cannula is only designed for single use in combination with the S⁴ Spinal System augmentation screw.

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!
► Do not reprocess the product.

The product is gamma-sterilized and supplied in sterile packaging.
The product must not be reused.
► Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
► Read, follow, and keep the instructions for use.
► Use the product only in accordance with its intended use, see Intended use.
► Do not use products from open or damaged sterile packaging.
► Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
► Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
► Do not use the product after its use-by date.

Safe operation

Risk of injury and/or malfunction!
► Always carry out a function check prior to using the product.

Risk of injury from cement leakages!
► Screw the injection cannula into the screw head without applying force.
► Make certain the axes of the pedicle screw and the injection cannula are aligned with each other.
► Cement application only under high-quality radiographic control (max. volume of cement per pedicle screw: 2 ml in the pedicle).

Contamination of the interlocking thread of the S⁴ Spinal System augmentation screw and/or insufficient fixing when removing the application system due to incorrect viscosity of the cement!
► Make certain the cement has the correct viscosity.

Risk of injury and/or malfunction!
► Screw the injection cannula into the thread of the pedicle screw, hand-tight.
► Screw the injection cannula into the pedicle screw and leave in place until the cement has completely set.

Risk of breaking the injection cannula by inducing excessive shear forces!
► Do not subject the injection cannula to transverse loads.
► Do not bend the injection cannula.

Injection cannula SR146SU for open approach

Risk of injury when using the product beyond the field of view!
► Apply the product only under visual control.

- Screw in the pedicle screw in accordance with the instructions. In doing so, keep the K-wire in the pedicle screw.
- Remove the injection cannula SR146SU 1 from the sterile packaging and rinse with sterile saline solution.
- Screw the injection cannula SR146SU 1 over the K-wire, hand-tight, in a line with the bone thread of the pedicle screw, into the screw head.
- Remove K-wire.
- Check that the injection cannula SR146SU 1 is positioned properly.
- Connect the cement applicator onto the the Luer lock connector.
- Apply the cement.
- Remove the cement applicator from the the Luer lock connector.
- Unscrew the injection cannula SR146SU 1, once the cement has completely set, out of the pedicle screw.
- Dispose of the injection cannula SR146SU 1 and do not reuse or resterilize it.
- Continue the operation and insert the rod.

Injection cannula SR148SU for percutaneous approach

Risk of injury and/or malfunction!
► Insert pestle 2b exclusively into the injection cannula SR148SU 2a.

Risk of injury and/or malfunction!
► Insert injection cannula SR148SU 2a over the K-wire or pestle 2b.

- Screw in the pedicle screw in accordance with the instructions. In doing so, keep the K-wire in the pedicle screw.
- Remove the injection cannula SR148SU 2a from the sterile packaging and rinse with sterile saline solution.
- Screw the injection cannula SR148SU 2a over the K-wire, hand-tight, in a line with the bone thread of the pedicle screw, into the screw head.
- or –
- Insert pestle into the injection cannula. Insert the cannula, over the pestle, into the screw head.
- In order to check that the injection cannula SR148SU 2a is correctly positioned, insert pestle 2b all the way into the injection cannula SR148SU 2a. If the position is correct, the marking on the pestle 2b disappears completely in the injection cannula SR148SU 2a.
- If the marking remains visible, check the alignment of the injection cannula SR148SU 2a. If necessary, screw in the injection cannula SR148SU 2a again.
- Remove the pestle 2b from the injection cannula SR148SU 2a.
- Connect the cement applicator onto the the Luer lock connector.
- Apply the cement.
- Remove the cement applicator from the the Luer lock connector. If necessary, press cement into the screw using the pestle 2b.
- Screw the injection cannula SR148SU 12a, once the cement has completely set, out of the pedicle screw.
- Dispose of the injection cannula SR148SU 2a and do not reuse or resterilize it.
- Continue the operation and insert the rod.

Validated reprocessing procedure

Single-use products

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product!

Designation	Art. no.
S ⁴ injection cannula for augmentation screws	SR146SU
S ⁴ percutaneous injection cannula for augmentation screws	SR148SU

Storage

- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Technical Service

Risk of injury and/or malfunction!

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems LLC

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042

Aesculap Repair Hotline

Phone: +1 (800) 214-3392

Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!
Incineration of the product, its components and its packaging produces harmless residues.
- Comply with national regulations for incineration under all circumstances!

Distributor in the US/Contact in Canada for product information and complaints

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