## Aesculap Spine

**(GB)** Instructions for use/Technical description

**(USA)** S<sup>4</sup> Instruments

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













26

### Legend

- 1 Set screw starter FW177R
- 2 Counter torque FW178R (with L-handle)
- 3 Rod persuader FW208R
- 4 Reduction lever (consisting of inner pipe FW229R, screwdriver FW228R, outer sleeve FW235R and spacer FW141P or FW143P)
- 5 Percutaneous reduction lever (consisting of lever inner pipe FW229R, screwdriver FW228R and outer sleeve with long slot FW353R and spacer FW141R or FW143P)
- 6 Counter torque FW236R
- 7 Wrench FW237R
- 8 Distractor (consisting of distractor FW238R and two distraction arms FW239R)
- 9 Rod insertion instrument FW240R
- 10 Distraction spindle FW241R (consisting of distractor and extension)
- 11 Rod length measuring instrument FW242R
- 12 Screw length measuring instrument FW351R
- 13 Slotted hammer FW243R
- 14 Tissue protection sleeve FW244P
- 15 Tissue protection sleeve FW355P
- 16 Kirschner wire FW247S
- 17 Trocar FW271M
- 18 Aiming device FW258M
- 19 Aiming device extraction tool FW274M
- 20 K-wire protection sleeve FW352R
- 21 Monoaxial screwdriver FW262R
- 22 Pedicle dissector FW263R
- 23 Tapper FW264R, FW265R, FW266R, FW267R and FW268R
- 24 Polyaxial screw driver FW270R
- 25 Dilation tube FW272R
- 26 Dilation tube FW354R

### Symbols on product and packages

Not for reuse in intended applications as defined by the manufacturer

Caution, general warning symbol  
Caution, see documentation supplied with the product

### Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bb.raun.com>

### Intended use

The S<sup>4</sup> instruments are used for the implantation of the S<sup>4</sup> implants with the following systems:

- S<sup>4</sup> Spinal System
- S<sup>4</sup> CS Cannulated System
- S<sup>4</sup> FRI/S<sup>4</sup> FRI percutaneous

The S<sup>4</sup> fracture reduction instruments are used for the reduction and stabilization of fractures with the cannulated S<sup>4</sup> implants of the S<sup>4</sup> spinal system.

## Safe handling and preparation

### CAUTION

Federal law restricts this device to sale by, or on order of a physician!

**Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!**

- ▶ **Do not re-use the product after use.**

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.
- ▶ To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).
- ▶ Read, observe and keep the operating manuals S<sup>4</sup> Spinal System Instrumentation, S<sup>4</sup> FRI and S<sup>4</sup> SRI.

## Safe operation

**Danger of death of the patient caused by pushing the guide wire into the aorta!**

- ▶ When inserting the S<sup>4</sup> implants or S<sup>4</sup> instruments with the guide wire, always observe the markings on the guide wire.
- ▶ Hold the guide wire with a forceps.
- ▶ Remove the guide wire at the right moment.

**Excessive trocar insertion depth due to missing end stop!**

- ▶ Use trocar FW271M only in combination with aiming device FW258M.

**Risk of injury and/or malfunction!**

- ▶ Always carry out a function check prior to using the product.

## Disassembly

- Dismantle the set screw starter with an SW6 open-ended spanner, see Fig. 1.


Fig. 1

### Note

The set screw starter does not need to be dismantled for cleaning.

- Disassembling reduction lever 4, see Fig. 2. Fully disassemble reduction lever and unscrew lock nut.


Fig. 2

- Disassembling percutaneous reduction lever 5, see Fig. 3. Fully disassemble reduction lever and unscrew lock nut.


Fig. 3

- Dismantle distractor 8 with the distraction blades, see Fig. 4: Activate slider B and remove distraction blade FW239R.


Fig. 4

### Legend

B Slider

- Dismantle rod insertion instrument 9, see Fig. 5: Unscrew the inner part.


Fig. 5

- Dismount distraction spindle 10, see Fig. 6.


Fig. 6

- Dismount rod length measuring instrument 11, see Fig. 7.


Fig. 7

## Assembly

### Assembling the set screw starter

- Screw together the set screw and the starter so that they are hand-tight, see Fig. 8.

Fig. 8

### Note

*The set screw starter is not intended for loosening set screws. If the working end becomes loose, it can be tightened until hand-tight, without tools. Use the revision screwdriver FW193R to loosen the set screws.*

### Mounting the reduction lever

- Mount reduction lever 4, see Fig. 9.


Fig. 9

## Mounting the reduction lever percutaneous

- Screw lock nut C onto the inner sleeve, see Fig. 10.


Fig. 10

### Legend

C Lock nut

- Insert screwdriver into inner sleeve, secure with spacer and apply set screw D, see Fig. 11.


Fig. 11

- Insert the subassembly into the FRI outer sleeve, see Fig. 12.


Fig. 12

## Mounting other product components

- Mount distractor 8 with the distraction blades, see Fig. 13.


Fig. 13

- Mount rod insertion instrument 9, see Fig. 14.


Fig. 14

- Mount distraction spindle 10, see Fig. 15.


Fig. 15

- Mount rod length measuring instrument 11, see Fig. 16.


Fig. 16

## Validated reprocessing procedure

### General safety instructions

#### Note

*Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.*

#### Note

*For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.*

#### Note

*Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.*

#### Note

*Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.*

*The recommended chemistry was used for validation.*

#### Note

*If there is no final sterilization, then a virucidal disinfectant must be used.*

#### Note

*For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>*

*The validated steam sterilization procedure was carried out in the Aesculap sterile container system.*

## Single-use products

K-wire FW247S

**Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!**

- **Do not re-use the product after use.**

- Do not reprocess single-use products that have been directly or indirectly contaminated with blood.
- Do not reprocess or re-use any used, damaged or contaminated single-use products.

## General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to Publications, Red Brochure – Proper maintenance of instruments.

## Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

## Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

## Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembly.

## Cleaning/disinfection

### Product-specific safety notes on the reprocessing procedure

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
    - are approved for use, for example, on aluminum, plastic materials, and high-grade steel,
    - do not attack softeners (e.g. in silicone).
  - ▶ Observe specifications regarding concentration, temperature, and exposure time.
- ▶ To prevent increased contamination of loaded instrument trays during use, please ensure that contaminated instruments are collected separately and not returned to the instrument tray.
  - ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

- ▶ Carry out ultrasound cleaning:
  - as an effective mechanical supplement to manual cleaning/disinfecting.
  - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
  - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
  - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- ▶ Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

# Aesculap®

## S4 Instruments

### Validated cleaning and disinfection procedure

| Validated procedure                                                                                                                                                                                                                                        | Specific requirements                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | Reference                                                                                                                                                                           |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Manual cleaning with immersion disinfection only:</p> <p>FW141P, FW143P, FW178R, FW228R, FW229R, FW235R, FW236R–FW239R, FW241R–FW243R, FW244P, FW247S, FW258M, FW258R, FW262R–FW268R, FW270R, FW271M, FW271R, FW272R, FW274M, FW351R–FW354R, FW355P</p> | <ul style="list-style-type: none"> <li>■ Cleaning brush: TA011994/TE654202, GK469200</li> <li>■ 20 ml disposable syringe</li> <li>■ Keep working ends open for cleaning.</li> <li>■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> <li>■ Check visible surfaces for residues after manual cleaning/disinfecting.</li> <li>■ Repeat the cleaning process if necessary.</li> </ul> | <p>Chapter Manual cleaning/disinfection and sub-chapter:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with immersion disinfection</li> </ul>                |
| <p>Manual cleaning with ultrasound and immersion disinfection only:</p> <p>FW177R, FW208R, FW240R</p>                                                                                                                                                      | <ul style="list-style-type: none"> <li>■ Cleaning brush: TA011994/TE654202, GK469200</li> <li>■ 20 ml disposable syringe</li> <li>■ Keep working ends open for cleaning.</li> <li>■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>                                                                                                                                         | <p>Chapter Manual cleaning/disinfection and sub-chapter:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual cleaning with ultrasound and immersion disinfection</li> </ul> |

| Validated procedure                                                                                                                                                                                                                                                                      | Specific requirements                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | Reference                                                                                                                                                                                                                                                       |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Mechanical alkaline cleaning and thermal disinfection only:<br/>FW141P, FW236R–FW238R, FW243R, FW247S, FW274M</p>                                                                                                                                                                     | <ul style="list-style-type: none"> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> <li>■ Keep working ends open for cleaning.</li> <li>■ Place products in the tray with their hinges open.</li> </ul>                                                                                                                    | <p>Chapter Mechanical cleaning/disinfecting and sub-chapter:</p> <ul style="list-style-type: none"> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>                                                                              |
| <p>Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection only:<br/>FW143P, FW177R, FW178R, FW208R, FW228R, FW229R, FW235R, FW238R, FW239R, FW240R–FW242R, FW258M, FW258R, FW262R–FW268R, FW270R, FW271R, FW271M, FW272R, FW351R–FW354R</p> | <ul style="list-style-type: none"> <li>■ 20 ml disposable syringe</li> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ FW178R, FW228R, FW235R, FW258M, FW262R–FW268R, FW270R, FW352R, FW272R, FW354R: Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> <li>■ Keep working ends open for cleaning.</li> <li>■ Place products in the tray with their hinges open.</li> </ul> | <p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with a brush</li> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul> |

# Aesculap®

## S4 Instruments

| Validated procedure                                                                                                                                 | Specific requirements                                                                                                                                                                                                                                                                                                                                                                                                                                       | Reference                                                                                                                                                                                                                                                                       |
|-----------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual pre-cleaning with brush and subsequent mechanical neutral or mild alkaline cleaning and thermal disinfection only:<br>FW355P, FW244P         | <ul style="list-style-type: none"> <li>■ 20 ml disposable syringe</li> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> <li>■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> <li>■ Keep working ends open for cleaning.</li> <li>■ Place the product on the tray with all product links and joints open.</li> </ul>                            | <p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with a brush</li> <li>■ Chapter Mechanical neutral or mild alkaline cleaning and thermal disinfecting</li> </ul> |
| Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection only:<br>FW208R, FW240R, FW242R | <ul style="list-style-type: none"> <li>■ 20 ml disposable syringe</li> <li>■ Mount jaws protection on the product.</li> <li>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).</li> <li>■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li> <li>■ Keep working ends open for cleaning.</li> <li>■ Place products in the tray with their hinges open.</li> </ul> | <p>Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:</p> <ul style="list-style-type: none"> <li>■ Chapter Manual pre-cleaning with ultrasound and brush</li> <li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li> </ul>     |

### Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

## Manual cleaning with immersion disinfection

| Phase | Step                     | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                             |
|-------|--------------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------------|
| I     | Disinfecting<br>cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| II    | Intermediate rinse       | RT (cold)    | 1          | -            | D-W              | -                                                                    |
| III   | Disinfection             | RT (cold)    | 15         | 2            | D-W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| IV    | Final rinse              | RT (cold)    | 1          | -            | FD-W             | -                                                                    |
| V     | Drying                   | RT           | -          | -            | -                | -                                                                    |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

### Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

### Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

# Aesculap®

## S4 Instruments

### Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

### Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

### Manual cleaning with ultrasound and immersion disinfection

| Phase | Step                | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|---------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Ultrasonic cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse  | RT (cold)    | 1          | -            | D-W              | -                                                              |
| III   | Disinfection        | RT (cold)    | 15         | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse         | RT (cold)    | 1          | -            | FD-W             | -                                                              |
| V     | Drying              | RT           | -          | -            | -                | -                                                              |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

#### **Phase I**

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### **Phase II**

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### **Phase III**

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### **Phase IV**

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

#### **Phase V**

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

## Mechanical cleaning/disinfecting

### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

## Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical/Note                                                                                                                                                                                                                                       |
|-------|----------------------|--------------|------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse             | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                                   |
| II    | Cleaning             | 55/131       | 10         | FD-W             | <div> <div>■</div> Concentrate, alkaline: <ul style="list-style-type: none"> <li>- pH = 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> <div>■</div> 0.5 % working solution <ul style="list-style-type: none"> <li>- pH = 11*</li> </ul> </div> |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                                   |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                                   |
| V     | Drying               | -            | -          | -                | According to the program for cleaning and disinfection device                                                                                                                                                                                       |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

## Mechanical cleaning/disinfection with manual pre-cleaning

### Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

### Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

## Manual pre-cleaning with a brush

| Phase | Step                     | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|--------------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfectant<br>cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing                  | RT (cold)    | 1          | -            | D-W              | -                                                              |

D-W: Drinking water

RT: Room temperature

\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

### Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual pre-cleaning with ultrasound and brush

| Phase | Step                | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|---------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Ultrasonic cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing             | RT (cold)    | 1          | -            | D-W              | -                                                              |

D-W:

Drinking water

RT:

Room temperature

\*Recommended:

BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

### Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical                                                                                                                                                                                                                                |
|-------|----------------------|--------------|------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse             | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                       |
| II    | Cleaning             | 55/131       | 10         | FD-W             | <div><div>■</div>Concentrate, alkaline:<ul style="list-style-type: none"><li>- pH = 13</li><li>- &lt;5 % anionic surfactant</li></ul><div>■</div>0.5 % working solution<ul style="list-style-type: none"><li>- pH = 11*</li></ul></div> |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                       |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                       |
| V     | Drying               | -            | -          | -                | According to the program for cleaning and disinfection device                                                                                                                                                                           |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

### Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------|----------------------|--------------|------------|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse             | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| II    | Cleaning             | 55/131       | 10         | FD-W             | <b>Neutral:</b> <ul style="list-style-type: none"> <li>■ Concentrate: <ul style="list-style-type: none"> <li>- pH neutral</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 %* working solution</li> </ul> <b>Mildly alkaline:</b> <ul style="list-style-type: none"> <li>■ Concentrate: <ul style="list-style-type: none"> <li>- pH = 9.5</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % solution</li> </ul> |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| V     | Drying               | -            | -          | -                | According to the program for cleaning and disinfection device                                                                                                                                                                                                                                                                                                                                                                                    |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner neutral

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning/disinfecting process if necessary.

## Inspection, maintenance and checks

CAUTION

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ **Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).**

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Assemble dismountable products, see Assembly.
- ▶ Check for compatibility with associated products.

## Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Store products with ratchet locks fully opened or locked no further than in the first notch.
- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

## Steam sterilization

### Note

*FW240R, FW270R: The product may only be sterilized when dismantled.*

### Note

*FW177R, FW208R, FW238R, FW239R, FW240R, FW242R, FW262R–FW268R, FW272R: The product can be sterilized either in disassembled or in assembled condition.*

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
  - Disassemble the product
  - Steam sterilization through fractionated vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

### Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of  $10^{-6}$ , Aesculap recommends the following parameters:

| Aesculap Orga Tray/Sterile container (perforated bottom)<br>Minimum cycle parameters* |                   |       |                     |
|---------------------------------------------------------------------------------------|-------------------|-------|---------------------|
| Sterilization method                                                                  | Temp.             | Time  | Minimum drying time |
| Prevacuum                                                                             | 270 °F/27<br>5 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

#### WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

### Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

### Technical Service

**Risk of injury and/or malfunction!**

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

#### Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 (7461) 95-1602  
Fax: +49 (7461) 16-5621  
E-Mail: [ats@aesculap.de](mailto:ats@aesculap.de)  
Or in the US:  
Aesculap Implant Systems LLC  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042  
Aesculap Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

**Accessories/Spare parts**

| Art. no. | Designation                                                            |
|----------|------------------------------------------------------------------------|
| FW233201 | Replacement nut (for distractor FW238R and distraction spindle FW241R) |

---

**Distributor in the US/Contact in Canada for product information and complaints**

Aesculap Implant Systems LLC  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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- S<sup>4</sup> Spinal System
- S<sup>4</sup> CS - Cannulated System
- S<sup>4</sup> FRI/S<sup>4</sup> FRI perkutan

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany