Aesculap Spine

(GB) Instructions for use/Technical description

(USA) S⁴® MIS Spyder Instrumentation

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


















Symbols on product and packages

Caution, general warning symbol
Caution, see documentation supplied with the product

Applicable to Spyder Instrumentation Set

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

The S⁴ MIS Spyder Instrumentation set is intended for use during posterior fusion of the lumbar spine with the S⁴ System to retract the soft tissue and provide access to the operative region. The instrumentation can be used for either the minimal invasive or the percutaneous approach to the pedicle.

Available sizes

Name	Sizes		
Tissue Dilator	6 mm 10 mm 16.5 mm 20 mm 24 mm		
Spyder Retractor C/C Blade	16.5x40 mm 16.5x50 mm 16.5x60 mm 16.5x70 mm 16.5x80 mm 16.5x90 mm 16.5x100 mm 16.5x110 mm	20x40 mm 20x50 mm 20x60 mm 20x70 mm 20x80 mm 20x90 mm 20x100 mm 20x110 mm	24x40 mm 24x50 mm 24x60 mm 24x70 mm 24x80 mm 24x90 mm 24x100 mm 24x110 mm
Spyder Retractor M/L Blade Cvd	7x40 mm 7x50 mm 7x60 mm 7x70 mm 7x80 mm 7x90 mm 7x100 mm	20x40 mm 20x50 mm 20x60 mm 20x70 mm 20x80 mm 20x90 mm 20x100 mm	
Delta Blade	7x17x40 mm 7x19x50 mm 7x22x60 mm 7x24x70 mm 7x27x80 mm 7x29x90 mm 7x31x100 mm 7x34x110 mm		

Aesculap®

S⁴® MIS Spyder Instrumentation

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.

Safe operation

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

Dilation

- ▶ Dilate the paraspinal muscles starting with the smallest diameter (6 mm) Tissue Dilator (ME389T) and work up to the diameter of the cranial-caudal (C/C) retractor blades you intend to use (16.5 mm, 20 mm, or 24 mm diameters).
- ▶ The dilator handle may be attached to each dilator if desired.

Lower Frame and Cranial Caudal (C/C) Blade Insertion

- ▶ After inserting the final dilation diameter, select the length of the cranial-caudal (C/C) retractor blades based on the depth indicator on the dilator.
- ▶ Insert the C/C blades into the assembled "lower" frame (ME302R, ME303R, ME304R) by loosening the grub screw, engaging the blades into the frame and then tightening the grub screw.
- ▶ Close the retractor, creating a tube from the blades.
- ▶ Slide blades over the final dilator, and remove the tissue dilator.
- ▶ Expand the lower arm of the retractor by rotating the thumb ratchets of the retractor.
- ▶ The blades of the lower retractor may be flared using the retractor rotator arms (ME380R). Insert the rotator arms into the retractor and rotate the arms inwards to flare blades.

Mini-open frame disassembles or jams due to unsecured blades!

- ▶ When placing blades into the frame arms always finger-tighten the grub screw to fix each blade in place to ensure blades do not pop out during use.
- ▶ Do not force blade into the frame arm as this will cause it to jam and make it difficult to remove.

Damage to the surrounding muscle and soft tissues by retractor blades!

- ▶ Do not over retract blades.

Upper Frame and Medial Lateral (M/L) Blade Insertion

- ▶ The M/L blades of the upper retractor may be inserted in a similar fashion using the upper frame of the retractor. The upper retractor frame has gold anodization on the collar for easy identification.
- ▶ Place the upper retractor arm over the lower retractor arm and expand and flare the upper retractor in a similar fashion if desired.

Frames may not nest correctly if not assembled in right order!

- ▶ For the mini-open procedure ensure that the lower frame (silver collar) is inserted first.
- ▶ If a second frame is to be used, insert the upper frame (gold collar) over and perpendicular to the lower frame, ensuring that the frames nest together correctly.

Blade Replacement

- ▶ Blade lengths and shapes may be changed in situ through the use of the blade holding instrument (ME381R).
- ▶ Before removing any blades, release the force of the soft tissue on the retractor blades by returning them to their parallel position.
- ▶ Insert the rotator arm into the retractor, pull back on the collar and rotate the arm to the upright position.
- ▶ Re-inserting the final dilator will hold the soft tissue while the retractor blade is exchanged.

- ▶ To remove a blade, insert the retractor blade handle into the blade notches and rotate the knob counter-clockwise. Loosen the grub screw and remove the blade.
- ▶ To insert a new blade, place the new blade in the retractor arm, tighten the grub screw and rotate the blade removal knob clockwise to disengage the tool.

Inability to remove retractor blades due to pressure from retracted tissues!

- ▶ Always set the retractor blades back to their vertical position prior to removing them.
- ▶ Pressure exerted by the tissues on the blades can make retraction of the collar difficult. Use the Blade Arm Rotator to relax the pressure on the teeth under the collar while gently retracting the collar.

Damage to Blade Holder Instrument or C/C blades due to over tightening knob!

- ▶ Over tightening the knob on the Blade Holder Instrument will put unnecessary load on the threads and jaws of the instrument and may cause expansion of the caudal cranial blades that makes removal or insertion of the blade in the Spyder frame arm impossible.

Aesculap®

S⁴® MIS Spyder Instrumentation

Lighting

- For additional visibility attach a light source to either the C/C Blade Light Rod Clip (ME387R) or the M/L Blade Light Rod Clip (ME388R).

For additional visibility, the medial blade and M/L light clip will provide the best visibility.

WARNING

Damage to the light clip due to overheating!

- Using a light source greater than 300 watts will lead to excessive heating of the light cable and light clip.

Securing the Lower Frame

If desired, the lower arms of the retractor can be attached to the table arm.

- Attach the table arm connection block (ME382R) to the retractor arm (ME304R).
- Attach table arm connection block (ME382R) to the table arm (FF280R) by sliding the (ME382R) into the collar and tighten the knob.

The Spyder Retractor is now attached to the table arm assembly.

Disassembling

Frame Arm (ME300R–ME303R):

- Remove any blades attached to the frame arm by loosening the grub screw.
- Open ratchet lock and slide the Frame Arm off the Mini-Open Frame Rack (ME304R).
- Remove the thumb ratchet from the Frame Arm.
- Unthread grub screw so that it remains in the holding arm but is completely disengaged from its threading.
- Retract and rotate the collar 180° so that it is held by the positive stop in the retracted position for cleaning.

Blade Handle (ME381R):

- Hold the blade handle head upright.
- Unscrew the nut on the shaft proximally towards the handle.
Nut disengages the threading on the shaft.
- Turn the knob on the bottom counter clockwise until the head/shaft detaches from the handle.
- Remove the head/shaft from the handle.
- Remove the knob from the handle.

Table Arm Pivot Connection Block (ME381R):

- Unscrew nut from the threading on the shaft away from the connector.
- Unscrew wheel from connection block.

Assembling

Frame Arm (ME300R–ME303R):

- Screw grub screws until they fully engage their threads.
- Insert thumb screws into the Frame Arm block.
- Engage Frame Arm onto the Mini-Open Frame Rack (ME304R):
 - The Left and Right Upper Frame arms (gold collar) should be paired with one another.
 - The Right and Left Lower Frame arms (silver collar) should be paired with one another.
 - The arm of the Frame Arm should face away from the teeth on the rack.
 - The thumb screws should be on the distant ends of the rack,
 - The frame arms are facing inward.
- Retract the collar and rotate the arms until the concave frame faces inward, and the grub screws point out and upward.

Blade Handle (ME381R):

- ▶ Insert the knob into the slot on the bottom of the handle.
- ▶ Insert the head/shaft until it engages the threads inside the knob.
- ▶ Hold the blade handle head upright.
- ▶ Turn the knob clockwise until the head/shaft fully engages with the handle.
- ▶ Orient the head so that the pegs on the head will mesh with the slots on the handle.
- ▶ Ensure that the head is fully engaged with the handle.
- ▶ Thread the nut on the handle distally towards the head until it is fully threaded.
- ▶ Ensure the Blade Handle operates correctly by turning the knob counter clockwise to spread open the arms and clockwise to close the arms.

Table Arm Pivot Connection Block (ME381R):

- ▶ Thread the wheel onto the block until it is fully engaged.
- ▶ Thread the nut on the shaft towards the connector until it is fully engaged.

Validated reprocessing procedure

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt–Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bb Braun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary.

Aesculap®

S⁴® MIS Spyder Instrumentation

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open up products with hinges.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which are approved for (e.g. titanium, high-grade steel).
- Observe specifications regarding concentration, temperature and exposure time.
- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ME389T–ME393T ME394T ME304R ME380R ME319T–ME326T ME327T–ME334T ME335T–ME341T ME350T–ME356T ME357T–ME364T ME387R ME388R	■ Cleaning brush: e.g. TA011944 ■ Disposable syringe 20 ml ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ME300R ME301R ME302R ME303R ME381R ME382R	■ Cleaning brush: TA011944, TE654202, or GK469200 ■ Disposable syringe 20 ml ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ME304R ME380R ME319T–ME326T ME387R ME388R	■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfecting and sub-chapter: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Aesculap®

S⁴® MIS Spyder Instrumentation

Validated procedure	Specific requirements	Reference
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ME389T–ME394T** ME327T–ME334T ME335T–ME341T ME350T–ME356T ME357T–ME364T	■ Cleaning brush: e.g. TA011944, TE654202, or GK469200 ■ Disposable syringe 20 ml ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ **Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ME300R ME301R ME302R ME303R ME381R ME382R	■ Cleaning brush: e.g. TA011944, TE654202, or GK469200 ■ Disposable syringe 20 ml ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

Manual cleaning and wipe disinfecting

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	-	-	D-W	-
II	Drying	RT	-	-	-	-
III	Wipe disinfection	-	>1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol
IV	Final rinse	RT (cold)	0.5	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

Phase I

- ▶ Clean the product under running faucet water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.

Phase II

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Phase III

- ▶ Wipe all surfaces of the product with a single-use disinfectant wipe.

Phase IV

- ▶ After the specified exposure time (at least 1 min), rinse the disinfected surfaces under running FD water.
- ▶ Drain any remaining water fully.

Aesculap®

S⁴® MIS Spyder Instrumentation

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Aesculap®

S⁴® MIS Spyder Instrumentation

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfecting

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Aesculap®

S⁴® MIS Spyder Instrumentation

Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Neutral: ■ Concentrate: - pH neutral - <5 % anionic surfactant ■ 0.5 %* working solution Mildly alkaline: ■ Concentrate: - pH = 9.5 - <5 % anionic surfactant ■ 0.5 % solution
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner neutral

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning/disinfecting process if necessary.

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Aesculap®

S⁴® MIS Spyder Instrumentation

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W:

Drinking water

RT:

Room temperature

*Recommended:

BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Aesculap®

S⁴® MIS Spyder Instrumentation

Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Neutral: ■ Concentrate: - pH neutral - <5 % anionic surfactant ■ 0.5 %* working solution Mildly alkaline: ■ Concentrate: - pH = 9.5 - <5 % anionic surfactant ■ 0.5 % solution
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner neutral

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning/disinfecting process if necessary.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- **Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).**

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note

The product can be sterilized either in disassembled or in assembled condition.

Note

To avoid breakage due to stress crack corrosion, sterilize the product with the lock fully open or locked no further than on the first ratchet tooth.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process:
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Aesculap®

S⁴® MIS Spyder Instrumentation

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
 - Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.
- To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom)			
Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/ 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Technical Service

Risk of injury and/or malfunction!

► Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
 Am Aesculap-Platz
 78532 Tuttlingen / Germany
 Phone: +49 (7461) 95-1602
 Fax: +49 (7461) 16-5621
 E-Mail: ats@aesculap.de
 Or in the US:
 Aesculap Implant Systems LLC
 Attn. Aesculap Technical Services
 615 Lambert Pointe Drive
 Hazelwood
 MO, 63042
 Aesculap Repair Hotline
 Phone: +1 (800) 214-3392
 Fax: +1 (314) 895-4420
 Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany