Aesculap Spine

(GB)

Instructions for use

(USA)

S⁴ percutaneous instruments

(D)

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S⁴ perkutan-Instrumente

(F)

Mode d'emploi

Instruments percutanés S⁴

(E)

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Instrumentos percutáneos S⁴

(I)

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S⁴ strumenti percutanei

(P)

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Instrumentos percutâneos S⁴

(NL)

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S⁴ percutaan-instrumenten

(RUS)

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
S⁴ percutaneous instruments

Legend

- 1 Cleaning instrument FW692R
- 2 Clamping sleeve FW693R
- 3 Inner sleeve FW694R
- 4 Screwdriver FW695R
(for polyaxial S⁴ screws)
- 5 Screwdriver FW696R
(for monoaxial S⁴ screws)
- 6 Axial spacer FW143P
- 7 Insertion instrument for locking screws FW697R

Symbols on product and packages

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture

Intended use

S⁴ percutaneous instruments are used for the implantation of cannulated S⁴ implants of the S⁴ Spinal System in percutaneous operating technique.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only for its intended purpose, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either manually or mechanically, prior to its initial sterilization.

- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- Read, observe and keep the operating manuals for S⁴ Spinal System Instrumentation, S⁴ FRI and, if applicable, S⁴ SRI.

Safe operation

DANGER

Danger of death of the patient caused by pushing the guide wire into the aorta!

- When inserting the S⁴ implants or S⁴ instruments with the guide wire, always observe the markings on the guide wire.
- Hold the guide wire with a forceps.
- Remove the guide wire at the right moment.

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

Cleaning instrument FW692R

The cleaning instrument is used for cleaning/clearing the cannulation, if a cannulated instrument becomes blocked intra-operatively with soft tissue or bone chips.

Risk of injury if cannulated instruments are punctured with the cleaning instrument!

- Remove the obstruction carefully with mild pressure.
- Always insert the cleaning instrument in the direction of the working end.
- Do not obstruct the outlet with your hand.

Damage to cannulated instruments if the cleaning instrument is incorrectly applied!

- Always insert the cleaning instrument in the direction of the working end.
- Check the cleaning instrument for bending, kinks and other damage before use.

- Disassemble the instrument as far as possible on a separate sterile table, if applicable, see Disassembling.
- Carefully insert the cleaning instrument into the cannulation from the handle side in the direction of the working end.
- Remove the obstruction by applying mild pressure with the handle part. Ensure that the cannulation does not become occluded.
- Assemble the instrument, if applicable, see Assembling.

Clamping sleeve FW693R and screwdriver FW695R/FW696R

Screwdrivers facilitate the percutaneous application of pedicle screws from the S⁴ Spinal System. Depending on the type of screw, the corresponding screwdriver for monoaxial or polyaxial screws must be chosen.

- Assemble the clamping sleeve, see Assembling.

Inserting the pedicle screws in the clamping sleeve

- Turn fixing screw **C** in the direction of the arrow until the thread in the direction of the working end becomes visible, see Fig. 1.


Fig. 1

- Pull back outer sleeve **A** as far as possible in the direction of fixing screw **C**, see Fig. 2.


Fig. 2

- Insert screwdriver **F** into clamping sleeve **2** from the rear end. Ensure that the lateral lugs slide into the grooves right to the bottom, see Fig. 3.


Fig. 3

S⁴ percutaneous instruments

- Insert the working ends of the pedicle screws into tension sleeve B, see Fig. 4.

When using polyaxial pedicle screws: Ensure that the hexagon of the screwdriver is fully inserted in the hexagon of the polyaxial pedicle screw.

When using monoaxial pedicle screws: Ensure that the dihedron sits securely in the screw head.


Fig. 4

- Slide outer sleeve A in the direction of the pedicle screw and with fixing screw C clamp pedicle screw in clamping sleeve 2, see Fig. 5.
- Check to ensure that the implant is properly seated in the clamping sleeve 2.


Fig. 5

- For implanting the pedicle screw and inserting the rod, see operating technique or operating manual for the S⁴ Spinal System, S⁴ FRI and S⁴ SRI.

Removing the clamping sleeve from the applied pedicle screw

- Turn fixing screw C in the direction of the arrow until the thread in the direction of the working end becomes visible, see Fig. 6.


Fig. 6

- Pull back outer sleeve A in the direction of fixing screw C.
- Completely remove clamping sleeve 2, see Fig. 7.


Fig. 7

Disassembling

Clamping sleeve FW693R

- Remove outer sleeve A over collet B in the direction of the working end.
- Unscrew fixing screw C from the rear end and remove from collet B, see Fig. 8.


Fig. 8

Inner sleeve FW694R

- Unscrew nut D from inner sleeve E, see Fig. 9.

Fig. 9

Assembling

Clamping sleeve FW693R

- Screw fixing screw **C** onto collet **B** until the thread of the collet is no longer visible.
- Push on outer sleeve **A** from the front end over collet **B**, see Fig. 10. Ensure that the slots are aligned with the slot nuts. The grooved end of the clamping sleeve should point in the direction of fixing screw **C** for improved usability.


Fig. 10

Inner sleeve FW694R

- Screw nut **D** onto inner sleeve **E**, see Fig. 11.

Fig. 11

Validated reprocessing procedure

Note

National laws, national and international standards and directives, and product-specific hygiene regulations for reprocessing must be observed.

Note

For patients with Creutzfeldt-Jakob Disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/sterile processing technician is responsible for this.

Note

Up-to-date information on reprocessing can be found on the Aesculap Extranet at www.aesculap-extra.net

S⁴ percutaneous instruments

General information

To prevent increased contamination of loaded instrument trays during use, please ensure that contaminated instruments are collected separately and not returned to the instrument tray.

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion of stainless steel. Consequently, the time interval between application and processing should not exceed 6 hours, pre-cleaning fixating temperatures >45 °C should not be used, and no fixating disinfectants (containing aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chloride – e.g. in surgical residues, drugs, saline solutions and water for cleaning, disinfection and sterilization – may cause corrosion damage to stainless steel (pinholing, stress corrosion), thus rendering the products unusable. These must be removed by rinsing thoroughly with demineralized water and then drying.

Only use process chemicals that have been tested and approved (with VAH, DGHM or FDA approval, for instance, or feature a CE mark) and are recommended by the chemical manufacturer as being compatible with the materials. All the chemical manufacturer's application specifications regarding temperature, concentration and contact time should be strictly observed. Failure to do so can result in the following problems:

- Visible changes to materials e.g. fading or changes in the color of titanium or aluminum. As regards aluminum, visible changes to the surface may already occur at a pH level of >8 in the application/working solution.
- Material damage such as corrosion, cracks, fractures, premature deterioration or swelling.
- Do not use process chemicals that cause stress cracks or brittleness in plastics.
- Clean the product immediately after use.

Please see www.a-k-i.org for more detailed information on hygienically safe reprocessing which is protective of materials and retains their value.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

Preparation at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Rinse non-visible surfaces such as those in instruments with concealed crevices, lumens or complex geometries, preferably with distilled water, using a disposable syringe for instance.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/disinfection

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for aluminum, plastics, and high-grade steel, for example, according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfection.

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning process if necessary.

S⁴ percutaneous instruments

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	>15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QUAT-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QUAT-free
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

Phase I

- Fully immerse the product in the cleaning/disinfecting solution for at least 15 minutes. Ensure that all accessible surfaces are moistened.
- Clean the instrument under running tap water with a suitable cleaning brush where necessary for as long as it takes to remove all discernable residues.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush (TA007747 or TE654202).
- Maneuver non-rigid components, such as adjustable screws, hinges etc during cleaning.

- Then flush these areas thoroughly at least five times with the cleaning solution using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Phase II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Maneuver non-rigid components, such as adjustable screws, hinges, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the instrument in the disinfectant solution.
- Maneuver non-rigid components, such as adjustable screws, hinges, etc. during disinfection.
- Rinse the lumens at least five times at the beginning of the contact period using a disposable syringe (20 ml) and a suitable irrigation adapter. Ensure that all accessible surfaces are moistened.
- Drain any remaining water fully.

Phase IV

- Rinse/flush the instrument thoroughly (all accessible surfaces).
- Maneuver non-rigid components, such as adjustable screws, hinges, etc. during final rinse.
- Rinse the lumens at least five times using a disposable syringe (20 ml) and a suitable irrigation adapter.
- Drain any remaining water fully.

Phase V

- Dry the instrument with a lint-free cloth or medical compressed air.

Mechanical cleaning/disinfecting

Cleaning instrument FW692R and insertion instrument for locking screws FW697R only

Note

The disinfectant must be of proven effectiveness (e.g. DGHM or FDA approval or CE mark).

Note

For thermal disinfection, always use fully desalinated water (demineralized, low microbial contamination: drinking water quality at least) and ensure that Ao is >3 000 for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

- Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots).
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Keep working ends open for cleaning.
- Place instruments in the tray with their hinges open.

S⁴ percutaneous instruments

Mechanical alkaline cleaning and thermal disinfection

Machine type: Single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER alkaline with tensides, application solution 0.5%
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfectant program

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

Mechanical cleaning/disinfection with manual pre-cleaning

Clamping sleeve FW693R, inner sleeve FW694R, screwdriver FW695R and FW696R and axial spacer FW143P

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note

For thermal disinfection, always use fully desalinated water (demineralized, low microbial contamination: drinking water quality at least) and ensure that Ao is >3 000 for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant Cleaning	RT (cold)	>15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QUAT-free
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

Phase I

- Fully immerse the product in the cleaning/disinfecting solution for at least 15 minutes. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all discernable residues have been removed.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush (TA007747 or TE654202).
- Maneuver non-rigid components, such as adjustable screws, hinges etc during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Phase II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Maneuver non-rigid components, such as adjustable screws, hinges, etc. during rinsing.

S⁴ percutaneous instruments

Mechanical alkaline cleaning and thermal disinfection

Machine type: Single-chamber cleaning/disinfection device without ultrasound

- Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots).

- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Keep working ends open for cleaning.
- Place instruments in the tray with their hinges open.

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER alkaline with tensides; application solution 0.5%
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfectant program

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).
- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable instruments, see Assembly.
- Check for compatibility with associated products.

Packaging

- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization

Note

The product may only be sterilized when dismantled.

- Ensure that the sterilant can reach all external and internal surfaces (e.g. by opening valves and stopcocks).
- Validated sterilization process
 - Disassemble the instrument.
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer specified in DIN EN 285 and validated in accordance with DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

S⁴ percutaneous instruments

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom)			
Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F— 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Aesculap Implant Systems LLC
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

**Distributor in the US/Contact in Canada
for product information and complaints**

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG