- GB

Instructions for use/Technical description
 S^{4®} Spinal System – supplied in sterile condition
- USA

Instructions for use/Technical description
 S^{4®} Spinal System – supplied in sterile condition
 Note for U.S. users
 This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapImplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
- D

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Aesculap® S⁴® Spinal System – supplied in sterile condition

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

Note

Generally, the S⁴ Spinal System – supplied in sterile condition is covered by the instructions for use S⁴ Spinal System – Lumbar/Deformity TA011187. The present information on the S⁴ sterile packed implants is intended to complement the respective information contained in the instructions for use S⁴ Spinal System – Lumbar/Deformity.

The S⁴ Spinal System implants are used for dorsal monosegmental and multisegmental stabilization of the lumbar and thoracic spine.

The parallel, domino and axial rod connectors are anchored to S⁴ Spinal System rods to join one rod in parallel to or in line with another rod.

The lateral offset rod connectors are anchored to the S⁴ Spinal System rods to offset the placement of a screw. The rod connectors therefore extend the rod to adjacent spinal segments.

The S⁴ Spinal System – supplied in sterile condition comprises:

- Rod connectors – parallel, domino, axial and lateral offset rod connectors

Note

There are special S⁴ instruments which must be used for implanting these components, as well as for the distraction, compression and reduction of the lumbar and thoracic spine.

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®: Titanium forged alloy Ti6Al4V acc. to ISO 5832-3

The titanium implants are anodized with a colored oxide layer. Slight changes in coloration may occur, but do not affect the implant quality.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Surgical implants serve to support normal healing processes. They are not intended to either replace normal structures of the human body or permanently bear the loads they would be subjected to in cases of incomplete healing.

Use for:

- Fractures
- Pseudarthrosis or delayed healing
- Dislocation
- Degenerative instability
- Post-discectomy syndrome
- Spondylolisthesis
- Kyphosis
- Scoliosis
- Stenosis
- Post-traumatic instabilities
- Tumors

Risk of screw breakage caused by pedicle screws used for spondylolisthesis!
► Support the stabilization measure with an intercorporeal fusion (ALIF, PLIF or TLIF).

Contraindications

Do not use in cases of:

- Fever
- Acute or chronic vertebral infections in the spine, local or systemic
- Pregnancy
- Severe osteoporosis or osteopenia
- Infection
 - Systemic
 - In the spine
 - Local
- Medical or surgical conditions that could preclude a successful implantation
- Mental illness
- Medicine abuse, drug abuse or alcoholism
- Adiposity
- Severe damage to bone structures, potentially precluding any stable implantation of the implant components
- Neuromuscular disorders or illnesses
- Bone tumors in the region of implant anchoring
- Wound healing disorders
- Lack of patient compliance
- Foreign body sensitivity to the implant materials
- Conditions not listed under indications

Side effects and interactions

The application or improper use of this system entails the following risks:

- Implant failure due to excessive load
 - Bending
 - Loosening
 - Breakage
- Inadequate fixation
- Failed or delayed fusion
- Infection
- Vertebral body fracture
- Injuries to
 - nerve roots
 - spinal cord
 - blood vessels
 - organs
- Changes of the normal spine lordosis
- Impairment of the gastrointestinal, urological and/or reproductive systems
- Pain or discomfort
- Bursitis
- Diminishing bone density due to load avoidance
- Bone atrophy/fracture above or below the spine section treated with the implant
- Limited physical capacity
- Persistence of symptoms treated by the implantation
- Death

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present instructions for use.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from wrong indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The instructions for use for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically removed components.
- Implants that have been used before must not be reused.
- Only combine Aesculap modular implant components with each other.

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

- Delayed healing can cause implant breakage due to material fatigue.
- The attending physician shall make any decision with regard to the removal of implant components that have been used.
- Damage to the load-bearing implant and/or bone structures can lead to loosening of the components, fracture of the bone or implant, dislocation and migration and other severe complications.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.

Damage to implants caused by processing and resterilization!
► Do not reprocess or resterilize the implants.

Application

Risk of injury caused by incorrect operation of the product!
► Attend appropriate product training before using the product.
► For information about product training, please contact your national B. Braun/Aesculap agency.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following points:

- Selection of the implant components and their sizes
- Positioning of the implant components in the bone
- Definition of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components available
- Highly aseptic operating conditions
- Implantation instruments, including the special Aesculap implant system instruments, complete and in working condition
- Operating surgeon and operating team familiar with information on the operating technique and available range of implants and instruments; information materials in this respect are complete and ready to hand.
- Familiarity with the rules governing surgical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors
- Implant manufacturer consulted in case of ambiguous preoperative situation and if implants were found in the area to be treated.

The surgical procedure was explained to the patient and the patient's understanding of the following information has been documented:

- Excessive loads can cause implant loosening or breakage in cases of delayed or incomplete fusion.
- The life-span of the implant depends on the patient's body weight.
- The implant components must not be overstrained by excessive loads, hard physical labor or sporting activities.
- Revision surgery can become necessary in cases of implant loosening, implant breakage or loss of correction.
- Smokers present an increased risk of fusion failure.
- The patient must undergo postoperative medical checks of the implant components at regular intervals.

Implantation of the S⁴ Spinal System rod connectors requires the following steps:

- Select the appropriate S⁴ Spinal System variant and anchoring implants according to the individual indication, preoperative planning and bone situation found intraoperatively.
- To prevent internal stresses on, and weakening of the implants: avoid scoring or scratching of the implant components.
- Position the rod connector on the rod at the location planned by the operating surgeon using the connector insertion/holding instrument intended for this purpose.
- Do not remove the clamping screws from the rod connector assembly.
- Do not reassemble the clamping screws with the other rod connector components.
- Do not implant the rod connector if the clamping screws are missing.

CAUTION

Rod connector is unusable due to the clamping screws being removed or missing!

- **Select another rod connector for implantation.**

- Make certain the rod connector is correctly positioned before tightening the clamping screws.
- For tightening and loosening of the clamping screws, always use the screwdriver and countering instrument provided for these purposes.
- Always tighten the clamping screws with the 4 Nm torque wrench FW207R intended for this purpose.

WARNING

Risk of injury due to insufficient clamping stability provided by the clamping screws!

- Check that no patient tissue is lodged within the rod connector assembly.
- Position the rod connector correctly.
- Check that the rods are fully inserted into the rod connector.
- Tighten the clamping screws with the 4 Nm torque wrench FW207R intended for this purpose.

WARNING

Risk of injury caused by over-/under-tightening of the clamping screws!

- Always apply exactly the specified torque of 4 Nm when tightening the clamping screws.

WARNING

Inadequate fixation due to incorrect rod position!

- Always position the rod connector so that the rod is fully inserted into the rod connector.
- Always position rods in such a way that the hexagonal-end or the tip is outside the rod connector clamping region.

WARNING

Inadequate fixation due to patient tissue lodged within the rod connector assembly!

- Check that no patient tissue is lodged within the rod connector assembly.
- Dislodge patient tissue from the rod connector assembly or select another rod connector for implantation.

- Only allow axial load on the clamping screws during tightening or loosening. Avoid induction of lateral forces when tightening or loosening the clamping screws.

WARNING

Damage to the hexagonal socket in the clamping screws due to incorrect application of the screwdriver or torque wrench!

- Make certain the hexagonal tip of the screwdriver or torque wrench is fully inserted in the hexagonal socket of the clamping screws during tightening or loosening.

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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