
Aesculap Sterile Technology


PrimeLine

Sterile container system
Edition 06/2005

-
- (GB)** **Instructions for use/Technical description**
Sterile container system
 - (D)** **Gebrauchsanweisung/Technische Beschreibung**
Sterilcontainer-System
 - (F)** **Mode d'emploi/Description technique**
Système de conteneurs de stérilisation
 - (E)** **Instrucciones de manejo/Descripción técnica**
Sistema de contenedores estériles
 - (I)** **Istruzioni per l'uso/Descrizione tecnica**
Sterile container system

Sterile container system

Legend

- 1 Outer lid
- 2 Outer lid lock
- 3 Inner lid
- 4 Inner lid lock
- 5 Bottom
- 6 Indicator seal retainer
- 7 Indicator seal
- 8 Handle
- 9 Slot (for the indicator seal on the inner lid lock)
- 10 Universal filter retainer
- 11 Cap
- 12 pushbuttons
- 13 Lid latch (for optional outer lid of BASIS container)
- 14 Plastic cover
- 15 Perforation field cover
- 16 Germ retention system
- 17 Adapter frame
- 18 Finger tray
- 19 Handling pin
- 20 Ribbed cover grid
- 21 Mounting handle
- 22 Plastic cover lock

Symbols on the product

	Attention, see instructions for use
---	-------------------------------------

Contents

1. Safe handling	3
2. Product description	3
2.1 Intended use	3
2.2 Mode of operation	3
3. Preparation and setup	3
3.1 First use	3
4. Working with the sterile container	4
4.1 Set-up system	4
4.2 Function check	5
4.3 Safe operation	5
5. Maintenance	7
6. Care	7
7. Troubleshooting list	8
8. Technical Service	9
9. Accessories/Spare parts	9
10. Technical specifications	9
11. Extracts from relevant standards	9
11.1 Standards cited	9
11.2 Storage periods for sterile materials – extract from DIN 58953-9, sterile materials supply	9

1. Safe handling

WARNING

Risk of contamination of sterile supplies due to inadequately sealed sterile container!

The sealing of the sterile container, and its function as a germ barrier, will be compromised if the sterile container is combined with components from other manufacturers.

➤ **Only combine Aesculap sterile container products with each other.**

- Prior to use, check for proper condition and functioning of the product.
- In order to avoid damaging the device through improper assembly/operation and mitigating the manufacturer's liability:
 - Use the product only according to these instructions for use.
 - Follow the safety instructions and maintenance advisories.
 - Never use damaged or faulty sterile containers. Replace any damaged components immediately with original spare parts.
 - If sterile container components are repaired in a way that could affect the container's germproof qualities: Inspect the sterile container thoroughly before use.
- Make sure that the product and accessories are operated and used only by qualified or trained and experienced personnel.
- Keep the instructions for use accessible to all staff.
- Adhere to applicable norms.
- Follow the general guidelines and aseptic principles when handling contaminated items that have undergone or are to undergo sterilization.

2. Product description

2.1 Intended use

The Aesculap sterile container system is intended for use as a sterile packaging for instruments and textiles due for steam sterilization according to a validated steam sterilization process (e.g. in a sterilizer complying with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated according to EN 554/ISO 13683).

After being sterilized, the sterile materials are stored in the sterile container until they are used.

Note

Please contact your Aesculap representative if your Aesculap sterile containers are to be used in any other steam sterilization process.

2.2 Mode of operation

The Aesculap sterile container system complies with the stipulations of DIN EN 868-1, DIN EN 868-8, DIN 58953-9 and ISO 11607.

Sterile containers with a perforated lid and a closed basin have been validated for steam sterilization in a sterilizer complying with EN 285, through:

- Fractional vacuum processes

Sterile containers with a perforated lid and a perforated basin are also suitable for steam sterilization in a sterilizer complying with EN 285, through:

- Gravitational procedures

Note

The suitability of any specific process must be proven through a validation at the respective site of application.

3. Preparation and setup

3.1 First use

- Thoroughly clean the new container prior to first use.
- After cleaning, use a suitable filter, see Changing filters.

Sterile container system PrimeLine:

The germ retention system **16** is integrated with the system.

Sterile container system

4. Working with the sterile container

4.1 Set-up system

Removing the outer lid/inner lid/plastic cover

If an outer lid is used, it can be removed to clean the sterile container and, if it is contaminated, it can be separated from the inner lid.

VARIO container (by default with outer and inner lid):

- Push down the outer lid lock **2** and remove the outer lid **1** and the inner lid **3**.

BASIS container (with retrofitted outer lid):

- Dismount the combined outer lid **1** and inner lid **3** from the bottom **5**.
- Loosen the lid latch **13** and remove the outer lid **1**.

Sterile container system PrimeLine:

- Remove the plastic cover **14** from the bottom **5**.

Cleaning and disinfecting the sterile container

WARNING

Risk of corrosion and stress cracks due to improper cleaning and disinfection!

- Follow the specifications provided by the cleaning and disinfecting agents' manufacturers when cleaning and disinfecting the product.

WARNING

Additional instruction for PrimeLine:

- Do not use solvents such as acetone or benzene.
- Final rinse with fully deionized water only.
- Make certain that there are no residues from the cleaning process left on the lid.

- Use fully deionized water of a quality according to EN 285 Annex B for reprocessing the product.

Note

Literature on cleaning and disinfecting:

"Test Series and Statements; Working Group Instrument Preparation; C 75911", available through Aesculap.

The recommendations of the Working Group Instrument Preparation offer further guidance for careful reprocessing of reusable medical products: "Proper Maintenance of Instruments", see "www.a-k-i.org".

Material-specific requirements for the cleaning and disinfection of sterile containers

	Disposal containers made of stainless steel	Standard containers made of anodized aluminum	PrimeLine container lids made of plastic
Ultrasound treatment, in manual and mechanical reprocessing	√	√	√
Cleaning/disinfecting agent*	acid, pH-neutral, alkaline	pH-neutral, mildly alkaline	pH-neutral, mildly alkaline
Thermal disinfection with deionized water, 5 min at 93 °C	√	√	√
Mechanical drying with hot air at up to 120 °C	√	√	√
*Mildly alkaline cleaning/disinfection may cause changes on surfaces or of the colors of the aluminum lids (fading, staining). Such changes do not affect the functionality of the product. For colored aluminum lids, use a process involving a neutral cleaning and disinfecting agent, and fully deionized water, if possible.			

Manual preparation

- Use suitable cleaning and disinfecting agents. Always adhere to the manufacturers' instructions.
- Wipe surfaces with a clean, lint-free cloth and a pH-neutral cleaning and disinfecting agent.
- Remove encrusted debris with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- Carry out the final rinse with fully deionized water.

Mechanical preparation

- Use suitable cleaning and disinfecting agents. Always adhere to the manufacturers' instructions.
- Use suitable storage racks according to the instructions of the manufacturers of the respective devices for mechanical cleaning.
- Carry out the final rinse with fully deionized water.

Changing filters

Change filters at the following intervals, depending upon filter type:

- Disposable filter: prior to each sterilization cycle
- Reusable filter: after max. 1000 sterilization cycles
- PrimeLine germ retention system: after max. 5000 sterilization cycles

VARIO container and BASIS container:

- Press simultaneously on both pushbuttons **12** on the universal filter retainer **10**.
- Remove the universal filter retainer **10**.
- Insert a new filter and remount the universal filter retainer **10**.
- Push down the cap **11** on the universal filter retainer **10** until you hear it click into place.

Sterile container system PrimeLine:

- Turn the ribbed cover grid **20** anticlockwise until it is unlocked.
- Remove the ribbed cover grid **20** of the germ retention system **16**.
- Turn the germ retention system **16** anticlockwise at the finger tray **18** or the mounting handle **21** until it disengages from the adapter frame **17**.
- Lift the germ retention system **16**, gripping it at the handling pin **19**, and remove it.
- For installing the germ retention system **16**, retrace the above steps.

4.2 Function check

- Visually inspect all components of the sterile container for damage and correct functioning.
 - Metal parts not deformed
 - Aluminum lid not warped
 - Seals intact
 - Plastic parts not cracked
 - Permanent filter intact
 - Lock fully functional (engages)
- Use sterile containers only if they are in immaculate condition. Replace any damaged components immediately with original spare parts or have the parts repaired.

4.3 Safe operation

Risk of contamination of sterile materials due to inadequately sealed sterile container!

The sealing of the sterile container, and its effectiveness as a germ barrier, will be compromised if the sterile container is combined with components from other manufacturers.

- **Only combine Aesculap sterile container products with each other.**

Loading sterile container

Instruments

According to DIN EN 868-8 and DIN 58953-9, the following maximum load of the container is recommended:

- Standard container: 10 kg
 - Half container: 5 kg
 - 3/4-container: 7 kg
- Store instruments in wire baskets in suitable supports. In doing this, lay hollow components, bowls, plates, etc. face down and on a slant.
 - For double-packaging of sterile materials in accordance with DIN 58953-9 wrap the loaded wire basket separately.

VARIO container and BASIS container:

- Load sterile containers in such a way that their universal filter retainers **10** remain unobstructed (VARIO and BASIS container).
- Use the inner lid lock **4** to closure the inner lid **3** to the container bottom **5**.

The inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile container system PrimeLine:

- Load the sterile container in such a way that the ribbed cover grid **20** in the plastic cover **14** remains unobstructed.
- Use the plastic cover lock **22** to lock the plastic cover **14** on the bottom **5**.
The plastic cover lock **22** has to engage perceptibly. If this is not the case: have the sterile container repaired, see Technical Service.

Sterile container system

Textiles

- Pack the textiles together in such a way that they are upright in the sterile container.
- Make certain that when the sterile container is fully loaded, it is still possible to insert an outstretched hand between the individual items without difficulty.
- Load sterile containers in such a way that their universal filter retainers **10** remain unobstructed (with VARIO and BASIS containers).
- Use the inner lid lock **4** to closure the inner lid **3** to the container bottom **5**.
The inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired see Technical Service.

Sterile container system PrimeLine:

- Load the sterile container in such a way that the ribbed cover grid **20** in the plastic cover **14** remains unobstructed.
- Use the plastic cover lock **22** to lock the plastic cover **14** on the bottom **5**.
The plastic cover lock **22** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Inserting the indicator seal

- Having loaded the sterile container, make notes of the following data on the indicator seal **7**: date of sterilization, sterilization number, expiry date, name and signature.
- Slide the indicator seal **7** from the outer side into the indicator seal retainer **6** in such a way that the indicator segment engages in the slot on the lid lock and seals the closure.

Loading the sterilizer

Prepare the sterile container and the sterilizer for loading of the latter in the following way:

Risk of vacuum damage to the sterile container due to insufficient pressure venting!

- **Do not use outer packaging for the sterile containers.**
- **Do not block the air flow through the perforation fields in the basin and inner lid under any circumstances (applicable to all container versions).**
- **Do not cover the sterile container with foil packaging.**
- Follow the sterilizer manufacturer's instructions.
- If possible, load the sterilizer with a mixture of textiles and instrument containers.
- Always place heavy sterile containers at the bottom of the sterilizer.
- Sterilize VARIO containers with their outer lids attached.
- Always carry the sterile containers by their handles.
- Stack sterile containers in the sterilizer if necessary.
- Transport stacks of sterile containers in such a way that the stacks do not topple over.

Sterilization method and parameters

- Sterilize with steam, taking note of the following:
The sterilization has to be done according to a validated steam sterilization procedure (e.g. in a sterilizer according to EN285/ANSI/AAMI/ISO11134-1993, ANSI/AAMI ST46-1993 and validated according to EN 554/ISO 13693).

Unloading the sterilizer

Risk of burns due to hot sterile container after sterilization!

- **Always wear protective gloves when handling a container.**
- Always carry the sterile containers by their handles.
- Transport stacks of sterile containers in such a way that the stacks do not topple over.
- To prevent formulation of condensate: Never allow sterile containers to cool down on a cold floor or in a drafty place.

Checking and providing sterile materials

The contents of a sterile container are only considered sterile if the container has been properly sterilized.

If this is not the case, the sterile materials must be reprocessed again.

Risk of contamination from improperly sterilized materials!

- **Before providing sterile materials, check to ensure that the sterilization was successful.**

- Before opening the sterile container, make certain that its indicator seal is not damaged.
- Make certain that the color of the indicator point has changed.
- Make certain that all container components, particularly the lid locks, are intact.
- Make certain that there is no condensate in the container.
If any condensate is found, eliminate the cause, see Troubleshooting list.

Storing sterile containers

Note

The sterile containers may be stored in stacks.

- Store sterile containers in a dry, clean and protected place.
- Comply with DIN 58953-9 standards for storage periods and conditions, see Storage periods for sterile materials – extract from DIN 58953-9, sterile materials supply.
Recommended storage periods according to DIN 58953-9 guidelines: 6 months.

5. Maintenance

Note

The sterile containers must only be tested and repaired by appropriately trained, knowledgeable or experienced persons.

The seals in the container lids have been tested, in conformance with DIN EN 868-8, Appendix G for a minimum service life of 5000 sterilization cycles.

- In case of visible damage, change the seal immediately.
- Have the lid sent to Aesculap, see Technical Service.
– or –
Order seal and sticker, see Accessories/Spare parts in brochure no. C 404 01 and replace seal.
- If sterile container components are repaired in a way that could affect the container's germproof qualities: Inspect the sterile container thoroughly before use.

6. Care

Note

Lightly lubricate moving parts such as latch hinges with a sterilizable and steam-permeable maintenance oil (e.g. Aesculap Sterilit spray JG 600 or maintenance oil JG 598).

Sterile container system

7. Troubleshooting list

Error	Cause	Remedy
Too much condensate inside the sterile container	Temperature of sterile materials too low prior to sterilization	Allow sterile materials to come to room temperature (approx. 20 °C) before sterilization
	Textiles too damp	Sterilize dry textiles only
	Sterile container too heavy	Standard container with instruments: max. load 10 kg with textiles: max. load 8 kg
		Half container with instruments: max. load 5 kg
		3/4-container with instruments: max. load 7 kg
	Sterile materials improperly positioned	Lay hollow components, bowls, plates etc. face down and on a slant
		Arrange textiles in loose vertical piles, do not press them together
	Sterile container incorrectly positioned in sterilizer	Always place heavy sterile containers at the bottom
	Sterile containers supplied for use immediately after sterilization	Allow sterile containers to cool down to room temperature prior to processing
	Sterile containers improperly positioned during cooling phase	Do not store sterile containers on a floor or in a drafty place Store sterile containers in air-conditioned areas according to DIN 58953-9 standards
	Sterilizer properties do not meet DIN EN 285 standards	Have sterilizer serviced regularly Check drying vacuum Check drying time
		Check steam quality and upgrade if necessary
	Empty-cycle and vacuum test not run daily before sterilization begins	Run empty-cycle and vacuum test daily before beginning sterilization
Condensate in groove of (inner) lid	Unsuitable sterilizer cycle selected	Select a program according to the load
	Sterilizer door left open too long, sterilizer cooled down	Load and unload sterilizer quickly
	Multiple sterile containers or heavy sterile materials stacked	Subdivide sterilizer with shelving Use outer lid
Indicator seal color unchanged	Sterilization improperly performed, sterilizer faulty	Have sterilizer repaired by manufacturer
BASIS/VARIO sterile container or PrimeLine sterile container system deformed	Perforations covered during sterilization	Never cover perforations

Error	Cause	Remedy
Lid lock faulty	Outer lid used to carry sterile container	Always carry the sterile containers by their handles

8. Technical Service

For service, maintenance or repairs, contact your national B. Braun/Aesculap representative.

All guarantees and warranty rights as well as licenses will be void if modifications on the device are performed.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.

Other service addresses can be obtained from the address indicated above.

9. Accessories/Spare parts

Accessories and consumables can be found in brochure no. C 404 01.

10. Technical specifications

The various models and sizes of the sterile containers can be found in brochure no. C 404 01.

11. Extracts from relevant standards

11.1 Standards cited

The following standards are cited in connection with the sterile containers:

- DIN 58953-9
- DIN 58946-1, EN 285, EN 866-3
- EN 868-1, EN 868-8
- EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, with applicable EN 554/ISO 13683
- ISO 11607

11.2 Storage periods for sterile materials – extract from DIN 58953-9, sterile materials supply

Loss of sterility is not so much connected to the storage periods as to outside influences and the effects of storage, transport and handling. Therefore, blanket statements cannot be made regarding appropriate storage periods. The appropriate length of storage periods is determined by the hygiene control committee. Responsibility for storage conditions and duration of storage periods lies with the hospital operator or its medical director.

Storage periods should be defined differently (where appropriate) for each department, and in writing. In doing this, the type of packaging and the storage conditions should be taken into account.

Recommended storage periods for sterile materials in sterile containers

according to EN 868-1 and EN 868-8: 6 months

These information are recommendations which, if followed, will decrease the risk of contamination when sterile containers are transported and opened. Storage for longer periods can lead to contamination of outer surfaces, depending upon storage conditions. Although such longer periods do not necessarily lead to recontamination of the contents of the package during storage, they do increase the risk of contamination when sterile containers are transported or opened. These storage recommendations are made on the assumption that the materials will be stored in a dry area that is protected against dust. Longer storage periods are feasible, depending upon conditions at the storage site.

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Aesculap AG & Co. KG