Aesculap Sterile Technology

(GB) Instructions for use/Technical description

(USA) Sterile Container System Edition 03/2016

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




13

Legend

- 1 Outer lid
- 2 Outer lid lock
- 3 Inner lid
- 4 Inner lid lock
- 5 Bottom
- 6 Indicator seal retainer
- 7 Indicator seal
- 8 Screw handle
- 9 Slot (for the indicator seal on the inner lid lock)
- 10 Universal filter retainer
- 11 Cap
- 12 Push buttons
- 13 Lid latch (for optional outer lid of BASIS container)
- 14 Plastic cover
- 15 Perforation field cover
- 16 Germ barrier system (permanent filter)
- 17 Adapter frame
- 18 Handling pin
- 19 Ribbed cover grid
- 20 Mounting handle
- 21 Plastic cover lock

Symbols on product and packages

Caution, general warning symbol
Caution, see documentation supplied with the product

Follow the instructions for use

Contents

1.	Safe handling	3
2.	Product description	3
2.1	Intended use.	3
2.2	Operating principle	3
3.	Preparation and setup	3
3.1	First use	3
4.	Working with the sterile container	4
4.1	System set-up	4
4.2	Function checks	4
4.3	Safe operation	5
5.	Validated reprocessing procedure	7
5.1	General safety notes	7
5.2	General information	7
5.3	Cleaning/disinfection.	7
5.4	Manual cleaning/disinfecting	8
5.5	Mechanical cleaning/disinfecting.	9
5.6	Inspection, maintenance and checks	9
6.	Troubleshooting list	10
7.	Technical Service	11
8.	Accessories/Spare parts	11
9.	Technical data	11
10.	Disposal	11
11.	Extracts from relevant standards	11
11.1	Standards cited	11
12.	Distributor in the US/Contact in Canada for product information and complaints	11

1. Safe handling

Risk of contamination of sterile goods from unsealed sterile containers!

The sealing of the sterile container and its germ barrier function will be compromised if the sterile container is combined with components from other manufacturers.

► **Only combine Aesculap sterile container products with each other.**

- Before use, check the product is in full working order and in good condition, see function check section and function check poster C63301.
- To prevent damage caused by improper setup or use and in order not to compromise the warranty and liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Never use damaged or faulty sterile containers. Replace any damaged components immediately with original spare parts.
 - If sterile container components are repaired in a way that could affect the container's germproof qualities: Inspect the sterile container thoroughly before use.
- Ensure that only persons with the requisite training, expertise or experience will handle the product and its accessories.
- Keep the instructions accessible for personnel.
- Follow general guidelines and aseptic principles when handling contaminated items that have undergone or are to undergo sterilization.

2. Product description

2.1 Intended use

The Aesculap sterile container system is a reusable sterile barrier system that preserves the sterility of medical products until they are used or reach their use-by date. Medical devices can be sterilized, stored and transported in the Aesculap optics container. The container can also be used for transporting medical devices back to the decontamination center after use. The Aesculap sterile container system is suitable for steam sterilization.

Note

Please contact your Aesculap representative if your Aesculap sterile containers are to be used in any other steam sterilization process.

2.2 Operating principle

The Aesculap sterile container system meets the requirements of DIN 58953-9 and EN ISO 11607.

Sterile containers with a perforated lid and a closed bottom have been validated for steam sterilization in a sterilizer acc. to EN 285 through the:

■ fractionated vacuum process.

Sterile containers with a perforated lid and a perforated bottom are also suitable for steam sterilization in a sterilizer acc. to EN 285 through:

■ gravitational procedures.

Note

The suitability of any specific process must be validated at the site of application.

3. Preparation and setup

3.1 First use

- Thoroughly clean the new sterile container prior to first use.
- After cleaning, use a suitable filter, see Filter change.

Sterile Container System PrimeLine:

The germ barrier system **16** (permanent filter) is integrated in the system.

4. Working with the sterile container

4.1 System set-up

Removing the outer lid / inner lid / plastic cover

If an outer lid is used, it can be removed to clean the sterile container and, if it is soiled, it can be separated from the inner lid.

VARIO container (by default with outer and inner lid):

- ▶ Push down the outer lid lock **2** and remove the outer lid **1** and the inner lid **3**.

BASIS container (with retrofitted outer lid):

- ▶ Dismount the combined outer lid **1** and inner lid **3** from the bottom **5**.
- ▶ Loosen the lid latch **13** and remove the outer lid **1**.

Sterile Container System PrimeLine:

- ▶ Remove the plastic cover **14** from the bottom **5**.

Filter change

Replace the filter at the following intervals, depending upon the filter type:

- Disposable filter: prior to each sterilization cycle
- PTFE permanent filter: after max. 1 000 processing cycles
- PrimeLine germ retention system: after max. 5 000 processing cycles

VARIO container and BASIS container:

- ▶ Press simultaneously both push buttons **12** on the universal filter retainer **10**.
- ▶ Remove the universal filter retainer **10**.
- ▶ Insert a new filter and remount the universal filter retainer **10**.
- ▶ Push down the cap **11** on the universal filter retainer **10** until you hear it click into place.

Sterile Container System PrimeLine:

- ▶ Turn the ribbed cover grid **19** anticlockwise until it is unlocked.
- ▶ Remove the ribbed cover grid **19** of the germ barrier system **16**.
- ▶ Turn the germ barrier system **16** anticlockwise with the mounting handle **20** until it is unlocked from the adapter frame **17**.
- ▶ Lift the germ barrier system **16**, gripping it at the handling pin **18**, and remove it.
- ▶ Install the germ barrier system **16** by the reverse sequence of steps.

4.2 Function checks

Note

*With the PrimeLine sterile container system, the ribbed cover grid **19** must be removed prior to checks!*

- ▶ Visually inspect all components of the sterile container before each use to ensure correct function and that there is no damage, see also function check poster C63301:
 - Metal parts are not deformed.
 - Aluminum lids are not warped.
 - The edges of the universal filter retainer **10** are in full-surface contact.
 - The seals are intact (no cracks).
 - The disposable single use filter has been changed.
 - The disposable single use filter is not damaged (holes, cracks).
 - Plastic parts are not damaged.
 - The permanent filter/PrimeLine sterile barrier system is undamaged (no kinks, holes, cracks or gaps).
 - The lock functions properly (engages).
 - Check both sides of the plastic lid for cracks.
 - ▶ Use sterile containers only if they have been visually inspected and are free from defects listed above. Replace any damaged components immediately with original spare parts or have the parts repaired.
- An illustrated description of the function check can be found on the function check poster C63301.

4.3 Safe operation

WARNING

Risk of contamination of sterile goods from unsealed sterile containers!

The sealing of any sterile container and its function as a germ barrier will be compromised if the sterile container is combined with components from other manufacturers.

- ▶ Only combine Aesculap sterile container products with each other.

CAUTION

Risk of non-sterility of container contents!

- ▶ Always carry the sterile containers by their handles.
- ▶ Never carry or lift the sterile container at the lid.
- ▶ Transport the sterile container in such a way that mechanical damage will not occur.

Loading the sterile containers

instruments

According to DIN EN 868-8 and DIN 58953-9 we recommend the following maximum load for the containers:

- Standard container: 10 kg
- Half container: 5 kg
- 3/4-container: 7.5 kg
- ▶ Store instruments in wire baskets in suitable supports. When doing so, lay hollow materials, dishes, plates, etc with the opening facing downwards at a slant.

VARIO container and BASIS container:

- ▶ Load the sterile containers in such a way that their universal filter retainers **10** remain unobstructed. Max. loading height: Up to about 2 cm below the edge of the container bottom.
- ▶ Use the inner lid lock **4** to fasten the inner lid **3** to the container bottom **5**.

The inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile Container System PrimeLine:

- ▶ Load the sterile container in such a way that the ribbed cover grid **19** in the plastic cover **14** remains unobstructed. Max. loading height: up to approx. 2 cm below the rim of the container base.
- ▶ Use the plastic cover lock **21** to lock the plastic cover **14** on the bottom **5**. The plastic cover lock **21** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Textiles

- ▶ Pack folded textiles in such a way that they fit vertically in the sterile container.
- ▶ Make certain that when the sterile container is fully loaded, it is still possible to insert a flat hand between the individual items without difficulty.
- ▶ Load sterile containers in such a way that their universal filter retainers **10** remain unobstructed.
- ▶ Use the inner lid lock **4** to fasten the inner lid **3** to the container bottom **5**.

The inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile Container System PrimeLine:

- ▶ Load the sterile container in such a way that the ribbed cover grid **19** in the plastic cover **14** remains unobstructed.
- ▶ Use the plastic cover lock **21** to lock the plastic cover **14** on the bottom **5**. The plastic cover lock **21** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Inserting the indicator seal

- ▶ After loading the sterile container, label the indicator seal **7** according to DIN 58953-7 (e.g. contents, batch no., expiry date etc.).
- ▶ Slide the indicator seal **7** from outside into the indicator seal retainer **6** so that the indicator segment engages in the slot on the lid lock and seals the lock.
 - or -
- ▶ Insert a plastic seal (e.g. JG739) at the lock of the closed container.

Loading the sterilizer

WARNING

Risk of vacuum damage to the sterile container due to inadequate pressure equalization!

- ▶ Do not use outer packaging for the sterile containers.
- ▶ Do not obstruct the air flow through the perforation fields in the bottom and in the inner lid under any circumstances (applicable to all container versions).
- ▶ Do not place foil packaging directly on the sterile container.

Note

VARIO containers and BASIS containers can be sterilized with the outer lid in place!

- ▶ Follow the sterilizer manufacturer's recommendations.
- ▶ Always place heavy sterile containers at the bottom of the sterilizer.
- ▶ Always carry the sterile containers by their handles.

Note

Sterile containers can be stacked in the sterilizer.

- ▶ Transport stacks of sterile containers in such a way that the stacks do not topple over.

Sterilization

CAUTION

Risk of sterilization failure!

- ▶ **Sterilize the containers only by approved and validated sterilizing processes.**

- ▶ Sterilize with steam, taking note of the following: Sterilization must be carried out through a validated steam sterilizing process (e.g. in a steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665).

Unloading the sterilizer

WARNING

Risk of burns due to a hot sterile container after sterilization!

- ▶ Always wear protective gloves when unloading the sterilizer.

Transporting the sterile container

CAUTION

Risk of non-sterility of container contents!

- ▶ Always carry the sterile containers by their handles.
- ▶ Never carry or lift the sterile container at the lid.
- ▶ Transport the sterile container in such a way that mechanical damage will not occur.

Storing the sterile containers

Note

The sterile containers may be stored in stacks.

- ▶ Store sterile containers in a dry, clean and protected place.
- ▶ Observe storage periods and storage conditions acc. to DIN 58953-9.

Checking and commissioning the sterile materials

The contents of a sterile container can be commissioned as sterile only if the container has been properly sterilized, stored and transported.

If this is not the case, the sterile materials must be processed again.

DANGER

Risk of contamination from improperly sterilized materials!

- ▶ Prior to commissioning the sterile materials, check to ensure that the sterilization was successful.

- ▶ Make certain that the color of the indicator point has changed.
- ▶ Make certain that all container components, particularly the lid locks, are intact.
- ▶ Check the plastic lid for cracks.
- ▶ Ensure the container seal is intact.

5. Validated reprocessing procedure

5.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

The recommended chemistry was used for validation.

Due to process tolerances, the manufacturer's specifications can only serve as an approximate guide for assessing the processing procedures applied by the individual operator/processors.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com.

5.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and reprocessing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehyde) should be used.

Excessive dosages of neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking on stainless steel becoming unreadable either visually or by machine.

Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Visual material changes, e.g. fading or discoloration of aluminum. For aluminum, the application/process solution only needs to be pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ Do not use process chemicals that cause stress cracking of plastics such as PPSU, or attack silicone softeners and cause brittleness of the material.
- ▶ For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

5.3 Cleaning/disinfection

Product-specific safety instructions for the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents approved for aluminum, plastics and high-grade steel according to manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and exposure time.

5.4 Manual cleaning/disinfecting

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfection process if necessary.

Manual cleaning with disinfecting cleaning by wipe disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals	
I	Cleaning	RT (cold)	-	-	D-W	-	
II	Dry	RT	-	-	-	-	
III	Wipe disinfection	-	>1	-	-	I Alcohol Denat. 70 % (Ethanol B. Braun)	II Aldehyde-free surface disinfectant (Stabimed)
IV	Final rinse	RT (cold)	0,5	-	FD-W	not required	Rinse off cleaning chemicals; do not leave residues
V	Dry	RT	-	-	-	-	

D-W: Drinking water

FD-W: Demineralized water

RT: Room temperature

Phase I

- Wipe off gross debris prior to cleaning.

Phase II

- Dry the instrument with a lint-free cloth or medical compressed air.

Phase III

- Carry out wipe disinfection.

Phase IV

- When using aldehyde-free disinfectants:
 - After the specified exposure time (at least 1 min), rinse the disinfected surfaces under running FD water.
- Drain any remaining water fully.

Phase V

- Dry the product with a lint-free cloth or in a heating cabinet.

Note

For PrimeLine use only Alcohol Denat. 70 % (e.g. ethanol B. Braun).

5.5 Mechanical cleaning/disinfecting

Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Neutral: B. Braun Helimatic Cleaner neutral, pH neutral, 0.5% working solution Mildly alkaline: ■ Concentrate: – pH = 9.5 – <5% anionic surfactant ■ 0.5 % solution Alkaline: Reprocessing is possible up to pH 10.5, provided the cleaning agent is manufacturer-approved for cleaning aluminum or plastic sterile containers. Note: Application of mild-alkaline or alkaline cleaning agents may cause discoloration of colored anodized aluminum, which however will not affect the functionality of the product.
III	Intermediate rinse	>10/50	1	FD-W	Especially for PrimeLine it must be ensured that the surface is rinsed without leaving any residues.
IV	Thermal disinfection	90/194	5	FD-W	Other process parameters may be feasible with agreement by the hospital hygienist.
V	Dry	-	-	-	According to mechanical cleaning program Do not use rinsing agent for PrimeLine.

D-W: Drinking water

FD-W: Demineralized water

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning/disinfection process if necessary.

Note

Ultrasonic treatment can be applied in manual processing as well as in mechanical processing. Temperatures of up to 120 °C are permitted for machine drying with hot air.

5.6 Inspection, maintenance and checks

Note

The sterile containers may be tested and repaired only by persons with the appropriate training, expertise or experience.

The seals in the container lids have been tested, in conformance with DIN EN 868-8, Appendix G for a minimum service life of 5 000 sterilization cycles.

- In case of visible damage, change the seal immediately.
- Send lid to Aesculap, address see Technical Service.
- If necessary, lightly lubricate moving metal parts, e.g. latch hinges, with sterilizable and steam-permeable maintenance oil (e.g. STERILIT® oil spray JG600 or STERILIT® I drip feed oiler JG598).

6. Troubleshooting list

Malfunction	Cause	Remedy
Excessive amounts of condensate inside the sterile container	Temperature of sterile materials too low prior to sterilization	Allow sterile materials to reach room temperature (approx. 20 °C).
	Textiles too damp	Sterilize dry textiles only.
	Sterile container too heavy	Standard container With instruments: max. load 10 kg With textiles: max. load 8 kg
		Half container with instruments: max. load 5 kg
		3/4-container with instruments: max. load 7 kg
	Sterile materials improperly positioned	Lay hollow materials, dishes, plates, etc with the opening facing downwards at a slant.
		Arrange textiles in loose vertical piles, do not press them together.
	Sterile container incorrectly positioned in sterilizer	Always place heavy sterile containers at the bottom.
	Sterile containers processed immediately after sterilization	Allow sterile containers to cool down to room temperature prior to processing.
	Sterile containers improperly positioned during cooling phase	Do not store sterile containers on a floor or in a drafty place. Store sterile containers in air-conditioned areas according to DIN 58953-9.
Condensate in groove of (inner) lid	Sterilizer properties do not comply with DIN EN 285	Have sterilizer serviced regularly. Check drying vacuum. Check drying time.
		Check steam quality and upgrade if necessary.
		Run empty-cycle and vacuum test daily before beginning sterilization.
	Empty-cycle and vacuum test not run daily before sterilization begins	
	Unsuitable sterilizer cycle selected	Select cycle in accordance with load.
Indicator seal color unchanged	Sterilization performed incorrectly Sterilizer defective	Have sterilizer repaired by manufacturer.
	Indicator seals stored improperly	Observe storage conditions as specified on the packaging of the indicator seals.
Basis/vario container or PrimeLine sterile container system deformed	Perforations covered during sterilization	Never cover perforation fields.
Inner or outer lid cannot be positioned or locked on the bottom component	Container lid or bottom deformed or damaged due to improper handling (e.g. carrying by the lid)	Replace container lid or bottom, or have components repaired by the manufacturer.

7. Technical Service

Risk of injury and/or malfunction!

► **Do not modify the product.**

► For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: ats@aesculap.de

Or in the US:

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042 USA

Aesculap Repair Hotline

Phone: +1 (800) 214 -3392

Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

8. Accessories/Spare parts

Accessories and supplies are listed in brochure no. C40401.

9. Technical data

The variants and sizes of the sterile containers are listed in the following brochures:

- C40401 (German)
- C40402 (English)
- C40404 (Spanish)

10. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

► Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

11. Extracts from relevant standards

11.1 Standards cited

The following standards are cited in connection with the sterile containers:

- DIN 58953-9
- EN 285
- EN 868-8
- EN 285 in combination with EN ISO 17665
- EN ISO 11607

12. Distributor in the US/Contact in Canada for product information and complaints

3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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