

Aesculap Sterile Technology

Sterile container system
Edition 07/2009

-
- (GB)** Instructions for use/Technical description
Sterile container system
 - (D)** Gebrauchsanweisung/Technische Beschreibung
Sterilcontainer-System
 - (F)** Mode d'emploi/Description technique
Système de containers de stérilisation
 - (E)** Instrucciones de manejo/Descripción técnica
Sistema de contenedores estériles
 - (I)** Istruzioni per l'uso/Descrizione tecnica
Sterile container system
 - (NL)** Gebruiksaanwijzing/Technische beschrijving
Steriele container-systeem
 - (DK)** Brugsanvisning/Teknisk beskrivelse
Sterilcontainer-system
 - (S)** Bruksanvisning/Teknisk beskrivning
Sterilcontainersystem

Sterile container system

Legend

- 1 Outer lid
- 2 Outer lid lock
- 3 Inner lid
- 4 Inner lid lock
- 5 Bottom
- 6 Indicator seal retainer
- 7 Indicator seal
- 8 Handle
- 9 Slot (for the indicator seal on the inner lid lock)
- 10 Universal filter retainer
- 11 Cap
- 12 Push buttons
- 13 Lid latch (for optional outer lid of BASIS container)
- 14 Plastic cover
- 15 Perforation field cover
- 16 Germ barrier system (permanent filter)
- 17 Adapter frame
- 18 Handling pin
- 19 Ribbed cover grid
- 20 Mounting handle
- 21 Plastic cover lock

Symbols on product and packages

Symbol	Explanation
	Caution: See documentation supplied with the product
	Follow the instructions for use

Contents

1. Safe handling	3
2. Product description	3
2.1 Intended use	3
2.2 Operating principle	3
3. Preparation and setup	4
3.1 First use	4
4. Working with the sterile container	4
4.1 System set-up	4
Removing the outer lid/inner lid/plastic cover	4
Filter change	4
4.2 Function checks	4
4.3 Safe operation	4
Loading the sterile containers	5
Inserting the indicator seal	5
Loading the sterilizer	5
Sterilization method and parameters	5
Unloading the sterilizer	6
Transporting the sterile container	6
Storing the sterile containers	6
Checking and commissioning the sterile materials	6
5. Validated processing procedure	6
5.1 General notes	6
5.2 Cleaning/Disinfecting	7
5.3 Manual cleaning/disinfecting	8
Manual cleaning and wipe disinfection	8
5.4 Mechanical neutral or mild-alkaline cleaning and thermal disinfection	9
5.5 Inspection, maintenance and checks	9
6. Troubleshooting list	10
7. Technical Service	11
8. Accessories and replacement parts	11
9. Technical specifications	11
10. Extracts from relevant standards	11
10.1 Standards cited	11

1. Safe handling

WARNING

Risk of contamination of sterile materials due to sealing failure of sterile container!

The sealing of the sterile container and its germ barrier function will be compromised if the sterile container is combined with components from other manufacturers.

➤ **Only combine Aesculap sterile container products with each other.**

- Prior to use, check for proper condition and functioning of the product.
- To avoid damaging the product through improper setup or operation and mitigating the warranty and liability on the part of the product's manufacturer:
 - Use the product only according to these instructions for use.
 - Follow the safety instructions and maintenance advisories.
 - Never use damaged or faulty sterile containers. Replace any damaged components immediately with original spare parts.
 - If sterile container components are repaired in a way that could affect the container's germproof qualities: Inspect the sterile container thoroughly before use.
- Ensure that only persons with the requisite training, expertise or experience will handle the product and its accessories.
- Keep the instructions accessible for personnel.
- Follow general guidelines and aseptic principles when handling contaminated items that have undergone or are to undergo sterilization.

2. Product description

2.1 Intended use

The Aesculap sterile container system is a reusable sterile barrier system that preserves the sterility of medical products until they are used or reach their use-by date. Medical products can be sterilized, stored and transported in the Aesculap optics container. The container can also be used for transporting medical products back to sterile processing after use. The Aesculap sterile container System is suitable for steam sterilization.

Note

Please contact your Aesculap representative if your Aesculap sterile containers are to be used in any other steam sterilization process.

2.2 Operating principle

The Aesculap sterile container system meets the requirements of DIN 58953-9 and EN ISO 11607.

Sterile containers with a perforated lid and a closed bottom have been validated for steam sterilization in a sterilizer acc. to EN 285 through the:

- fractionated vacuum process.

Sterile containers with a perforated lid and a perforated bottom are also suitable for steam sterilization in a sterilizer acc. to EN 285 through:

- gravitational procedures

Note

The suitability of any specific process must be validated at the site of application.

Sterile container system

3. Preparation and setup

3.1 First use

- Thoroughly clean the new sterile container prior to first use.
- After cleaning, use a suitable filter, see Filter change.

Sterile container system PrimeLine:

The germ barrier system **16** (permanent filter) is integrated in the system.

4. Working with the sterile container

4.1 System set-up

Removing the outer lid/inner lid/plastic cover

If an outer lid is used, it can be removed to clean the sterile container and, if it is soiled, it can be separated from the inner lid.

VARIO container (by default with outer and inner lid):

- Push down the outer lid lock **2** and remove the outer lid **1** and the inner lid **3**.

BASIS container (with retrofitted outer lid):

- Dismount the combined outer lid **1** and inner lid **3** from the bottom **5**.
- Loosen the lid latch **13** and remove the outer lid **1**.

Sterile container system PrimeLine:

- Remove the plastic cover **14** from the bottom **5**.

Filter change

Replace the filter at the following intervals, depending upon the filter type:

- Disposable filter: prior to each sterilization cycle
- PTFE permanent filter: after max. 1 000 processing cycles
- PrimeLine germ barrier system: after max. 5 000 processing cycles

VARIO container and BASIS container:

- Press simultaneously both push buttons **12** on the universal filter retainer **10**.
- Remove the universal filter retainer **10**.
- Insert a new filter and remount the universal filter retainer **10**.
- Push down the cap **11** on the universal filter retainer **10** until you hear it click into place.

Sterile container system PrimeLine:

- Turn the ribbed cover grid **19** anticlockwise until it is unlocked.
- Remove the ribbed cover grid **19** of the germ barrier system **16**.
- Turn the germ barrier anticlockwise system **16** with the mounting handle **20** until it is unlocked from the adapter frame **17**.

- Lift the germ barrier system **16**, gripping it at the handling pin **18**, and remove it.
- Install the germ barrier system **16** by the reverse sequence of steps.

4.2 Function checks

Note

*With the PrimeLine sterile container system, the ribbed cover grid **19** must be removed prior to checks!*

- Visually inspect all components of the sterile containers for damage and correct functioning prior to each use. Check that:
 - metal parts are not deformed
 - aluminum lids are not warped
 - the edges of the universal filter retainer **10** are seated in full-surface contact
 - the seals are intact (no cracks, ...)
 - plastic parts are not cracked
 - the permanent filter/PrimeLine sterile barrier system is undamaged (no kinks, holes, cracks or gaps)
 - the lock functions properly (engages)
- Use sterile containers only if they are in mint condition. Replace any damaged components immediately with original spare parts or have the parts repaired.

4.3 Safe operation

WARNING

Risk of contamination of sterile materials due to sealing failure of sterile container!

The sealing of any sterile container and its function as a germ barrier will be compromised if the sterile container is combined with components from other manufacturers.

- **Only combine Aesculap sterile container products with each other.**

CAUTION

Risk of non-sterility of container contents!

- **Always carry the sterile containers by their handles.**
- **Never carry or lift the sterile container at the lid.**
- **Transport the sterile container in such a way that mechanical damage will not occur.**

Loading the sterile containers

Instruments

According to DIN EN 868-8 and DIN 58953-9 we recommend the following maximum load for the containers:

- Standard container: 10 kg
 - Half container: 5 kg
 - 3/4-container: 7 kg
- Store instruments in appropriate storage devices in a tray, with hollow components, bowls, dishes, etc. with their openings at a downward angle.

VARIO container and BASIS container:

- Load the sterile containers in such a way that their universal filter retainers **10** remain unobstructed. Max. loading height: Up to about 2 cm below the edge of the container bottom.
- Use the inner lid lock **4** to fasten the inner lid **3** to the container bottom **5**.
- The inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile container system PrimeLine:

- Load the sterile container in such a way that the ribbed cover grid **19** in the plastic cover **14** remains unobstructed. Max. loading height: up to approx. 2 cm below the rim of the container base.
- Lock the plastic lid **14** with the plastic lid lock **21** on the bottom **5**. The plastic lid lock **21** must engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Textiles

- Pack folded textiles in such a way that they fit vertically in the sterile container.
- Make certain that when the sterile container is fully loaded, it is still possible to insert a flat hand between the individual items without difficulty.
- Load sterile containers in such a way that their universal filter retainers **10** remain unobstructed.
- Use the inner lid lock **4** to fasten the inner lid **3** to the container bottom **5**.
- The inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile container system PrimeLine:

- Load the sterile container in such a way that the ribbed cover grid **19** in the plastic cover **14** remains unobstructed.
- Lock the plastic lid **14** with the plastic lid lock **21** on the bottom **5**. The plastic lid lock **21** must engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Inserting the indicator seal

- After loading of the sterile container, enter the following data on the indicator seal **7**: Date of sterilization, sterilization number, expiry date, name and signature.
- Slide the indicator seal **7** from outside into the indicator seal retainer **6** so that the indicator segment engages in the slot on the lid lock and seals the lock.
- or -
- Insert a plastic seal (e.g. JG739) at the lock of the closed container.

Loading the sterilizer

WARNING

Risk of vacuum damage to the sterile container due to inadequate pressure equalization!

- **Do not use outer packaging for the sterile containers.**
- **Do not obstruct the air flow through the perforation fields in the bottom and in the inner lid under any circumstances (applicable to all container versions).**
- **Do not put foil packaging directly on the sterile container.**

Note

VARIO containers and BASIS containers can be sterilized with the outer lid in place!

- Follow the sterilizer manufacturer's recommendations.
- Always place heavy sterile containers at the bottom of the sterilizer.
- Always carry the sterile containers by their handles.

Note

Sterile containers can be stacked in the sterilizer.

- Transport stacks of sterile containers in such a way that the stacks do not topple over.

Sterilization method and parameters

CAUTION

Risk of sterilization failure!

- **Sterilize the containers only by approved and validated sterilizing processes.**
- Sterilize with steam, taking note of the following: Sterilization must be carried out through a validated steam sterilizing process (e.g. in a steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665).

Sterile container system

Unloading the sterilizer

WARNING

Risk of burns due to a hot sterile container after sterilization!

- Always wear protective gloves when unloading the sterilizer.

Transporting the sterile container

CAUTION

Risk of non-sterility of container contents!

- Always carry the sterile containers by their handles.
- Never carry or lift the sterile container at the lid.
- Transport the sterile container in such a way that mechanical damage will not occur.

Storing the sterile containers

Note

The sterile containers may be stored in stacks.

- Store sterile containers in a dry, clean and protected place.
- Observe storage periods and storage conditions acc. to DIN 58953-9.

Checking and commissioning the sterile materials

The contents of a sterile container can be commissioned as sterile only if the container has been properly sterilized, stored and transported.

If this is not the case, the sterile materials must be processed again.

DANGER

Risk of contamination from improperly sterilized materials!

- Prior to commissioning the sterile materials, check to ensure that the sterilization was successful.

- Make certain that the color of the indicator point has changed.
- Make certain that all container components, particularly the lid locks, are intact.
- Ensure the container seal is intact.

5. Validated processing procedure

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, the manufacturer's specifications can only serve as an approximate guide for assessing the processing procedures applied by the individual operator/processors.

5.1 General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. To avoid this, the time interval between application and processing should not exceed 6 h, and neither fixating pre-cleaning temperatures >45 °C nor any fixating disinfecting agents (active ingredients: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only process chemicals that have been tested and approved (e.g. VAH/ DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- optical changes of materials, e.g. fading or discoloration of aluminum. For aluminum, pH >8 in the application/process solution is sufficient to cause visible surface changes.
 - material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use process chemicals that cause stress cracking of plastics such as PPSU, or attack softeners e.g. in silicone and cause brittleness of the material.

5.2 Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions. The cleaning and disinfecting agent must
 - be approved for plastics and high-grade steel
 - not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.

Sterile container system

5.3 Manual cleaning/disinfecting

- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning and wipe disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical	
I	Cleaning	RT (cold)	-	-	D-W	-	
II	Drying	RT	-	-	-	-	
III	Wipe disinfection	-	>1	-	-	I Alcohol Denat. 70 % (ethanol B BRAUN)	II Aldehyde-free surface disinfectant (Stabimed)
IV	Final rinse	RT (cold)	-	-	FD-W	not required	Rinse off cleaning chemicals; do not leave residues
V	Drying	RT	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Wipe off gross debris prior to cleaning.

Stage II

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

Stage III

- Carry out wipe disinfection.

Stage IV

- When using aldehyde-free disinfectants:
After the specified exposure time (at least 1 minute), rinse the disinfected surfaces under running FD water.
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with lint-free tissue or in a heating cabinet.

5.4 Mechanical neutral or mild-alkaline cleaning and thermal disinfection

Machine type: Single-chamber washer/disinfector without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Neutral: BBRAUN HELIMATIC CLEANER neutral, pH-neutral, process solution 0.5 % Mild-alkaline: - Concentrate: pH = 9.5 < 5 % anionic tensides - 0.5 % solution Alkaline: Processing possible at pH up to 10.5, provided the cleaning agent is manufacturer-approved for container cleaning. Note Application of mild-alkaline or alkaline cleaning agents may cause discoloration of colored anodized aluminum, which however will not affect the functionality of the product.
III	Intermediate rinse	>10/50	1	FD-W	Especially for PrimeLine it must be ensured that the surface is rinsed without leaving any residues.
IV	Thermal disinfecting	90/194	5	FD-W	Other process parameters may be feasible with agreement by the hospital hygienist.
V	Drying	-	-	-	According to mechanical cleaning program

D-W: Drinking water

FD-W: Fully desalinated (demineralized) water

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Note

Ultrasonic treatment can be applied in manual processing as well as in mechanical processing. Temperatures of up to 120 °C are permissible for mechanical drying with hot air.

5.5 Inspection, maintenance and checks

Note

The sterile containers may be tested and repaired only by persons with the appropriate training, expertise or experience.

The seals in the container lids were tested according to DIN EN 868-8, Annex G, for a minimum service life of 5 000 sterilization cycles.

- In case of visible damage, change the seal immediately.
- Have the lid sent to Aesculap, see Technical Service.
- If necessary, lightly lubricate moving metal parts, e.g. latch hinges, with sterilizable and steam-permeable maintenance oil (e.g. STERILIT® oil spray JG600 or STERILIT® I drip feed oiler JG598).

Sterile container system

6. Troubleshooting list

Malfunction	Cause	Remedy
Excessive amounts of condensate inside the sterile container	Temperature of sterile materials too low prior to sterilization	Allow sterile materials to come to room temperature (appr. 20 °C) before sterilization
	Textiles too damp	Sterilize dry textiles only
	Sterile container too heavy	Normal container with instruments: max. load 10 kg with textiles: max. load 8 kg
		Half container with instruments: max. load 5 kg
		3/4-container with instruments: max. load 7 kg
	Sterile materials improperly positioned	Lay hollow components, bowls, plates etc. face down and on a slant
		Arrange textiles in loose vertical piles, do not press them together
	Sterile container incorrectly positioned in sterilizer	Always place heavy sterile containers at the bottom
	Sterile containers processed immediately after sterilization	Allow sterile containers to cool down to room temperature prior to processing
	Sterile containers improperly positioned during cooling phase	Do not store sterile containers on a floor or in a drafty place Store sterile containers in air-conditioned areas according to DIN 58953-9
	Sterilizer properties do not meet DIN EN 285 standards	Have sterilizer serviced regularly Check drying vacuum Check drying time
		Check steam quality and upgrade if necessary
	Empty-cycle and vacuum test not run daily before sterilization begins	Run empty-cycle and vacuum test daily before beginning sterilization
Condensate in groove of (inner) lid	Unsuitable sterilizer cycle selected	Select a program according to the load
	Sterilizer door left open too long, sterilizer cooled down	Load and unload sterilizer quickly
Indicator seal color unchanged	Multiple sterile containers or heavy sterile materials stacked	Subdivide sterilizer with shelving Use outer lid
	Sterilization performed incorrectly Sterilizer defective	Have sterilizer repaired by manufacturer
	Indicator seals stored improperly	Observe storage conditions as specified on the packaging of the indicator seals

Malfunction	Cause	Remedy
Basis/Vario sterile container or PrimeLine sterile container system deformed	Perforations covered during sterilization	Never cover perforations
Inner or outer lid cannot be positioned or locked on the bottom component	Container lid or bottom deformed or damaged due to improper handling (e.g. carrying by the lid)	Replace container lid or bottom, or have components repaired by the manufacturer

7. Technical Service

Risk of injury and/or malfunction!

➤ Do not modify the product.

WARNING

➤ For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

8. Accessories and replacement parts

Accessories and supplies are listed in brochure no. C40401.

9. Technical specifications

The variants and sizes of sterile containers are listed in brochure no. C40401.

10. Extracts from relevant standards

10.1 Standards cited

The following standards are cited in connection with the sterile containers:

- DIN 58953-9
- EN 285
- EN 868-8
- EN 285 in combination with EN ISO 17665
- EN ISO 11607

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Aesculap AG