

# Aesculap® T-Space® PEEK

Transforaminal Lumbar Interbody Fusion System



Aesculap Spine

New Generation  
T-Space PEEK Cages

# Aesculap® T-Space® PEEK



## Content

|                        |    |
|------------------------|----|
| A Foreword             | 3  |
| B Implant material     | 4  |
| C Implant features     | 5  |
| D Surgical technique   | 6  |
| E Ordering information |    |
| E1 Implants            | 16 |
| E2 Instruments         | 17 |

## Foreword

The high incidence of spinal disorders and consecutive symptoms calls for optimized diagnostics and therapies. Minimally invasive surgical procedures, being generally cost-efficient, significantly less invasive with fewer complications, and producing better outcomes, are of particular interest. Minimally invasive spine surgery relies on various retractor systems to create small ventral and dorsal access channels to the spine. At the same time innovative implants reducing tissue trauma and new percutaneous surgical techniques are coming up. Accordingly, the innovative S<sup>4</sup> Spinal System, the Spine Classics retractor system and the T-Space cage for intercorporeal fusion form an excellent treatment concept. In this way, minimally invasive mono- and bisegmental fusion surgeries at the lumbar spine can be successfully performed after a relatively short learning curve.

The new generation of T-Space PEEK cages is a further development of the 2006 launched and clinically proven Aesculap transforaminal interbody fusion system. The changes mainly regard the implant design – the new bulletted nose facilitates the implantation of the cage especially in very degenerated discs. Furthermore the connection of the implant with the inserter is improved by means of a screw thread. Besides, the system offers an extended size range presenting the right implant to fit the patient. The instrumentation stays reduced, clearly arranged and simple in handling. Moreover the surgeon can choose between three inserters as the preference or situation in-situ may be.

Spine Classics



S<sup>4</sup> Spinal System



# Aesculap® T-Space® PEEK

# B

## Implant material T-Space PEEK

The material used is biocompatible PEEK-Optima®, which was introduced by Invibio in 1999. PEEK stands for PolyEtherEtherKetone. PEEK-Optima® polymer comply with ISO10993-1, USP Class VI and ASTM F2026 for use as a medical implant material.

The use of PEEK-Optima® as an orthopedic device material enjoys increased popularity in recent years due to the material's unique combination of characteristics. It's properties include radiolucency, high mechanical strength, biocompatibility and compatibility with standard sterilization methods.

The intrinsic radiosopic transparency of the material on X-rays and CT scans makes it possible to view bone growth adjacent to the implant. This allows quick and simple assessment of the bone

structure and progress towards bone fusion. To verify the position of the PEEK implant on radiosopic images, non-radiolucent tantalum marker serve as location marker.


Of particular interest is the modulus of elasticity of PEEK-Optima® of 3.6 GPa, which is similar to that of cortical bone. This specific stiffness encourages load sharing between implant material and natural bone, thereby stimulating bone healing activity. The material provides excellent strength and rigidity. PEEK-Optima® also exhibits high fatigue resistance and low wear factor.

Extensive investigations into the biocompatibility of PEEK-Optima® have proven that the material is suitable for the use as a long-term implant.

Lateral view



AP view



## Implant features T-Space PEEK




### Position verification despite X-ray transparency



- PEEK-Optima® is radiolucent and allows therefore a quick and simple assessment of the bone structure and progress towards bone fusion
- Tantalum marker allow for easy and exact implant localization

### Intelligent implant design



- Bulleted nose for easier implantation especially in very degenerated discs
- Screw thread connection for safe and easy connection with the inserter
- 1 mm increments in height
- Optimized ratio between contact area and opening

### Thought-out instruments





- Straight or curved inserter as the preference or situation in-situ may be
- Offset inserter for comfortable use with microscope or C-arm

# Aesculap® T-Space® PEEK

Surgical technique

D









#### Positioning of the patient and incision marking (Fig. 1-2)

- A minimally invasive approach requires the patient to be placed on a radiolucent table which allows for AP views of the various anatomic structures.
- The appropriate position of the longitudinal incision (4–5 cm in length) is determined by using a C-arm. The intended skin incision is marked paraspinally on the right respectively on the left side.



# Aesculap® T-Space® PEEK

## Surgical technique



### Fascial incision (Fig. 3-4)

- A slightly arcuate fascial incision 1.5 cm from the midline is performed. This allows a firm hold of the speculum and counter retractor, facilitating the exposure of the individual segment.


### Exposure and blunt dissection of the paraspinal muscles (Fig. 5-6)

- After splitting of the thoracolumbar fascia a blunt dissection of the paraspinal muscles is performed with the fingertip. In accordance with the palpatory finding, a correction of the skin incision is still possible, as the muscle retractor should be introduced as vertically as possible and in the direction of the interlaminar space. The length of the retractor is selected by using the index finger.





- *Spine Classics retractor system*
- ▶ *See brochure 011402*




#### Introduction of the Spine Classics retractor system (Fig. 7-8)


- The muscle retractor is introduced with closed blades and with the handle in the longitudinal direction of the body. It is then turned 90° with the handle towards the assistant and afterwards expanded.

# Aesculap® T-Space® PEEK

## Surgical technique



- *S<sup>4</sup> Spinal System*
- ▶ *See surgical technique 026702*



- *Osteotome FJ658R*

### Insertion of S<sup>4</sup> screws (Fig. 9)

- Using the standard technique the S<sup>4</sup> Spinal System pedicle screws are inserted.

### Removal of facet joint (Fig. 10–11)

- A complete unilateral facetectomy should be considered on the side targeted for the implant insertion. The inferior articular process of the facet joint is resected first, then the subjacent superior articular process is resected.




Fig. 12




Fig. 13




Fig. 14

- *T-handle SJ033R*
- *Distractors FJ647R-FJ657R*

#### Opening of the disc and removal of disc material (Fig. 12-13)


- To open the disc a small window is cut into the annulus.
- Rongeurs are used to remove the opened annulus.
- Posterior osteophytes are removed by using Kerrisons.

#### Restoration of disc height (Fig. 14)

- The desired restoration of the natural disc height can be set using the distractors. They are available in heights from 7-17 mm in 1 mm increments.
- The distractor must be inserted horizontally and then rotated. Rotating clockwise the distractors are blunt. A special designed sharp rim allows removal of disc material. If so, the distractor has to be rotated counterclockwise.

# Aesculap® T-Space® PEEK

## Surgical technique




- Bone curettes, angled FJ679R-FJ680R or FJ698R-FJ699R
- Box curette, straight FJ681R
- Box curettes, angled FJ682R-FJ683R or FJ702R-FJ703R
- Bone rasps, angled FJ685R-FJ686R or FJ704R-FJ705R

### Cleaning of the intervertebral space (Fig. 15)

- The disc space is cleared using rongeurs, bone curettes and box curettes.
- The bone rasps are used to refresh the cartilaginous endplates. Alternatively the box curettes can be used.

Excessive preparation of the endplates may weaken the construct and cause subsidence of the interbody device.



- T-handle SJ033R
- Trials FJ667R-FJ677R or FJ619R-FJ629R

### Determination of implant size using trial implants (Fig. 16)

- The trial implants are available in heights from 7-17 mm in 1 mm increments.
- Using the T-handle the desired trial implant is inserted.

The trials are essential to ensure the correct implant size to be used.




Fig. 17a



Fig. 17b

- Packing block FJ615R
- Inserter FJ604R/FJ605R/FJ603R


### Implant insertion (Fig. 17a-b)

- The T-Space PEEK implant should be filled with bone or bone substitute by using the packing block.
- It is recommended to place bone graft in the anterior part.
- Partially the T-Space PEEK implant is inserted into the disc space. Three inserter options are available: a straight and a curved one as well as an offset inserter. The implant is connected with the inserter by means of a screw thread.




# Aesculap® T-Space® PEEK

## Surgical technique



■ *Impactor FJ613R*



### Final implant positioning (Fig. 18a-b and Fig. 19)

- Using the impactor the implant is rotated 90° to achieve the final positioning.
- X-ray control to verify the implant positioning.
- It is recommended to put bone material harvested from the facet joint around the T-Space implant.






Fig. 20




Fig. 21

- *S<sup>4</sup> Spinal System*
- ▶ *See surgical technique 026702*

#### Application of rod and set screw (Fig. 20)

- Final assembly of the S<sup>4</sup> Spinal System.
- Compression is applied to the pedicle screws to support the contact area between the T-Space implant and the endplates.
- Final tightening of the S<sup>4</sup> pedicle screws and removal of the tabs.

#### S<sup>4</sup> screw positioning on the contra-lateral side (Fig. 21)

- S<sup>4</sup> Spinal System is applied on the contra-lateral side.



# Aesculap® T-Space® PEEK

E<sub>1</sub>

## Ordering information – T-Space PEEK implants









| Art. no. | Description  | Height | Width   | Length | Angle |
|----------|--------------|--------|---------|--------|-------|
| SJ907P   | T-Space PEEK | 7 mm   | 11.5 mm | 26 mm  | 5°    |
| SJ908P   | T-Space PEEK | 8 mm   | 11.5 mm | 26 mm  | 5°    |
| SJ909P   | T-Space PEEK | 9 mm   | 11.5 mm | 26 mm  | 5°    |
| SJ910P   | T-Space PEEK | 10 mm  | 11.5 mm | 26 mm  | 5°    |
| SJ911P   | T-Space PEEK | 11 mm  | 11.5 mm | 26 mm  | 5°    |
| SJ912P   | T-Space PEEK | 12 mm  | 11.5 mm | 26 mm  | 5°    |
| SJ913P   | T-Space PEEK | 13 mm  | 11.5 mm | 26 mm  | 5°    |
| SJ915P   | T-Space PEEK | 15 mm  | 11.5 mm | 26 mm  | 5°    |
| SJ917P   | T-Space PEEK | 17 mm  | 11.5 mm | 26 mm  | 5°    |
| SJ937P   | T-Space PEEK | 7 mm   | 11.5 mm | 30 mm  | 5°    |
| SJ938P   | T-Space PEEK | 8 mm   | 11.5 mm | 30 mm  | 5°    |
| SJ939P   | T-Space PEEK | 9 mm   | 11.5 mm | 30 mm  | 5°    |
| SJ940P   | T-Space PEEK | 10 mm  | 11.5 mm | 30 mm  | 5°    |
| SJ941P   | T-Space PEEK | 11 mm  | 11.5 mm | 30 mm  | 5°    |
| SJ942P   | T-Space PEEK | 12 mm  | 11.5 mm | 30 mm  | 5°    |
| SJ943P   | T-Space PEEK | 13 mm  | 11.5 mm | 30 mm  | 5°    |
| SJ945P   | T-Space PEEK | 15 mm  | 11.5 mm | 30 mm  | 5°    |
| SJ947P   | T-Space PEEK | 17 mm  | 11.5 mm | 30 mm  | 5°    |
| SJ967P   | T-Space PEEK | 7 mm   | 11.5 mm | 34 mm  | 5°    |
| SJ968P   | T-Space PEEK | 8 mm   | 11.5 mm | 34 mm  | 5°    |
| SJ969P   | T-Space PEEK | 9 mm   | 11.5 mm | 34 mm  | 5°    |
| SJ970P   | T-Space PEEK | 10 mm  | 11.5 mm | 34 mm  | 5°    |
| SJ971P   | T-Space PEEK | 11 mm  | 11.5 mm | 34 mm  | 5°    |
| SJ972P   | T-Space PEEK | 12 mm  | 11.5 mm | 34 mm  | 5°    |
| SJ973P   | T-Space PEEK | 13 mm  | 11.5 mm | 34 mm  | 5°    |
| SJ975P   | T-Space PEEK | 15 mm  | 11.5 mm | 34 mm  | 5°    |
| SJ977P   | T-Space PEEK | 17 mm  | 11.5 mm | 34 mm  | 5°    |

## Ordering information – Preparation instruments



**FJ610 T-Space PEEK instrumentation complete  
(FJ633R and FJ611R)**




consisting of:

|   | Art. no. | Description                      | Recommended | Optional |
|---|----------|----------------------------------|-------------|----------|
|   | FJ679R   | Left angled bone curette, 45°    |             | 1        |
|   | FJ680R   | Right angled bone curette, 45°   |             | 1        |
|  | FJ698R   | Left angled bone curette, 20°    | 1           |          |
|   | FJ699R   | Right angled bone curette, 20°   | 1           |          |
|  | FJ681R   | Straight curette                 | 1           |          |
|   | FJ682R   | Left angled curette, 45°         |             | 1        |
|  | FJ683R   | Right angled curette, 45°        |             | 1        |
|   | FJ702R   | Left angled curette, 20°         | 1           |          |
|  | FJ703R   | Right angled curette, 20°        | 1           |          |
|   | FJ658R   | Straight osteotome, 8 mm         | 1           |          |
|  | FJ685R   | Left angled bone rasp, 45°       |             | 1        |
|   | FJ686R   | Right angled bone rasp, 45°      |             | 1        |
|   | FJ704R   | Left angled bone rasp, 20°       | 1           |          |
|   | FJ705R   | Right angled bone rasp, 20°      | 1           |          |
|   | FJ633R   | Tray for preparation instruments | 1           |          |
|   | JH217R   | Wide perforated basket lid       | 1           |          |
|   | TE989    | Graphic template for FJ633R      | 1           |          |

# Aesculap® T-Space® PEEK


## Ordering Information – implantation instruments



|   | Art. no. | Description                 | Recommended | Optional |
|---|----------|-----------------------------|-------------|----------|
|  | FJ647R   | Distractor, 7 mm            | 1           |          |
|   | FJ648R   | Distractor, 8 mm            | 1           |          |
|   | FJ649R   | Distractor, 9 mm            | 1           |          |
|   | FJ650R   | Distractor, 10 mm           | 1           |          |
|   | FJ651R   | Distractor, 11 mm           | 1           |          |
|   | FJ652R   | Distractor, 12 mm           | 1           |          |
|   | FJ653R   | Distractor, 13 mm           | 1           |          |
|   | FJ655R   | Distractor, 15 mm           | 1           |          |
|   | FJ657R   | Distractor, 17 mm           | 1           |          |
|  | FJ667R   | T-Space trial, 7 mm         | 1           |          |
|   | FJ668R   | T-Space trial, 8 mm         | 1           |          |
|   | FJ669R   | T-Space trial, 9 mm         | 1           |          |
|   | FJ670R   | T-Space trial, 10 mm        | 1           |          |
|   | FJ671R   | T-Space trial, 11 mm        | 1           |          |
|   | FJ672R   | T-Space trial, 12 mm        | 1           |          |
|   | FJ673R   | T-Space trial, 13 mm        | 1           |          |
|   | FJ675R   | T-Space trial, 15 mm        | 1           |          |
|   | FJ677R   | T-Space trial, 17 mm        | 1           |          |
|  | FJ619R   | T-Space curved trial, 7 mm  |             | 1        |
|   | FJ620R   | T-Space curved trial, 8 mm  |             | 1        |
|   | FJ621R   | T-Space curved trial, 9 mm  |             | 1        |
|   | FJ622R   | T-Space curved trial, 10 mm |             | 1        |
|   | FJ623R   | T-Space curved trial, 11 mm |             | 1        |
|   | FJ624R   | T-Space curved trial, 12 mm |             | 1        |
|   | FJ625R   | T-Space curved trial, 13 mm |             | 1        |
|   | FJ627R   | T-Space curved trial, 15 mm |             | 1        |
|   | FJ629R   | T-Space curved trial, 17 mm |             | 1        |



| Art. no. | Description                         | Recommended | Optional |
|----------|-------------------------------------|-------------|----------|
| SJ033R   | T-handle for distractors and trials | 2           |          |
| FJ051R   | Retractor S                         | 1           |          |
| FJ052R   | Retractor M                         | 1           |          |
| FJ053R   | Retractor L                         | 1           |          |
| FJ054R   | Retractor XL                        | 1           |          |
| FJ615R   | Packing block                       | 1           |          |
| FJ603R   | Offset inserter                     | 1           |          |
| FJ604R   | Straight inserter                   | 1           |          |
| FJ605R   | Curved inserter                     | 1           |          |
| FJ613R   | Impactor                            | 1           |          |
| FJ611R   | Tray for implantation instruments   | 1           |          |
| JH217R   | Wide perforated basket lid          | 1           |          |
| TF023    | Graphic template for FJ611R         | 1           |          |



Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 7461 95-0 | Fax +49 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

The main product trademark 'Aesculap' and the product trademarks 'T-Space' and 'S+' are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. 042502

1211/1/1