- GB** Instructions for use/Technical description
- USA** Targon® femur and tibia nail system
- D** Gebrauchsanweisung/Technische Beschreibung
- Targon®-Femur- und -Tibia-Nagelsystem
- F** Mode d'emploi/Description technique
- Système de clous Targon® pour fémur et tibia
- E** Instrucciones de manejo/Descripción técnica
- Sistema de clavos para fémur y tibia Targon®
- I** Istruzioni per l'uso/Descrizione tecnica
- Sistema di chiodi femorali e tibiali Targon®
- P** Instruções de utilização/Descrição técnica
- Sistema de pregos femorais e tibiais Targon®
- NL** Gebruiksaanwijzing/Technische beschrijving
- Targon®-femur- en -tibia-nagelsysteem
- S** Bruksanvisning/Teknisk beskrivning
- Targon® lår- och skenbensstiftsystem
- RUS** Инструкция по применению/Техническое описание
- Система гвоздей Targon® для бедренной и большеберцовой кости
- CZ** Návod k použití/Technický popis
- Systém femorálních a tibialních hřebů Targon®
- PL** Instrukcja użytkowania/Opis techniczny
- System gwoździ udowych i piszczelowych Targon®
- SK** Návod na používanie/Technický opis
- Femorálny a tibialný klincový systém Targon®
- TR** Kullanım Kılavuzu/Teknik açıklama
- Targon® femur ve tibia çivi sistemi

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 010451 01/14 V6 Änd.-Nr. 48295

Technical alterations reserved

Aesculap®

Targon® femur and tibia nail system

Intended use

The Targon® femur and tibia nail system is used as a splint for fractures of long tubular bones, the femur and the tibia.

System configurations

- Targon® universal femur nail: various lengths and diameters
- Targon® universal tibia nail: various lengths and diameters
- Targon® solid titanium tibia nail: various lengths and diameters
- Targon® solid titanium femur nail: various lengths and diameters
- For further information regarding diameters and lengths, see operating manual O113 and brochure O111.
- Targon® retrograde femur nail: various lengths and diameters
- For further information regarding diameters and lengths, see brochure O190.

Associated implants:

- Fixation screws in various diameters and lengths
- Threaded sleeves for fixation screws
- Fixation screws for the end of the nail on the connecting side

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_F Titanium forged alloy Ti6Al4V acc. to ISO5832-3
- Implant steel acc. to ISO 5832-1

The titanium implants are coated with a colored oxide. Their color may change slightly, but this has no bearing on the quality of the implant.

ISOTAN®_F is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use Targon® universal femur nails, Targon® universal tibia nails, Targon® solid titanium tibia nails and Targon® solid titanium femur nails for the following indications:

- Open or closed fractures of the femur shaft or tibia shaft
- Repositioning osteotomies in the shaft region
- Reconstruction following tumor resection in the shaft region
- Osteotomies in the shaft region to lengthen or shorten the limb
- Pathological fractures
- Pseudarthroses
- Knee joint arthrodesis (Targon® F special length)

Use Targon® retrograde femur nails for the following indications:

- Femur shaft fractures
- Extra- and intra-articular fractures of the distal femur (these fractures also with ipsilateral femur shaft fractures)
- Patients with ipsilateral knee prosthesis, hip prosthesis or other proximal implants, e.g. dynamic hip screw
- Shaft fractures necessitating additional proximal implants, e.g. for ipsilateral femoral neck fractures
- Shaft fracture and additional nailing of the ipsilateral tibia (identical, extended approach incision)
- Pathological fractures
- Pseudarthroses
- Ankle arthrodesis

The operating surgeon is responsible for establishing the above indications and indications beyond those, taking into account the specific clinical, biological, and biomechanical situation.

Contraindications

Do not use Targon® universal femur nails, Targon® universal tibia nails, Targon® solid titanium tibia nails and Targon® solid titanium femur nails for the following:

- Open or closed fractures in joint regions
- Trochanteric fractures of the femur

Do not use Targon® retrograde femur nails for the following:

- Septic knee arthritis

Generally, do not use in cases of:

- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Anticipated overloading of the implant
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Side effects and interactions

- Dislocation, loosening and fracture of implant components
- Delayed healing, or non-union and development of pseudarthrosis
- Primary and secondary infections
- Venous thrombosis, lung embolism, cardiac arrest
- Tissue reaction to implant materials
- Trauma to nerves, tendons, or vessels
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Compartment syndrome
- Pain in the region of the nail insertion point and in the region of locking components

Safety notes

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

Implant steel

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- In case of long steel implants, MRI can cause critical, localized heat development.
- Implants produce clear MRI artifacts.

Titanium/titanium alloy

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- Implants produce moderate MRI artifacts.

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Mobilization, dynamic locking or use of compression involve the risk of migration of the nail and/or the locking component into the joint region or out of the bone. Therefore the migration paths of the implants and those of the fragments must be considered in the implantation and taken into account when selecting the implant length.
- Considerably higher strain due to increased leverage is to be expected if there are locking holes at the height or in the region of the fracture line. This can lead to implant failure. Accordingly, postoperative strain has to be reduced, or full load may be envisaged only after complete callus formation.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudarthrosis develops, or if the implant is overloaded for an excessively long period. One way to achieve this is through mobilization.
- Always insert the Targon® implants in such a way that forces to be transmitted are slight and can be transmitted by the bone at an early stage.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Do not apply strong force when implanting or explanting the implants. Should a problem arise, check the position and orientation of the implants, instruments, and fragments and analyze the origin of the problem. If necessary, repeat previous operating steps and check the instruments (e.g. for blocked drill helix).
- Targon® implants may only be implanted using the Aesculap implantation instruments specially designed for this purpose.
- Use implant components only in the nail holes intended for this purpose.
- Only use cutting instruments, e.g. drill bits, guide pins, etc., that are sharp.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- The attending physician shall make any decision with regard to the removal of implant components that have been used,
- Increase the load on the fracture and implant in accordance with the progress of healing.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- At the time of follow-up exams, the scope of postoperative load and aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. Also, the dimensions of the implant should be taken into account.
- Recommended follow-up intervals
 - prior to the patient's discharge from hospital
 - 10-12 weeks post operationem
 - 6 months post operationem
 - 12 months post operationem

For implanting the nail via a guide pin:

- Use the appropriate guide pin
- Avoid jamming
- Monitor the tip of the guide pin under X-ray control
- Remove the guide pin prior to locking

Aesculap and the Aesculap Academy regularly hold general and product-specific training events, offering users helpful application advice, especially for complicated or critical situations.

Sterility

- The implant components are supplied in an unsterile condition.
 - The implant components are packaged individually.
 - Store the implant components in their original packaging and only remove them from their original and protective packaging immediately prior to application.
 - Use the implant system storage devices for processing, sterilization and sterile setup.
 - Use a suitable tray for cleaning/disinfection.
 - Use the system storage device only for sterilization and sterile provision.
 - Ensure that the implant components in their implant system storage devices do not come into contact with each other or with instruments.
 - Ensure that the implant components are not damaged in any way.
- Prior to initial sterilization and subsequent resterilization, the implant components must be cleaned using the following validated reprocessing procedure:

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

The recommended chemistry was used for validation.

Note

Sterilization should always be favored over disinfection under all circumstances. Only if final sterilization is not possible should disinfection with a virucidal agent be preferred.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com.

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Validated reprocessing procedure

Mechanical alkaline cleaning and thermal disinfection

- ▶ Process the implant in its system storage device.
 - ▶ Place the implants on a tray that is suitable for cleaning (avoiding rinsing blind spots and damage).
- Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemistry/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH ~ 13 - <5% anionic surfactant ■ 0.5% working solution - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Dry	-	-	-	In accordance with the program for the cleaning and disinfecting machine

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

For implant components that are to be resterilized:

WARNING

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- ▶ Handle the implants with new gloves only.
- ▶ Keep the implant system storage devices covered or closed.
- ▶ Process implant system storage devices separately from instrument trays.
- ▶ Clean implants must not be processed together with contaminated implants.
- ▶ Process the implant components individually and separately if no implant system storage devices are available, ensuring that the implant components are not damaged in the process.
- ▶ Mechanically clean and disinfect the implant components.
- ▶ Do not reuse surgically contaminated implants!

WARNING

Direct or indirect contamination may render implants unsuitable for resterilization!

- ▶ Do not reprocess implants that have been directly or indirectly contaminated with blood.

Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Immediately sort out damaged or inoperative products.

Packaging

- ▶ Place the product in its holder or on a suitable tray.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Make sure that the packaging will prevent a recontamination of the product whilst in storage.

Sterilization

- ▶ Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Application

WARNING

Risk of injury caused by incorrect operation of the product!

- ▶ Attend appropriate product training before using the product.
- ▶ For information about product training, please contact your national B. Braun/ Aesculap agency.

A detailed operation manual for the Targon® nail system is available upon request from Aesculap or from your local Aesculap representative.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- Information has been obtained from the manufacturer in the event of an ambiguous preoperative situation and if implants are present in the area to be treated

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- Under certain circumstances it can happen that, when a fracture is treated with an implant, the original anatomy of the extremity will not be fully restored.
- Following treatment of a fracture, shoulder and elbow joint function may be impaired.
- Pains may occur following the fracture treatment.
- The implants must not be exposed to overload due to excessive strains, hard physical work or sports activities. Overtraining entails the risk of implant loosening or breakage.
- The patient must be informed about the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The risks of transgressing these rules must be explained to the patient.
- Revision surgery may become necessary in cases of implant loosening or pseudarthrosis.
- The patient must undergo regular medical follow-up examinations.
- The attending physician decides when to remove, in part or completely, the osteosynthesis implants through another surgical procedure.

WARNING

Complications due to incorrect use and/or incorrect assessment/treatment of the fracture!

- ▶ Only apply nails of appropriate length and diameter.
- ▶ Make certain that implants are aligned correctly.

The implantation and implantation site are prepared in the following way:

The application steps for Targon® implant components are described in detail in the operating manuals. The following notes must be observed for all Targon® implantations:

Check for correct mounting of the nail on targeting device: Insert the drill through the drill sleeve and the targeting hole and guide it through the nail hole.

WARNING

Improper use can cause bones to shatter and nails to jam or bend!

- ▶ Do not use excessive force, under any circumstances, when handling/applying the instruments and implants.
- ▶ Only apply mild hammer taps to Targon® F/T nails and, if necessary, RF nails to insert them into the medullary cavity.
- ▶ Check the nail insertion site, and correct if necessary.
- ▶ Should a nail be difficult to insert, drill open the medullary cavity to a smaller (for solid nails) or larger (for drilled nailing) extent.

The operating surgeon determines the following, depending upon the indications:

- if implant locking is necessary
- if the implant is to be locked, whether it will be locked statically or dynamically.

Note

For locking an implant without a targeting attachment, Aesculap recommends free-hand locking supported by a C-bow and laser pilot beam, or with a radiolucent angle transmission

WARNING

Avoid drill mistargeting, drill breakage, guide pin breakage and premature implant failure due to crack formation as a result of damaged surfaces in the drill hole regions!

- ▶ Use a drill bit of the correct diameter.
- ▶ Level the drill insertion site with a flat counterbore or use a center punch.
- ▶ Only use sharp drill bits and guide pins and drill carefully and without applying excessive force. Avoid changes of the drilling direction.
- ▶ During center punching, flat countersinking and drilling, do not apply force to the targeting device; do not bend or twist the targeting device.

Implant removal requires the following steps:

- ▶ Prior to removing the nail, remove all screws.
- ▶ Remove the nail with quick, light taps, and with the aid of the impactor, extraction pin and slotted hammer.

Note

Implant removal can give rise to complications due to implants that are stuck, ingrown bone tissue, etc. Such complications can lead to damage to implants and/or instruments. For such cases we recommend using a set of special instruments, which can be ordered from Aesculap. This instrument set should be at hand for every explantation. Other special instruments are available for dealing with broken implants. Note the instructions for use for such case!

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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