Targon® FN

GB	Instructions for use Targon® FN
D	Gebrauchsanweisung Targon® FN
F	Mode d'emploi Targon® FN
E	Instrucciones de manejo Targon® FN
I	Istruzioni per l'uso Targon® FN
P	Instruções de utilização Targon® FN
RUS	Инструкция по применению Targon® FN
CZ	Návod k použití Targon® FN
PL	Instrukcja użytkowania Targon® FN
SK	Návod na používanie Targon® FN
CN	使用说明书 带锁髓内钉系统（植入物）

B | BRAUN

SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 012010 05/11 V5 Änd.-Nr. 43219

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Simbolo CE, em conformidade com a Directiva 93/42/CEE
Маркировка CE согласно директиве 93/42/ЕЭС
CE-označení podle směrnice 93/42/EHS
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
Označenie CE podľa smernice 93/42/EHS
CE 标记符合 93/42/EEC 指令

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Мы оставляем за собой право вносить технические изменения
Technické zmeny vyhradeny
Zmiany techniczne zastrzeżone
Technické zmeny vyhradené
保留技术更改的权利

Intended use

The Targon® FN system is used for splinting, stabilization, and fixation of fractures of the proximal femur.

System variants:

- Targon® FN femoral plate
- Targon® FN Telescrew support screw system in various lengths

Materials

The materials used in the implant are listed on the packaging.

- ISOTAN®_F titanium forged alloy Ti6Al4V according to ISO 5832-3

The titanium implants are coated with a colored oxide. Their color may change slightly, but this has no bearing on the quality of the implant.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for:

- Medial femoral neck fractures

The operating surgeon is responsible for establishing the above indications and indications beyond those, taking into account the specific clinical, biological, and biomechanical situation.

Contraindications

Do not use in the presence of:

- Trochanteric femur fractures
- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Anticipated overloading of the implant
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Side effects and adverse interactions

- Changes in position, loosening, wear and tear or fracture of implant components
- Position change and loosening of fragments
- Delayed healing, or non-union and development of pseudarthrosis
- Primary and secondary infections
- Venous thrombosis, lung embolism, cardiac arrest
- Tissue reaction to implant materials
- Trauma to nerves, tendons, or vessels
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Pain in the implant entry region and in the region of locking components

Safety notes

- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- Always observe the instructions for use of the individual Aesculap implant components.
- Follow the operating manual.
- Use implant components only in the appropriate plate drill holes.
- The operating surgeon must be thoroughly familiar with the bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Prior to the operation, the surgeon must make themselves familiar with the Targon® implants and the operating technique to be used.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- The Targon® implants may only be implanted with the Aesculap implantation instruments intended for this purpose.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Implants that have been used before must not be reused.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Considerably higher strain due to increased leverage is to be expected if there are locking holes at the height or in the region of the fracture line. This can lead to implant failure. Accordingly, postoperative strain has to be reduced, or full load may be envisaged only after complete callus formation.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudarthrosis develops, or if the implant is overloaded for an excessively long period. One way to achieve this is through mobilization.
- Always insert the Targon® implants in such a way that forces to be transmitted are slight and that any such forces will be borne by the bone at an early stage.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- Never use damaged implant components or implant components that have been removed surgically.
- Do not apply strong force when implanting or explanting the implants. Should a problem arise, check the position and orientation of the implants, instruments, and fragments and analyze the origin of the problem. If necessary, repeat previous operating steps and check the instruments (e.g. for blocked drill helix).

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.**
- MRI causes non-critical, localised heat development.**
- Implants produce moderate MRI artifacts.**

- Increase the load on the fracture and implant in accordance with the progress of healing.
- The timing of follow-up examinations, the type of postoperative load and the aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. Also, the dimensions of the implant should be taken into account.

- Recommended follow-up intervals:
 - prior to the patient's discharge from hospital
 - 10-12 weeks post operationem
 - 6 months post operationem
 - 12 months post operationem

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.

Risk of contamination of the patient!

- **Never resterilize or reuse implants contaminated e.g. with blood or secretions.**

- Ensure that the implant components in their implant system storage devices do not come into contact with each other or with instruments.

- Ensure that the implant components are not damaged in any way.

Prior to resterilization, the implant components must be cleaned using the following validated reprocessing procedure:

Note

National laws, national and international standards and directives, and product-specific hygiene regulations for reprocessing must be observed.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/sterile processing technician is responsible for this.

Note

Up-to-date information on reprocessing can be found on the Aesculap Extranet at www.aesculap-extra.net

Validated reprocessing procedure

Mechanical alkaline cleaning and thermal disinfection

- Place the implants on a tray that is suitable for cleaning (avoiding rinsing blind spots).

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C°/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER alkaline with tensides, application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD-W	–
IV	Thermal disinfection	90/194	5	FD-W	–
V	Drying	–	–	–	According to disinfectant program

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- **Handle the implants with new gloves only.**
- **Keep the implant system storage devices covered or closed.**
- **Process implant system storage devices separately from instrument trays.**
- **Clean implants must not be processed together with contaminated implants.**
- **Process the implant components individually and separately if no implant system storage devices are available, ensuring that the implant components are not damaged in the process.**
- **Mechanically clean and disinfect the implant components.**
- **Do not reuse intraoperatively contaminated implants!**

Direct or indirect contamination may render implants unsuitable for resterilization!

- **Do not reprocess implants that have been directly or indirectly contaminated with blood.**

Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational and undamaged.
- Immediately sort out damaged or inoperative products.

Packaging

- Place the product in its holder or on a suitable tray.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization

- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134°C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- Implantation instruments, including the special Aesculapimplantsystem instruments, must be complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure has been explained to the patient and their agreement concerning the following information has been recorded:

- Treatment of a fracture with implants may not always be successful in completely restoring the original bone anatomy.
 - Following the fracture treatment, the functionality of adjacent joints may be restricted.
 - Pains may occur following the fracture treatment.
 - The implants must not be subjected to overload through extreme strains, hard physical work or exercise. The implant may loosen or suffer material fracture due to excessive strain e.g. suffered in a fall.
 - The patient must be informed about the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The implants will not sustain the load and functionality of a healthy bone. The risks of transgressing these rules must be explained to the patient.
 - Corrective surgery may become necessary in cases of implant loosening, fragment migration, pseudarthroses or other complications.
 - The patient must undergo regular medical follow-up examinations.
 - The attending physician decides the appropriate timing of the surgical removal, in part or completely, of the osteosynthesis implants.
- Check for correct mounting of the plate at the aiming device: Insert the drill through the drill sleeve and the aiming hole and guide it through the plate hole.

WARNING

Complications due to incorrect use and/or incorrect assessment/treatment of the fracture!

- Use support screws of the appropriate length.
- Make certain that implants are aligned correctly.
- Handle the implants with appropriate care.

Explanation

Note
Implant removal can give rise to complications due to implants that are stuck, ingrown bone tissue, etc. Such complications can lead to damage to implants and/or instruments. For such cases a set of special instruments is recommended, which can be ordered from Aesculap. This instrument set should be at hand for every explantation. Other special instruments are available for dealing with broken implants. Note the instructions for use for such case!

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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KA351S, KA354S, KA356S, KA358S, KA360S, KA362S, KA364S, KA366S, KA368S, KA370S, KA372S, KA374S, KA458S, KA460S, KA462S, KA464S, KA466S, KA468S, KA470S, KA472S, KA474S, KA500S, KA558S, KA560S, KA562S, KA564S, KA566S, KA568S, KA570S, KA572S, KA574S, KA576S, KA600S, KA662S, KA664S, KA666S, KA668S, KA670S, KA672S, KA674S, KA676S, KA764S, KA766S, KA768S, KA770S, KA772S, KA774S, KA776S, KA864S, KA866S, KA868S, KA870S, KA872S, KA874S, KA876S	TARGON 股骨通用髓内钉
KA351SS, KA354SS, KA356SS, KA358SS, KA360SS, KA362SS, KA364SS, KA366SS, KA368SS, KA370SS, KA372SS, KA374SS, KA458SS, KA460SS, KA462SS, KA464SS, KA466SS, KA468SS, KA470SS, KA472SS, KA474SS, KA562SS, KA664SS, KA666SS, KA668SS, KA670SS, KA672SS, KA674SS, KA676SS, KA764SS, KA766SS, KA768SS, KA770SS, KA772SS, KA774SS, KA776SS	TARGON 股骨通用髓内钉 消毒包装
KA558SS, KA560SS, KA562SS, KA564SS, KA566SS, KA568SS, KA570SS, KA572SS, KA574SS, KA576SS,	TARGON 股骨通用髓内钉消毒包装
KB336TS	TARGON FN 锁定螺钉 Ø 4.5 x 36 MM
KB340T	TARGON 锁定螺钉钛 4.5/40MM
KB340TS	TARGON FN 锁定螺钉 Ø 4.5 x 40 MM
KB344T	TARGON 锁定螺钉钛 4.5 /44MM
KB344TS	TARGON FN 锁定螺钉 Ø 4.5 x 44 MM
KB348T	TARGON 锁定螺钉钛 4.5/48MM
KB348TS	TARGON FN 锁定螺钉 Ø 4.5 x 48 MM
KB352T	TARGON 锁定螺钉钛 4.5/52MM
KB352TS	TARGON FN 锁定螺钉 Ø 4.5 x 52 MM
KB356T	TARGON 锁定螺钉钛 4.5/56MM
KB356TS	TARGON FN 锁定螺钉 Ø 4.5 x 56 MM
KB360T	TARGON 锁定螺钉钛 4.5/60MM
KB360TS	TARGON FN 锁定螺钉 Ø 4.5 x 60 MM

TA-Nr.: 012010 05/11 Änd.-Nr.: 43219