- GB

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0123 – DIR 93/42/EEC

Intended use

The implant system is used for the intramedullary splinting, stabilization, and fixation of fractures in long tubular bones.

System configurations:

- Targon® PH/H
 - Proximal humerus nail
 - Humerus shaft nail
- Targon® PF
 - Proximal femur nail
- Targon® F
 - Femur nail
- Targon® T
 - Tibia nail
- Targon® RF
 - Retrograde femur nail

Each system consists of locking nails, locking components (e.g., screws, sleeves, pins) and associated special implantation instruments.

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®: titanium forged alloy Ti6Al4V acc. to ISO 5832-3
- Implant steel acc. to ISO 5832-1

The titanium implants are coated with a colored oxide. Slight changes in coloration may occur, but do not affect the implant quality.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

PEEK-OPTIMA® is a registered trademark of Invivio, Ltd Lancashire FY5 4QD / UK.

Indications

Use Targon® PH/H for:

- Multi-fragment fractures of the humerus head, up to 4-fragment fracture (Targon® PH)
- Combined humerus head and shaft fractures (Targon® PH long)
- Stable or unstable humeral shaft fractures (Targon® H)
- Pathological fractures
- Pseudarthroses

Use Targon® PF for:

- Pertrochanteric, intertrochanteric and subtrochanteric femur fractures
- Above-mentioned fractures combined with ipsilateral femoral shaft fractures (Targon® PF long)
- Lateral femoral neck fractures combined with ipsilateral femoral shaft fractures (Targon® PF long)
- Lateral femoral neck fractures extending into the trochanter zones
- Pathological fractures
- Pseudarthroses

Use Targon® F/T for:

- Open or closed fractures of the femur shaft or tibia shaft
- Reconstruction following tumor resection in the shaft region
- Lengthening, shortening and corrective osteotomies in the shaft region
- Pathological fractures
- Pseudarthroses

Use Targon® RF for:

- Extra- and intra-articular fractures of the distal femur (these fractures also with ipsilateral femur shaft fractures)
- Patients with ipsilateral knee prosthesis, hip prosthesis or other proximal implants, e.g. a dynamic hip screw
- Shaft fractures necessitating additional proximal implants, e.g. for ipsilateral femoral neck fractures
- Shaft fractures and additional nailing of the ipsilateral tibia (identical, extended approach incision)
- Pathological fractures
- Pseudarthroses

Absolute contraindications

Do not use Targon® PH/H for:

- Fractures up to 7 cm proximal of the olecranon
- Crushed head cap of the humerus head (C3 fractures)
- Isolated avulsion fracture of a tuberosity

Do not use Targon® PF for:

- Medial femoral neck fractures

Do not use Targon® F/T for:

- Open or closed fractures in joint regions
- Trochanteric fractures of the femur

Do not use Targon® RF for:

- Septic knee arthritis

Generally, do not use in cases of:

- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation

Relative contraindications

The following conditions, individual or combined, can lead to delayed healing or compromise the success of the operation:

- Anticipated overload on the implant
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Expected lack of patient cooperation
- Foreign body sensitivity to the implant materials

Use of the products requires individual, critical assessment by the surgeon.

Side effects and interactions

- Changes in position, loosening, wear and tear on, or fracture of implant components
- Position change and loosening of bone fragments
- Periimplant fractures
- Delayed healing, or non-union and development of pseudarthrosis
- Primary and secondary infections
- Thromboses, embolisms
- Tissue reaction to implant materials
- Injury to the surrounding tissue, including nerve and vascular damage
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Compartment syndrome
- Pain in the region of the nail insertion point and in the region of locking components

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- Use implant components only in the nail holes intended for this purpose.
- Prior to the operation, the operating surgeon must familiarize himself/herself with the Targon® implants and the associated operating technique.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- The Targon® implants may only be implanted with the Aesculap implantation instruments intended for this purpose.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants which have already been used must not be used again.
- Mobilization, dynamic locking or use of compression involve the risk of migration of the nail and/or the locking component into the joint region or out of the bone. Therefore the migration paths of the implants and those of the fragments must be considered in the implantation and taken into account when selecting the implant length.
- Considerably higher strain due to increased leverage is to be expected if there are locking holes at the height or in the region of the fracture line. This can lead to implant failure. Accordingly, postoperative strain has to be reduced, or full load may be envisaged only after complete bone formation.
- The functional lifespan of the system is limited to the duration until bone consolidation. Bone consolidation normally occurs within three months.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudarthrosis develops, or if the implant is overloaded for an excessively long period. One way to achieve this is through mobilization.
- There are currently no findings regarding the behavior of the implant after bone consolidation. Experiences with the implant show no uncharacteristic accumulation of complications for this case. After consolidation has occurred, leaving the implant in the body is associated with risks that the user must weigh against the risks of removal by operation. The indication for implant removal is relative, and rarely mandatory. It is always made on an individual basis. There is no evidence-based recommendation. Otherwise, refer to the corresponding applicable guidelines of the Deutsche Gesellschaft für Unfallchirurgie e.V. (DGU).
- Once the nail is locked, remove transmedullary bone and fragment correction screws again, or take into account, when determining the allowable patient load, that adverse lever effects can expose the nail to increased load.
- Always insert the Targon® implants in such a way that forces to be transmitted are low and can be absorbed from the bone at an early stage.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Always use cutting instruments, e.g., drill bits, guide pins, that are sharp.
- Do not apply strong force when implanting or explanting the implants. Should a problem arise, check the position and orientation of the implants, bone fragments and instruments and analyze the origin of the problem. If necessary, repeat previous operating steps and check the instruments (e.g., for a blocked drill helix).

The implant was not examined for safety and compatibility in MR environments. It has not been tested for heat, movements or image artifacts during MR examinations. MR-scanning of a patient with this implant can lead to injuries to the patient.

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Increase the load on the fracture and implant in accordance with the progress of healing.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- At the time of follow-up exams, the scope of postoperative load and aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. Also, the dimensions of the implant should be taken into account.
- Recommended follow-up intervals
 - prior to the patient's discharge from hospital
 - 10-12 weeks post operationem
 - 6 months post operationem
 - 12 months post operationem
- Implanting the nail via a guide pin
 - Use an appropriate guide pin
 - Avoid jamming the implant or guide pin
 - Monitor the guide pin end by X-ray control
 - Remove the guide pin before locking the implant

Sterility

- The implant components are supplied in unsterile condition.
- The implant components are packaged individually.
- ▶ Store the implant components in their original packaging and only remove them from their original and protective packaging immediately prior to application.
- ▶ Use the implant system storage devices for processing, sterilization and sterile setup.
- ▶ Use a suitable tray for cleaning/disinfection.
- ▶ Use the system storage device only for sterilization and sterile provision.
- ▶ Ensure that the implant components do not come into contact with each other or with instruments in their implant system storage devices.
- ▶ Ensure that the implant components are not damaged in any way.

Validated reprocessing procedure

General safety instructions

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.
The recommended chemistry was used for validation.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.
Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.
Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.
Additional drying, if necessary.
Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Cleaning/Disinfection

Device-specific safety information on the reprocessing method

Targon® implants not falling under the above restrictions can be resterilized after proper processing.

DANGER

Danger to the patient!
The product must not be manually prepared.

- Only mechanically reprocess the product.

WARNING

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- Handle the implants with new gloves only.
- Keep the implant system storage devices covered or closed.
- Process implant system storage devices separately from instrument trays.
- Clean implants must not be cleaned together with contaminated implants.
- Reprocess the implant components individually and separately if no implant system storage devices are available. When doing this, make certain that the implant components are not damaged.
- Mechanically clean and disinfect the implant components.
- Do not reuse surgically contaminated implants.

WARNING

Direct or indirect contamination may render implants unsuitable for resterilization!

- Do not reprocess implants that have been directly or indirectly contaminated with blood.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Mechanical alkaline cleaning and thermal disinfection	■ 20 ml disposable syringe ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).	Chapter Mechanical cleaning/disinfection and subsection: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Mechanical cleaning/disinfection

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	■ Concentrate, alkaline: – pH = 13 – <5 % anionic surfactant ■ 0.5 % working solution – pH = 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functional and undamaged.
- Immediately sort out damaged or inoperative products.

Packaging

- Place the product in its holder or on a suitable tray.
- Adequately wrap the trays in accordance with the sterilization procedure used (e. g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against contamination of the product during storage.

Sterilization method and parameters

- Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Intraoperative orientation points

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- Treatment of a fracture with implants may not always be successful in completely restoring the original bone anatomy.
- Following the fracture treatment, the functionality of adjacent joints may be restricted.
- Pains may occur following the fracture treatment.
- The implants must not be exposed to overload due to excessive strains, hard physical work or sports activities. Overtraining entails the risk of implant loosening or breakage.
- The patient must be informed about the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The implants will not sustain the load and functionality of a healthy bone. The risks of transgressing these rules must be explained to the patient.
- If the implant loosens, fragments migrate, or pseudoarthrosis develops, a revision operation may be required.
- The patient must undergo regular medical follow-up examinations.

WARNING

Complications due to incorrect use and/or incorrect assessment/treatment of the fracture!

- Make certain that nails and screws of a suitable length and diameter are used.
- Make certain that implants are aligned correctly.

- Check correct mounting of the nail on the target device: insert the drill bit through the drill sleeve and the targeting hole and guide it through the nail hole.

WARNING

Risk of nail seizure, nail bending, or bone rupture!

- Do not apply excessive force to instruments or implants.
- Apply only gentle hammer taps to insert Targon® F/T/RF nails and, if necessary, PF nails into the medullary cavity.
- To avoid bone rupture, do not hammer down Targon® PH/H nails under any circumstances.
- Check and, if necessary, correct the nail insertion point.
- If a nail proves difficult to insert, slightly drill open the medullary cavity.

The operating surgeon determines the following, depending upon the indications:

- whether implant locking is necessary
- if the implant is to be locked, whether it will be locked statically or dynamically.

Note
For locking an implant without a targeting attachment, Aesculap recommends free-hand locking supported by a C-bow and laser pilot beam, or with a radiolucent angle transmission!

Drill targets missed, drill bit breakage and guide pin breakage as well as early implant failure due to crack formation caused by surface damage in the drill hole regions!

- Use a drill bit of the correct diameter.
- Level the drill insertion site with a flat counterbore or use a center punch.
- Only use sharp drill bits and guide pins, and drill carefully and without applying excessive force. Avoid direction deviations.
- During center punching, flat countersinking and drilling, do not apply force to the targeting device; do not bend or twist the targeting device.

Implant removal

The attending physician decides when to remove, in part or completely, the osteosynthesis implants through another surgical procedure.

- Prior to removing the nail, remove all screws.
- Remove the nail with quick, light taps, and with the aid of the impactor, extraction pin and slotted hammer.

Note

During implant removal, complications can arise due to jammed implants, grown-in bone tissue, etc. Under certain circumstances, the implants and/or the instruments could get damaged. For such cases we recommend using a set of special instruments, which can be ordered from Aesculap.

This instrument set should be at hand for every explantation. Other special instruments are available for dealing with broken implants. Note the instructions for use for such cases!

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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