## Aesculap Spine


- |              |                                                    |              |                                                      |
|--------------|----------------------------------------------------|--------------|------------------------------------------------------|
| <b>(GB)</b>  | <b>Instructions for use/Technical description</b>  | <b>(S)</b>   | <b>Bruksanvisning/Teknisk beskrivning</b>            |
| <b>(USA)</b> | TSPACE PEEK/TSPACE <sup>XP</sup> instruments       |              | TSPACE PEEK/TSPACE <sup>XP</sup> instrument          |
| <b>(D)</b>   | <b>Gebrauchsanweisung/Technische Beschreibung</b>  | <b>(RUS)</b> | <b>Инструкция по применению/Техническое описание</b> |
|              | TSPACE PEEK/TSPACE <sup>XP</sup> Instrumente       |              | Инструменты TSPACE PEEK/TSPACE <sup>XP</sup>         |
| <b>(F)</b>   | <b>Mode d'emploi/Description technique</b>         | <b>(CZ)</b>  | <b>Návod k použití/Technický popis</b>               |
|              | Instruments TSPACE PEEK/TSPACE <sup>XP</sup>       |              | Nástroje TSPACE PEEK/TSPACE <sup>XP</sup>            |
| <b>(E)</b>   | <b>Instrucciones de manejo/Descripción técnica</b> | <b>(PL)</b>  | <b>Instrukcja użytkowania/Opis techniczny</b>        |
|              | Instrumental TSPACE PEEK/TSPACE <sup>XP</sup>      |              | Instrumenty TSPACE PEEK/TSPACE <sup>XP</sup>         |
| <b>(I)</b>   | <b>Istruzioni per l'uso/Descrizione tecnica</b>    | <b>(SK)</b>  | <b>Návod na používanie/Technický opis</b>            |
|              | Strumenti TSPACE PEEK/TSPACE <sup>XP</sup>         |              | Nástroje TSPACE PEEK/TSPACE <sup>XP</sup>            |
| <b>(P)</b>   | <b>Instruções de utilização/Descrição técnica</b>  | <b>(TR)</b>  | <b>Kullanım Kılavuzu/Teknik açıklama</b>             |
|              | Instrumental TSPACE PEEK/TSPACE <sup>XP</sup>      |              | TSPACE PEEK/TSPACE <sup>XP</sup> el aletleri         |
| <b>(NL)</b>  | <b>Gebruiksaanwijzing/Technische beschrijving</b>  |              |                                                      |
|              | TSPACE PEEK/TSPACE <sup>XP</sup> Instrumenten      |              |                                                      |


A

B


C

D

E

## Legend

- A** Insertion instrument SN305R (assembled)
- B** Insertion instrument SN305R (disassembled):
  - 1** Tube
  - 2** Switch
  - 3** Rotary nut
  - 4** Irrigation connector
  - 5** Pin
  - 6** Rod
  - 7** Dihedron
- C** Impact cylinder SN320R
- D** Trial implant
- E** Filling holder SN304R

## Symbols on product and packages

|                                                                                     |                                                                                         |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
|    | Caution, general warning symbol<br>Caution, see documentation supplied with the product |
|    | Date of manufacture                                                                     |
|  | Implant coupled                                                                         |
|  | Implant uncoupled                                                                       |

## Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bb Braun.com>

## Intended use

The TSPACE PEEK/TSPACE<sup>XP</sup> instruments are used to insert TSPACE PEEK (third generation) and TSPACE<sup>XP</sup> trial implants and implants into the thoracolumbar spine.

The insertion instrument SN305R is used for implanting implants and for determining the implant size by inserting trial implants into the intervertebral space.

The impact cylinder SN320R is coupled to the insertion instrument. The impact cylinder is used to expel the trial implants from the intervertebral space or to revise the implants.

The filling holder SN304R is used to hold the implants on the operating table while the inner window is filled with bone material.

## Available sizes

| Art. no. | Designation                                           |
|----------|-------------------------------------------------------|
| SN304R   | TSPACE PEEK/TSPACE <sup>XP</sup> filling holder       |
| SN305R   | TSPACE PEEK/TSPACE <sup>XP</sup> insertion instrument |
| SN320R   | TSPACE PEEK/TSPACE <sup>XP</sup> impact cylinder      |

## TSPACE PEEK/TSPACE<sup>XP</sup> trial implants 26 mm

| Art. no. | Designation                        |
|----------|------------------------------------|
| SN322R   | Trial implant 5° 26 x 11.5 x 7 mm  |
| SN323R   | Trial implant 5° 26 x 11.5 x 8 mm  |
| SN324R   | Trial implant 5° 26 x 11.5 x 9 mm  |
| SN325R   | Trial implant 5° 26 x 11.5 x 10 mm |
| SN326R   | Trial implant 5° 26 x 11.5 x 11 mm |
| SN327R   | Trial implant 5° 26 x 11.5 x 12 mm |
| SN328R   | Trial implant 5° 26 x 11.5 x 13 mm |
| SN329R   | Trial implant 5° 26 x 11.5 x 14 mm |
| SN330R   | Trial implant 5° 26 x 11.5 x 15 mm |
| SN332R   | Trial implant 5° 26 x 11.5 x 17 mm |

## TSPACE PEEK/TSPACE<sup>XP</sup> trial implants 30 mm

| Art. no. | Designation                        |
|----------|------------------------------------|
| SN352R   | Trial implant 5° 30 x 11.5 x 7 mm  |
| SN353R   | Trial implant 5° 30 x 11.5 x 8 mm  |
| SN354R   | Trial implant 5° 30 x 11.5 x 9 mm  |
| SN355R   | Trial implant 5° 30 x 11.5 x 10 mm |
| SN356R   | Trial implant 5° 30 x 11.5 x 11 mm |
| SN357R   | Trial implant 5° 30 x 11.5 x 12 mm |
| SN358R   | Trial implant 5° 30 x 11.5 x 13 mm |
| SN359R   | Trial implant 5° 30 x 11.5 x 14 mm |
| SN360R   | Trial implant 5° 30 x 11.5 x 15 mm |
| SN362R   | Trial implant 5° 30 x 11.5 x 17 mm |

## TSPACE PEEK/TSPACE<sup>XP</sup> trial implants 34 mm

| Art. no. | Designation                        |
|----------|------------------------------------|
| SN382R   | Trial implant 5° 34 x 11.5 x 7 mm  |
| SN383R   | Trial implant 5° 34 x 11.5 x 8 mm  |
| SN384R   | Trial implant 5° 34 x 11.5 x 9 mm  |
| SN385R   | Trial implant 5° 34 x 11.5 x 10 mm |
| SN386R   | Trial implant 5° 34 x 11.5 x 11 mm |
| SN387R   | Trial implant 5° 34 x 11.5 x 12 mm |
| SN388R   | Trial implant 5° 34 x 11.5 x 13 mm |
| SN389R   | Trial implant 5° 34 x 11.5 x 14 mm |
| SN390R   | Trial implant 5° 34 x 11.5 x 15 mm |
| SN392R   | Trial implant 5° 34 x 11.5 x 17 mm |

## Safe handling and preparation

### CAUTION

**Federal law restricts this device to sale by, or on order of a physician!**

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.

- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.

## Safe operation

**WARNING**

**Risk of injury and/or malfunction!**

- ▶ Always carry out a function check prior to using the product.

**WARNING**

**Risk of injury when using the product beyond the field of view!**

- ▶ Apply the product only under visual control.

## Preparing the implant bed

**WARNING**

**Risk of injury to nerves, blood vessels, and tissue caused by sharp instruments!**

- ▶ To avoid injuries, proceed with extreme care during the operation.

- ▶ Resect the facet joint.
- ▶ Open the intervertebral disk space.
- ▶ Restore the height of the intervertebral disk using distractors.
- ▶ Clear the intervertebral disk space using rongeurs, sharp spoons and curettes.
- ▶ If necessary, apply bone material.

## Determining the implant height with the trial implant

### Pick up the trial implant with the insertion instrument

- ▶ Turn the rotary nut **3** on the insertion instrument **A** counterclockwise ("loosen" direction) to push the rod **6** forward, see Fig. 1.


Fig. 1 Close-up view of the insertion instrument, switch open

- ▶ Couple the trial implant **D** to the dihedron **7** on the insertion instrument **A**, see Fig. 2. When doing so, pay attention to the orientation of the trial implant (orientation according to the marking on the insertion instrument).


Fig. 2 Coupling the trial implant to the insertion instrument

#### Note

*The switch 2 can only be turned when the rod is in the front position and the pin 5 is fully visible in the switch.*

- Turn the switch 2 to the right (↻ direction), see Fig. 3.


Fig. 3 Close-up view of the insertion instrument, switch closed

- Hold the trial implant D at 0° position on the insertion instrument A, see Fig. 8.

- Turn the rotary nut 3 on the insertion instrument A clockwise ("fix" direction) until the trial implant is clamped, see Fig. 4.


Fig. 4 Close-up view of the insertion instrument, trial implant/implant clamped

- Check manually that the trial implant is securely in place.

#### Inserting the trial implant

- Insert the trial implant D with the insertion instrument 1 straight (0° position) into the prepared implant bed. To do so, carefully hit with a hammer on the impact surface of the rotary nut 3.
- To turn the trial implant D:
  - Gently turn the rotary nut 3 counterclockwise ("loosen" direction). When doing so, make certain that the pin 5 remains outside of the switch 2.
  - Impact the trial implant D further with the hammer until the required end position is reached.
- Check the position of the trial implant D using intraoperative X-ray control.

### Removing the trial implant

- Mount the impact cylinder **C** onto the insertion instrument **A**, see Fig. 5.


Fig. 5 Mounting the impact cylinder onto the insertion instrument

- Carefully expel the trial implant **D** using the impact cylinder **C** out of the intervertebral space, see Fig. 6.


Fig. 6 Expelling the trial implant/implant with the impact cylinder

- Uncouple the impact cylinder **C** from the insertion instrument **A**.

### Uncoupling the trial implant from the insertion instrument

- Turn the rotary nut **3** counterclockwise ("loosen" direction), until the pin **5** is fully in the switch **2**, see Fig. 3. When doing so, ensure that the pin **5** does not jam the switch **2**.
- Turn the switch **2** to the left (↶ direction), see Fig. 1.
- Remove the trial implant.
- If necessary, repeat the process with the next biggest trial implant, until the correct implant size is determined.

### Inserting the implant in the implant bed

#### Picking up the implant with the insertion instrument

- Select an implant according to the size of the trial implant.
- Remove the implant with protective packaging from the sterile packaging.
- Mount the implant **a** in the protective packaging onto the dihedron **7** of the insertion instrument **A**, see Fig. 7. When doing so, pay attention to the orientation of the implant (orientation according to the marking on the insertion instrument).


Fig. 7 Coupling the implant to the insertion instrument

#### Legend

**a** Implant

- Turn the switch **2** to the right (↷ direction), see Fig. 3.
- Hold the implant at 0° position, see Fig. 8.


Fig. 8 Implant at 0° position

- Turn the rotary nut **3** on the insertion instrument **A** clockwise (direction: "fix"), until the implant is clamped, see Fig. 4.

- Remove the protective packaging from the implant.
- Check manually that the implant is securely in place.

### Filling the implant with bone material (optional)

- Place the implant **a** in the filling holder **E**, see Fig. 9. When doing so, observe the implant length (marking on the filling holder: L26/L30/L34).


Fig. 9 Implant in filling holder

- Compress the bone grafts or bone replacement material in the inner window of the implant **a** using the pestle provided in the set. When doing so, ensure that the pestle does not damage the surface of the implant.

### Inserting the implant

- Insert the implant **a** with the insertion instrument **A** straight (0° position) into the prepared implant bed. To do so, carefully hit with a hammer on the impact surface of the rotary nut **3**.


Fig. 10 Inserting the implant into the implant bed

- To turn the implant **a**:
  - Gently turn the rotary nut **3** counterclockwise ("loosen" direction). When doing so, make certain that the pin **5** remains outside of the switch **2**.
  - Impact the implant **a** further with the hammer until the required end position is reached.
- Check the position of the implant **a** using intraoperative X-ray control.
- As soon as the end position (see Fig. 11) is reached: Turn the rotary nut **3** counterclockwise ("loosen" direction), until the pin **5** is fully in the switch **2**, see Fig. 3. When doing so, ensure that the pin **5** does not jam the switch **2**.


Fig. 11 Implant in end position

- Turn the switch **2** to the left (↶ direction), see Fig. 1.

- ▶ Carefully remove the insertion instrument **A** from the implant **a**. Avoid any tilting or canting of the instrument when doing so.
- ▶ Place bone or bone replacement material around the implant **a**.

### Removing the implant intraoperatively

If, in the opinion of the surgeon, it becomes necessary to remove the implant intraoperatively, carry out the following steps:

- ▶ Check the position for coupling the insertion instrument **A**, see Fig. 1.
- ▶ Couple insertion instrument **A** to implant **a** in situ, see Fig. 7.
- ▶ Turn the switch **2** to the right (↻ direction), see Fig. 3.
- ▶ Mount the impact cylinder **B** onto the insertion instrument **A**, see Fig. 5.
- ▶ Carefully expel the implant **a** using the impact cylinder **C** out of the intervertebral space, see Fig. 6.

## Disassembling

### Insertion instrument SN305R

- ▶ Turn the rotary nut **3** counterclockwise ("loosen" direction).
- ▶ Turn the switch **2** to the left (↻ direction), see Fig. 1.
- ▶ Turn the rotary nut **3** clockwise ("fix" direction), until the rod **6** is out of the thread.
- ▶ Remove the rod **6** from the tube **1**, see Fig. 12.


Fig. 12 Insertion instrument disassembled

## Assembling

- ▶ Turn the switch **2** to the left (↻ direction).
- ▶ Insert the rod **6** into the tube **1**, see Fig. 12.
- ▶ Turn the rotary nut **3** counterclockwise ("loosen" direction) to push the rod **6** forward, see Fig. 1.

## Validated reprocessing procedure

### General safety instructions

#### Note

*Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.*

#### Note

*For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.*

#### Note

*Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.*

#### Note

*Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.*

*The recommended chemistry was used for validation.*

#### Note

*If there is no final sterilization, then a virucidal disinfectant must be used.*

#### Note

*For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>.*

*The validated steam sterilization procedure was carried out in the Aesculap sterile container system.*

### General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to Publications, Red Brochure – Proper maintenance of instruments.

### Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.
- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- ▶ In case of bone or tissue residues on the product: Pre-clean the product manually with a brush.

### Preparation before cleaning

- ▶ Dismantle the product prior to cleaning, see Disassembling.
- ▶ For mechanical cleaning: Mount the irrigation connector 4 to the rear end of the rotary nut 3.

### Cleaning/disinfection

#### Product-specific safety notes on the reprocessing procedure

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
  - are approved for plastic material and high-grade steel,
  - do not attack softeners (e.g. in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 96 °C.

## Validated cleaning and disinfection procedure

| Validated procedure                                                                                                                 | Specific requirements                                                                                                                                                                                                | Reference                                                                                                                                                                                                        |
|-------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection<br>■ SN304R<br>■ Trial implants                                                         | <ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>                             | Chapter Manual cleaning/disinfection and sub-chapter:<br>■ Chapter Manual cleaning with immersion disinfection                                                                                                   |
| Manual cleaning with ultrasound and immersion disinfection<br>■ SN305R<br>■ SN320R                                                  | <ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Drying phase: Use a lint-free cloth or medical compressed air</li> </ul>                             | Chapter Manual cleaning/disinfection and sub-chapter:<br>■ Chapter Manual cleaning with ultrasound and immersion disinfection                                                                                    |
| Mechanical alkaline cleaning and thermal disinfection<br>■ SN304R                                                                   | <ul style="list-style-type: none"> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> </ul>                                                                        | Chapter Mechanical cleaning/disinfecting and sub-chapter:<br>■ Chapter Mechanical alkaline cleaning and thermal disinfecting                                                                                     |
| Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection<br>■ SN305R<br>■ Trial implants | <ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> </ul> | Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter:<br>■ Chapter Manual pre-cleaning with a brush<br>■ Chapter Mechanical alkaline cleaning and thermal disinfecting              |
| Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection<br>■ SN320R     | <ul style="list-style-type: none"> <li>■ Suitable cleaning brush</li> <li>■ 20 ml disposable syringe</li> <li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li> </ul> | Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter:<br>■ Chapter Manual pre-cleaning with ultrasound and brush<br>■ Chapter Mechanical alkaline cleaning and thermal disinfecting |

# Aesculap®

## TSPACE PEEK/TSPACE<sup>XP</sup> instruments

### Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

### Manual cleaning with immersion disinfection

| Phase | Step                     | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                             |
|-------|--------------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------------|
| I     | Disinfecting<br>cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| II    | Intermediate rinse       | RT (cold)    | 1          | -            | D-W              | -                                                                    |
| III   | Disinfection             | RT (cold)    | 15         | 2            | D-W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| IV    | Final rinse              | RT (cold)    | 1          | -            | FD-W             | -                                                                    |
| V     | Drying                   | RT           | -          | -            | -                | -                                                                    |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.

#### Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.

- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

#### Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

#### Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

#### Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

#### Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

### Manual cleaning with ultrasound and immersion disinfection

| Phase | Step                | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|---------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Ultrasonic cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse  | RT (cold)    | 1          | -            | D-W              | -                                                              |
| III   | Disinfection        | RT (cold)    | 15         | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse         | RT (cold)    | 1          | -            | FD-W             | -                                                              |
| V     | Drying              | RT           | -          | -            | -                | -                                                              |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\*Recommended: BBraun Stabimed

# Aesculap®

## TSPACE PEEK/TSPACE<sup>XP</sup> instruments

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

### Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

### Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

### Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

### Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

### Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

## Mechanical cleaning/disinfecting

*Note*

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

*Note*

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

## Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical/Note                                                                                                                                                                                                                           |
|-------|----------------------|--------------|------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse             | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                       |
| II    | Cleaning             | 55/131       | 10         | FD-W             | <div><div>■</div>Concentrate, alkaline:<ul style="list-style-type: none"><li>- pH = 13</li><li>- &lt;5 % anionic surfactant</li></ul><div>■</div>0.5 % working solution<ul style="list-style-type: none"><li>- pH = 11*</li></ul></div> |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                       |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                       |
| V     | Drying               | -            | -          | -                | According to the program for cleaning and disinfection device                                                                                                                                                                           |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Aesculap®

TSPACE PEEK/TSPACE<sup>XP</sup> instruments

Mechanical cleaning/disinfection with manual pre-cleaning

*Note*

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

*Note*

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfectant cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing               | RT (cold)    | 1          | -            | D-W              | -                                                              |

D-W:            Drinking water

RT:            Room temperature

\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual pre-cleaning with ultrasound and brush

| Phase | Step                | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|---------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Ultrasonic cleaning | RT (cold)    | >15        | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing             | RT (cold)    | 1          | -            | D-W              | -                                                              |

D-W:            Drinking water  
RT:            Room temperature  
\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

# Aesculap®

## TSPACE PEEK/TSPACE<sup>XP</sup> instruments

### Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical                                                                                                                                                                                                                                            |
|-------|----------------------|--------------|------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse             | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                                   |
| II    | Cleaning             | 55/131       | 10         | FD-W             | <div> <div>■</div> Concentrate, alkaline: <ul style="list-style-type: none"> <li>- pH = 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> <div>■</div> 0.5 % working solution <ul style="list-style-type: none"> <li>- pH = 11*</li> </ul> </div> |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                                   |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                                   |
| V     | Drying               | -            | -          | -                | According to the program for cleaning and disinfection device                                                                                                                                                                                       |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

## Inspection, maintenance and checks

CAUTION

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ Prior to function checks, lubricate the switch and rotary nut on the insertion instrument SN305R with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Assemble dismountable products, see Assembling.
- ▶ Check for compatibility with associated products.

## Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

## Steam sterilization

### Note

*Insertion instrument SN305R may only be sterilized in disassembled condition.*

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
  - Disassemble the product
  - Steam sterilization through fractionated vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

# Aesculap®

## TSPACE PEEK/TSPACE<sup>XP</sup> instruments

### Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10<sup>-6</sup>, Aesculap recommends the following parameters:

| Aesculap Orga Tray/Sterile container (perforated bottom) |               |       |                     |
|----------------------------------------------------------|---------------|-------|---------------------|
| Minimum cycle parameters*                                |               |       |                     |
| Sterilization method                                     | Temp.         | Time  | Minimum drying time |
| Prevacuum                                                | 270 °F/275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

### WARNING for the US market

**If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.**

### Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

### Technical Service

**Risk of injury and/or malfunction!**

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

### Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 (7461) 95-1602  
Fax: +49 (7461) 16-5621  
E-Mail: [ats@aesculap.de](mailto:ats@aesculap.de)

Or in the US:  
Aesculap Implant Systems LLC  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042

Aesculap Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

## **Disposal**

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

## **Distributor in the US/Contact in Canada for product information and complaints**

Aesculap Implant Systems LLC  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany