- (GB)** Instructions for use/Technical description
- (USA)** YASARGIL aneurysm clip removal forceps
- (D)** Gebrauchsanweisung/Technische Beschreibung  
YASARGIL Aneurysmen-Clip-Abnehmezange
- (F)** Mode d'emploi/Description technique  
Pince de retrait de clips d'anévrisme YASARGIL
- (E)** Instrucciones de manejo/Descripción técnica  
Pinza de extracción de clips de aneurismas YASARGIL
- (I)** Istruzioni per l'uso/Descrizione tecnica  
Pinze per rimozione di clip da aneurisma YASARGIL
- (P)** Instruções de utilização/Descrição técnica  
Pinça de remoção de clips para aneurisma YASARGIL
- (NL)** Gebruiksaanwijzing/Technische beschrijving  
YASARGIL aneurysmaclip-verwijderingstang
- (S)** Bruksanvisning/Teknisk beskrivning  
YASARGIL aneurysm-clip-borttagningstång
- (RUS)** Инструкция по применению/Техническое описание  
Зажимы для снятия аневризматических клипсов YASARGIL
- (CZ)** Návod k použití/Technický popis  
Kleště na sejmutí aneurysmatické cévní svorky YASARGIL
- (PL)** Instrukcja użytkowania/Opis techniczny  
YASARGIL kleszczyki do wyjmowania zacisków tętniakowych
- (SK)** Návod na používanie/Technický opis  
Kliešte na odstránenie svoriek pre aneurizmy YASARGIL
- (TR)** Kullanım Kılavuzu/Teknik açıklama  
YASARGIL anevrizma klipsi çıkarma penseti

## B | BRAUN

SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

TA-Nr. 013147 11/14 V6 Änd.-Nr. 51256

CE 0123 – DIR 93/42/EEC

Technical alterations reserved

# Aesculap®

## YASARGIL aneurysm clip removal forceps

### Legend

- 1 Jaw part
- 2 Branches left/right
- 3 Clip removal forceps

### Symbols on product and packages

Caution, general warning symbol  
Caution, see documentation supplied with the product

### Intended use

Aesculap YASARGIL aneurysm clip removal forceps are used to remove Aesculap YASARGIL titanium aneurysm clips in neurosurgical procedures.

*Note*  
In the following, the Aesculap YASARGIL aneurysm clips will be referred to simply as "clips".

### Available sizes

Detailed information on available sizes can be found in the main Neurosurgery Catalog and in the YASARGIL aneurysm Clips brochures.  
All removal forceps are marked according to size (Mini or Standard) and clip material (titanium) so that their correct application with clips of the appropriate size and material is ensured.

### Safe handling and preparation

- CAUTION**  
**Federal law restricts this device to sale by, or on order of a physician!**
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
  - ▶ Read, follow, and keep the instructions for use.
  - ▶ Use the product only in accordance with its intended use, see Intended use.
  - ▶ Clean the new product mechanically after removing its transport packaging and prior to its initial sterilization.
  - ▶ Store any new or unused products in a dry, clean, and safe place.
  - ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
  - ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
  - ▶ Replace any damaged components immediately with original spare parts.

### Safe operation

**Damage, imprecise function and incorrect closing force of the clips due to incorrect choice of removal forceps!**  
▶ Only use Aesculap clip removal forceps.  
▶ Use clips only with the appropriate Aesculap clip removal forceps (see labels).

- ▶ Use the clip removal forceps only in the following combinations:
  - Titanium Standard clips with the clip removal forceps for titanium Standard clips
  - Titanium Mini clips with the clip removal forceps for titanium Mini clips
- ▶ Adapt the applied clip in jaw piece 1.
- ▶ Actuate branches left/right 2.
- ▶ Take clip removal forceps 3 with the clip out of the operating field.

### Validated reprocessing procedure

#### General safety notes

*Note*  
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

*Note*  
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

*Note*  
It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.  
The specified chemistry was used for validation.

*Note*  
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at [www.extranet.bbraun.com](http://www.extranet.bbraun.com)  
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

#### General information

- Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.  
Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.  
Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result damage to metallic products. These must be removed by rinsing thoroughly with demineralized water and then drying.  
Perform additional drying, if necessary.  
Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:
- Optical changes in materials, e.g., fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
  - Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.

- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see [www.a-k-l.org](http://www.a-k-l.org), link to Publications, Red Brochure – Proper maintenance of instruments.

#### Disassembling the product before carrying out the reprocessing procedure

- ▶ Open up products with hinges.

#### Preparations at the place of use

- ▶ If applicable, rinse surfaces that are not accessible to visible inspection (preferably with demineralized water), using a disposable syringe, for example.
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

#### Cleaning/disinfection

##### Product-specific safety notes on the reprocessing procedure

**DANGER**

**Danger to the patient!**  
▶ Reprocess the product only with manual pre-cleaning followed by mechanical cleaning.

**DANGER**

**Risk to patient due to cross contamination!**  
▶ Do not clean contaminated products together with uncontaminated products in a tray.

**CAUTION**

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**  
▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which

- are approved for use, for example, on aluminum, plastic materials, and high-grade steel,
- do not attack softeners (e.g. in silicone).

▶ Observe specifications regarding concentration, temperature and exposure time.  
▶ Do not exceed the maximum permitted cleaning temperature of 55 °C.

- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- ▶ Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

#### Mechanical cleaning/disinfection with manual pre-cleaning

*Note*  
The cleaning and disinfecting machine must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

*Note*  
The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

##### Manual pre-cleaning with a brush

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemicals                          |
|-------|-----------------------|--------------|------------|--------------|------------------|------------------------------------|
| I     | Disinfectant cleaning | RT<br>(cold) | >15        | 1            | D–W              | Dr. Weigert neodisher® SeptoClean* |
| II    | Rinsing               | RT<br>(cold) | 1          | –            | D–W              | –                                  |

D–W: Drinking water  
RT: Room temperature  
\*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

- Phase I**
- ▶ Fully immerse the product in the cleaning/disinfecting solution for at least 15 minutes. Ensure that all accessible surfaces are moistened.
  - ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed.
  - ▶ If applicable, brush non-visible surfaces for at least 1 minute with a suitable cleaning brush.
  - ▶ Mobilize non-rigid components, such as set screws and hinges, during cleaning.
  - ▶ Then flush these areas thoroughly at least five times with the cleaning disinfectant solution using a disposable syringe (20 ml).

- Phase II**
- ▶ Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
  - ▶ Mobilize non-rigid components, such as set screws and hinges, during rinsing.

##### Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemicals                                                                |
|-------|----------------------|--------------|------------|------------------|--------------------------------------------------------------------------|
| I     | Pre-rinse            | <25/77       | 3          | D–W              | –                                                                        |
| II    | Cleaning             | 55/131       | 10         | FD–W             | Dr. Weigert neodisher® SeptoClean 1 %* working solution                  |
| III   | Intermediate rinse   | >10/50       | 1          | FD–W             | –                                                                        |
| IV    | Thermal disinfection | 90/194       | 5          | FD–W             | –                                                                        |
| V     | Drying               | –            | –          | –                | In accordance with the program for the cleaning and disinfecting machine |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

- ▶ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots)
- ▶ Place instruments in the tray with their hinges open.
- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Check for compatibility with associated products.

Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Store products with ratchet locks fully opened or locked no further than in the first notch.
- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
  - Steam sterilization using fractional vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractional vacuum process at 134 °C, holding time 18 minutes for prion inactivation
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10<sup>-6</sup>, Aesculap recommends the following parameters:

| Aesculap Orga Tray/Sterile container (perforated bottom) |               |       |                     |
|----------------------------------------------------------|---------------|-------|---------------------|
| Minimum cycle parameters*                                |               |       |                     |
| Sterilization method                                     | Temp.         | Time  | Minimum drying time |
| Prevacuum                                                | 270 °F/275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

**Risk of injury and/or malfunction!**

- ▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 (7461) 95-1602  
Fax: +49 (7461) 16-5621  
E-Mail: ats@aesculap.de  
Or in the US:  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042  
Aesculap Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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