





Aesculap® Acculan 4

- (GB)** Instructions for use/Technical description
(USA) Dermatome GA340 / Dermatome 0.1 mm GA341
- (D)** Gebrauchsanweisung/Technische Beschreibung
Dermatom GA340 / Dermatom 0,1 mm GA341
- (F)** Mode d'emploi/Description technique
Dermatome GA340 / Dermatome 0,1 mm GA341
- (E)** Instrucciones de manejo/Descripción técnica
Dermátomo GA340 / Dermátomo 0,1 mm GA341
- (I)** Istruzioni per l'uso/Descrizione tecnica
Dermatomo GA340 / Dermatomo 0,1 mm GA341
- (P)** Instruções de utilização/Descrição técnica
Dermátomo GA340 / Dermátomo 0,1 mm GA341
- (NL)** Gebruiksaanwijzing/Technische beschrijving
Dermatoom GA340 / dermatoom 0,1 mm GA341

- (S)** Bruksanvisning/Teknisk beskrivning
Dermatom GA340 / Dermatom 0,1 mm GA341
- (RUS)** Инструкция по применению/Техническое описание
Дерматом GA340/дерматом 0,1 мм GA341
- (CZ)** Návod k použití/Technický popis
Dermatom GA340 / Dermatom 0,1 mm GA341
- (PL)** Instrukcja użytkowania/Opis techniczny
Dermatom GA340 / dermatom 0,1 mm GA341
- (SK)** Návod na použitie/Technický opis
Dermatóm GA340/Dermatóm 0,1 mm GA341
- (TR)** Kullanım Kılavuzu/Teknik açıklama
Dermatom GA340 / Dermatom 0,1 mm GA341
- (KR)** 사용 설명서 / 기술 설명
피부 절편기 GA340 / 피부 절편기 0.1 mm GA341

Aesculap® Acculan 4

Dermatome GA340 / Dermatome 0.1 mm GA341

Legend

- 1 Product (Dermatome)
- 2 Trigger (for oscillation frequency control)
- 3 Safety catch
- 4 Battery bay
- 5 Locking pin
- 6 Sterile funnel
- 7 Battery
- 8 Lid
- 9 Lid release element
- 10 Battery removal device
- 11 Dermatome blade
- 12 Sliding surfaces
- 13 Blade cover
- 14 Roll pin
- 15 Clamping lever
- 16 Blade guide peg
- 17 Pusher dog
- 18 Setting lever
- 19 Lock (of setting lever)
- 20 Dial wheel
- 21 Surface
- 22 Guide slots
- 23 Flap rod
- 24 Flaps
- 25 Flap tab
- 26 Nut
- 27 Symbol for alignment of flaps (on the flap rod)

The depictions are schematic only.

Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
 YYYY-MM	Maintenance label Indication of the next maintenance appointment (Date: Year-Month)
	Machine-readable, two-dimensional code The code contains a unique serial number which can be used for electronic tracking of the individual instrument. The serial number is based on the global standard sGTIN (GS1).
	Date of manufacture
	Follow the instructions for use
	Labeling of electrical and electronic devices pursuant to directive 2002/96/EG (WEEE)
	Classification Type BF
	Motor speed control
	Rotational direction for loosening the nut
	Direction of rotation to tighten the nut
	Alignment of flaps on the flap rod
	Manufacturer's batch designation

	Manufacturer's serial number
	Manufacturer's article number
	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage
S9	Operating mode

Contents

1.	Applicable to.....	3
2.	General information.....	4
2.1	Intended use.....	4
2.2	Main functions and design characteristics	4
2.3	Indications	4
2.4	Absolute contraindications	4
2.5	Relative contraindications.....	4
3.	Safe handling	4
4.	Product description	5
4.1	Scope of supply.....	5
4.2	Components required for operation	5
4.3	Operating principle.....	5
5.	Preparation	5
6.	Working with the device	5
6.1	System set-up.....	6
6.1.1	Connecting the accessories.....	6
6.1.2	Inserting the rechargeable battery.....	6
6.1.3	Intraoperative battery change.....	6
6.1.4	Removing the rechargeable battery	6
6.1.5	Protection against inadvertent activation	7
6.1.6	To mount the flap rod and flaps	7
6.1.7	Inserting the dermatome blade.....	7
6.1.8	Removing the dermatome blade.....	7
6.1.9	Intraoperative storage	7
6.2	Function checks	7
6.3	Safe operation	8

6.3.1	Adjusting the cutting thickness	8
6.3.2	Adjusting the cutting width	8
6.3.3	Operating the product.....	8
6.3.4	Taking skin grafts.....	8
7.	Validated reprocessing procedure	9
7.1	General safety notes	9
7.2	General information	9
7.3	Preparations at the place of use.....	10
7.4	Disassembling the product before carrying out the reprocessing procedure	10
7.4.1	Removing the flap rod of the dermatome	10
7.5	Preparation before cleaning	10
7.6	Single-use products.....	10
7.7	Cleaning/disinfection.....	10
7.7.1	Product-specific safety instructions for the reprocessing procedure	10
7.8	Manual cleaning with wipe disinfection	11
7.9	Automatic cleaning/disinfection with manual pre-cleaning ..	12
7.9.1	Manual pre-cleaning with a brush.....	12
7.9.2	Mechanical alkaline cleaning and thermal disinfection.....	12
7.10	Inspection, maintenance and checks.....	13
7.11	Packaging	13
7.12	Steam sterilization.....	13
7.13	Sterilization for the US market.....	13
7.14	Storage	13
8.	Maintenance	14
9.	Troubleshooting list.....	14
10.	Technical Service	16
11.	Accessories/Spare parts.....	16
12.	Technical data	17
12.1	Classification acc. to Directive 93/42/EEC.....	17
12.2	Performance data, information about standards	17
12.3	Operating mode.....	17
12.4	Environmental conditions	17
13.	Disposal	17
14.	Distributor in the US/Contact in Canada for product information and complaints	17

1. Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

Aesculap® Acculan 4

Dermatome GA340 / Dermatome 0.1 mm GA341

2. General information

2.1 Intended use

Task/Function

The Dermatome GA340 / Dermatome 0.1 mm GA341, combined with dermatome blade, is used for split skin removal in adjustable depth.

Application Environment

The product fulfills the requirements for type BF pursuant to IEC/DIN EN 60601-1 and is used in operating rooms in sterile environments of explosion risk areas (such as areas with pure oxygen or anaesthesia gases).

2.2 Main functions and design characteristics

Oscillation frequency	min. 0 min ⁻¹ to max. 6 500 min ⁻¹
Operating mode	Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)
	■ 60 second application, 60 second pause ■ 10 repetitions ■ 30 min cooling time ■ Max. temperature 48 °C

Electrical systems generally heat up during continual operation. It is advised to give the system a break after use to cool down, as listed in the table on operating mode.

Heating depends on the tool used and the load. After a certain number of repetitions, the system should cool down. This procedure prevents the system overheating as well as possible injury to the patient or user.

The user is responsible for the use and adherence to the pause sequence described.

2.3 Indications

The type and area of application depend on the tool selected.

2.4 Absolute contraindications

The product is not licensed for use on the central nervous system or central circulatory system.

2.5 Relative contraindications

The safe and effective use of the product greatly depends on influences which can only be controlled by the user. Therefore the specifications provided represent framework conditions only.

Clinically successful use of the product is dependent on the knowledge and experience of the surgeon. The surgeon must decide which structures it is sensible to treat and take into account the safety and warning information contained in these instructions for use.

3. Safe handling

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

WARNING

Risk of injury and material damage if this product is not used as intended!

► Use the product only for its intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

► Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Prior to use, check that the product is in good working order.
- Observe TA022130 "Notes on Electromagnetic Compatibility (EMC)".
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.
- Ensure that the electrical installation of the room is consistent with the requirements of IEC/DIN EN.

- ▶ Do not operate the product in explosion-hazard areas.
- ▶ Sterilize product before use.
- ▶ When using the ECCOSfixation system, adhere to the appropriate instructions for use for TA009721, see Aesculap Extranet at <https://extranet.bbraun.com>

4. Product description

4.1 Scope of supply

Art. no.	Designation
GA340	Dermatome
GA341	– or – Dermatome 0.1 mm
GA678	Sterile funnel
TA014549	Instructions for use of Dermatome GA340/GA341 (flyer)

4.2 Components required for operation

- Li-Ion battery short GA346 (charged)
- Sterile funnel GA678
- Lid GA675
- Dermatome blade GB228R

4.3 Operating principle

The product **1** contains an electric motor, which is powered by a replaceable battery **7**.

The charged non-sterile battery **7** is inserted into the product **1** via the sterile funnel **6** and sealed with the lid **8** so that it is sterile.

The product moves the dermatome blade **11** in oscillation. The oscillation frequency is controlled electronically and can be continuously regulated with the trigger **2**.

5. Preparation

Aesculap assumes no liability if the following rules are not followed:

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Prior to use, inspect the product and its accessories for any visible damage.
- ▶ Use the products and their accessories only if they are in perfect technical condition.

6. Working with the device

WARNING

Risk of infection and contamination!

Product is delivered unsterilized!

- ▶ Sterilize the product before use pursuant to the operating instructions.

WARNING

Risk of injury and material damage due to accidental activation of the product!

- ▶ Products which are not being actively used must be secured against accidental activation (position OFF).

WARNING

Risk of injury and material damage due to inappropriate use of tools!

- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ When coupling/uncoupling, handle tools with cutting edges with care.

WARNING

Damage to the product if dropped!

- ▶ Use the products only if they are in perfect technical condition, see Function check

WARNING

Risk of injury and/or malfunction due to unrecognized change to cutting setting!

The basic adjustment of the cutting setting can change if the product falls on the floor, for example.

- ▶ Do not use product.
- ▶ Have the manufacturer inspect the product.

WARNING

Risk of burns to skin and tissue through blunt tools or if product has not been maintained properly!

- ▶ Use tools only if they are in perfect condition.
- ▶ Replace blunt tools.
- ▶ Maintain the product properly, see maintenance guide.

Aesculap® Acculan 4

Dermatome GA340 / Dermatome 0.1 mm GA341

6.1 System set-up

6.1.1 Connecting the accessories

DANGER

Risk of injury due to unapproved configuration using additional components!

- Ensure that the classification of all components used is consistent with the classification of the product (such as type BF or type CF).

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

All configurations must fulfill the fundamental standard IEC/DIN EN 60601-1. The person connecting the devices with each other is responsible for the configuration and must ensure that the fundamental standard IEC/DIN EN 60601-1 or relevant national standards.

- Follow the instructions for use of individual accessories.
- Please address your B. Braun/Aesculap partner or Aesculap Technical Service with any inquiries in this respect; for a contact address, see Technical Service.

6.1.2 Inserting the rechargeable battery

- Product 1 with the battery shaft 4 must be turned upwards and the sterile funnel 6 (sterile) attached, see Fig. A.
- Battery 7 (non-sterile) is inserted into the battery shaft 4 by a second (non-sterile) person, see Fig. A.

Note

After the battery has been inserted, multiple signals will sound, which indicate that the product is ready to use.

- After insertion of the battery, have the sterile funnel 6 (non-sterile) removed by a second person.
- Attach the lid 8 (sterile) in such a manner that it clicks into place with both release catches 9.

Note

The sterility of the product is only guaranteed with a correctly placed lid.

6.1.3 Intraoperative battery change

WARNING

Risk of injury from hot battery!

Battery may be hot after use in the machine.

- Remove the battery with the battery removal device and let it cool.

The battery removal device is used for intraoperative battery changes ensuring that sterile conditions are safely maintained.

- Turn the product 1 with the battery shaft 4 upwards.
- Press both safety catches 9 on the lid 8 at the same time and remove the lid 8.
- Attach the sterile battery removal device 10, see Fig. B.
- Shake the product 1 with the attached battery removal device 10 with battery shaft 4 slightly downwards.
The battery 7 will slide easily into the battery removal device 10.
- Give the battery removal device 10 including dead battery 7 to the non-sterile person.
- Insert a charged battery 7, see Inserting the rechargeable battery.

6.1.4 Removing the rechargeable battery

WARNING

Risk of injury from hot battery!

Battery may be hot after use in the machine.

- Let the battery cool in the machine and then remove.
– or –
- Remove the battery with the battery removal device and let it cool.

CAUTION

Damage to the battery caused by knocking it against hard objects!

- Only remove battery by tapping the product on the palm of the hand.

CAUTION

Damage to, or destruction of the batteries due to processing!

- Do not sterilize batteries.

Rechargeable battery must be removed after each surgical procedure and prior to processing of the device.

Note

The battery removal device can be used for easier battery removal, see Fig. B.

- ▶ Turn the product **1** with the battery shaft **4** upwards.
- ▶ Press both safety catches **9** on the lid **8** at the same time and remove the lid **8**.
- ▶ Tightly grasp the battery shaft **4** at the lower end.
- ▶ Tap the battery shaft **4** on the palm of the hand until the battery **7** slides out of the battery shaft **4** and can be removed.

6.1.5 Protection against inadvertent activation

To prevent inadvertent activation of the product during blade change, for example, the trigger can be locked.

To lock trigger:

- ▶ Twist the safety catch **3** to position OFF.
The trigger **2** is blocked and the product **1** cannot be operated.

To unlock trigger:

- ▶ Twist the safety catch **3** to position ON.
The trigger **2** is unlocked and the product **1** can be operated.

6.1.6 To mount the flap rod and flaps

- ▶ Install the flaps **24** on flap rod **23** to the edge of the rear surface. Observe the symbol **27** on the face of the thread, see Fig. C.
- ▶ Screw the nut **26** counterclockwise onto the thread of flap rod **23** (left-handed thread).
- ▶ Turn the nut **26** up to the edge of the visible surface.
- ▶ Insert and twist the complete flap rod **23** in its guide slots **22** on the side surfaces.
- ▶ Push the flap rod **23** laterally all the way so that the cross pin of the flap rod **23** comes to rest in the guide slot **22**.
- ▶ Tighten nut **26** by turning it counterclockwise.

6.1.7 Inserting the dermatome blade

DANGER

Risk of infection and contamination due to torn or cut surgical gloves!

- ▶ Avoid contact with the cutting edge when inserting the dermatome blade.

DANGER

Risk of injury and material damage due to inadvertent activation of the dermatome while inserting the dermatome blade!

- ▶ Prior to inserting the dermatome blade, engage the safety catch, see Protection against inadvertent activation.

- ▶ Engage safety catch **3**.
- ▶ Push both clamping levers **15** in the direction of the arrows, see Fig. C. Blade cover **13** is now unlocked.

- ▶ Open blade cover **13** in the direction of the arrow.
- ▶ Safely insert dermatome blade **11** in pusher dog **17** and blade guide peg **16**, see Fig. D.
- ▶ Close blade cover **13**.
- ▶ Tighten both clamping levers **15** so that roll pin **14** engages in the recess of clamping lever **15**.

6.1.8 Removing the dermatome blade

- ▶ Release both clamping levers **15**.
- ▶ Open blade cover **13**.
- ▶ Remove the dermatome blade **11** from pusher dog **17** and blade guide peg **16**.

6.1.9 Intraoperative storage

Note

The ECCOS fixation can be used for intraoperative storage. The dermatome can not be processed in this position, however. For the correct processing position, see Validated reprocessing procedure.

- ▶ Engage safety catch **3**.
- ▶ Fold back the clamp of the ECCOS fixation.
- ▶ Insert product **1** in the ECCOS fixation, see Fig. E

6.2 Function checks

The function checks must be carried out prior to each use and after each intraoperative battery change.

- ▶ Make certain there is a battery in the battery compartment.
- ▶ Ensure that the lid has completely clicked into place.
- ▶ Ensure that the dermatome blades are not mechanically damaged.
- ▶ Ensure that a dermatome blade is correctly seated in the dermatome.
- ▶ Check that the dermatome blade is securely seated.
- ▶ Check that the roll pin is securely seated.
- ▶ Release product for use (ON position).
- ▶ Briefly operate product at maximum oscillation frequency.
- ▶ Make certain that the rotational direction is correct in each case.
- ▶ Do not use the product if it is damaged or defective.
- ▶ Set aside the product if it is damaged.

Aesculap® Acculan 4

Dermatome GA340 / Dermatome 0.1 mm GA341

6.3 Safe operation

WARNING

Risk of infection or injury due to aerosol formation or particles!

- Use suitable protection (such as waterproof protective clothing, face mask, safety glasses, suction).

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

WARNING

Risk of injury when using the product beyond the field of view!

- Apply the product only under visual control.

Note

The product drive motor is operated using a magnetic sensor system. In order to prevent inadvertent activation of the motor, the product may not be exposed to any magnetic fields (such as magnetic instrument pads).

Note

The slight whistling noise when the product is starting up is due to the construction.

6.3.1 Adjusting the cutting thickness

Note

The skin condition (e.g. due to age) of the patient must be considered when adjusting the cutting depth (graft thickness).

- Set the cutting depth with setting lever **18**. The lock **19** of the setting lever **18** prevents any unintentional change to the setting of the cutting depth, see Fig. C.
- Read the cutting depth on dial wheel **20** (1/10 mm gradation).
- Follow the troubleshooting information, see Troubleshooting list.

6.3.2 Adjusting the cutting width

WARNING

Risk of injury/cuts by the dermatome blade if the flaps are opened incorrectly!

- Always use the flap tabs for lifting the flaps.

Note

The width of the skin graft will be approx. 2mm less than the cutting width set by the flaps.

- Set the required cutting width by opening the appropriate number of flaps **24**.
- Safely lift flaps **24** by pressing flap tabs **25**, see Fig. F.

6.3.3 Operating the product

Note

The slight whistling noise when the product is starting up is due to the construction.

- Activate trigger **2**.
The saw runs at the continuously controlled oscillation frequency.

6.3.4 Taking skin grafts

Note

To familiarize with the operation of the Dermatome prior to its first surgical application, the user should carry out some test cuts on a specimen.

- Stretch flabby skin.
- With the sliding surface in level contact with the skin, push product **1** forward in a smooth movement, applying gentle pressure. Be careful not to cant the product **1**.
- For very thin grafts, use a forceps, from time to time, to lay up the resected portion of the graft on surface **21**, see Fig. C.
- Turn off the motor and pull the freshly cut skin graft out of the device
- or -
- Lower the handle and cut off the skin graft with the product **1**.

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this. The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

*For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.*

7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

On stainless steel, residues containing chlorine or chloride (such as surgical residues, drugs, saline solutions in water for cleaning, disinfection and sterilization) may lead to corrosion (pitting corrosion, tensile corrosion) and thus to the destruction of the product. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Visual material changes (such as fading or color changes in titanium or aluminum). For aluminum, the application/process solution only needs to be pH >8 to cause visible surface changes.
- Material damage (such as corrosion, cracks, breaks, premature aging or swelling).
- Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Aesculap® Acculan 4

Dermatome GA340 / Dermatome 0.1 mm GA341

7.3 Preparations at the place of use

- ▶ Remove all attached components from the product (dermatome blade, battery, accessories).
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

7.4 Disassembling the product before carrying out the reprocessing procedure

- ▶ Remove battery 7 from the product 1, see Removing the rechargeable battery.
- ▶ Open the clamping levers 15 to remove the dermatome blade 11, see Removing the dermatome blade.
- ▶ Disassemble the product 1 immediately after use, according to instructions.

7.4.1 Removing the flap rod of the dermatome

- ▶ Loosen nut 26 by turning it clockwise (left-handed thread).
- ▶ Screw the nut 26 back to the end of the visible surface.
- ▶ Press on nut 26 and push flap rod 23 approx. 4 mm to one side.
- ▶ Twist flap rod 23 until it can be removed.
- ▶ Remove flap rod 23.
- ▶ Slide off the flaps 24 from flap rod 23.
The flap rod 23 has been dismantled.

7.5 Preparation before cleaning

- ▶ Prior to the first mechanical cleaning/disinfection process, mount the ECCOS fixations in a suitable tray (e.g. JC254R) or use the ECCOS tray GB256R with fixations equipped.
- ▶ Insert the products in the correct position into the ECCOS fixations, see Fig. G.

7.6 Single-use products

Dermatome blade GB228R

WARNING

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- ▶ Do not reprocess the dermatome blade GB228R.

7.7 Cleaning/disinfection

7.7.1 Product-specific safety instructions for the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for plastic material and high-grade steel,
 - do not attack softeners (e.g. in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C during chemical cleaning and/or disinfection.
- ▶ Do not exceed the maximum temperature of 96 °C during thermal disinfection with FD water.
- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

CAUTION

Damage or destruction of the batteries due to processing!

- ▶ Protect battery from moisture.

Note

The indicated drying time is a guide time only. It must be checked taking into account the specific conditions (e.g. load) and if applicable adjusted.

7.8 Manual cleaning with wipe disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Pre-cleaning	RT (cold)	≥2	-	D-W	Until visually clean
II	Cleaning with enzyme solution	RT (cold)	≥2	0.8	D-W	pH-neutral*
III	Intermediate rinse	RT	≥5	-	D-W	-
IV	Drying	RT	-	-	-	-
V	Wipe disinfection	-	>1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol
VI	Final rinse	RT (cold)	0.5	-	FD-W	-
VII	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

* Validated with enzyme cleaner "Cidezyme Johnson & Johnson"

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

Phase I

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Clean the product under running water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.
- ▶ Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.

Note

For details on difficult to access surfaces, see *Acculan pre-cleaning and care information TA016000* (available in *Aesculap Extranet* at <https://extranet.bbraun.com>).

Phase II

- ▶ Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- ▶ Spray products with a pH neutral enzyme solution, let soak in for at least 2 minutes and then wipe off.

Phase III

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Rinse product under running tap water for at least 5 minutes.
- ▶ Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- ▶ Contamination should be removed with a lint-free cloth or soft brush moistened with enzyme cleaner.

- ▶ Flexible components should be rinsed for 20 seconds with the water pistol (cold water, at least 2.5 bar).
- ▶ After manual cleaning, check visible surfaces and areas of flexible components for residues.
- ▶ If necessary, repeat the cleaning process (phase 1 to 3).

Phase IV

- ▶ Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

Phase V

- ▶ Wipe all surfaces of the product with a single-use disinfecting wipe.

Phase VI

- ▶ Rinse disinfected surfaces after the prescribed reaction time for at least 1 minute under running demineralized water.
- ▶ Drain any remaining water fully.

Stage VII

- ▶ Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

Aesculap® Acculan 4

Dermatome GA340 / Dermatome 0.1 mm GA341

7.9 Automatic cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must fundamentally have a tested efficacy (such as FDA approval or CE label pursuant to DIN EN ISO 15883).

Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

7.9.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemistry/Note
I	Rinsing	RT (cold)	-	-	D-W	Until visually clean
II	Brushes	RT (cold)	-	-	D-W	Until visually clean

D-W: Drinking water
RT: Room temperature

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

Phase I

- ▶ Dismantle the flap rod for cleaning. Store flap rods, flaps and nut separately in tray.
- ▶ Mobilize non-rigid components during cleaning.
- ▶ Thoroughly clean the product under running water.

Phase II

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.
- ▶ After manual preparation, check visible surfaces for residue and repeat the pre-cleaning process as needed.

Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available in Aesculap Extranet at <https://extranet.bbraun.com>).

7.9.2 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemistry/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH ~ 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	mind. 10 min at max. 120 °C

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Insert the product in its proper position in the ECCOS fixation. Insert the battery compartment over the clamp and swivel device into the ECCOS fixation so that the trigger 2 points upwards.
- ▶ Ensure that the blade cover 13 is kept open through the entire processing procedure.
- ▶ Dismantle the flap rod 23 for cleaning. Store flap rod 23, flaps 24 and nut 26 separately in tray.
- ▶ After automatic cleaning/disinfection, check visible surfaces for residues and repeat the cleaning/disinfection process as needed.

7.10 Inspection, maintenance and checks

Note

Aesculap recommends lubricating the flap rod, flaps, sliding surfaces and other moving parts, (e.g. trigger, cover lids) with STERILIT-Power-Systems oil spray GB600 or the STERILIT-Power Systems drip lubricator GA059 prior to every sterilization. Remove excess oil with a lint-free wipe.

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Check the product for any damage, abnormal running noise, overheating or excessive vibration.
- ▶ Set aside the product if it is damaged.

7.11 Packaging

- ▶ Always follow the instructions for use of the respective packaging and storage devices (e.g. instructions for use TA009721 for Aesculap-ECCOS holder system).
- ▶ Insert the products in the correct position into the ECCOS fixations, see Fig. G.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

7.12 Steam sterilization

CAUTION

Damage or destruction of the batteries due to processing!

- ▶ Do not sterilize the rechargeable battery.

Note

The product may only be sterilized with the blade, battery and lid as well as opened blade cover all demounted.

Note

The product can be sterilized with a mounted flap tab and flaps in the open position.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Use a validated sterilization method:
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer DIN EN 285 and validated pursuant to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min

When sterilizing multiple products in one steam sterilizer:

- ▶ Ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

7.13 Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

7.14 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- ▶ Store sterile single-use GB258R product in germ-proof packaging in a dry, dark and temperature-controlled room.

Aesculap® Acculan 4

Dermatome GA340 / Dermatome 0.1 mm GA341

8. Maintenance

To ensure reliable operation, the product must be maintained in accordance with the maintenance labeling or at least once a year.

YYYY-MM

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

9. Troubleshooting list

► Have defective products repaired by Aesculap Technical Service, see Technical Service.

Issue	Cause	Detection	Remedy
Product not running	No battery	No battery in the battery shaft	Insert battery.
	Battery not charged	No signal tone when inserting the battery	Charge battery in charger.
	Battery defective	No signal tone when inserting the battery	Have the manufacturer repair the battery.
	Product in locked OFF position	Safety catch is in OFF position	Turn the safety catch into the ON position.
	Product defective	Product not running	Have the manufacturer repair the product.
Product becomes too hot	Excessive use	Product heating	Follow operating instructions (operating mode).
	Reprocessing/maintenance carried out incorrectly	Product heating	Follow operating instructions (preparation, care). Preventative: lubricate product before every sterilization.
	Fall damage, product defective	Product heating	Have the manufacturer repair the product.
	Use of a defective attachment	Product heating	Have the manufacturer repair the attachment.
	Blunt tool	Heating of tool and product	Change tool.
Insufficient power	Product defective	Insufficient product performance	Follow operating instructions (preparation, care). Preventative: lubricate product before every sterilization. Have the manufacturer repair the product.
		Large amount of heat generated after a short time	Follow operating instructions (operating mode). Have the manufacturer repair the product
	Product is operating in left rotation mode	Serrated tool run in counterclockwise mode	Operate serrated tool in right rotation.
	Use of a defective attachment	Attachment very hot	Have the manufacturer repair the attachment.
	Blunt tool	Tool blades worn out	Change tool.

Issue	Cause	Detection	Remedy
Loud running noise	Defective product drive/ball bearings	Loud, perceptible noise during operation	Follow operating instructions (preparation, care). Preventative: lubricate product before every sterilization. Have the manufacturer repair the product
	Defective attachment drive/ball bearings	Loud, perceptible noise during operation	Follow operating instructions (preparation, care). Preventative: lubricate attachment before every sterilization. Have the attachment repaired by the manufacturer
Lid cannot be attached/removed	Lid not compatible	Lid does not click into place	Use the compatible lid for GA340/GA341.
	Lid deformed/defective	Lid is difficult or cannot be mounted/removed	Have the manufacturer repair the lid.
	Connection to product deformed/defective	Lid is difficult or cannot be mounted/removed	Have the manufacturer repair the product.
Battery cannot be inserted/removed	Battery not compatible	Battery cannot be completely inserted into the battery shaft	Use the compatible battery for GA340/GA341.
	Battery deformed/defective	Battery is difficult or cannot be mounted/removed	Have the manufacturer repair the battery.
	Battery shaft to product deformed/defective	Battery is difficult or cannot be mounted/removed	Have the manufacturer repair the product
Sterile funnel or battery removal device not attachable	Sterile funnel or battery removal device not compatible	Sterile funnel or battery removal device not attachable to the battery shaft	Use compatible sterile funnel or battery removal device for GA340/GA341.
	Sterile funnel or battery removal device deformed/defective	Sterile funnel or battery removal device difficult or not attachable	Change sterile funnel or battery device.
	Battery shaft to product deformed/defective	Sterile funnel or battery removal device difficult or not attachable	Have the manufacturer repair the product.
Speed control knob not working	Product in locked OFF position	Safety catch is in OFF position	Turn the safety catch into the ON position.
	Speed control knob sticks/is defective	Speed control knob not working	Have the manufacturer repair the product.
Dermatome blade not moving	Defective gearings	Product runs loudly	Have the manufacturer repair the product
Insufficient cutting performance of dermatome blade	Dermatome blade blunted	Worn cutting edge of dermatome blade	Replace dermatome blade
	Adjusted dermatome	Product has insufficient power	Have the manufacturer repair the product
Cutting depth uneven/deviating from setting	Skin grafts uneven / too thick / too thin	Mechanical influence (e.g. fall on floor)	Have the manufacturer repair the product

Aesculap® Acculan 4

Dermatome GA340 / Dermatome 0.1 mm GA341

10. Technical Service

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- ▶ **Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.**
- ▶ **Do not modify the product.**

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ For service and repairs, please contact your national B. Braun/Aesculapagency.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95 -1601
Fax: +49 (7461) 14 -939
E-Mail: ats@aesculap.de
Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042 USA
Aesculap Repair Hotline
Phone: +1 (800) 214 -3392
Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

11. Accessories/Spare parts

Art. no.	Designation
GA059	STERILIT-Power-Systems drip lubricator (50 ml)
GA346	Li-Ion battery short
GA340345	Nut
GA643417	Flap rod
GA670210	Flap
GA675	Lid
GA678	Sterile funnel
GA679	Battery removal device
GB228R	Dermatome blade
GB256R	ECCOS set for dermatome
GB487R	ECCOS fixation for lid
GB488R	ECCOS fixation for sterile funnel
GB489R	ECCOS fixation for battery removal device
GB498R	ECCOS fixation for straight machine
GB600	STERILIT Power Systems oil spray
TA014548	Instructions for use of Dermatome GA340/GA341 (A4 for ring binder)
TA014549	Instructions for use of Dermatome GA340/GA341 (flyer)

12. Technical data

12.1 Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
GA340	Dermatome	Ila
GA341	0.1 mm Dermatome	Ila

12.2 Performance data, information about standards

Max. power	250 W
Max. oscillation frequency	6 500 min ⁻¹
Blade travel	approx. 3.1 mm
Cutting depth setting	0.2 mm to 1.2 mm (with GA340) 0.1 mm to 1.1 mm (with GA341) Scale gradations of 1/10 mm
Cutting width	max. approx. 78 mm, set by 10 flaps of width 8 mm
Weight (operational)	1.65 kg ± 10 %
Dimensions (L x W x H, ready for operation)	277 mm x 113 mm x 83 mm ± 5 %
Applied part	Type BF
EMC	IEC/DIN EN 60601-1-2
Conforming to standard	IEC/DIN EN 60601-1

The product was subject to 350 preparation cycles as a test by the manufacturer and has passed.

12.3 Operating mode

Operating mode	Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)
	■ 60 second application, 60 second pause ■ 10 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C

12.4 Environmental conditions

	Operation	Storage and transport
Temperature	10 °C to 27 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

13. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

Adhere to national regulations when disposing of or recycling the product, its components and its packaging! The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.) Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

14. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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