Aesculap Spine

GB	Instructions for use
USA	activ C intervertebral disk prosthesis
D	Gebruuchsanweisung
	activ C-Bandscheibenprothese
F	Mode d'emploi
	Prothèse de disque intervertébral activ C
E	Instrucciones de manejo
	Disco intervertebral protésico activ C
I	Istruzioni per l'uso
	Protesi discale activ C
P	Instruções de utilização
	Prótese de disco intervertebral activ C
NL	Gebruiksaanwijzing
	activ C-tussenwervelschijfprothese
CZ	Návod k použití
	activ C-protéza meziobratlové ploténky
PL	Instrukcja użytkowania
	Proteza dysku międzykręgowego activ C
CN	使用说明书
	颈椎人工椎间盘（商品名：activ C）
KR	사용방법및사용시주의사항
	추간채고장보형재 (activ C, 형명개별기재)

B|BRAUN
SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 011995 02/10 Änd.-Nr. 31985/32607

CE
0123

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Símbolo CE, em conformidade com a Directiva 93/42/CEE
CE-certificering conform richtlijn 93/42/EEG
CE-označení podle směrnice 93/42/EHS
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
符合指令 93/42/EEC 的 CE 标志

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Technische wijzigingen voorbehouden
Technické změny vyhrazeny
Zmiany techniczne zastrzeżone
保留技术更改的权利

activ C intervertebral disk prosthesis

CAUTION
Federal law restricts this device to sale by or on order of a physician!

Intended use

The activ C intervertebral disk prosthesis is used for replacing intervertebral disks in the cervical spine (between C3 and C7). It serves to restore the disk height and segmental mobility.

The activ C intervertebral disk prosthesis consists of three components:

- superior prosthesis plate with spikes for anchoring at the vertebral body
- inferior prosthesis plate with central anchoring fin for anchoring at the vertebral body
- Polyethylene inlay

The prosthesis plates and the polyethylene inlay together form a ball joint. The polyethylene inlay is anchored in form-fit in the inferior prosthesis plate. activ C intervertebral disk prostheses are available in six different sizes (XS, S, M, L, XL and XXL) and in up to three different heights (5 mm, 6 mm and 7 mm). The activ C intervertebral disk prostheses are supplied completely pre-assembled.

Materials

The materials used in the implant are listed on the packaging:

Materials
ISODUR®F Cobalt-Chromium forged alloy CoCr29Mo according to ISO 5832-12
PLASMAPORE® surface coating, pure titanium according to ISO 5832-2
Ultra-high molecular low-pressure polyethylene according to ISO5834-2

ISODUR® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for the following indications:

- Symptomatic cervical diskopathy
- Fresh disk prolapse

Contraindications

Do not use in the presence of:

- Fever
- Infection
 - acute
 - systemic
 - in the spine
 - local
- Systemic diseases and metabolic disorders
- Pregnancy
- Osteoporosis, osteochondrosis, severe osteopenia
- Medical conditions that could preclude a successful implantation (e.g. pre-existing kyphosis)
- Spondylarthrosis
- Chronic degenerative spinal stenosis
- Segmental instability
- Spinal deformity
- Chronic spondylotic myelopathy
- Ossification of the anterior longitudinal ligament (ALL) and/or the posterior longitudinal ligament (PLL)
- Distinct osteophytic and sclerotic changes of the vertebral bodies
- Preoperative immobile segment (mobility less than 2° in flexion/extension)
- Preoperative disk height of the segment to be operated on less than 3 mm
- Narrow, hard disk symptoms
- Poor quality of the articular complex
- Facet joint degeneration
- Vertebral body fracture of the cervical spine
- Tumors of the cervical spine
- Pharmaceutical or other drug abuse; alcoholism
- Foreign body sensitivity to the implant materials
- Inadequate patient compliance
- Cases not listed under indications

Side-effects and adverse interactions

- Bent, loose, worn or broken implant components
- Loss of fixation, dislocation and migration
- Infections
- Allergic reactions to the implant materials
- Tissue reaction to implant materials
- Venous thrombosis, lung embolism, and cardiac arrest
- Hematomas and wound healing disorders
- Periarticular calcifications and fusions
- Injuries to
 - nerve roots
 - spinal cord
 - blood vessels
 - organs

Preoperatively existing segmental kyphosis and/or so-called "straight neck" can lead to postoperative kyphotic position of the prosthesis and/or immobility of the segment.

Acute myelopathies can have various reasons. Not all of them can be treated successfully with a prosthesis.

Safety notes

- The implant components were tested and approved for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- Always observe the instructions for use of the individual Aesculap implant components.
- Only combine Aesculap modular implant components with each other.
- Do not use damaged or surgically excised components under any circumstances.
- Implants that have already been used must not be reused.
- The physician in charge decides whether an implant component implanted should be removed. In this case he must be conscious of the potential risks of a further operation as well as of the difficulties of an implant removal.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Damage to the load-bearing structures of the implant may cause the loosening of components, their dislocation and migration, and other severe complications.
- To allow the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals through appropriate examination procedures.
- Only use instruments intended for Aesculap activ C intervertebral disk prostheses.
- Follow the instructions for use of trial implants, handle and reamer guide TAO11996 and insertion instrument and spacer TAO11997.
- Follow the operating manual (no. 031301).

Sterility

- The implant components are supplied pre-assembled in labeled protective packaging.
- The implant components have been sterilized by irradiation (min. dose 25 kGy).
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiry date or whose packaging is damaged.

WARNING

Implants supplied in sterile condition must not be resterilized under any circumstances. If the present instructions for use do not offer explicit information in this respect, the product user carries all responsibility for any resterilization!

WARNING

- Danger to the patient and possible loss of implant functionality due to resterilization!
- Never resterilize and reuse surface-coated implants (e.g. PLASMAPORE®), implants made of ceramics or polyethylene, or implants containing plastic components.
 - Never resterilize used, explanted or damaged implants.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating team are aware of information concerning the operating technique and the range of implants and associated instruments; this information is complete and ready to hand.
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific publications by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.
- Select implant components according to the indication, preoperative planning and the bone situation found intraoperatively

The operating procedure has been explained to the patient, whose understanding of the following information has been documented:

- The patient is aware of the risks involved in neurosurgery, general surgery, orthopedic surgery and general anesthesia.
- The patient has been informed of the advantages and drawbacks of an activ C intervertebral disk prosthesis, and has been made aware of possible alternative therapeutic methods.
- Excessive load, wear or infection can lead to malfunctions of the activ C intervertebral disk prosthesis.
- The life-span of the implant depends on the patient's body weight.
- The implant components must not be subjected to overload due to extreme strains or through work-related or sports activities.
- Implant failure may necessitate revision surgery.

- If revision surgery is carried out, it may not be possible to preserve the segmental movement.
- The patient must undergo regular check-ups of the implant components, performed by a physician.

The height indicated on the implant packaging is the total posterior implant height.

Implantation of activ C intervertebral disk prostheses requires the following steps:

- Surgical access
- Midline marking
- Disectomy and preparation of the vertebral body end plates
- Application of the distractor and distraction
- Verification of implant size
- Fin bed reaming
- Insertion of the activ C intervertebral disk prosthesis

Surgical access

WARNING

Risk of damage to blood vessels, nerves and tissue caused by sharp instruments!

- Work with extreme care throughout the operation.

- Establish the surgical approach to the affected segment.

Midline marking

WARNING

Mispositioning of the activ C intervertebral disk prosthesis caused by inaccurate marking of the midline!

- Always mark the midline under X-ray control.

- Determine the center of the intervertebral disk under X-ray control, using the midline marker.
- Apply markings on the superior/ inferior vertebral body.

Disectomy and preparation of the vertebral body end plates

WARNING

Insufficient anchoring stability of the prosthesis plates due to incomplete removal of cartilage!

- Remove all cartilage from the vertebral body end plates.

WARNING

Increased risk of migration due to over-preparation of the vertebral body endplates!

- Work carefully when preparing the vertebral body end plates.
- Do not remove the vertebral body end plates.

- Completely remove the disk, using standard instruments.
- Resect cartilage from the vertebral body end plates.

Application of the distractor and distraction

WARNING

Compression of the spinal canal and of other posterior elements caused by the distraction forceps inserted too deeply!

- Always use X-ray control when inserting the distraction forceps.

WARNING

Irritation of the facet joints/ strain on spinal cord and nerve roots due to excessive distraction!

- Extend the distraction forceps slowly and carefully under X-ray control.
- Observe the tactile response.
- Use the heights of adjacent (healthy) disks as a measure for the degree of distraction.

WARNING

Risk of loosening/tear-out of the distraction screws during distraction with the distractor and the distraction screws!

- Use the distraction forceps for the distraction step.
- Do not distract with the distractor and the distraction screws.
- Lock the distractor when it is extended.

- Place the distractor.
- Under X-ray control, insert the closed distraction forceps into the intervertebral disk space and extend the distraction forceps.
- Lock the distractor.
- Close the distraction forceps and remove it from the intervertebral disk space.

Verification of implant size

WARNING

Compression of the spinal canal and of other posterior elements caused by trial implants inserted too deeply!

- Use the adjusting wheel to set the depth stop of the trial implant to minimum insertion depth.
- Always use X-ray control when inserting the trial implants in the intervertebral disk space.

WARNING

Increased risk of migration caused by selecting a prosthesis plate of insufficient size!

- Select the size of the trial implant/implant to achieve maximum coverage of the vertebral body end plates.
- On the X-ray image, in AP and lateral view, check to ensure that the size, position and alignment of the trial implant are correct.

WARNING

Irritation of the facet joints/strain on spinal cord and nerve roots due to excessive/insufficient distraction!

- Select the correct trial implant/implant size.
- Inspect the height of the trial implant in the lateral-view X-ray image.

- Mount the trial implant of the selected size on the handle, observing the "CRANIAL" and "CAUDAL" markings on the trial implant/handle.
- Use the adjusting wheel to set the depth stop of the trial implant to minimum insertion depth.
- Under X-ray control, insert the trial implant with the handle in the intervertebral disk space.
- Align the markings of the trial implant center to the midline marking.
- If necessary, correct the depth stop position using the adjusting wheel until the intended position is reached. Ensure that depth stop abuts the anterior edge of the vertebral body.
- On the intraoperative X-ray image, in AP and lateral view, check the size/height and position of the trial implant.
- Dismount the handle from the trial implant.

Fin bed reaming

WARNING

Mispositioning of the activ C intervertebral disk prosthesis caused by disregard of the midline!

- Check the position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- Ensure secure positioning of the trial implant during the reaming procedure.

WARNING

Risk of damage to blood vessels, nerves and tissue caused by sharp reamer!

- Carefully insert the handpiece with the coupled reamer.
- Do not activate the motor before the handpiece with the coupled reamer is properly positioned in the appropriate guide hole and the respective guide slot of the reamer guide.

WARNING

Risk of damage to blood vessels, nerves and spinal cord caused by reamer applied too deeply!

- For reaming the holes, apply the reamer only in combination with a correctly selected and positioned trial implant and with the reamer guide coupled.
- Ensure that depth stop of the trial implant abuts the anterior edge of the vertebral body.
- Only use activ C instruments.
- Use instruments according to their intended use.
- Ream the holes under lateral X-ray control.

- Mount the reamer guide on the trial implant, observing the "CRANIAL" and "CAUDAL" markings on the reamer guide.
- Clamp in the activ C reamer in the handpiece.
- Thread the handpiece with the coupled reamer into the appropriate guide hole and the first guide slot of the reamer guide.
- Ream the first hole under lateral X-ray control.
The reaming depth will be limited by a stop at the reamer guide.
- Thread the handpiece with the coupled reamer into the appropriate guide hole and the second guide slot of the reamer guide.
- Ream the second hole under lateral X-ray control.
The reaming depth will be limited by a stop at the reamer guide.
- Remove the handpiece with the coupled reamer, the reamer guide and the trial implant from the intervertebral disk space.
- Clean the fin bed.

Insertion of the activ C intervertebral disk prosthesis

WARNING

Damage to vertebral body end plates when impelling the activ C intervertebral disk prosthesis!

- Apply proper care when impelling the activ C intervertebral disk prosthesis.
- Do not remove the vertebral body end plates.

WARNING

Compression of the spinal canal and of other posterior elements caused by the activ C intervertebral disk prosthesis inserted too deeply!

- Prior to insertion, adjust the activ C intervertebral disk prosthesis to minimum insertion depth.
- Always maintain X-ray control while inserting the activ C intervertebral disk prosthesis in the intervertebral disk space.
- On the X-ray image, in AP and lateral view, check for correct position/alignment of the activ C intervertebral disk prosthesis.

CAUTION

Damage to the activ C intervertebral disk prosthesis due to improper handling!

- Mount the implant on the insertion instrument immediately after removing it from the implant packaging.
- Handle with caution the friction surfaces of the prosthesis plates/polyethylene inlay.

- Put out the implant ready for use.
- Connect the spacer corresponding to the selected implant height to the insertion instrument, observing the "CRANIAL" and "CAUDAL" markings on the spacer/insertion instrument.
- Use the adjusting wheel to set the depth stop of the insertion instrument to minimum insertion depth.
- Mount the intervertebral disk implant on the insertion instrument immediately after removing it from the implant packaging, observing the "CRANIAL" and "CAUDAL" markings.
The superior prosthesis plate must abut at the part of the insertion instrument that is marked "CRANIAL", while the inferior prosthesis plate must be at the "CAUDAL" part.
- Turn the clamping sleeve to clamp in the intervertebral disk prosthesis on the instrument.
- Remove the disk prosthesis from its packaging, holding it with the insertion instrument.
- Under X-ray control, carefully insert the disk prosthesis into the disk space.
When doing this, the superior prosthesis plate must be oriented towards cranial, the inferior prosthesis plate towards caudal.
- If necessary, correct the depth stop position using the adjusting wheel until the intended position is reached.
- Release the clamping sleeve on the insertion instrument and remove the insertion instrument.
- On the intraoperative X-ray image, in AP and lateral view, check the size/height and position of the intervertebral disk implant.

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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