- (GB)** Instructions for use/Technical description
activ C-insertion instrument
- (USA)** Gebrauchsanweisung/Technische Beschreibung
activ C-Einsetzinstrument
- (D)** Mode d'emploi/Description technique
Instrument d'insertion activ C
- (F)** Instrucciones de manejo/Descripción técnica
Instrumento de inserción activ C
- (E)** Istruzioni per l'uso/Descrizione tecnica
Strumento inseritore activ C
- (I)** Instruções de utilização/Descrição técnica
Instrumento introdutor activ C
- (P)** Gebruiksaanwijzing/Technische beschrijving
activ C-inbrenginstrument
- (NL)** Bruksanvisning/Teknisk beskrivning
activ C-insättningsinstrument
- (S)** Инструкция по применению/Техническое описание
Инструмент для установки актив С
- (RUS)** Návod k použití/Technický popis
Zaváděcí nástroj aktiv C
- (CZ)** Instrukcja użytkowania/Opis techniczny
Instrument do wprowadzania aktiv C
- (PL)** Návod na používanie/Technický opis
Zavádzací nástroj aktiv C
- (SK)** Kullanım Kılavuzu/Teknik açıklama
Yerleştirme ekipmanı aktiv C
- (TR)**


B | BRAUN
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 011997 09/15 V6 Änd.-Nr. 52765

CE 0123 – DIR 93/42/EEC

Aesculap®

activ C-insertion instrument

Legend

- A Insertion instrument and accessories**
- 1 Key (for insertion instrument)
 - 2 Insertion instrument
 - 3 Repositioning instrument
 - 4 Irrigation connector
 - 5 Button
 - 6 Adjusting wheel
 - 7 Clamping sleeve
 - 8 Connector piece
 - 9 Spacer
- B Revision instrument**

Symbols on product and packages

Caution, general warning symbol

Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

Intended use

The insertion instrument is used for picking up the activ C intervertebral disc prosthesis from its sterile packaging, inserting it in the prepared intervertebral space, and carrying out position corrections of the prosthesis.

Available sizes

Art. no.	Designation
FW857R	Insertion instrument with irrigation connector H6
FW863R	Spacer for implant height 5 mm
FW864R	Spacer for implant height 6 mm
FW866R	Insertion instrument with irrigation connector H5
FW867R	Repositioning instrument
FW868R	Revision instrument
FW945R	Key for insertion instrument

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).

Safe operation

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

- Follow OP-Manual no. Q31302 and instructions for use "activ C disk prosthesis", see TA011995.
- Mounting the insertion instrument, see Assembling.

Coupling the implant

Damage to the implant due to improper handling!

- Follow the instructions for use of activ C disk prosthesis TA011995.
- Mount the implant on the insertion instrument immediately after removing it from the implant packaging.
- Handle the friction surfaces of the prosthesis plates and of the polyethylene inlay with care.

- Open the implant packaging.

- Note*
- Do not remove the implant from the packaging!
- Mount insertion instrument **2** on the disc prosthesis with the latter still in its packaging, observing the "CRANIAL" and "CAUDAL" markings. The superior prosthesis plate (the one with the spikes) must abut the part of insertion instrument **2** marked "CRANIAL", while the inferior prosthesis plate (with the fin) has to abut the part of insertion instrument **2** marked "CAUDAL".
 - Turn clamping sleeve **7** clockwise to fixate the intervertebral disc prosthesis on insertion instrument **2**.
 - Remove the disc prosthesis from its packaging, holding it with the insertion instrument **2**.

Inserting the implant

Damage to vertebral body end plates when the intervertebral disk prosthesis is tapped in!

- Apply proper care when tapping in the intervertebral disk prosthesis.

- Carefully, under X-ray control, insert the disk prosthesis into the intervertebral disk space. The superior prosthesis plate (plate with the spikes) must be oriented cranially, the inferior prosthesis plate (the plate with the fin) caudally.
- If necessary impel the implant with light hammer taps on the proximal end of insertion instrument **2**.

Compression of the spinal canal and of other posterior elements caused by the intervertebral disk prosthesis inserted too deeply!

- Prior to inserting the intervertebral disk prosthesis, adjust the depth stop to minimum insertion depth.
- Introduce the disk prosthesis into the intervertebral space under X-ray control.

- If necessary, correct the depth stop position using the adjusting wheel **6** until the intended position is reached.
- Should repositioning become necessary, screw repositioning instrument **3** onto the proximal end of insertion instrument **2** and carefully, with a slotted hammer, move the implant to its intended position.

Decoupling the implant

- Note*
- To avoid damage to the implant, do not couple an insertion instrument that already has been decoupled once.
- Loosen clamping sleeve **7** on insertion instrument **2** by turning it counterclockwise, and remove insertion instrument **2**.
 - If clamping sleeve **7** is difficult to loosen, insert key **1** in one of the lateral drill holes of clamping sleeve **7** and loosen clamping sleeve **7** in this way.
 - On the intraoperative X-ray image, in AP and lateral view, check the size, height and position of the intervertebral disk implant.

Disassembling

- Uncouple spacer **9** from insertion instrument **2** by pressing "push to clean" button **5** at the handle of insertion instrument **2** and simultaneously pulling it out, axially, towards the working end.
- If applicable, unscrew repositioning instrument **3** from insertion instrument **2** by turning it counterclockwise.
- Before reprocessing, ensure that the clamping sleeve **7** is sitting loosely on the instrument shaft.

Assembling

- Note*
- Account for the size of the implant when selecting the insertion instrument!
- Remove irrigation connector **4** from the handle end of insertion instrument **2**.
 - Verify that clamping sleeve **7** on insertion instrument **2** has been loosened.
 - Turn clamping sleeve **7** counterclockwise until it is no longer in connection with connecting piece **8** of insertion instrument **2**.
 - Connect the spacer **9** corresponding to the intraoperatively selected implant height to the insertion instrument **2**.
 - Insert spacer **9** in the hole of the connecting piece, observing the "CRANIAL" and "CAUDAL" markings on the spacer **9** and insertion instrument **2**. At the same time keep button **5**, "push to clean", at the handle of insertion instrument **2** pressed down.
 - Insert spacer **9** down to the stop in insertion instrument **2**.
 - Release "push to clean" button **5**. By axially pulling at spacer **9**, check that the spacer is properly locked in insertion instrument **2**.
 - Turn adjusting wheel **6** counterclockwise to set the depth stop to minimum insertion depth.

Validated reprocessing procedure

General safety instructions

- Note*
- Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.
- Note*
- For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

- Note*
- Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.
- Note*
- Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.
- The recommended chemistry was used for validation.

- Note*
- If there is no final sterilization, then a virucidal disinfectant must be used.

- Note*
- For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>
- The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

- Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixing pre-cleaning temperatures >45 °C nor fixing disinfecting agents (active ingredient: aldehydes/alcohols) should be used.
- Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading of the laser marking becoming unreadable visually or by machine for stainless steel.
- Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.
- Additional drying, if necessary.
- Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:
- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
 - Material damage such as corrosion, cracks, fracturing, premature aging or swelling.

- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- ▶ Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- ▶ Dismantle the product prior to cleaning, see Disassembling.

Note
During cleaning the clamping sleeve 7 must sit loosely on the instrument shaft.

Note
To ensure proper cleaning of spacers 9, the PEEK forks must sit loosely (so that they can be turned in both directions) on the thread.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for the material in question e.g., aluminum, plastics, high-grade steel,
 - do not attack softeners (e.g. in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 60 °C.

- ▶ Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfection.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FW867R ■ FW945R	■ Cleaning brush: 250 mm/Ø 3.7 mm, e.g. GK469200 ■ 20 ml disposable syringe ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-section: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ FW868R	■ Cleaning brush: 250 mm/Ø 3.7 mm, e.g. GK469200 ■ 20 ml disposable syringe ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ FW867R ■ FW945R	■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection and sub-section: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW868R	■ Cleaning brush: 250 mm/Ø 3.7 mm, e.g. GK469200 ■ 20 ml disposable syringe ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Disinfection	RT (cold)	15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	–	FD–W	–
V	Drying	RT	–	–	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Disinfection	RT (cold)	15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	–	FD–W	–
V	Drying	RT	–	–	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II**
- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - ▶ Drain any remaining water fully.
- Phase III**
- ▶ Fully immerse the product in the disinfectant solution.
 - ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.
- Phase IV**
- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
 - ▶ Rinse lumens with an appropriate disposable syringe at least five times.
 - ▶ Drain any remaining water fully.
- Phase V**
- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfection

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	Concentrate, alkaline: – pH = 13 – <5 % anionic surfactant 0.5 % working solution – pH = 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disin- fecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concen- trate, pH ~ 9*
II	Rinsing	RT (cold)	1	–	D–W	–

D–W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I**
- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
 - ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
 - ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II**
- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	Concentrate, alkaline: – pH = 13 – <5 % anionic surfactant 0.5 % working solution – pH = 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disin- fecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Assemble dismountable products, see Assembling.
- ▶ Check for compatibility with associated products.

Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Store products with ratchet locks fully opened or locked no further than in the first notch.
- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note
The product may only be sterilized when dismantled.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Disassemble the instrument
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{–6}, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

Risk of injury and/or malfunction!
► Do not modify the product.

► For service and repairs, please contact your national B. Braun/Aesculap agency.
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Aesculap Implant Systems LLC
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Disposal

► Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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