- GB** Instructions for use
- USA** activ L Intervertebral disc prosthesis
- D** **Gebrauchsanweisung**  
activ L-Bandscheibenprothese
- F** **Mode d'emploi**  
Prothèse de disque intervertébral activ L
- E** **Instrucciones de manejo**  
Prótesis de disco intervertebral activ L
- I** **Istruzioni per l'uso**  
Protesi discale activ L
- P** **Instruções de utilização**  
Protéze de disco intervertebral activ L
- NL** **Gebruiksaanwijzing**  
activ L-tussenwervelschijfprothese
- S** **Bruksanvisning**  
activ L-diskprotes
- RUS** **Инструкция по применению**  
Протез межпозвоночного диска activ L
- CZ** **Návod k použití**  
Protéza meziobratlové ploténky activ L
- PL** **Instrukcja użytkowania**  
Proteza krążka międzykręgowego activ L
- SK** **Návod na používanie**  
Náhrada medzistavcovej platničky activ L
- TR** **Kullanım Kılavuzu**  
activ L-Protez disk
- KR** 사용방법 및 사용시 주의사항  
추간체고정보형제 (activ L, 형명개별기재)

**B | BRAUN**  
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

TA-Nr. 011430 06/14 V6 Änd.-Nr. 50611

**CE** 0123 - DIR 93/42/EEC

Technical alterations reserved

# Aesculap®

## activ L Intervertebral disc prosthesis

### Intended use

The activ L intervertebral disc prosthesis is used for replacing intervertebral discs in the lumbar spine. They restore the disc height and the segmental mobility.

They comprise:

- inferior and superior prosthesis plates for anchoring in the vertebral bodies
  - in sizes S, M, L and XL
  - with a central anchoring fin (can only be fitted with an anterior approach) or without central anchoring fins (can be fitted with an anterior or lateral approach, at approximately 45°)
- a polyethylene inlay, which is inserted into the inferior prosthesis plate
  - for total heights of 8.5 mm, 10 mm, 12 mm and 14 mm
  - to be combined with any plate size

The prosthesis plates and the polyethylene inlay together form a ball joint.

The operating surgeon chooses, according to the anatomic conditions, between using prosthesis plates with anchoring fin and prosthesis plates without such a fin.

### Materials

The materials used in the implant are listed on the packaging:

- ISODUR®: Cobalt-chromium forged alloy CoCrMo according to ISO 5832-12
- PLASMAPORE®:  $\mu$ -CaP with calcium phosphate surface
- Ultra-high molecular low-pressure polyethylene according to ISO 5834-2

ISODUR® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttingen / Germany.

### Indications

Lumbar discopathy

Defined in detail as mono- or bisegmental intervertebral disc generation between L3 and S1. There is no significant degeneration in the facet joints. There is no degeneration to the discs adjacent to a segment to be treated and there is no, or only slight, posterior instability. The preoperative intervertebral disk height is 3 mm to 4 mm or more.

### Contraindications

Do not use in the presence of:

- Fever
- Infection: Infection: systemic, in the spinal column, local
- Systemic diseases and metabolic disorders
- Pregnancy
- Osteoporosis, osteochondrosis, severe osteopenia
- Medical or surgical conditions that could preclude the success of the implantation
- Spinal stenosis, radiculopathy
- Degenerated facets
- Increased segmental instability
- Fractured vertebral body or bodies
- Spinal deformity
- Spondylolisthesis over 25 %
- Pharmaceutical or other drug abuse; alcoholism
- Obese or overweight patients
- Foreign body sensitivity to the implant materials
- Inadequate patient compliance
- All cases that are not listed under indications

### Side effects and interactions

- Bent, loose, worn or broken implant components
- Loss of fixation, dislocation and migration
- infections.
- Allergic reactions to implant materials
- Tissue reaction to implant materials
- Venous thrombosis, lung embolism, cardiac arrest
- Hematomas and wound healing disorders
- Periarticular calcifications and fusions
- Trauma to: nerve roots, spinal cord, blood vessels, organs

### Safety notes

#### Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

#### CAUTION

Federal law restricts this device to purchase by, or on instruction by a physician!

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- The attending physician shall make any decision with regard to the removal of implant components that have been used, taking into account the potential risks of another operation and the difficulties involved in the removal of an implant.
- Damage to the load-bearing structures of the implant can lead to loosening of components, dislocation, migration, and other severe complications.

- To ensure the earliest possible detection of such factors favoring implant malfunction, the activ L intervertebral disc prosthesis must be inspected postoperatively, on a regular basis, by suitable procedures.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

### Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are sterilized by irradiation.
- ▶ Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- ▶ Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- ▶ Do not use implant components that are past their expiration date or whose packaging is damaged.

#### Damage to implants caused by processing and reesterilization!

- ▶ Do not reprocess or reesterilize the implants.

### Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- Implantation instruments including the special Aesculap implant system instruments are complete and in working condition
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The patient is aware of the risks involved in neurosurgery, general surgery, orthopedic surgery and general anesthesia.
- The patient has been informed of the advantages and drawbacks of an activ L intervertebral disc prosthesis, and has been made aware of possible alternative therapeutic methods.
- Excessive load, wear or infection can lead to malfunctions of the activ L-intervertebral disc prosthesis.
- The life-span of the implant depends on the patient's body weight.
- The activ L intervertebral disc prosthesis must not be subjected to excessive strains through extreme loads, hard physical labour or sports activities.
- Implant failure may necessitate revision surgery.
- If revision surgery is carried out, it may not be possible to preserve the segmental movement.
- The patient must undergo regular medical follow-up examinations of the activ L-intervertebral disc prosthesis

The implantation involves the following steps:

- ▶ Only use instruments for activ L intervertebral disc prostheses from Aesculap.
- ▶ Observe the instructions for use of insertion instrument TA011458, the distractor TA011540, the chisel/chisel guide TA013425 and the distraction forceps revision instrument TA011904.
- ▶ Follow the O27901 operating manual.
- ▶ Select implant components according to the indication, preoperative planning and the bone situation found intraoperatively.

#### Note

The height printed on the packaging of the polyethylene inlay is equivalent to the total posterior height of the disc prosthesis, i.e. of the inferior prosthesis plate, the polyethylene inlay and the superior prosthesis plate. The final height of the polyethylene inlay is determined intraoperatively, in most cases.

#### Preoperative planning

- Risk of incorrect size selection due to using an X-ray template of the wrong scale!
- ▶ Make certain to choose an X-ray template of the correct scale.

- ▶ Carry out preoperative planning, using X-ray templates.

- Increased risk of migration caused by selecting a prosthesis plate that is too small!
- ▶ The prosthesis plate must completely cover the end plate of the vertebral body.

- ▶ Select the size of the prosthesis plate so that maximum coverage of the vertebral bodies is achieved.

#### Surgical access

- Risk of damage to blood vessels and nerves caused by sharp instruments!
- ▶ Operate with extreme caution to avoid damage to blood vessels, nerves or tissue.
- ▶ Have a vascular surgeon assisting or on call.

- ▶ Establish the surgical access to the affected segment.

#### Note

Prosthesis plates with an anchoring fin can only be implanted through anterior access. Prosthesis plates without an anchoring fin can also be implanted through lateral access (approx. 45°).

#### Marking the midline

- Mispositioning of the vertebral disc prosthesis caused by inaccurate marking of the midline!
- ▶ Always mark the midline under X-ray control.
- ▶ Make certain to use a midline marker of the correct size when installing the lateral access.

- ▶ Determine the center of the vertebral disc, using the midline marker, under X-ray control.
- ▶ Apply markings on the superior and inferior vertebral bodies.
- ▶ When operating through an anterior access, use the "anterior midline marker".
- ▶ When operating through a lateral access, use the lateral midline marker of the appropriate size.

## Discectomy and preparation of the vertebral body end plates

WARNING

**Insufficient anchoring stability due to incomplete removal of cartilage!**  
▶ Remove the cartilage layer from the vertebral body end plates.

WARNING

**Increased risk of migration due to over-preparation of the vertebral body end-plates!**  
▶ Avoid any aggressive preparation of the vertebral body end plates.  
▶ Do not remove the vertebral body end plates completely.

- ▶ Completely remove the intervertebral disc, using standard instruments.
- ▶ Remove all cartilage from the vertebral body end plates.

### Verifying the size and distraction

- ▶ Mount the trial plates of the selected size on the distractor in such a way that the marking for the lateral/anterior access corresponds to the marking on the working tip of the distractor.

WARNING

**Compression of the spinal canal and of other posterior elements caused by the distractor inserted too deeply!**  
▶ Introduce the distractor into the intervertebral space under X-ray control.  
▶ On the X-ray image, under AP and lateral irradiation, check to ensure that the size, position and alignment of the trial plates are correct.

- ▶ Cautiously introduce the undistracted distractor into the intervertebral space, making sure that it is positioned centrally.
- ▶ Check the size of the vertebral body end plates. A margin of 1 mm to 2 mm width is allowed between the lateral and posterior edges of the trial plates and the edge of the vertebral body.

WARNING

**Irritation of the facet joints and strain on dura and nerve roots due to excessive or insufficient distraction!**  
▶ Distract only to such an extent that the instrument firmly fits in the intervertebral space.  
▶ On the X-ray image, under AP and lateral irradiation, check to ensure that the size, position and alignment of the trial plates are correct.

- ▶ Open the distractor until it is firmly seated in the intervertebral space.
- ▶ Read the distraction height on the distractor scale.
- ▶ If the mark is between to height values on the scale, choose the lower height.
- ▶ Release the distraction and remove the distractor from the intervertebral space.
- ▶ If necessary, use a small spacer for keeping the vertebral bodies apart.

### Chiseling out the anchoring for prosthesis plates with anchoring fin

The chiseling procedure will determine the central position and the alignment of intervertebral disc prosthesis in the intervertebral disc plane.

The chisel depth is defined by the bedstop of the chisel on the chisel guide.

- ▶ Mount the chisel guide of the selected height on the handle.
- ▶ Turn the adjusting wheel of the depth bedstop to adjust it to the minimum insertion depth.

WARNING

**Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!**  
▶ Align the chisel guide with the midline mark.  
▶ Make certain that the position of the chisel guide does not change during the chiseling procedure.

WARNING

**Compression of the spinal canal and of other posterior elements caused by the chisel guide inserted too deeply!**  
▶ Prior to insertion, adjust the chisel guide to the minimum insertions depth.  
▶ Introduce the chisel guide into the intervertebral space under X-ray control.

- ▶ Insert the chisel guide into the intervertebral space under X-ray control. Align the groove on the chisel guide with the midline marks.
- ▶ If necessary, readjust the depth bedstop until the posterior edge of the chisel guide is at the height of the posterior edge of the vertebral body, or 1 mm to 2 mm in front of it.
- ▶ Check to ensure that the chisel guide is firmly seated. Choose a different height or angle variant, if necessary.
- ▶ Dismount the handle from the chisel guide.

WARNING

**Risk of injury caused by sharp chisel!**  
▶ Apply the chisel with extreme caution to avoid damage to critical blood vessels, nerves or other tissue, as well as injury to yourself.

- ▶ Carefully slide the chisel over the shaft of the chisel guide, until it engages in the guide groove.
- ▶ Drive the chisel into the vertebral body, but only down to the bedstop on the chisel guide, not any further.
- ▶ Carefully knock the chisel out of the vertebral body, using a slotted hammer.
- ▶ Remove the chisel guide from the intervertebral space.
- ▶ If necessary, use a narrow spacer for keeping the vertebral bodies apart.

### Inserting the activ L intervertebral disc prosthesis

- ▶ Select the implant components of the size that was determined intraoperatively.

CAUTION

**Damage to the friction surfaces due to improper handling!**  
▶ Handle with caution the friction surfaces of the prosthesis plates and the polyethylene inlay.

- ▶ Insert the polyethylene inlay into the depression in the inferior prosthesis plate, in such a way that the two thinner protrusions of the polyethylene inlay engage in the undercut of the prosthesis plate. The polyethylene inlay can be mounted in only one position, without having to exert any force.
- ▶ Mount the complete intervertebral disc assembly on the insertion instrument, according to the surgical access, observing the markings.  
The inferior prosthesis plate must be at the part of the insertion instrument clamp that is marked "CAUDAL", while the superior prosthesis plate must be at the "CRANIAL" part
- ▶ Turn the clamping sleeve to clamp in the intervertebral disc prosthesis on the instrument.

WARNING

**Compression of the spinal canal and of other posterior elements caused by the intervertebral disk prosthesis inserted too deeply!**  
▶ Introduce the disk prosthesis into the intervertebral space under X-ray control.  
▶ On the X-ray image, under AP and lateral irradiation, check to ensure that the position and alignment of the intervertebral disc prosthesis are correct.

## Note

For prosthesis plates with anchoring fin, the central position is determined by the previously applied slits in the vertebral bodies.

For prosthesis plates without anchoring fin, the intervertebral disc prosthesis has to be aligned with midline marks. The central anchoring tooth on the superior and inferior prosthesis plates must be in line with the midline markings.

- ▶ Carefully introduce the intervertebral disc prosthesis into the intervertebral space.  
When doing this, the superior prosthesis plate must be oriented towards cranial, the inferior prosthesis plate towards caudal.

WARNING

**Damage to vertebral body end plates when the intervertebral disk prosthesis is tapped in!**  
▶ Apply proper care when tapping in the intervertebral disk prosthesis.

WARNING

**Vascular damage caused by protruding implants!**  
▶ Check to make certain that the intervertebral disc prosthesis is positioned centrally.

- ▶ Under X-ray control, drive the intervertebral disc prosthesis into the intervertebral space until the posterior edge of the prosthesis is at the posterior edge of the vertebral body, or 1 mm to 2 mm in front of it. While doing this, check the insertion depth on the lateral X-ray image, for this depth is not limited by the instrument.
- ▶ Release the clamping sleeve on the insertion instrument and remove the insertion instrument.
- ▶ Check the position of the intervertebral disc prosthesis on the intraoperative X-ray image, in AP and lateral view. Further information on Aesculap implant systems is available at any time from B. Braun/Aesculap or from the responsible B. Braun/Aesculap-branch.

## Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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