activ L

GB	Instructions for use
USA	activ L intervertebral disc prosthesis
D	Gebrauchsanweisung
	activ L-Bandscheibenprothese
F	Mode d'emploi
	Prothèse de disque intervertébral activ L
E	Instrucciones de manejo
	Prótesis de disco intervertebral activ L
I	Istruzioni per l'uso
	Protesi discale activ L
P	Instruções de utilização
	Próteses de disco intervertebral activ L
NL	Gebruiksaanwijzing
	activ L-tussenwervelschijfprothese
CZ	Návod k použití
	Protézy meziobratlové ploténky activ L
PL	Instrukcja użytkowania
	Proteza dysku międzykręgowego activ L
SK	Návod na používanie
	Protéza medzistavcovej platničky activ L
CN	使用说明
	activ L 人工椎间盘 椎间盘假体
KR	사용방법및사용시주의사항 추간체고정보형체 (activ L, 형명개별기재)

B | BRAUN

SHARING EXPERTISE

Aesculap AG

Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 011430 02/11 Änd.-Nr. 32568/32607

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Símbolo CE, em conformidade com a Directiva 93/42/CEE
CE-certificering conform richtlijn 93/42/EEG
CE-označení podle směrnice 93/42/EHS
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
Označenie CE podľa smernice 93/42/EHS
CE 标志符合 93/42/EEC 指令

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Technische wijzigingen voorbehouden
Technické změny vyhrazeny
Zmiany techniczne zastrzeżone
Technické zmeny vyhradené
保留技术变更权利

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Intended use

The activ L intervertebral disc prosthesis is used for replacing intervertebral discs in the lumbar spine. They restore the disc height and the segmental mobility.

They comprise:

- inferior and superior prosthesis plates for anchoring in the vertebral bodies
 - in sizes S, M, L and XL
 - with a central anchoring fin (can only be implanted through anterior access)
 - or -
 - without a central anchoring fin (for implantation through anterior or lateral (approx. 45°) access)
- a polyethylene inlay, which is inserted into the inferior prosthesis plate
 - for total heights of 8.5 mm, 10 mm, 12 mm and 14 mm
 - to be combined with any plate size

The prosthesis plates and the polyethylene inlay together form a ball joint.

The operating surgeon chooses, according to the anatomic conditions, between using prosthesis plates with anchoring fin and prosthesis plates without such a fin.

Materials

The materials used in the implant are listed on the packaging:

- ISODUR®: Cobalt-Chromium wrought alloy CoCr29Mo acc. to ISO 5832-12
 - PLASMAPORE®_{U-CaP} surface coating made of pure titanium acc. to ISO 5832-2, with an additional calcium phosphate coating
 - Ultra-high molecular low-pressure polyethylene acc. to ISO 5834-2
- ISODUR® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Lumbar discopathy

Detailed definition: Monosegmental or bisegmental disc degeneration between L3 and S1. No existing significant degeneration of the facet joints. No existing degeneration, or only minor degeneration, of the discs adjacent to the segments to be treated, and no existing posterior instability. Preoperative disc height of at least 3-4 mm.

Contraindications

Do not apply in the presence of:

- Fever
- Infection: acute, systemic, in the spine, local
- Systemic diseases and metabolic disorders
- Pregnancy
- Osteoporosis, osteochondrosis, severe osteopenia
- Any medical condition that could preclude the success of the implantation
- Spinal stenosis, radiculopathy
- Degenerated facets
- Increased segmental instability
- Vertebral body fracture
- Spinal deformity
- Spondylolisthesis above 25 %
- Pharmaceutical or other drug abuse; alcoholism
- Obese or overweight patients
- Foreign body sensitivity to the implant materials
- Inadequate patient compliance
- All cases that are not listed under indications

Side-effects or adverse interactions

- Bent, loose, worn or broken implant components
- Loss of fixation, dislocation and migration
- Infections
- Allergic reactions to implant materials
- Tissue reaction to implant materials
- Venous thrombosis, lung embolism, and cardiac arrest
- Hematomas and wound healing disorders
- Periarticular calcifications and fusions
- Injuries to: nerve roots, spinal cord, blood vessels, organs

Safety information

- The implant components were tested and approved for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the operative procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Under no circumstances may modular implant components from different suppliers be combined.
- Under no circumstances should damaged components or surgically excised components be used.
- Implants that have already been used must not be reused.

- The physician in charge decides whether an implant component implanted should be removed, taking into account the potential risks of another operation and the difficulties involved in the removal of an implant.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keep the patient well informed.
- Damage to the load-bearing structures of the implant may cause the loosening of components, their dislocation and migration, and other severe complications.
- To ensure the earliest possible detection of such factors favoring implant malfunction, the activ L intervertebral disc prosthesis must be inspected postoperatively, on a regular basis, by suitable procedures.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation (min. dose 25 kGy).
- The implant components must not be re-sterilized.
- Store the implant components in their original packaging. Do not remove them from their original and protective packaging until immediately before using them.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiry date or whose packaging is damaged.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bones
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components must be ready to hand
 - Operating conditions must be highly aseptic
 - Implantation instruments including the special Aesculap implant system instruments are complete and in working condition
 - The operating surgeon and operating room team must be thoroughly conversant with the operating technique, as well as the range of implants and instruments to be applied; complete information on these subjects must be readily available at the workplace.
 - The operating surgeon must be familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant articles by medical specialists from professional literature.
 - The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.
- The operative procedure has been explained to the patient, and the latter's understanding of the following information has been documented:
- The patient is aware of the risks involved in neurosurgery, general surgery, orthopedics and general anesthesia.
 - The patient has been informed of the advantages and drawbacks of an activ L intervertebral disc prosthesis, and has been made aware of possible alternative therapeutic methods.
 - Excessive load, wear or infection can lead to malfunctions of the activ L intervertebral disc prosthesis.
 - The life-span of the implant depends on the patient's body weight.
 - The activ L intervertebral disc prosthesis must not be subjected to excessive strains through extreme loads, hard physical labor or sports activities.
 - Implant failure may necessitate revision surgery.
 - If revision surgery is carried out, it may not be possible to preserve the segmental movement.
 - The patient must undergo regular medical check-ups of the activ L intervertebral disc prosthesis.
 - Only use instruments intended for Aesculap activ L intervertebral disc prostheses.

- Observe the instructions for use of insertion instrument TA011458, distractor TA011540 and distraction forceps revision instrument TA011904.
- Observe the OR manual 027901.
- Select implant components according to the indication, preoperative planning and the bone situation found intraoperatively.
- The height printed on the packaging of the polyethylene inlay is equivalent to the total posterior height of the disc prosthesis, i. e. of the inferior prosthesis plate, the polyethylene inlay and the superior prosthesis plate. The final height of the polyethylene inlay is determined intraoperatively, in most cases.

Preoperative planning

WARNING

Risk of incorrect size selection due to using an X-ray template of the wrong scale!

- Make certain to choose an X-ray template of the correct scale.

- Carry out preoperative planning, using X-ray templates.

WARNING

Increased risk of migration caused by selecting a prosthesis plate that is too small!

- The prosthesis plate must completely cover the end plate of the vertebral body.

- Select the size of the prosthesis plate so that maximum coverage of the vertebral bodies is achieved.

Surgical access

WARNING

Risk of damage to blood vessels and nerves caused by sharp instruments!

- Operate with extreme caution to avoid damage to blood vessels, nerves or tissue.
- Have a vascular surgeon assisting or on call.

- Establish the surgical access to the affected segment.

Note

Prosthesis plates with an anchoring fin can only be implanted through anterior access. Prosthesis plates without an anchoring fin can also be implanted through lateral access (approx. 45°).

Marking the midline

WARNING

Mispositioning of the vertebral disc prosthesis caused by inaccurate marking of the midline!

- Always mark the midline under X-ray control.
- Make certain to use a midline marker of the correct size when installing the lateral access.

- Determine the center of the vertebral disc, using the midline marker, under X-ray control.
- Apply markings on the superior and inferior vertebral bodies.
- When operating through an anterior access, use the anterior midline marker.
- When operating through a lateral access, use the lateral midline marker of the appropriate size.

Discectomy and preparation of the vertebral body end plates

WARNING

Insufficient anchoring stability due to incomplete removal of cartilage!

- Remove the cartilage layer from the vertebral body end plates.

WARNING

Increased risk of migration due to over-preparation of the vertebral body endplates!

- Avoid any aggressive preparation of the vertebral body end plates.
- Do not remove the vertebral body end plates completely.

- Completely remove the intervertebral disc, using standard instruments.
- Remove all cartilage from the vertebral body end plates.

Verifying the size and distraction

- Mount the trial plates of the selected size on the distractor in such a way that the marking for the lateral/anterior access corresponds to the marking on the working tip of the distractor.

WARNING

Compression of the spinal canal and of other posterior elements caused by the distractor inserted too deeply!

- Introduce the distractor into the intervertebral space under X-ray control.
- On the X-ray image, under AP and lateral irradiation, check to ensure that the size, position and alignment of the trial plates are correct.

- Cautiously introduce the undistracted distractor into the intervertebral space, making sure that it is positioned centrally.

- Check the size of the vertebral body end plates. A margin of 1 mm to 2 mm width is allowed between the lateral and posterior edges of the trial plates and the edge of the vertebral body.

WARNING

Irritation of the facet joints and strain on dura and nerve roots due to excessive or insufficient distraction!

- Distract only to such an extent that the instrument firmly fits in the intervertebral space.
- On the X-ray image, under AP and lateral irradiation, check to ensure that the size, position and alignment of the trial plates are correct.

- Open the distractor until it is firmly seated in the intervertebral space.
- Read the distraction height on the distractor scale.
- If the mark is between to height values on the scale, choose the lower height.
- Release the distraction and remove the distractor from the intervertebral space.
- If necessary, use a small spacer for keeping the vertebral bodies apart.

Chiseling out the anchoring for prosthesis plates with anchoring fin
The chiseling procedure will determine the central position and the alignment of intervertebral disc prosthesis in the intervertebral disc plane.

The chisel depth is defined by the bedstop of the chisel on the chisel guide.

- Mount the chisel guide of the selected height on the handle.
- Turn the adjusting wheel of the depth bedstop to adjust it to the minimum insertion depth.

WARNING

Mispositioning of the intervertebral disc prosthesis due to failure to align the chisel guide with the midline!

- Align the chisel guide with the midline mark.
- Make certain that the position of the chisel guide does not change during the chiseling procedure.

WARNING

Compression of the spinal canal and of other posterior elements caused by the chisel guide inserted too deeply!

- Prior to insertion, adjust the chisel guide to the minimum insertions depth.
- Introduce the chisel guide into the intervertebral space under X-ray control.

- Insert the chisel guide into the intervertebral space under X-ray control. Align the groove on the chisel guide with the midline marks.

- If necessary, readjust the depth bedstop until the posterior edge of the chisel guide is at the height of the posterior edge of the vertebral body, or 1 mm to 2 mm in front of it.

- Check to ensure that the chisel guide is firmly seated. Choose a different height or angle variant, if necessary.
- Dismount the handle from the chisel guide.

CAUTION

Risk of injury caused by sharp chisel!

- Apply the chisel with extreme caution to avoid damage to critical blood vessels, nerves or other tissue, as well as injury to yourself.

- Carefully slide the chisel over the shaft of the chisel guide, until it engages in the guide groove.

- Drive the chisel into the vertebral body, but only down to the bedstop on the chisel guide, not any further.
- Carefully knock the chisel out of the vertebral body, using a slotted hammer.
- Remove the chisel guide from the intervertebral space.
- If necessary, use a narrow spacer for keeping the vertebral bodies apart.

Inserting the activ L vertebral disc prosthesis

- Select the implant components of the size that was determined intraoperatively.

WARNING

Damage to the friction surfaces due to improper handling!

- Handle with caution the friction surfaces of the prosthesis plates and the polyethylene inlay.

- Insert the polyethylene inlay into the depression in the inferior prosthesis plate, in such a way that the two thinner protrusions of the polyethylene inlay engage in the undercut of the prosthesis plate. The polyethylene inlay can be mounted in only one position, without having to exert any force.
- Mount the complete intervertebral disc assembly on the insertion instrument, according to the surgical access, observing the markings. The inferior prosthesis plate must be at the part of the insertion instrument clamp that is marked "CAUDAL", while the superior prosthesis plate must be at the "CRANIAL" part.
- Turn the clamping sleeve to clamp in the intervertebral disc prosthesis on the instrument.

WARNING

Compression of the spinal canal and of other posterior elements caused by the intervertebral disc prosthesis inserted too deeply!

- Introduce the distractor into the intervertebral space under X-ray control.
- On the X-ray image, under AP and lateral irradiation, check to ensure that the position and alignment of the intervertebral disc prosthesis are correct.

Note

For prosthesis plates with anchoring fin, the central position is determined the previously applied slits in the vertebral bodies.

For prosthesis plates without anchoring fin, the intervertebral disc prosthesis has to be aligned with midline marks. The central anchoring tooth on the superior and inferior prosthesis plates must be in line with the midline markings.

- Carefully introduce the intervertebral disc prosthesis into the intervertebral space. When doing this, the superior prosthesis plate must be oriented towards cranial, and the inferior prosthesis plate towards caudal.

WARNING

Damage to vertebral body end plates when driving in the intervertebral disc prosthesis!

- Apply proper care when driving in the intervertebral disc prosthesis.

WARNING

Vascular damage caused by protruding implants!

- Check to make certain that the intervertebral disc prosthesis is positioned centrally.

- Under X-ray control, drive the intervertebral disc prosthesis into the intervertebral space until the posterior edge of the prosthesis is at the posterior edge of the vertebral body, or 1 mm to 2 mm in front of it. While doing this, check the insertion depth on the lateral X-ray image, for this depth is not limited by the instrument.
- Release the clamping sleeve on the insertion instrument and remove the insertion instrument.
- Check the position of the intervertebral disc prosthesis on the intraoperative X-ray image, in AP and lateral view.

Further information about B. Braun/Aesculap implant systems is available from B. Braun/Aesculap or the B. Braun/Aesculap office responsible.

Handling of explanted devices

- Place explanted device in appropriate biohazard or designated container for transport to decontamination area (Central Services, etc.). Ensure that the appropriate documentation accompanies the explant separately to avoid contamination.

Decontaminate explanted device according to facility procedures or reference the following:

- Ensure that a pair of rubber, chemical resistant gloves, goggles, and a facemask (or face shield) is being worn.
- Place contaminated instrumentation into a container with Metricide or equivalent disinfecting solution.
- Soak the instruments for a minimum of five minutes at room temperature (25 °C).
- Using a soft toothbrush below the surface of the solution removes all contamination from the instruments or follows the "Manufactures Recommended Procedure".
- Continue to soak the instruments for an additional 15 minutes.
- After a total of 20 minutes, remove instrumentation, rinse thoroughly and dry.
- Discard residual solution in the drain and flush thoroughly with water.
- Rinse the container thoroughly and dry.
- Once the explanted device is decontaminated, package it appropriately along with documentation and contact your Aesculap representative for further instructions on return of the device.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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