GB RECYCLING PASSPORT
D RECYCLINGPASS

LAPAROSCOPIC SURGERY

AESCULAP® EinsteinVision® 3.0

EV3.0 Camera Control Unit / EV3.0 Kamerakontrolleinheit PV630

EV3.0 Camera Head 0° / EV3.0 Kamerakopf 0° PV631

EV3.0 Camera Head 30° / EV3.0 Kamerakopf 30° PV632

STATEMENT

We, Aesculap AG

hereby confirm that the products

EV3.0 Camera Control Unit PV630

EV3.0 Camera Head 0° PV631

EV3.0 Camera Head 30° PV632

comply with the following guidelines:

Waste electrical and electronic equipment (WEEE) 2012/19/EC (dated 04/07/2012)

Category 8 product according to Annex 11a

Electrical and Electronic Equipment Act - ElektroG

as well as the statutory requirements of the EEA member states

November, 2017

Aesculap AG, Am Aesculap-Platz, 78532 Tuttlingen / Germany

INFORMATION FOR RE-USE AND HANDLING PV630, PV631, PV632

POS.	MATERIAL	WEIGHT	INFORMATION
	Re-use		
	-	-	-
	Scrap		
010	Circuit board a) Camera motherboard	240 g	Remove: Circuit board > 10 cm ² x 2
020	Circuit board b) Ethernet board	12 g	Remove: Circuit board > 10 cm ²
030	Circuit board c) Back board	12 g	Remove: Circuit board > 10 cm ²
040	Circuit board d) Front board	40 g	Remove: Circuit board > 10 cm ²
050	Circuit board e) Additional board	100 g	Remove: Circuit board > 10 cm ²
060	Circuit board f) Power switch board	10 g	Remove: Circuit board > 10 cm ²
070	Circuit board d) Front board, small	10 g	Remove: Circuit board > 10 cm ²
080	Mains filter	100 g	Remove: Mixed scrap
090	Fan	103 g	Remove: Mixed scrap
100	Front membrane and front board	210 g	Remove: Mixed scrap
110	Power supply	320 g	Remove: Circuit board > 10 cm ²
120	Connector socket	100 g	Remove: Mixed scrap
130	Switch membrane incl. carrier plate	5 g	Remove: Mixed scrap
140	Front membrane with product name	5 g	Remove: Mixed scrap
150	Heat pipe	6 g	Remove: Mixed scrap
160	Camera head with cable and plug	1230 g	Remove: Mixed scrap
	Recycling		
	Metal		
200	Aluminum, pressure die-casting	1150 g	Housing cover
210	Aluminum, pressure die-casting	1850 g	Housing base
220	Aluminum, pressure die-casting	245 g	Front frame
230	Stainless steel	120 g	Power supply cover

INFORMATION FOR RE-USE AND HANDLING PV630, PV631, PV632

POS.	MATERIAL	WEIGHT	INFORMATION
	Recycling		
	Metal		
240	Stainless steel	626 g	Housing spacer plate
250	Stainless steel	409 g	Rear wall
260	Aluminum	180 g	Heat sink, x 2
270	Aluminum	1777 g	Heat sink, large
280	Stainless steel	10 g	Mounting bracket, x 4
290	Stainless steel	90 g	Slider, x 2
300	Stainless steel	94 g	Retaining plate
310	Stainless steel	580 g	Side part, x 2
	Plastic		
400	Polyurethane	5 g	Housing feet
410	PPSU	5 g	Socket frame
420	PA 66 GF20	5 g	Sleeves
430	PA 66-GF25FR	5 g	Ethernet socket holder
440	TGF-W2000-SI [silicone ceramic]	10 g	Heat conductor pads
450	Silicone	150 g	Contoured foam
	Composite		
500	PCV-copper	200 g	Cable
	Total	approx. 11.43 kg	

GB REMOVAL OF POLLUTANTS PV630, PV631, PV632

Device overview

Fig. 1: Camera controller

Fig. 2: Endoscope

GB REMOVAL OF POLLUTANTS PV630, PV631, PV632

Disassembly

Fig. 3: Housing

- Unscrew the screws (x2) on the back
- Lift the housing hood upwards
- Unplug the protective earth cable and remove the cover
- Unscrew the screws (x8) on the side parts and remove the side parts

REMOVAL OF POLLUTANTS PV630, PV631, PV632


Fig. 4: Camera control unit opened

- Disconnect all connectors and remove cables
- Unscrew the nuts (x2) on the power supply bracket and remove the power supply with cover
- Unscrew the nuts (x4) on the power supply and remove the power supply
- Unscrew the screws (x2) on the mains filter and remove the mains filter
- Unscrew the screws (x3) on the retaining plate and remove the retaining plate
- Unscrew the nuts (x3) and spacer bolts (x3) on circuit board e) and remove circuit board e)
- Unscrew the screws (x4) on circuit board c) and remove circuit board c)
- Unscrew the nut on the Ethernet socket mount and remove the Ethernet socket mount
- Unscrew the screws (x2) on circuit board b) and remove circuit board b)
- Unscrew the nuts (x4) on the mounting bracket and remove the mounting bracket (x4)
- Unscrew the screws (x5) and spacer bolts (x6) as well as the spacer bolts (x6) and nuts (x2) on the back on circuit board a) and remove circuit board a)
- Remove the heat sinks, heat conductive pad and heat pipe
- Unscrew the spacer bolts (x10) as well as the spacer bolts (x4) and nuts (x2) on the back on circuit board a) and remove circuit board a)
- Unscrew the screws (x3) on the rear wall and remove the rear wall
- Unscrew the screws (x7) on the front frame and remove the front

REMOVAL OF POLLUTANTS PV630, PV631, PV632

Disassembly


Fig. 5: Housing base

- Unscrew the nuts (x4) on the fan and remove the fan
- Remove the heat sinks, heat conductive pad and heat pipe
- Unscrew the screws (x4) on the large heat sink and remove the large heat sink
- Remove contoured foam
- Unscrew the screws (x6) on the housing spacer plate and remove the housing spacer plate
- Unscrew the screws (x6) on the sleeves and remove the sleeves with slider
- Remove the housing feet from the housing base

REMOVAL OF POLLUTANTS PV630, PV631, PV632


Fig. 6: Front

- Unscrew the nuts (x2) on circuit board f) and remove circuit board f)
- Unscrew the spacer bolts (x2) on the carrier plate and remove the switch membrane incl. carrier plate
- Remove the front membrane with the product name from the front frame
- Unscrew the nuts (x6) on the front plate and remove the front plate with front membrane
- Unscrew the nuts (x5) on circuit board d) and remove circuit board d)
- Unscrew the screws (x8) on the connector socket and remove the connector socket
- Unscrew the nuts (x4) on the socket frame and remove the socket frame

REMOVAL OF POLLUTANTS PV630, PV631, PV632

Disassembly

Fig. 7: Camera head with cable and plug (Item 160)

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Waste electrical and electronic equipment (WEEE) 2012/19/EC (vom 04.07.2012)

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The main product mark 'Aesculap' is a registered mark of Aesculap AG. Subject to technical changes. All rights reserved.
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Die Hauptproduktmarke 'Aesculap' ist eine eingetragene Marke der Aesculap AG. Technische Änderungen vorbehalten.