Aesculap Surgical Instruments

-
- | | |
|-----|--|
| GB | Instructions for use/Technical description |
| USA | CO ₂ cartridge adapter |
| D | Gebrauchsanweisung/Technische Beschreibung |
| | CO ₂ -Kartuschenadapter |
| F | Mode d'emploi/Description technique |
| | Adaptateur de cartouche de CO ₂ |
| E | Instrucciones de manejo/Descripción técnica |
| | Adaptador para cartuchos de CO ₂ |
| I | Istruzioni per l'uso/Descrizione tecnica |
| | Adattatore per cartucce di CO ₂ |





Fig. 2

Fig. 3

Fig. 4

CO₂ cartridge adapter

Legend

- 1 CO₂ cartridge adapter
- 2 Cartridge sleeve
- 3 Connector piece
- 4 Cartridge (unsterile)
- 5 Plastic seal
- 6 Sliding sleeve
- 7 Pneumatic holding arm

Symbols on product and packages

Symbol	Explanation
	Not for reuse in intended applications as defined by the manufacturer
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The CO₂ cartridge adapter is used for the provision of compressed gas for operating the UNITRAC holding system in unsterile areas.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

CAUTION

Loss of sterility of close-by instruments caused by use of the cartridge in sterile area or by unsterile exhaust air!

- Use the cartridge and the CO₂ cartridge adapter only in unsterile areas of the O.R.

Note

Optionally, for joint storage and processing the CO₂ cartridge adapter can be sterilized with the UNITRAC holding system. The CO₂ cartridge adapter is intended for use in unsterile areas only.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

DANGER

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not process the cartridge.

- Do not use the product after its use-by date.

Assembling

CAUTION

Impact damage to the plastic seal!

- **Remove the protective foam cover only immediately before use.**

- Remove protective foam cover from cartridge 4.
- Insert cartridge 4 into cartridge sleeve 2, see Fig. 2.
- Align connector piece 3 on the top end of cartridge sleeve 2 and screw it down to the positive stop, see Fig. 3.

Cartridge 4 is broached.

Note

The cartridge sleeve can be difficult to release if tightened too forcefully.

Installing the cartridge adapter at the pneumatic holding arm

- Pull back sliding sleeve 6 and couple CO₂ cartridge adapter 1 at pneumatic holding arm 7, see Fig. 4.
- Release sliding sleeve 6.
Sliding sleeve 6 engages audibly, with a click.

Note

Heavy or sluggish movement of the holding arm indicates that the cartridge is spent.

Note

Use a new cartridge for each operation.

Replacing the cartridge

- Uninstall the pressure unit and remove cartridge 4, see Disassembling.
- Insert new cartridge 4 and couple CO₂ cartridge adapter 1, see Assembling.

Disassembling

WARNING

Risk of injury caused by very cold cartridge, due to expansion of residual gas!

- **Touch the cartridge only at the plastic seal.**

CAUTION

Risk of alarm caused by residual gas hissing out of the partially emptied cartridge at the moment of removal!

- **Inform O.R. personnel well in time of the impending cartridge change.**
- Pull back sliding sleeve 6 and decouple CO₂ cartridge adapter 1 from pneumatic holding arm 7.
- Unscrew cartridge sleeve 2 from connector piece 3.
As soon as residual gas escapes audibly, pause unscrewing until all residual gas has escaped.
- Remove cartridge 4 from cartridge sleeve 2.

CO₂ cartridge adapter

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Therefore the time interval between application and cleaning should not exceed 6 h, and neither potentially fixating pre-cleaning temperatures >45 °C nor any fixating disinfectants (active ingredient: aldehyde) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfecting and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use process chemicals that would cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.

For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.

Manual cleaning/disinfecting

- Disassemble the product before cleaning.
- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Preparation prior to cleaning

- Ensure that the cartridge has been removed from the CO₂ cartridge adapter, see Disassembling.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions. The cleaning and disinfecting agents must
 - be approved for plastic materials and high-grade steel,
 - not attack softeners.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55 °C.

CAUTION

Contamination or blocking of the CO₂ cartridge adapter with cleaning and disinfecting agents!

- Do not immerse/clean the CO₂ cartridge adapter in an ultrasonic cleaning bath and do not apply immersion disinfection.

CO₂ cartridge adapter

Manual cleaning and wipe disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	-	-	D-W	-
II	Drying	RT	-	-	-	-
III	Wipe disinfection	-	≥1	-	-	Meliseptol HBV wipes 50 % propan-1-ol
IV	Final rinse	RT (cold)	0.5	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Clean the product under running tap water, using a suitable cleaning brush if necessary, until all visible residues have been removed from the surfaces.
- Make certain the product is held in such a way that water will not enter e.g. through coupling interfaces. (Immediately remove any fluid that entered the product inadvertently.)
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Dry the product with lint-free tissue or medical-quality filtered compressed air.
- Dry the product with a lint-free tissue.

Stage III

- Wipe all surfaces of the product with a single-use disinfecting wipe.

Stage IV

- After the specified exposure time (at least 1 minute), rinse the disinfected surfaces under running FD water.
- Hold the product in such a way that water cannot enter e.g. through couplings or adapters.
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting

Note

The disinfectant must be of proven effectiveness (e.g. DGHM or FDA approval or CE mark).

Note

For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that A_0 is $>3\,000$ for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfectant without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	–
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER with tensides, application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD-W	–
IV	Thermal disinfecting	90/194	5	FD-W	–
V	Drying	–	–	–	According to disinfectant program

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

CO₂ cartridge adapter

Inspection, maintenance and checks

CAUTION

Damage to sealing components caused by oil spray on the CO₂ cartridge adapter!

- Do not lubricate the CO₂ cartridge adapter.

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Ensure that the cartridge has been removed from the CO₂ cartridge adapter, see Disassembling.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

Packaging

- Sort the product into its appropriate storage device or put it on a suitable tray.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

Note

The product may only be sterilized in disassembled condition.

CAUTION

Loss of sterility/damage to the CO₂ cartridge adapter, close-by instruments and the sterilizer caused by sterilization of the CO₂ cartridge adapter with cartridge!

- **Make certain that no cartridge is ever put into the sterilizer.**
- Make certain that all external and internal surfaces of the product will be exposed to the sterilizing agent.
- Validated sterilization process:
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom)			
Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F– 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

Technical Service

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Designation	Art. no.
UNITRAC holding arm	RT040R
UNITRAC CO ₂ cartridge (unsterile)	RT044SU
Cartridge sleeve for CO ₂ cartridge adapter	RT043201

Disposal

- Adhere to national regulations when disposing of or recycling the product or its components and their packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG