


# CONTAINERS: GUARDIANS OF STERILITY

## SCIENTIFIC INFORMATION


Rigid sterile containers deliver surgical instruments in a sterile, functional condition to the right place at the right time. To achieve this, they have to be strong, stable, maintain sterility during storage and have a long shelf life – all whilst withstanding thousands of cycles of sterilization and usage. **The following tests demonstrate that the AESCULAP Aicon® Sterile Containers fulfill all these requirements.**



### TESTS

### RESULTS

#### STABILITY / ENDURANCE TESTING




Container components (face plates, identification tags, filter holders, protection covers, enhanced drying system EDS, closing latches) were cleaned, sterilized **(zero, 250, 500 cycles)**

and underwent thousands of cycles of **assembly + disassembly (2,000-10,000)**.



All components fulfilled success criteria, **no** tearing, breaking or loss of function.

#### MECHANICAL LOAD TESTING



##### Stacking (without slipping or toppling over):

A load of up to **5,700 N or roughly 581.24 kg** was applied for **10 min**; a tensile force of **50 N** was applied in **longitudinal, transversal and diagonal direction for 5 s** (according to DIN<sup>1</sup>).

##### Stability after pressure change with and without stacking:

Container was **loaded with 64 kg** and then underwent sterilization.

##### Stability of container grips:

A tensile machine using **4x** the maximum container load of **11.34 kg** was lifted for **10 min**.



**No** slipping or sliding

**No** optical transformation, measurements before and after loading differed **< 1mm**

**No** deformations

#### STERILIZATION PERFORMANCE



##### Sterilization performance

Containers were loaded and placed in sterilizer<sup>2</sup>.

Three different times were noted:

##### Equilibration time

Time until desired temperature is established within chamber.


ISO norm: < 15 s with a chamber volume of < 800 l

##### Holding time

ISO norm: sterilization temperature of 134–137°C for ≥ 180 s

##### Plateau time (holding time + equilibration time)

ISO norm: 1–60 s: delta T ≤ 5 K; > 60 s: delta T ≤ 2 K



##### Sterilization validation after accelerated ageing

Sterilized, "aged" containers (**1,000 sterilization cycles, 500 cleaning cycles**) were transferred to a test chamber, challenged with a test germ and incubated for **30 min**<sup>3</sup>.



##### Sterilization efficacy and thermal profile

Two containers were stacked in a steam cycle of **132°C for 3 min**.

Instruments in the container were inoculated with spores in the most difficult to sterilize locations<sup>4</sup>.



**Performance:** All normative criteria were met for all test samples<sup>5</sup>.

**Validation:** No entrance of microorganisms was observed.

##### Sterilization efficacy and thermal profile:

All indicator test samples were negative for growth after incubating. Steady state thermal conditions were confirmed. Positive controls were performed.

#### SHELF LIFE



The containers were fully loaded, sterilized and stored for **365 days**, followed by sterility testing in a microbiological culture medium<sup>6</sup>.



**No** growth of microorganisms after sterilization, storage and incubation. Containers can maintain sterility after **365 days**.

#### STERILE UPPER EDGE



**Upper edges of the container bottom** were sterilized and then inoculated with **bacterial spores** along the outer rim; test swabs were incubated for **7 days at 55–60°C**.



Aseptic presentation of internal contents was considered to be at a **safe distance from potential contamination** of the interior lip of the container and was determined to be at the **+2 mm** mark from the top rim of the container.

#### STERILE BARRIER INTEGRITY AFTER TRANSPORT AND VERTICAL MOVEMENT IN AN ELEVATOR



Containers were sterilized, **transported** via a route with physical obstacles and changes in environmental conditions, taken up and down an **elevator (minimum 10 floors)** and then tested for sterility again<sup>6</sup>.



Containers **maintained** whole package integrity and sterility after transportation and elevator transit.

#### Sources

- (1) Packaging for terminally sterilized medical devices, Part 8: re-usable sterilization containers for steam sterilizers conforming to EN 285 – requirements and test methods; German version EN 868-8:200.
- (2) DIN EN 868-8:2009-09, Packaging for terminally sterilized medical devices – Part 8: re-usable sterilization containers for steam sterilizers conforming to EN 285-Requirements and test methods.
- (3) DIN EN ISO 11607-1:2017-10, Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems.
- (4) Association for the Advancement of Medical Instrumentation (AAMI), International Organization for Standardization (ISO) Guidelines 4161:2009(R2014) Sterilization of Health Care Products.
- (5) DIN EN 285-2016-05, Sterilization – Steam Sterilization – Large Sterilizers.
- (6) Premarket Notification [501(k)] Submissions for Medical Sterilization Packaging Systems in Health Care Facilities; Draft Guidance for Industry and FDA; <https://federalregister.gov/documents/2002/03/07/02-5489/medical-devices-draft-guidance-for-industry-and-fda-on-premarket-notification-submissions-for->