CRANIOFIX2

- GB

USA

D

F

E

I

P

NL

DK

S

GR

Instructions for use

CRANIOFIX2 instrument set

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CRANIOFIX2-Instrumenten-Set

Mode d'emploi

Kit d'instruments CRANIOFIX2

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Jogo de instrumentos CRANIOFIX2

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CRANIOFIX2-instrumentenset

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
- 1

CRANIOFIX2 applier FF494R
- 2

Lubrication points of the CRANIOFIX2 applier
- 3

CRANIOFIX2-cutting forceps FF103R
- 4

CRANIOFIX2 extraction forceps FF104R
- 5

CRANIOFIX2 holding forceps FF105R
- 6

CRANIOFIX2 titanium clamp FF490T, Ø 11 mm
- 7

CRANIOFIX2 titanium clamp FF491T, Ø 16 mm
- 8

CRANIOFIX2 titanium clamp FF492T, Ø 20 mm

Symbols on product

	Attention, see instructions for use
--	-------------------------------------

Intended use

The Aesculap CRANIOFIX2 instrument set is used for implanting CRANIOFIX2 titanium clamps.

System components

The CRANIOFIX2 instrument set includes the following components:

- CRANIOFIX2 applier
- CRANIOFIX2 cutting forceps
- CRANIOFIX2 extraction forceps
- CRANIOFIX2 holding forceps
- CRANIOFIX storage tray

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

The CRANIOFIX2 instrument set is supplied in unsterile condition and must be cleaned and sterilized prior to use.

- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to the initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside damaged products immediately.

Safe operation

CRANIOFIX2 applier

Damage to/destruction of the titanium clamps caused by using a non-compatible applier!

- Use CRANIOFIX2 applier FF494R only for CRANIOFIX2 titanium clamps FF490T–FF492T.
- Use CRANIOFIX2 appliers FF102R, FF106R and FF107R only for CRANIOFIX2 titanium clamps FF099T–FF101T.

The CRANIOFIX2 applier 1 is not to be disassembled. The CRANIOFIX2 applier is compatible with CRANIOFIX storage tray FF094P. CRANIOFIX2 applier 1 is fitted with an automatic release mechanism to ensure that the maximum force sustainable by the pin is not exceeded when a CRANIOFIX2 titanium clamp is applied.

CRANIOFIX2 holding forceps

- Hold the CRANIOFIX2 titanium clamp by applying CRANIOFIX2 holding forceps 5 in the holding groove. This serves to avoid damage to the pin indentation.

CRANIOFIX2 cutting forceps

The cutting edges of CRANIOFIX2 cutting forceps 3 have been finished (rounded) in such a way that the forceps are self-centering in the depression of the CRANIOFIX2 titanium clamp. This design also ensures that the cutting plane does not protrude, or protrudes only very slightly, over the upper edge of the CRANIOFIX2 titanium clamp.

The CRANIOFIX2 cutting forceps must be replaced as soon as it exerts excessive forces on the CRANIOFIX2 titanium clamp during the cutting process and/or when the cutting burr on the pin becomes too large. Both are signs that the cutting edges of the CRANIOFIX2 cutting forceps are worn out.

- To prevent uncontrolled snapping off of the protruding pin end into the operating field while cutting with CRANIOFIX2 cutting forceps 3, keep hold of the pin of the CRANIOFIX2 titanium clamp.

CRANIOFIX2 extraction forceps

- Apply CRANIOFIX2 extraction forceps 4 with their teeth in the 4 notches milled into the CRANIOFIX2 titanium clamp for this purpose.

Care and handling

Note

Observe all relevant national regulations and standards concerning process ing.

Note

For patients with Creutzfeldt-Jakob Disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net.

Preparations

- Pre-clean the product, if necessary.
- Put away the product in dry condition.
- Open jointed products.
- Reprocess the product immediately after use.

Cleaning/Disinfecting

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Only use a cleaning/disinfecting agent that has been approved for surgical steel; follow the agent manufacturer's instructions.

- Use cleaning/disinfecting agents that are suitable for the product. Always follow the manufacturer's instructions with regard to concentration, temperature and exposure time.
- Avoid encrustation of residues/proteins (e.g. caused by aldehyde/ alcohol).
- Only use bactericidal, fungicidal and virucidal disinfecting agents.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/ disinfection.
 - to prepare products with encrusted debris for mechanical cleaning/ disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - as an aftertreatment for products that are still dirty after mechanical cleaning/disinfection.
- Preferably use thermal disinfecting processes.

Mechanical cleaning/disinfecting

Pre-cleaning

- Thoroughly pre-rinse under running water.
- Carry out ultrasound treatment.
- Insert the special cleaning brush (art. no. FF494801) through the distal opening of the applier and clean the inner surfaces of the tensioning forceps.
- Carry out the final rinse under running water.

Cleaning

- Place the product on a tray that is suitable for cleaning (to avoid rinsing blind spots).
- Store any products with hinges or joints on the tray in such a way that the joints are kept open.
- Ensure that water can flow out of openings.
- Process the product in a cleaning/disinfecting device. Follow the instructions provided by the manufacturer of the machine.
- Processing cycle:
 - Use an appropriate alkaline cleaning/disinfecting agent according to manufacturer's instructions.
 - Do not exceed the maximum cleaning temperature of 55 °C.
 - Clean for at least 10 min.
 - Neutralize.
 - Apply intermediate rinse for 1 min.
 - Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
 - Carry out thermal disinfection: Rinse for 10 min at 93 °C with distilled, demineralized or fully desalinated water.
 - Complete the program with a drying phase of at least 20 min at a temperature of not more than 110 °C.
- After completion of the mechanical cleaning/disinfecting cycle, inspect surfaces, cavities, lumens, and openings for visible debris.
- Clean manually if necessary.

Manual cleaning/disinfecting

- Use an appropriate alkaline cleaning/disinfecting agent according to manufacturer's instructions.
- Immerse the product in the cleaning/disinfecting agent in such a way that all surfaces, cavities, lumens, and openings are covered.
- After the end of the disinfection period, thoroughly rinse the product under running water, making sure that water flows through every lumen and channel and all blind holes are repeatedly filled and drained.
- Clean hinged or jointed products in open and closed positions.
- Remove encrusted debris with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- Clean lumens, channels and blind holes with soft round plastic brushes of fitting diameter.
- Carry out an intensive final rinse with distilled, demineralized or fully desalinated water.
- Inspect surfaces, cavities, lumens and openings for visible debris. If necessary, repeat the cleaning/disinfection process.
- Use a lint-free cloth or a compressed-air gun for drying the product.
- Make certain that lumens, channels and blind holes are dried, too.

Control, care and inspection

- Allow the product to cool down to room temperature.
- Lightly lubricate moving parts such as hinges and latches with a steriliz able, steam-permeable and tissue-compatible maintenance oil (e.g. Aesculap Sterilit spray JG600 or maintenance oil JG598).
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Set aside damaged products immediately.

Packaging

- Secure product with retractor in the first notch.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected. Observe the weight limit for each tray/container.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Apply steam sterilization, observing the following rules: Sterilization must be carried out through a validated steam sterilization process (e.g. in a sterilizer conforming to EN 285/ANSI/AAMI/ISO 11134–1993, ANSI/AAMI ST46–1993, validated acc. to EN 554/ISO 13683). For the fractionated vacuum process, the sterilization has to be carried out for a minimum of 5 minutes at 134 °C and at 2 bar pressure.
- When sterilizing several products at the same time in one steam steril- ization container (sterile container) for processing under generally accepted hospital in-use conditions.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid steril- ization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile con- tainer system
Pre-vacuum	270— 275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

Storage

- Store processed products under conditions as germ–free as possible, in a dry, dark, cool and dust-protected room.

Technical Service

For service, maintenance or repairs, contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
AESCULAP Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420
Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Designation	Art. no.
CRANIOFIX2 applier	FF494R
CRANIOFIX2 holding forceps	FF105R
CRANIOFIX2 cutting forceps	FF103R
CRANIOFIX2 extraction forceps	FF104R
CRANIOFIX2 storage tray	FF094P
Cleaning brush for CRANIOFIX2 applier	FF494801
CRANIOFIX2 titanium clamp, Ø 11 mm	FF490T *
CRANIOFIX2 titanium clamp, Ø 16 mm	FF491T *
CRANIOFIX2 titanium clamp, Ø 20 mm	FF492T *

* Regarding the operation and function of this product, always observe instructions for use TA011200.

Distributor in the US/Contact in Canada for product information and complaints

AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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