CRANIOFIX2

- GB

Instructions for use
CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm
- USA

Instructions for use
CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm
- D

Gebruiksaanweisung
CRANIOFIX2-Titanklammern Ø 11 mm, Ø 16 mm, Ø 20 mm
- F

Mode d'emploi
Pincen en titane CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- E

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
Symbols on product

	Sterilization using irradiation
	This product is not for reuse in terms of its proper applications as defined by the manufacturer
	Use by
	Attention, see instructions for use

Product description

The Aesculap **CRANIOFIX2** titanium clamp system comprises sterile **CRANIOFIX2** titanium clamps for single use and diverse reusable neurosurgical instruments for applying titanium clamps.

Intended use

The Aesculap **CRANIOFIX2** titanium clamp system is used for fixating craniotomized bone flaps and fractures at the brain skull (neurocranium). The **CRANIOFIX2** titanium clamp FF490T/FF490T-UNI, Ø 11 mm, is applied in the craniotomy gap (see Fig. 2) or in the fracture gap. The **CRANIOFIX2** titanium clamp FF491T/FF491T-UNI, Ø 16 mm, is applied in a trepanation hole of Ø 12 mm (see Fig. 3) or for bridging the defect in case of a fracture. The **CRANIOFIX2** titanium clamp FF492T/FF492T-UNI, Ø 20 mm, is applied in a trepanation hole of Ø 15 mm (see Fig. 3) or for bridging the defect in case of a fracture.

Materials

The **CRANIOFIX2** titanium clamps are made of titanium alloy (Ti6Al4V) acc. to ISO 5832-3. This titanium alloy is a biocompatible material. The **CRANIOFIX2** titanium clamps are non-ferromagnetic and MR-compatible up to 3.0 Tesla. Consequently, the magnetic fields applied during NMR tomography do not present any risk. Stronger magnetic fields or an enlargement of the imaging field can lead to a dramatically increased incidence of positioning errors and artifacts. The size of the artifacts can vary considerably, depending on the NMR-pulse frequency, and compromise the significance of the NMR image if the area concerned is in the immediate vicinity (within a distance of some millimeters) of the clamp position.

Contraindications

- The following list shows some, but not all conceivable contraindications.
- Patients with systemic illnesses and metabolic disorders, oversensitivity to metals or allergies against the implant materials, inflammations in the region of the implant site, and patients that show insufficient compliance and co-operation
 - Bone conditions that rule out the application of **CRANIOFIX2** titanium clamps
 - Artificial cranial bone flaps
 - Bone tumors in the area on which the implant is supported
 - Degenerative bone diseases
 - Missing dura mater
 - Application in the facial skull (viscerocranium) and in the orbita or skull-base regions

Notes

- Improper fixation or combination with other fixation methods (sutures, mini plates or clamp systems from other manufacturers) can lead to step formation, position shifts of the skull calotte and loosening or breakage of implant components. For maximum stability of a craniotomized bone flap, we recommend using three **CRANIOFIX2** titanium clamps in a triangular arrangement. For larger bone flaps or multiple fragments of a skull fracture, additional **CRANIOFIX2** titanium clamps may be required to ensure adequate stability.

- During application of **CRANIOFIX2** titanium clamps, friction between the upper titanium disc and the pin indentation can lead to spark formation. Appropriate safety measures must be taken in order to prevent this. Because of the fire hazard, surgical interventions of this kind must not be carried out in rooms where flammable anesthetics or other flammable gases, liquids, objects or oxidizing agents are stored.
- After the operation, exceptionally high mechanical loads on the fixated bone flap must be avoided.

Safety measures

The operating surgeon should be thoroughly familiar with the surgical techniques required for using this product. The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. The operating surgeon is also responsible for complications resulting from misdiagnosis and flawed operation techniques.

Possible side effects

Infections, damage to nervous tissue, hematomas and wound dehiscence are among the general perioperative risks. The patient should be made aware of these risks as well as other risks of neurosurgery, general surgery and the administration of anesthetics.

Delivery format

Aesculap **CRANIOFIX2** titanium clamps are delivered in sterile condition and are designed for single use only. The **CRANIOFIX2** titanium clamps must not be resterilized.

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product.

Provision and handling of the implants

Damage to/destruction of the titanium clamps caused by using a non-compatible applier!

- Use the **CRANIOFIX2** applier (FF494R) only for **CRANIOFIX2**-titanium clamps (FF490T/FF490T-UNI-FF492T/FF492T-UNI).
- Use the **CRANIOFIX** applier (FF102R, FF106R and FF107R) only for the **CRANIOFIX** titanium clamps (FF099T-FF101T).

Prior to the surgical intervention, the operating surgeon must prepare an operation plan in order to ensure that a sufficient number of sterile **CRANIOFIX2** titanium clamps are available, aseptic operation conditions are in place, the **CRANIOFIX2** applier required and other **CRANIOFIX** instruments are complete and ready for operation, and the operating surgeon and the O.R. team are familiar with the instruments and operation technique. The craniotomy incisions must be performed in an osteoplastic fashion. Correct handling of the **CRANIOFIX2** titanium clamps before and during the operation is essential for successful bone flap fixation. The **CRANIOFIX2** titanium clamps are individually wrapped in protective packaging and gamma-ray sterilized.

- The **CRANIOFIX2** titanium clamps have to be kept in their original packaging until immediately before use. Prior to use, check the expiry date and make sure that the sterile packaging is intact.
- Make certain that the surfaces of the **CRANIOFIX2** titanium clamps are not bent or damaged.
- The **CRANIOFIX2** titanium clamps have been designed as a disposable product and must not be reused.
- Damaged or surgically excised **CRANIOFIX2** titanium clamps must not be used again.
- To avoid spark formation, immerse the **CRANIOFIX2** titanium clamps in a sterile liquid prior to application.

Application of CRANIOFIX2 titanium clamps

Caution

The Aesculap **CRANIOFIX2** titanium clamp of Ø 16 mm may only be applied in a trepanation hole of Ø 12 mm (see Fig. 3) or for bridging the defect in case of a fracture. The Aesculap **CRANIOFIX2** titanium clamp of Ø 20 mm may only be applied in a trepanation hole of Ø 15 mm (see Fig. 3) or for bridging the defect in case of a fracture.

The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. For maximum stability of one bone flap approx. 20 cm² in size, we recommend applying three **CRANIOFIX2** titanium clamps in the craniotomy gap or the trepanation hole, respectively, at equal distances to each other. Figure 2 shows a typical application of a craniotomized bone flap.

The following is a brief description of the surgical procedure:

- 1 Position the three **CRANIOFIX2** titanium clamps in equal distances to each other at the bone flap.
- 2 The lower discs are pushed in between the dura and the calvarium.
- 3 Replace the bone flap in its original position.
- 4 Apply the holding forceps on the first **CRANIOFIX2** titanium clamp. To ensure a safe fixation and avoid damage to the pin indentation, the pin of the **CRANIOFIX2** titanium clamp must be held in the recess intended for this purpose.
- 5 Position the **CRANIOFIX2** applier over the pin of the **CRANIOFIX2** titanium clamps.
- 6 Compress the handle of the **CRANIOFIX2** applier so that the upper disc approaches the lower disc along the pin indentation. The spring reeds engage in the pin indentation of the **CRANIOFIX2** titanium clamp, with an audible "click".
- 7 Do not remove the holding forceps until immediately before closing the **CRANIOFIX2** titanium clamp.
- 8 The **CRANIOFIX2** applier is equipped with an automatic strain relief mechanism, which ensures that the **CRANIOFIX2** applier is disengaged as soon as the maximum allowable force is applied. Once the **CRANIOFIX2** applier has disengaged, the handle should be fully compressed several times until there is no more audible clicking of the spring reeds as they engage in the indentation of the pin of the **CRANIOFIX2** titanium clamp.

The top disk can loosen after application if the pin is not cut off properly!

- After tightening the top disk, cut the excess pin off at a right angle to the pin's axis using the **CRANIOFIX** cutting forceps.
- Be sure that the indentation above the top disk remains completely intact.

- 9 Use the cutting forceps to remove the remaining, outward projecting section of the pin. When doing so, the protruding pin of the **CRANIOFIX2** titanium clamp must be held securely to prevent the pin from snapping off into the operation field in an uncontrolled manner.
- 10 The same procedure is carried out with the other **CRANIOFIX2** titanium clamps.

The surgical procedure for applying **CRANIOFIX2** titanium clamps is described in the **CRANIOFIX2** brochure and in the literature list (www.aesculap.de). Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact us directly:

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-0
Fax: +49 7461 95-2600
Website: www.aesculap.de

Patient ID

The Patient ID serves to keep a written record of important information on the implanted **CRANIOFIX2** titanium clamps and the surgical intervention that was carried out:

- Patient data
 - Hospital data
 - Date of provision
 - Physician in charge and hospital
 - Art. no. and lot no. of the **CRANIOFIX2** titanium clamps implanted
- To facilitate the postoperative radiological examination, each patient should be issued with a Patient ID, which can be ordered separately from Aesculap. Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact Aesculap directly.

Removing the CRANIOFIX2 titanium clamps in a revision operation

In the rare case that the operation has to be repeated, the operating surgeon can leave the **CRANIOFIX2** titanium clamps in place on the bone flap. Aesculap recommends applying a second incision around the original clamp position. If this is not feasible (e.g. for anatomic reasons), Aesculap recommends a second incision within the original clamp site. Under no circumstances must the incision cut through a **CRANIOFIX2** titanium clamp. The operating surgeon is responsible for any decision about how to approach the revision intervention (see Fig. 4).

Accessories

Designation	Art. no.
CRANIOFIX2 titanium clamps Ø 11 mm	FF490T
CRANIOFIX2 titanium clamps Ø 11 mm	FF490T-UNI
CRANIOFIX2 titanium clamps Ø 16 mm	FF491T
CRANIOFIX2 titanium clamps Ø 16 mm	FF491T-UNI
CRANIOFIX2 titanium clamps Ø 20 mm	FF492T
CRANIOFIX2 titanium clamps Ø 20 mm	FF492T-UNI
CRANIOFIX2 applier, non-dismountable	FF494R*
CRANIOFIX holding forceps	FF105R*
CRANIOFIX cutting forceps	FF103R*
CRANIOFIX extractor forceps	FF104R*
CRANIOFIX storage tray	FF094P

* Concerning the handling, functioning and reprocessing of this product, the instructions for use TA011161 must be observed.

TA-Nr.: 011200 09/12 Änd.-Nr.: 45901

USA CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm

Figure captions

- Fig. 1**
- CRANIOFIX2** titanium clamp FF490T/FF490T-UNI, Ø 11 mm
 - CRANIOFIX2** titanium clamp FF491T/FF491T-UNI, Ø 16 mm
 - CRANIOFIX2** titanium clamp FF492T/FF492T-UNI, Ø 20 mm
- Fig. 2**
- CRANIOFIX2** titanium clamp FF491T/FF491T-UNI, y 16 mm/**CRANIOFIX2** titanium clamp FF492T/FF492T-UNI, y 20 mm
 - CRANIOFIX2** titanium clamp FF490T/FF490T-UNI, Ø 11 mm
- Fig. 3**
- Correct position of **CRANIOFIX2** titanium clamps FF491T/FF491T-UNI/FF492T/FF492T-UNI: Trepanation hole y 12/15 mm
 - Incorrect position of **CRANIOFIX2** titanium clamps FF491T/FF491T-UNI/FF492T/FF492T-UNI: Craniotomy gap
- Fig. 4**
- CRANIOFIX2** titanium clamp
 - New craniotomy line
 - Old craniotomy line

Symbols on product

	Sterilization using irradiation
	This product is not for reuse in terms of its proper applications as defined by the manufacturer
	Use by
	Attention, see instructions for use

Product description

The Aesculap **CRANIOFIX2** titanium clamp system comprises sterile **CRANIOFIX2** titanium clamps for single use and diverse reusable neurosurgical instruments for applying titanium clamps.

Intended use

The Aesculap **CRANIOFIX2** titanium clamp system is used for fixing craniotomized bone flaps and fractures at the brain skull (neurocranium). The **CRANIOFIX2** titanium clamp FF490T/FF490T-UNI, Ø 11 mm, is applied in the craniotomy gap (see **Fig. 2**) or in the fracture gap. The **CRANIOFIX2** titanium clamp FF491T/FF491T-UNI, Ø 16 mm, is applied in a trepanation hole of Ø 12 mm (see **Fig. 3**) or for bridging the defect in case of a fracture. The **CRANIOFIX2** titanium clamp FF492T/FF492T-UNI, Ø 20 mm, is applied in a trepanation hole of Ø 15 mm (see **Fig. 3**) or for bridging the defect in case of a fracture.

Materials

The **CRANIOFIX2** titanium clamps are made of titanium alloy (Ti6Al4V) acc. to ISO 5832-3. This titanium alloy is a biocompatible material. The **CRANIOFIX2** titanium clamps are non-ferromagnetic and MR-Conditional up to 3.0 Tesla. Consequently, the magnetic fields applied during NMR tomography do not present any risk. Stronger magnetic fields or an enlargement of the imaging field can lead to a dramatically increased incidence of positioning errors and artifacts. The size of the artifacts can vary considerably, depending on the NMR-pulse frequency, and compromise the significance of the NMR image if the area concerned is in the immediate vicinity (within a distance of some millimeters) of the clamp position.

Contraindications

- The following list shows some, but not all conceivable contraindications.
- Patients with systemic illnesses and metabolic disorders, oversensitivity to metals or allergies against the implant materials, inflammations in the region of the implant site, and patients that show insufficient compliance and co-operation
 - Bone conditions that rule out the application of **CRANIOFIX2** titanium clamps
 - Artificial cranial bone flaps
 - Bone tumors in the area on which the implant is supported
 - Degenerative bone diseases
 - Missing dura mater
 - Application in the facial skull (viscerocranium) and in the orbita or skull-base regions

Notes

- Improper fixation or combination with other fixation methods (sutures, mini plates or clamp systems from other manufacturers) can lead to step formation, position shifts of the skull calotte and loosening or breakage of implant components. For maximum stability of a craniotomized bone flap, we recommend using three **CRANIOFIX2** titanium clamps in a triangular arrangement. For larger bone flaps or multiple fragments of a skull fracture, additional **CRANIOFIX2** titanium clamps may be required to ensure adequate stability.

- During application of **CRANIOFIX2** titanium clamps, friction between the upper titanium disc and the pin indentation can lead to spark formation. Appropriate safety measures must be taken in order to prevent this. Because of the fire hazard, surgical interventions of this kind must not be carried out in rooms where flammable anesthetics or other flammable gases, liquids, objects or oxidizing agents are stored.
- After the operation, exceptionally high mechanical loads on the fixated bone flap must be avoided.

Safety measures

The operating surgeon should be thoroughly familiar with the surgical techniques required for using this product. The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. The operating surgeon is also responsible for complications resulting from misdiagnosis and flawed operation techniques.

Possible side effects

Infections, damage to nervous tissue, hematomas and wound dehiscence are among the general perioperative risks. The patient should be made aware of these risks as well as other risks of neurosurgery, general surgery and the administration of anesthetics.

MRI Information

MR-Conditional.

Non-clinical testing demonstrated that the **CRANIOFIX2** System is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 3,000 Gauss/cm or less
- First Level Controlled Operating Mode, with a maximum whole body averaged SAR of 4-W/kg

MRI-Related Heating

In non-clinical testing, the **CRANIOFIX2** System produced the following temperature rises during MRI performed for 15 min of scanning (i.e., per pulse sequence) in 1.5 Tesla/64 MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3 Tesla (3 Tesla/128 MHz, Excite, HDx, Software 14X, M5, General Electric Healthcare, Milwaukee, WI) MR Systems:

	1.5 Tesla	3 Tesla
First Level Controlled Operating Mode, scaled to whole body averaged SAR of 4 W/kg, highest temperature change	+2.6 °C	+2.8 °C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the **CRANIOFIX2** System. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 5 mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	503 mm²	196 mm²	701 mm²	563 mm²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Delivery format

Aesculap **CRANIOFIX2** titanium clamps are delivered in sterile condition and are designed for single use only. The **CRANIOFIX2** titanium clamps must not be resterilized.

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product.

Provision and handling of the implants

CAUTION

Federal law restricts this device to sale by or on order of a physician

Damage to/destruction of the titanium clamps caused by using a non-compatible applier!

- Use the **CRANIOFIX2** applier (FF494R) only for **CRANIOFIX2**-titanium clamps (FF490T/FF490T-UNI-FF492T/FF492T-UNI).
- Use the **CRANIOFIX** applier (FF102R, FF106R and FF107R) only for the **CRANIOFIX** titanium clamps (FF099T-FF101T).

Prior to the surgical intervention, the operating surgeon must prepare an operation plan in order to ensure that a sufficient number of sterile **CRANIOFIX2** titanium clamps are available, aseptic operation conditions are in place, the **CRANIOFIX2** applier required and other **CRANIOFIX** instruments are complete and ready for operation, and the operating surgeon and the O.R. team are familiar with the instruments and operation technique. The craniotomy incisions must be performed in an osteoplastic fashion. Correct handling of the **CRANIOFIX2** titanium clamps before and during the operation is essential for successful bone flap fixation. The **CRANIOFIX2** titanium clamps are individually wrapped in protective packaging and gamma-ray sterilized.

- The **CRANIOFIX2** titanium clamps have to be kept in their original packaging until immediately before use. Prior to use, check the expiry date and make sure that the sterile packaging is intact.
- Make certain that the surfaces of the **CRANIOFIX2** titanium clamps are not bent or damaged.
- The **CRANIOFIX2** titanium clamps have been designed as a disposable product and must not be reused.
- Damaged or surgically excised **CRANIOFIX2** titanium clamps must not be used again.
- To avoid spark formation, immerse the **CRANIOFIX2** titanium clamps in a sterile liquid prior to application.

Application of CRANIOFIX2 titanium clamps

Caution

The Aesculap **CRANIOFIX2** titanium clamp of Ø 16 mm may only be applied in a trepanation hole of Ø 12 mm (see **Fig. 3**) or for bridging the defect in case of a fracture. The Aesculap **CRANIOFIX2** titanium clamp of Ø 20 mm may only be applied in a trepanation hole of Ø 15 mm (see **Fig. 3**) or for bridging the defect in case of a fracture.

The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. For maximum stability of one bone flap approx. 20 cm² in size, we recommend applying three **CRANIOFIX2** titanium clamps in the craniotomy gap or the trepanation hole, respectively, at equal distances to each other. Figure 2 shows a typical application of a craniotomized bone flap.

The following is a brief description of the surgical procedure:

- Position the three **CRANIOFIX2** titanium clamps in equal distances to each other at the bone flap.
- The lower discs are pushed in between the dura and the calvarium.
- Replace the bone flap in its original position.
- Apply the holding forceps on the first **CRANIOFIX2** titanium clamp. To ensure a safe fixation and avoid damage to the pin indentation, the pin of the **CRANIOFIX2** titanium clamp must be held in the recess intended for this purpose.
- Position the **CRANIOFIX2** applier over the pin of the **CRANIOFIX2** titanium clamps.
- Compress the handle of the **CRANIOFIX2** applier so that the upper disk approaches the lower disc along the pin indentation. The spring reeds engage in the pin indentation of the **CRANIOFIX2** titanium clamp, with an audible "click".
- Do not remove the holding forceps until immediately before closing the **CRANIOFIX2** titanium clamp.
- The **CRANIOFIX2** applier is equipped with an automatic strain relief mechanism, which ensures that the **CRANIOFIX2** applier is disengaged as soon as the maximum allowable force is applied. Once the **CRANIOFIX2** applier has disengaged, the handle should be fully compressed several times until there is no more audible clicking of the spring reeds as they engage in the indentation of the pin of the **CRANIOFIX2** titanium clamp.

The top disk can loosen after application if the pin is not cut off properly!

- After tightening the top disk, cut the excess pin off at a right angle to the pin's axis using the **CRANIOFIX** cutting forceps.
- Be sure that the indentation above the top disk remains completely intact.

- Use the cutting forceps to remove the remaining, outward projecting section of the pin. When doing so, the protruding pin of the **CRANIOFIX2** titanium clamp must be held securely to prevent the pin from snapping off into the operation field in an uncontrolled manner.
- The same procedure is carried out with the other **CRANIOFIX2** titanium clamps.

The surgical procedure for applying **CRANIOFIX2** titanium clamps is described in the **CRANIOFIX2** brochure and in the literature list (www.aesculap.de). Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact us directly:

Aesculap AG

Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-0
Fax: +49 7461 95-2600
Website: www.aesculap.de

Or in the US:
Aesculap Inc.
Attn: Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Patient ID

The Patient ID serves to keep a written record of important information on the implanted **CRANIOFIX2** titanium clamps and the surgical intervention that was carried out:

- Patient data
- Hospital data
- Date of provision
- Physician in charge and hospital
- Art. no. and lot no. of the **CRANIOFIX2** titanium clamps implanted

To facilitate the postoperative radiological examination, each patient should be issued with a Patient ID, which can be ordered separately from Aesculap. Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact Aesculap directly.

Removing the CRANIOFIX2 titanium clamps in a revision operation

In the rare case that the operation has to be repeated, the operating surgeon can leave the **CRANIOFIX2** titanium clamps in place on the bone flap. Aesculap recommends applying a second incision around the original clamp position. If this is not feasible (e.g. for anatomic reasons), Aesculap recommends a second incision within the original clamp site. Under no circumstances must the incision cut through a **CRANIOFIX2** titanium clamp. The operating surgeon is responsible for any decision about how to approach the revision intervention (see **Fig. 4**).

Accessories

Designation	Art. no.
CRANIOFIX2 titanium clamps Ø 11 mm	FF490T
CRANIOFIX2 titanium clamps Ø 11 mm	FF490T-UNI
CRANIOFIX2 titanium clamps Ø 16 mm	FF491T
CRANIOFIX2 titanium clamps Ø 16 mm	FF491T-UNI
CRANIOFIX2 titanium clamps Ø 20 mm	FF492T
CRANIOFIX2 titanium clamps Ø 20 mm	FF492T-UNI
CRANIOFIX2 applier, non-dismountable	FF494R*
CRANIOFIX holding forceps	FF105R*
CRANIOFIX cutting forceps	FF103R*
CRANIOFIX extractor forceps	FF104R*
CRANIOFIX storage tray	FF094P

* Concerning the handling, functioning and reprocessing of this product, the instructions for use TA011161 must be observed.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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